


Short Instructions

Eppendorf ThermoMixer F0.5/F1.5/F2.0/FP

English (EN)

1 Operating controls



1 Temperature keys with control LEDs

Select a temperature or switch off temperature control.

2 Key short

Short Mix runs as long as the **short** key is pressed.

3 Key start/stop

Start or stop mixing/temperature control.

4 Arrow keys speed

Set the mixing frequency.
Keep the arrow key pressed: quick setting.

5 Arrow keys temp

Set the temperature.
Keep the arrow key pressed: quick setting.
As soon as the set temperature is changed, the device begins to perform temperature control.

6 Display


7 Menu arrow key

Navigate in the menu: set key lock or volume.

8 Key menu/enter

Open the menu.
Confirm your selection.

2 Display ThermoMixer F0.5/F1.5/F2.0/FP



- | | |
|---|--|
| <p>1 Device status
🔄 Device is performing mixing/temperature control.</p> <p>2 ThermoTop
🔊 ThermoTop has been attached. In order to prevent the formation of any condensate in a reliable manner, the device first heats the ThermoTop until it reaches the set temperature. The tempering of the thermoblock occurs with a delay. After inserting samples into a pre-heated thermoblock, the displayed actual temperature may fall temporarily.</p> <p>3 Speaker
🔊 Speaker switched on.
🔇 Speaker switched off.</p> | <p>4 Key lock
🔒 Key lock activated: parameters cannot be changed.
🔓 No key lock.</p> <p>5 Mixing frequency</p> <p>6 Actual temperature
When the actual temperature flashes on the display, the device does not operate.</p> <p>7 Set the temperature
When the set temperature has been reached, only one value is displayed.</p> |
|---|--|