

Pipette Performance Data Sheet



Eppendorf Xplorer® plus Move It®

Are you frustrated with the challenges of pipetting samples across different vessel formats, which prevent you from using a standard multi-channel pipette? Fret not – Eppendorf has your solution with the Eppendorf Xplorer plus Move It electronic multi-channel pipettes featuring adjustable tip spacing. Like the mechanical Research® plus Move It models, these pipettes accelerate and simplify the synchronous pipetting of a series of samples. The electronic Xplorer plus Move It models not only provide precision and durability but also introduce enhanced ergonomics and versatility with additional pipetting modes, such as multi-dispensing.

Features

- > 4-, 6-, 8- and 12-channel air-cushion pipettes with variable volume and adjustable tip spacing option
- > Intuitive handling via multi-function rocker with »up is up« (aspiration) and »down is down« (dispensing) logic and selection dial to access all pipetting modes
- > Manual adjustment knob for smooth vibration-free format switches and liquid transfers without dripping
- > Format limiter for quick switches backwards and forwards between source and destination formats
- > Individual programming, fixed volume settings, password protection, display of service intervals
- > Autoclavable lower parts for increased user and sample safety
- > 360° rotatable pipette head for comfortable readability of display and fast identification of parameters
- > Also available as mechanical Eppendorf Research® plus Move It pipettes with adjustable tip spacing

Technical Specifications

| | |
|---------------------------------------|---|
| Intended use | Electronic multi-channel pipette for general laboratory use. Designed to quickly transfer several samples between different vessel formats. |
| Pipetting type | Air-cushion pipette |
| Channel options | 4, 6, 8, 12 with adjustable cone spacing |
| Volume range | 1 µL – 1,200 µL |
| Volume increments | 0.02 µL – 1 µL |
| Operation | > Electronic (pipetting & settings) > Manual (tip spacing adjustment) |
| Pipetting modes | 9 (see Applications) |
| Pipette Manager enabled | Yes (sold separately) |
| Aspiration & dispensing speed | 8 levels |
| Display languages | Chinese, Dutch, English, French, German, Italian, Japanese, Portuguese, Spanish |
| Autoclavable | Lower parts (121 °C, 1 bar, 20 min) |
| UV resistance | Yes (254 nm) |
| Consumables | epT.I.P.S.® pipette tips |
| Open system (3rd party consumables) | Yes |
| Conforms with | ISO 8655 |
| Temporary adjustment | > Pre-programmed adjustments for liquid types, altitude or elongated pipette tips > Individual adjustments in 1, 2 or 3 points |
| Power supply | Power supply/mains with plug adapters included (100-240 V ± 10 %, 50-60 Hz) |
| Battery type | Rechargeable lithium polymer (3.7 V) |
| Charging ports | > Connector socket for plug > Charging contacts for Carousel or Stand |
| Interface | Micro USB (for software updates) |
| Storage accessories (sold separately) | > Charger Carousel for 6 pipettes > Charger Stand > Pipette Holder for Carousel, Stand or wall mounting |
| Storage & transport conditions | > Temperature: -5 °C – 45 °C > Humidity: 10 % – 95 % |
| Operating conditions | > Temperature: 5 °C – 40 °C > Humidity: 10 % – 95 % |
| Others | > Easy and fast tip spacing between 4.5 mm (min) and 33 mm (max) depending on the model > No tube connections between cone and piston-cylinder system for increased reliability and safety |

Applications

- > Efficiency gains in all applications with frequent vessel format changes, i.a. ELISA, cell culture, PCR, genotyping, blood analysis, FACS analysis, media transfer, dilution series, screening
- > Automated (forward) pipetting of liquids (**Pip**)
- > Reverse pipetting (**Spc**) for non-aqueous solutions
- > Dispensing (**Dis**) to quickly fill plates, gels, and reaction vessels
- > Automatic dispensing (**Ads**) for automated filling of plates, gels, and reaction vessels
- > Pipetting and mixing (**P/M**) for automated mixing of samples
- > Manual pipetting (**Man**) for volume determination
- > Multi aspiration (**Spc**): Pooling supernatants from different vessels
- > Sequential dispensing and pipetting (**Spc**) for user-defined volume steps
- > Diluting (**Spc**) for preparation of dilution series
- > Fix volume (**Fix**): Save up to 10 volumes for ease of operation

Dimensions

| | |
|-----------------------|-----------------------------------|
| Length | 290 mm (8-channel, 100 µL) |
| Weight | 307 g (8-channel, 100 µL) |
| Packaging (W x H x D) | 32 cm x 16 cm x 12 cm |
| Shipping weight | 913 g – 1 kg (depending on model) |

Measurement Deviations

Eppendorf Xplorer® plus Move It®, adjustable multi-channel, incl. mains/power supply with plug adapters

| Volume range/channel | Color code | Volume | Systematic error ¹⁾ | | Random error ¹⁾ | |
|---------------------------------|---|---------|--------------------------------|----------|----------------------------|----------|
| 1 – 20 µL 8/12-channel | ■ light pink (for epT.I.P.S.® 384, 20 µL) | 2 µL | ±8 % | ±0.16 µL | ±5 % | ±0.1 µL |
| | | 10 µL | ±4 % | ±0.4 µL | ±2 % | ±0.2 µL |
| | | 20 µL | ±2 % | ±0.4 µL | ±1 % | ±0.2 µL |
| 5 – 100 µL 8/12-channel | ■ light yellow (for epT.I.P.S.® 384, 100 µL) | 10 µL | ±3 % | ±0.3 µL | ±2 % | ±0.2 µL |
| | | 50 µL | ±1.2 % | ±0.6 µL | ±0.8 % | ±0.4 µL |
| | | 100 µL | ±1 % | ±1 µL | ±0.6 % | ±0.6 µL |
| 30 – 300 µL 4/6/8-channel | ■ orange (for epT.I.P.S.®, 300 µL) | 30 µL | ±3 % | ±0.9 µL | ±1 % | ±0.3 µL |
| | | 150 µL | ±1 % | ±1.5 µL | ±0.5 % | ±0.75 µL |
| | | 300 µL | ±0.6 % | ±1.8 µL | ±0.25 % | ±0.75 µL |
| 120 – 1,200 µL 4/6/8-channel | ■ green (for epT.I.P.S.®, 1,250 µL) | 120 µL | ±6 % | ±7.2 µL | ±0.9 % | ±1.08 µL |
| | | 600 µL | ±2.7 % | ±16.2 µL | ±0.4 % | ±2.4 µL |
| | | 1200 µL | ±1.2 % | ±14.4 µL | ±0.3 % | ±3.6 µL |

¹⁾ The error data, according to EN ISO 8655, only apply if original Eppendorf tips are used. Technical specifications are subject to change. Errors and omissions excepted.

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

www.eppendorf.com

For more information and article numbers, visit

www.eppendorf.com/move-it-pipette

