

CHO Cell Culture with Single-Use Eppendorf BioBLU® Packed-Bed Fibra-Cel® Basket

Taylor Hatton, Shaun Barnett, and Kamal Rashid*

*Associate Director, Center for Integrated BioSystems, Utah State University, Logan, UT, USA.

Contact: becken.u@eppendorf.com

Abstract

We used Eppendorf packed-bed bioreactors filled with Fibra-Cel disks for the cultivation of Chinese Hamster Ovary (CHO) cells. The objective of this study was to compare process performance in BioBLU single-use and traditional glass packed-bed bioreactors. Alkaline phosphatase (ALKP)-secreting CHO cells were used to measure protein production in each bioreactor. Overall,

the results from these comparisons suggest that there is no significant difference between the reusable and single-use FibraCel basket systems for bench-scale production of recombinant proteins. Productivity of cells and collection of secreted proteins will not be hindered by the implementation of single-use bioreactor systems.

Introduction

The packed-bed basket technology, developed by Eppendorf, provides a shear-free environment for production of animal cells. At present, little information is available on the utility of the Eppendorf BioBLU 5p Single-Use Vessel for the production of secreted proteins, especially in perfusion mode of operation. Thus, this study was conducted to measure the growth and productivity of alkaline phosphatase (ALKP)-secreting rCHO. Two packed-bed bioreactor types were used: an Eppendorf BioBLU 5p Single-Use Vessel (3.75 L working volume) and a 3 L autoclavable glass vessel (1.25 - 3.75 L working volume),

both operated by an Eppendorf bench scale bioprocess control station in perfusion mode. The perfusion process provides a homeostatic environment for optimal cell growth similar to that experienced by cells *in vivo*, where waste products are constantly removed and fresh nutrients are replenished. Cells cultured in packed-bed bioreactors are not exposed to hydrodynamic forces, thus, allowing for maximum cell growth and protein expression (1). The objective of this study was to compare the two types of bioreactors to determine if any differences are observed between the productivity of the reusable and the single-use system.

Materials and Methods

Culture procedures

In order to evaluate the impact of these bioreactor systems on protein production, we utilized a recombinant alkaline phosphatase-secreting CHO cell line (rCHO), a proprietary

cell line provided by CDI Bioscience, Inc. (Madison, WI). The rCHO cells were engineered with the IPTG-regulated RP Shift vector so that the rCHO cells stop replicating and shift to protein production when induced with IPTG. Serum-



Fig. 1:
BioBLU 5p Single-Use Vessel
preloaded with Fibracel disks.

free CD-CHO medium (Life Technologies®, USA) was used throughout these experiments. The media contains 6.3 g/L glucose and was supplemented with 8 mM L-glutamine and 100 µg/ml of an antibiotic/antimycotic solution (Life Technologies, USA). Frozen rCHO cells were thawed and transferred to T-75 flasks with CD-CHO medium and allowed to expand. Once a sufficient number of cells were achieved, sterile disposable spinner flasks were utilized to further expand the cells. Subculture of the cells continued until a sufficient number of viable cells was achieved for use as a seed culture at the density of 5×10^5 cells/mL. Both the glass bioreactor and the BioBLU 5p Single-Use Vessel were controlled with an Eppendorf bioprocess control station.

Packed-bed basket impeller operated in perfusion mode

Two experimental trials were performed using the packed-

Table 1: Comparison of perfusion volumes.

*Perfusion occurred every other day.

Perfusion	Glass	BioBLU
Day 1	0.5 L	1 L
Day 2	1 L	2 L
Days 3 - 15*	2 L	4 L

bed vessels in perfusion mode: 3 L autoclavable vessel (1.25 - 3.75 L working volume) and a BioBLU 5p Single-Use Vessel (3.75 L working volume, pre-loaded with 150 g of Fibracel disks). The perfusion process was initiated once the cells reached the exponential growth phase as shown in table 1.

Both experimental trials had the following parameters shown in Table 2.

Table 2: Bioprocess setpoints

Parameter	Glass	BioBLU
Temperature	37° C (± 0.1°C)	37° C (± 0.1°C)
Agitation	120 rpm (± 5 rpm)	120 rpm (± 5 rpm)
DO	35 % (± 1 %)	35 % (± 1 %)
pH	7.1 (± 0.01)	7.1 (± 0.01)
Gas flow	0.5 SLPM	1.5 SLPM

Biomarkers of cell growth and productivity

Cell productivity was assessed by measuring activity of the secreted ALKP protein using an enzyme assay (AnaSpec®, USA) according to the manufacturer’s protocol. For simplicity unit measurements were used in this study. A unit (U) of ALKP activity was defined as the amount of enzyme that hydrolyzes 1µmol of p-nitrophenylphosphate to p-nitrophenol in a total reaction volume of 1 ml in 1 minute at 37°C. The YSI® 2700 Select Biochemistry Analyzer (YSI, Inc., USA) was utilized to monitor the glucose and lactate levels in the culture media every 24 hours for the duration of each trial.

Results

Glucose utilization and lactate production

Glucose is the main energy source for cell proliferation and ALKP production. Thus, glucose levels were expected to directly correlate with ALKP production in each experiment. Because lactate is a secondary energy source, lactate levels were expected to decline following this initial increase and the utilization of glucose in the media. Lactate metabolism is beneficial to the system by reducing a major metabolic by-product from the system (2,3). Glucose levels measured

at the time of induction (day 3) were nearly 0 g/L in both experiments (Fig. 2). Media lactate concentrations increased in response to decreasing glucose availability. The use of lactate as a secondary energy source can also be observed as lactate levels decrease at each 2 L perfusion.

Comparison of bioreactor systems for ALKP production

The average total ALKP production per experiment trial is shown in Figure 3. Overall, there is not a significant

difference in ALKP production between the two bioreactor systems. The total amount of ALKP measured after five

media exchanges in the reusable vessel was 17.44 U/mL and 16.22 U/mL in the single-use vessel.




Fig. 2: Glucose consumption and lactate production by rCHO cells cultured in two packed-bed bioreactor system. Values shown are the amounts of glucose and lactate measured in the culture media at each media exchange. The time and volume of the media exchange is indicated at each dashed line. Induction of ALKP activity by IPTG began on culture day 5 and continued every two days throughout the remainder of the experiment. Results of two experimental trials are shown (A, reusable; B, single-use).




Fig. 3: ALKP production by rCHO cells cultured in two packed-bed bioreactor system. (A) ALKP concentration in culture media measured each day of each experimental trial. IPTG induction of ALKP began on culture day five and continued every two days for the remainder of the experiment. (B) Stacked bar charts show the cumulative production of ALKP throughout the experiment, with each bar representing a perfusion.

Conclusion

In summary, these results demonstrated comparable yields in ALKP production between the two packed-bed bioreactor systems when operated in perfusion. Given the greater productivity of cells cultured in the packed-bed bioreactor and the advantages of this system operated in perfusion

mode, researchers desiring to scale up mammalian cell culture for protein production should consider utilization of the Eppendorf BioBLU packed-bed, single-use bioreactor system.

Literature

1. Merten, O.W., P.E. Cruz, C. Rochette, C. Geny-Fiamma, C. Bouquet, D. Goncalves, O. Danos, and M.J.T. Carrondo (2001) Comparison of different bioreactor systems for the production of high titer retroviral vectors. *Biotechnology Progress* 17(2): p. 326-335.
2. Ma, N.N., J. Ellet, C. Okediadi, P. Hermes, E. McCormick, and S. Casnocha (2009) A single nutrient feed supports both chemically defined NS0 and CHO fed-batch processes: improved productivity and lactate metabolism. *Biotechnology Progress* 25(5): p. 1353-1363.
3. Hu, W.S., T.C. Dodge, K.K. Frame, and V.B. Himes (1987) Effect of glucose on the cultivation of mammalian cells. *Dev Biol Stand* 66: p. 279-90.

Ordering information

Description	Order no.
BioFlo® 120, Advanced	
Plug type B (USA, Canada, Mexico, Japan)	B120ACS000
Plug type CEE 7/7 (EU (except UK, Ireland, Switzerland), Russia, Korea)	B120ACS001
Plug type I (Australia, New Zealand, China, Argentina)	B120ACS002
Plug type J (Switzerland)	B120ACS003
Plug type G (UK, Ireland)	B120ACS004
Plug type N (Brazil)	B120ACS005
Plug type D (India)	B120ACS006
BioFlo® 320	
Base control station: All configured units include the same base control station.	1379963011
BioBLU® 5p Packed-Bed Vessel, microsparge	M1363-0119
BioBLU® 5p Packed-Bed Vessel, macrosparge	M1363-0133
Fibra-Cel® disks, a solid support growth matrix for mammalian, animal, and insect cells	
50 g	M1292-9984
250 g	M1292-9988

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
 eppendorf@eppendorf.com

www.eppendorf.com