

Comparing Culture Methods in *Escherichia coli* Fermentation: Batch, Fed-Batch, and Continuous

Ying Yang and Ma Sha Eppendorf, Inc., Enfield, CT, USA Presenting author: yang.y@eppendorf.com

Abstract

Scope

In this study, three different culture methods, batch, fed-batch, and continuous fermentation, were carried out to grow E. coli (ATCC[®] 25922GFPTM) in the Eppendorf BioBLU[®] 3f Single-Use Vessel controlled by the BioFlo[®] 320 bioprocess control station. The objectives of this study were (1) to evaluate bioprocess parameter setup, vessel preparation, and compare how different culture methods affect microbial growth, nutrient consumption, productivity, and yield; (2) to provide a detailed cost analysis for growing E. coli under different fermentation modes and further compare their pros and cons; and (3) to demonstrate the strong performance of BioBLU 3f Single-Use Vessel in high-density microbial applications. Batch fermentation is a relatively easy operation with competitive biomass yield on glucose, and it can be applied in early development for process optimization including medium selection. With a feasible feeding strategy, fed-batch fermentation gave the highest cell dry weight at 85 g/L among all culture methods in this study, and it was very cost-effective per unit biomass. In continuous fermentation, feeding and harvest take place simultaneously at the same rate. It has the potential to achieve steady state to greatly reduce the downtime, and it is scale-up friendly since the working volume is kept constant. Therefore, depending on the experimental needs and the laboratory settings, and with a brief estimate of the process budget and scheduling, this study can help fermentation scientists to choose the ideal culture method to meet their unique needs.

How to choose the ideal fermentation method.

Bioprocess system			Cultivation strategies				Terrific Broth supplemented with
<image/>	Bioprocess system BioBLU 3f Single-Use Vessel controlled by BioFlo 320 bioprocess control station. Process parameters common to all runs Parameter Setpoint/control Vessel BioBLU 3f Inoculation density 5 % (v/v) Dissolved oxygen (DO) 30 % Agitation Magnetic drive; controlled by D0 cascade; maximum 1,200 rpm Sparger Gassing control Automatic gas flow and mix, controlled by D0 cascade Temperature 37 °C, controlled with stainless-steel cooling finger (cooling water: 18 °C) Impeller mpeller 3x Rushton-type impellers pH 7.0 ± 0.1; controlled by 25 % (v/v) amonium hydroxide	<image/>	1. Batch/ Complex medium	2. Batch/ Chemically defined medium	3. Fed-batch/ Chemically defined medium	4. Continous/ Chemically defined medium	glycerol. Used as the inoculation medium in all runs and as the fermentation medium in basic batch fermentation.
							Chemically-defined medium Citrate phosphate buffer + 90 g/L glucose + MgSO ₄ + thiamine + trace elements. Used in glucose-enriched batch fermentation.
			Batch fermentation	Fed-batch fermentation	Continous fermentation	Chemically-defined medium Citrate phosphate buffer + 15 g/L glu- cose + MgSO ₄ + thiamine + trace ele- ments. Used in fed-batch and continu-	
	Unit 1 - Cascade		With microbial growth, n ally consumed, by-produ	nutrients are gradu- licts accumulate, and	After inoculation, the microorganisms are grown under batch regime until the	Fresh medium is continuously added to the fermentor, while used medium and cells	ous fermentation.
BioBLU 3f Single-Use Vessel	From To DO Agit S-Flo S-O2 + +		the culture environment	gradually changes.	nutrients are consumed. Fresh nutrients	are harvested at the same time. Consumed	
Specifically designed for microbial applications:	Out % SP SP SP # ## \$#00		The growth curve is usua	ally divided into lag	are added in increments throughout	nutrients are replaced, and toxic metabolites	Feeding medium
> Reduced turnaround time between runs, because no	0.0 300 X 1000 1000 50.0 1200 0.0 X 965 6.0 900 65.0 3.0 0.0 X 1000 1000 1000		phase, exponential phase	e, stationary phase,	the remaining duration of termentation	are removed. When addition and removal	614 g/L glucose MgSO thiamine and
cleaning and sterilization is required	100.0 100.0 X 730 45 - 60.0		and death phase.		stationary growth phase can be extended	stavs constant.	trace elements. Used in fed-batch
> Scalability through industrial design					with extensive biomass accumulation.		fermentation.
> Rigid-wall design reduces the risk of vessel damage	25 00 0.0 20.0 40.0 Out % 60.0 80.0 100.0						
> Autoclavable material, allowing medium sterilization			E.coli strain:			Inoculum preparation	
within the vessel	EFT(HH:MM) 00:00 FT(HH:MM) FT(HH:M) F	BioBLU 3f Single-Use Vessel	ATCC 25922	GFP		Mini cell bank in 2 mL cryovial C	vernight culture in 500 mL baffled shake flask
	Cascade for DO control		30,00				

Results

1. Batch/Complex medium



> Cascade for DO control: No O₂ enrichment throughout the fermentation

2. Batch/Chemically defined medium



> Cascade for DO control: O_2 got enriched up to 40 % (data not shown)

3. Fed batch/Chemically defined medium



> Cascade for DO control: O_2 got enriched up to 100 % (data not shown) > Feeding was initiated before glucose depletion and triggered robust exponential growth

4. Continuous/Chemically defined medium



> Feeding was initiated before glucose depletion > Broth harvesting was initiated simultaneously with feeding at the

> Standard growth curve with distinct lag, exponential growth, stationary, and death phases > Simple operation with low risk of contamination

> A relatively extended lag phase > The exponential growth of *E. coli* was accompanied by drastic glucose consumption > After glucose depletion, bacteria growth stopped

> With the exponential increment of feeding rate, a significant accumulation of biomass was observed corresponding to an OD₆₀₀ increase of more than 150 within 2 hours > Glucose was depleted or almost depleted when the maximum OD₆₀₀ was reached, then started to accumulate during death phase > More complex and time-consuming handling before and during fermentation

same speed ensuring a constant culture volume throughout the run

> Glucose depletion occured at t = 7 h and remained depleted > Oxygen demand increased to 100 % at t = 10.5 h (data not shown)

> The OD_{600} of the harvested broth was 95 > Growth control via control of nutrient supply is possible

Summa	ary					
Operation mode	CDW ¹ per vessel at the end of the run (g)	Total volume (L)	Volumetric biomass productivity (g/(L*h))	Total glucose consumption (g)	Biomass yield on glucose (g CDW/g glucose consumed)	Maximum OD ₆₀₀
Batch, complex medium	9.87	3	0.37	N/A		11
Batch, chemi- cally defined medium	77.16	3	2.2	273.6	0.282	77
Fed-batch	225.6	3	6.27	713.7	0.316	240

Conclusion

> BioBLU 3f Single-Use Vessels are lightweight, easy-to-use, labor-saving, and cost-effective compared to the traditional glass vessels.

> Batch fermentation is the easiest and cheapest fermentation mode, but considering the unit cost per CDW, fed-batch fermentation was the most economic mode, for this case study.

> For batch fermentation, a switch of the medium composition increased the maximum OD₆₀₀ from 11 to 77.

> With the very cost-efficient fed-batch fermentation, we successfully reached the highest OD₄₀₀ of 240 with an expo-

nentially increasing medium feeding rate.

> A maximum OD₆₀₀ of 159 was reached with the economically competitive and scale-up friendly continuous fermentation method.

> Depending on the experimental needs, the laboratory settings, and a brief estimate of the

process budget and scheduling, the ideal fermentation method should be chosen.

> Medium composition, cultivation mode, and process control contribute to shaping the

growth curve in fermentation.

¹CDW, cell dry weight ²In continuous fermentation, the biomass is collected from the culture in the vessel (3 L) and harvest (1.3 L)

4.44

4.3²

210.0

Continous

Your local distributor: www.eppendorf.com/contact · Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany · www.eppendorf.com

846.0

0.248

ATCC[®] is a registered trademark of the American Type Culture Collection, USA. 25922GFPTM is a trademark of the American Type Culture Collection, USA. Eppendorf[®], the Eppendorf[®], the Eppendorf Brand Design, and BioBLU[®] are registered trademarks of Eppendorf AG, Germany. BioFlo[®] is a registered trademark of Eppendorf, Inc., USA. All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

Maximum Cost per unit CDW (g/L) CDW (\$/g)

3.9

27.3

85.0

56.3

159

159.9

25.5

9.3

9.7

www.eppendorf.com

0.2

0.4

1.2

1.2

1.6

2.0

2.6

3.1

5.1

8.0