

# Kluyveromyces lactis

Multiporator/Eppendorf Eporator®

## Transformation Protocol

Protocol No. 4308 915.516 – 12/2001

<b>Microorganism</b>	<i>Kluyveromyces lactis</i> AWJ137
<b>Cell type</b>	Yeast
<b>Molecules injected</b>	DNA
<b>Growth medium</b>	Not given
<b>Washing solution</b>	YPD, 10 mM DTT, 40 mM HEPES (pH 8.0)
<b>Electroporation solution</b>	2.5 mM Tris/HCl, 1 mM MgCl <sub>2</sub> , (pH 8.0)
<b>Outgrowth medium</b>	YPD with 1 M sorbitol
<b>Cuvette</b>	4 mm gap width
<b>Reference</b>	Prof. Friedhelm Meinhardt • Universität Münster • Corrensstr. 3 • D-48149 Münster e-mail: meinhar@uni-muenster.de

### Making electrocompetent cells:

1. Grow a 50 ml cell culture of *Kluyveromyces lactis* to a density of  $5 \times 10^7$  cells/ml.
2. Harvest by centrifugation (5 min., at 4,000 x g)
3. Resuspend cells in 5 ml washing solution and incubate for 5 min.
4. Wash twice with electroporation solution and resuspend in 0.25 ml of the same solution. Incubate for 1 hour on ice.

### Electroporation of cells:

1. Add 5 µl DNA (1 µg/µl) to 100 µl of electrocompetent cells. Homogenize by gently mixing with pipette several times. Transfer mixture into a prechilled cuvette.
2. Wipe moisture from the cuvette and insert the cuvette into the device.
3. Electroporation:

<b>Mode</b>	Prokaryotes "O"
<b>Voltage (V)</b>	2,500 V
<b>Time constant (τ)</b>	5 ms

4. Add 1 ml YPD medium with 1 M sorbitol, incubate for 1 hour at 30 °C.
5. Centrifuge cells 30 sec., resuspend in 0.4 ml sterile water, plate in aliquots on selective plates.

### Expected results:

Transformation efficiency up to 200 transformants/µg of DNA.

**eppendorf**