















Master Any Type of Liquid

| Type of Liquid | Potential problems | Workaround | Recommendations | | | |
|--|---|---|---|---|---|---|
| | | | Observations | Air-cushion pipettes | Positive displacement dispenser | Positive displacement pipettes |
| Water |  | <ul style="list-style-type: none"> > Air-cushion pipettes are optimized to the physical properties of water | <ul style="list-style-type: none"> > Optimally suitable for the use of water > No adaptation necessary | <ul style="list-style-type: none"> > Serial pipetting for multiple samples and vessel formats | <ul style="list-style-type: none"> > Varitip S^{*3,4} system allows accurate pipetting from large bottles and narrow vessels | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles |
| Viscous e.g. glycerol, oil |  | <ul style="list-style-type: none"> > High resistance to flow > Liquid residues stay attached to inside tip wall > Imprecise results | <ul style="list-style-type: none"> > Work slowly > Reverse pipetting > Adjust to liquid type^{*1} | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing > No adjustment to liquid type needed | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles (with Varispenser[®] 2/2x up to a viscosity of 500 mm^{2/s}) |
| Dense e.g. sulfuric acid, caesium chloride |  | <ul style="list-style-type: none"> > Influence on size of air-cushion > Dispensed volume too low or too high | <ul style="list-style-type: none"> > Adjust pipette to liquid density > Adjust to liquid type^{*1} | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing > No adjustment to liquid type needed | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles up to a density of 2.2 g/cm³ |
| Volatile e.g. acetone, ethanol |  | <ul style="list-style-type: none"> > Air-cushion expands > Liquid drips out of the tip > Imprecise results | <ul style="list-style-type: none"> > Prewet at least 5 times > Reverse pipetting > Adjust to liquid type^{*1} | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing > No adjustment to liquid type needed | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers > Varitip S system and valve for drip-free dispensing | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles up to a vapor pressure of 500 mbar |
| Infectious / radioactive e.g. biohazard material |  | <ul style="list-style-type: none"> > Aerosols contaminate pipette > Threat to human health and sample safety | <ul style="list-style-type: none"> > Use filter tips > Automated systems protect user and sample | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles |
| Detergent / detergent-containing e.g. Tween 20, Triton™ X-100 |  | <ul style="list-style-type: none"> > Reduced surface tension > Liquid residues stick to the inner wall of the tip > Imprecise results | <ul style="list-style-type: none"> > Use tips with low retention effect > Adjust to liquid type^{*1} | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles (with Varispenser[®] 2/2x up to a viscosity of 500 mm^{2/s}) |
| Foaming e.g. protein-containing liquids |  | <ul style="list-style-type: none"> > Foam is created > Liquid residues remain in the tip > Imprecise results | <ul style="list-style-type: none"> > Reverse pipetting | <ul style="list-style-type: none"> > Higher precision regardless of physical properties of liquid > Serial dispensing | <ul style="list-style-type: none"> > Varitip P^{*2} allows accurate pipetting, for example from beakers | <ul style="list-style-type: none"> > Liquid dispensing directly from supply bottles |

*1 This option is only available on automated systems and electric pipettes
*2,3,4 See Varipette[®] 4720 for corresponding Eppendorf Varitips[®]

Eppendorf Solutions

| | | | | | |
|---------------------------|--|--|--|---|--|
| Mechanical systems | Advantages <ul style="list-style-type: none"> > Easy to clean > Economical > Lightweight | <ul style="list-style-type: none"> > Eppendorf Research[®] plus > Eppendorf Reference[®] 2 > Research plus Move It[®] > Pipet Helper[®]  | <ul style="list-style-type: none"> > Multipette[®] M4  | <ul style="list-style-type: none"> > Varipette[®] 4720  | <ul style="list-style-type: none"> > Varispenser[®] 2/2x for dispensing large volumes  |
| Electronic systems | Advantages <ul style="list-style-type: none"> > High reproducibility > Ergonomic working > Multifunctionality | <ul style="list-style-type: none"> > Eppendorf Xplorer[®] (plus) > VisioNize[®] pipette manager > Xplorer plus Move It[®] > Easypet[®] 3 > epMotion[®]  | <ul style="list-style-type: none"> > Multipette[®] E3/E3x  | | <ul style="list-style-type: none"> > Eppendorf Top Buret[™] for titration  |