Keep this instruction manual and “ROTORS, TUBES, BOTTLES AND CAPS” in the file.

Tube sealer for Ultracentrifuge

STF3

Before using this equipment, please carefully read this instruction manual and “ROTORS, TUBES, BOTTLES AND CAPS” for its efficient operation and for your safety. Keep this instruction manual for your reference and refer to it as required.

Contents

| | |
|---------------------------------------|-----|
| Safety Reminder _____ | (i) |
| 1. Introduction _____ | 1 |
| 2. External View _____ | 1 |
| 3. Specifications _____ | 2 |
| 4. Installation _____ | 3 |
| 5. Summary of Sealing Procedure _____ | 4 |
| 6. Operation _____ | 5 |
| 7. Seal Tube and Tube Rack _____ | 8 |
| 8. Maintenance After Operation _____ | 9 |
| 9. Troubleshooting _____ | 10 |
| 10. Packing List _____ | 12 |
| 11. Decontamination _____ | 13 |

Safety Reminder

Common safety precautions

The equipment is designed to seal the resin tube which contains liquid-suspended materials having different densities and particle size.

Carefully read the following safety precautions for a thorough understanding.

- Follow the instructions and procedures described in this manual to operate this equipment safely.
- Carefully read all safety messages in this manual and safety instructions on your equipment.
- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- Safety messages are labeled as indicated below. They are in combination with signal words of “WARNING” and “CAUTION” with the safety alert symbol to call your attention to items or operations that could be dangerous to you or other persons using this equipment. The definitions of signal words are as follows:

 **WARNING** : Warning notes indicate any condition or practice, which if not strictly observed, could result in serious injury or possible death.

 **CAUTION** : Caution notes indicate any condition or practice, which if not strictly observed or remedied, could result in personal injury, damage or destruction of the equipment.

NOTE : Notes indicate an area or subject of special merit, emphasizing either the product’s capabilities or common errors in operation or maintenance.

- Do not operate this equipment in any manner not described in this Instruction Manual. Should you have any troubles with this equipment, contact an authorized sales or service representative.
- The precautions described in this Instruction Manual are carefully developed in an attempt to cover all the possible risks. However, it is also important that you are alert for unexpected incidents. Be careful when operating this equipment.

Safety Reminder

WARNING

1. This equipment is not explosion proof equipment. Never use explosive or flammable samples or materials that chemically react vigorously.
2. Do not install the equipment in or near places where inflammable gases are generated or chemicals are stored.
3. Take all necessary safety measures before using samples that are toxic or radioactive, or blood samples that are pathogenic or infectious. You use such samples at your own responsibility.
 - 1) If the centrifuge, rotor, or an accessory (this equipment etc.) is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
 - 2) If there is a possibility that the centrifuge, rotor, or an accessory (this equipment etc.) is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory (this equipment etc.) properly before requesting repairs from an authorized sales or service representative.
 - 3) It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, an accessory (this equipment etc.), or parts properly before returning them to an authorized sales or service representative.
4. Your tube sealer must be grounded properly.
5. To avoid electrical shocks, do not detach the cover from this tube sealer.
6. To prevent electrical shocks, turn off the POWER switch and disconnect the power cord from the receptacle before replacing the fuse with the new one.
7. To avoid electrical shocks, do not handle the power cord or turn on or off the POWER switch with wet hands.
8. Repairs, disassembly, and other modifications to the equipment are strictly prohibited unless performed by the authorized service representative.





Safety Reminder

CAUTION

1. Do not touch the heater unit while this equipment is in operation. The heater unit is very hot and your hands may be burnt if touched.
2. If you observe some abnormality in this product, stop using it immediately and contact an authorized sales or service representative.
3. Remove the power cord from the receptacle when this equipment is not in use for a long time.
4. Use the provided power cord that is designated for this equipment by us.
Do not use the provided power cord for a machine that is not designated by us.
5. Depending on the magnitude, an earthquake might damage the equipment. If you observe some abnormality, contact an authorized sales or service representative.

Safety Reminder


 Warning Label and  Caution Label

1. Introduction

The STF3 tube sealer is a heat-welding machine to seal the seal tubes used for angle rotors vertical rotors, and swing rotors. For safe and proper use of the STF3 tube sealer, carefully read and follow the instructions below.

2. External View

*Tube rack is optionally available.

Figure 1. External View of STF3 Tube Sealer

3. Specifications

Sealing system: Direct welding by means of the molding heater

Power requirements:

110-120 V equipment: 99~132 V AC, 2 A, 50/60 Hz

200-240 V equipment: 180~264 V AC, 2 A, 50/60 Hz

Dimensions: 132 (W) x 255 (D) x 260 (H) mm (when raising the handle: 320 (H) mm)

Weight: 8 kg

 WARNING:

- This tube sealer can be damaged if connected to a wrong voltage.

CE

The STF3 satisfies CE marking requirements. The CE marking is an international symbol, which shows that the product conforms to EC directives.

Standards concerning these directives are as follows:

- Low Voltage Directive (2014/35/EU)
EN 61010-1:2010
Environment requirement: Altitude 2000M max, relative humidity \leq 80% between 5°C to 31°C and relative humidity \leq 50% between 31°C to 40°C
Installation category II ,Pollution degree 2
- EMC Directive (2014/30/EU)
EN 61326-1: 2013 Class A
- RoHS Directive (2011/65/EU)
EN 50581: 2012

4. Installation

⚠ CAUTION:

- When moving this equipment, support the bottom of this equipment.

This section describes the instructions that you should abide by when installing the equipment to ensure your safety and the maximum performance.

(1) Power source

110-120 V equipment: 99~132 V AC, 2 A, 50/60 Hz

200-240 V equipment: 180~264 V AC, 2 A, 50/60 Hz

- When abnormalities arise on a main part, please draw out a power plug from an electric socket.
- Do not position an object so that it is difficult to disconnect the power cord from the receptacle. If you do so, you cannot disconnect the power cord from the receptacle when you observe some abnormality.

(2) Location

- Position this equipment to a level surface, such a table or laboratory bench that can support the weight of this equipment and resist vibration.
- Ambient temperature for performance: 15°C to 25°C.
- Ambient temperature for operation: 6°C to 32°C
- Do not locate it under the direct sunlight.
- Do not install near by a heat-generating device or waterworks from where water may drip or splash. Such location may incur temperature rise in the precious samples or malfunction of this equipment.

(3) Connection of the power cord and grounding


⚠ WARNING:

- Do not touch the power cord with wet hands to avoid electrical shocks.
- Do not hold the cord when disconnecting the power cord from the receptacle. Hold the plug when disconnecting it.
- This equipment must be grounded properly.


⚠ CAUTION:

- Use the provided power cord that is designated for this equipment by us. Do not use the provided power cord for a machine that is not designated by us.


The power cord plug has a three-prong plug. One of the three prongs is provided for grounding. This equipment can be grounded just by putting the three-prong plug into a three-prong receptacle. Do not connect the three-prong plug to a two-prong outlet using an adapter. Do not hold the cord when disconnecting the power cord from the receptacle.


5. Summary of Sealing Procedure

Position the tube inlet just under the heater.

Welding ~ Cooling

Contour after welding

- (1) Adjust the position of the tube rack as shown in the figure so that the tube inlet is positioned just under the heater.

NOTE:

- Inject the sample into the tube up to the base of the inlet.
- Be sure to wipe away sample adhered on the surface of the inlet.

- (2) Pull the handle down until the heater contacts the tube inlet and keep at the position for 1 to 2 seconds. Then pull the handle down to the end and wait until the HEAT lamp goes out (about 40 seconds).

NOTE:

- Keep contacting the heater with the tube inlet for 1 to 2 seconds.
- If the sample adheres to the heater, remove the cord from the receptacle. Clean the heater using cotton swabs moistened with water and dry it well.

- (3) Check the height of the welded part.

Points

Check that the height of the welded part meets the following specifications.

| Tube size | \varnothing / L dimension (mm) |
|-----------------|----------------------------------|
| 1.5PA seal tube | L = 27.7 ~ 28.5 |
| 2PA seal tube | L = 34.5 ~ 35.3 |
| 40PA seal tube | \varnothing = 5.5 ~ 6 |
| 6PA seal tube | \varnothing = 5 ~ 6 |
| 10PA seal tube | |
| 15PA seal tube | |
| 20PA seal tube | |
| 30PA seal tube | |
| 33PA seal tube | \varnothing = 5 ~ 5.5 |
| Other tubes | |

- (4) Set the tube in the rotor and set the space cap on it. For some tube and rotor combinations, a float spacer must be set before setting the space cap.

NOTE:





- Use the correct space cap and the correct float spacer according to the seal tube and the rotor in use.



If sample leakage is found after centrifugation, carefully read this manual again and check the [⚠ WARNING] information and [⚠ CAUTION] information. For safety, do not perform centrifugation if sample leakage occurs.




6. Operation

⚠ CAUTION:

- If you observe some abnormality in this product, stop using it immediately and contact an authorized sales or service representative.
- Never use any tube or tube rack that is not designated for this equipment by us.
- Do not touch the heater unit while this equipment is in operation. The heater unit is very hot and your hands may be burnt if touched. Do not put the power cord near the heater unit while this equipment is in operation.

| Step | Operation and precaution | Illustration |
|------|---|--|
| (1) | <p>Inject the sample into the tube and check the balance (refer to the rotor instruction manual for the allowable imbalance) Inject the sample up to the base of the inlet.</p> <p>⚠ CAUTION: If the sample in the tube is not enough, the tube may be deformed or damaged during centrifugation.</p> |  |
| (2) | <p>Set the seal tubes in the proper tube rack (see Section 7). Before setting, check that no liquid such as sample is adhered on both inside and outside surfaces of the tube inlet. If adhered, wipe the liquid off the surface.</p> <p>⚠ CAUTION: Be sure to wipe the liquid off the surface. Failure to do so may cause imperfect welding and the heater failure may occur.</p> | <p>*Be sure to wipe the liquid off the surface.</p>  |
| (3) | <p>Set the tube rack on the guide of the STF3 tube sealer.</p> | |
| (4) | <p>Turn on the POWER switch. (POWER lamp lights up.)</p> <p>NOTE: The fuse is blown if the POWER lamp does not light up. Replace it with the 5A fuse (accessory).</p> | |
| (5) | <p>Preheat the heater before operation. Move the tube rack to the position where the tube does not contact the heater. Pull the handle forward. Then the HEAT lamp lights up. Wait until the HEAT lamp goes out (about 40 seconds) and put the handle back to the original position.</p> |  |
| (6) | <p>Move the tube rack from side to side then move the seal tube to and fro so that the tube inlet is positioned just under the heater.</p> <p>⚠ CAUTION: If the positioning before welding is improper, the seal tube may be improperly welded and deformed. Do not put such a tube in the centrifugal operation. Otherwise sample leakage may occur.</p> |  |

| Step | Operation and precaution | Illustration |
|------|---|--|
| (7) | <p>Pull the handle forward. Then the HEAT lamp lights and the heater heats up. Gradually pull the handle down until the heater contacts the inlet of the seal tube. Keep at this position for <u>1 to 2 seconds</u>.</p> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p>⚠ CAUTION:</p> <p>① As a guide, count slowly from 1 to 2 by yourself.</p> <p>② If the contacting time is too long, the heater temperature becomes high and it may cause imperfect welding.</p> </div> |  <p>Keep contacting the heater with the tube inlet for 1 to 2 seconds.</p> |
| (8) | <p>Pull the handle down to the end. The heater moves downward melting the inlet of the seal tube.</p> |  |
| (9) | <p>Keep pulling the handle down. The seal tube is welded and the heater automatically turns off upon completion of welding. The Heat lamp keeps lighting for about 40 seconds while the welded part is cooled.</p> | |
| (10) | <p>Put the handle back to the original position <u>when the HEAT lamp goes out</u>.</p> <div style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p>NOTE:</p> <p>When you pull the handle up and down time after time, although the HEAT lamp lights up, the seal tube may not be welded. If the seal tube is not welded, keep pulling the handle up for 40 seconds and then repeat the above procedure in Step (7), (8), and (9).</p> </div> | |
| (11) | <p>Slide the tube rack and seal the next tube in the same manner as above.</p> | |
| (12) | <p>When all the tubes are sealed, turn off the POWER switch.</p> | |

| Step | Operation and precaution | Illustration | | | | | | | | | | | | | | | | | |
|-----------------|---|--|-----------------------------------|-----------------|-------------|---------------|-------------|----------------|----------------------|---------------|--------------------|----------------|----------------|----------------|----------------|----------------|-------------|----------------------|--|
| (13) | <p>Check the height of the welded part. Checkpoints differ depending on the tube size. Check the tube size according to the following tables.</p> <table border="1" data-bbox="261 387 858 913"> <thead> <tr> <th>Tube Size</th> <th>\varnothing / L dimensions (mm)</th> </tr> </thead> <tbody> <tr> <td>1.5PA seal tube</td> <td>L=27.7~28.5</td> </tr> <tr> <td>2PA seal tube</td> <td>L=34.5~35.3</td> </tr> <tr> <td>40PA seal tube</td> <td>\varnothing =5.5~6</td> </tr> <tr> <td>6PA seal tube</td> <td rowspan="6">\varnothing =5~6</td> </tr> <tr> <td>10PA seal tube</td> </tr> <tr> <td>15PA seal tube</td> </tr> <tr> <td>20PA seal tube</td> </tr> <tr> <td>30PA seal tube</td> </tr> <tr> <td>33PA seal tube</td> </tr> <tr> <td>Other tubes</td> <td>\varnothing =5~5.5</td> </tr> </tbody> </table> | Tube Size | \varnothing / L dimensions (mm) | 1.5PA seal tube | L=27.7~28.5 | 2PA seal tube | L=34.5~35.3 | 40PA seal tube | \varnothing =5.5~6 | 6PA seal tube | \varnothing =5~6 | 10PA seal tube | 15PA seal tube | 20PA seal tube | 30PA seal tube | 33PA seal tube | Other tubes | \varnothing =5~5.5 |  |
| Tube Size | \varnothing / L dimensions (mm) | | | | | | | | | | | | | | | | | | |
| 1.5PA seal tube | L=27.7~28.5 | | | | | | | | | | | | | | | | | | |
| 2PA seal tube | L=34.5~35.3 | | | | | | | | | | | | | | | | | | |
| 40PA seal tube | \varnothing =5.5~6 | | | | | | | | | | | | | | | | | | |
| 6PA seal tube | \varnothing =5~6 | | | | | | | | | | | | | | | | | | |
| 10PA seal tube | | | | | | | | | | | | | | | | | | | |
| 15PA seal tube | | | | | | | | | | | | | | | | | | | |
| 20PA seal tube | | | | | | | | | | | | | | | | | | | |
| 30PA seal tube | | | | | | | | | | | | | | | | | | | |
| 33PA seal tube | | | | | | | | | | | | | | | | | | | |
| Other tubes | \varnothing =5~5.5 | | | | | | | | | | | | | | | | | | |
| (14) | <p>Set the sealed tubes in the rotor and the space caps on them. For some sealed tube and rotor combinations, a float spacer must be set before setting the space cap. Perform centrifugation by following the rotor instruction manual.</p> <div data-bbox="240 1263 874 1534" style="border: 1px solid black; padding: 5px;"> <p>⚠ CAUTION: Applicable space caps and applicable float spacers are different depending on the types of seal tubes and the rotor model in use. Use proper space caps and proper float spacers according to the “ROTOR, TUBES, BOTTLES AND CAPS Instruction Manual (Part No.S999204)”.</p> </div> <div data-bbox="240 1601 874 1765" style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p>NOTE: The tube might be partially deformed, especially the center side of the rotation. There is no problem with this deformation.</p> </div> | <ul style="list-style-type: none"> ● When using a space cap  <ul style="list-style-type: none"> ● When using a space cap and float spacer  | | | | | | | | | | | | | | | | | |

7. Seal Tubes and Tube Rack

⚠ CAUTION:

- Use of a tube rack that is not suitable for your seal tubes can cause imperfect welding or sample leakage because the length of the welded part is out of the specification.

Use a proper tube rack that is suitable for your seal tubes. Refer to Table 1 for applicable tube rack for each seal tubes.

Table 1. Correspondence of seal tubes and tube racks

| Applicable seal tubes | Tube racks | Part No. | Tube sealers | | |
|--|----------------|----------|--------------|------|------|
| | | | STF-1 | STF2 | STF3 |
| 1.5PA seal tube | Tube rack (G2) | S201778G | | | |
| 2PA seal tube * | Tube rack (G) | S201778F | | | |
| 3.5PA seal tube | Tube rack (B2) | S201778E | ○ | ○ | ○ |
| 4PA seal tube | Tube rack (B3) | S201778H | | | |
| 5PA seal tube | Tube rack (B) | S201778A | | | |
| 6PA seal tube | Tube rack (C4) | S206017A | × | ○ | ○ |
| 6.5PA seal tube * | Tube rack (B4) | S201778J | ○ | ○ | ○ |
| 8PA seal tube | Tube rack (C2) | S201778L | | | |
| 10PA seal tube | Tube rack (D) | S206017E | × | ○ | ○ |
| 12PA seal tube | Tube rack (C) | S201778B | ○ | ○ | ○ |
| 15PA seal tube | Tube rack (C3) | S206017F | | | |
| 20PA seal tube | Tube rack (E4) | S206017C | | | |
| 30PA seal tube | Tube rack (E3) | S206017B | × | ○ | ○ |
| 33PA seal tube | Tube rack (E2) | S206017D | | | |
| 40PA seal tube | Tube rack (E) | S201778C | ○ | ○ | ○ |
| 94PA seal tube | Tube rack (F) | S201778D | ○ | × | × |
| | Tube rack (F2) | S201778M | × | ○ | ○ |
| 2PA,4PA,5PA,6.5PA,12PA, 40PA seal tubes | *Multi-rack | S201778K | ○ | ○ | ○ |

○: Compatible ×: Not compatible

*Multi-rack can be used with 2PA, 4PA, 5PA, 6.5PA, 12PA, and 40PA seal tubes only.

8. Maintenance After Operation

 CAUTION:

- Using cleaning or sterilization methods other than recommended in this instruction manual may cause in corrosion or deterioration of this equipment
- Disconnect the power cord from the receptacle before cleaning or sterilizing the equipment.
- If the equipment is exposed to ultraviolet rays for a long time, the color of the covers may be changed or the coating may be peeled off. After use, cover the equipment with a cloth to protect it from direct exposure.
- Do not pour water, neutral detergent or disinfectant solution directly over the equipment. Otherwise fluids may leak into the inside of the equipment and cause corrosion or deterioration of the equipment.

- If the equipment is heavily stained, clean it with a cloth or sponge moistened with a dilute solution of neutral detergent.

9. Troubleshooting

(1) The melted tube inlet is stuck to the heater and the tube is raised when pulling up the handle (see figure 2).

- The cooling time is insufficient. Wait until the HEAT lamp goes out then pull up the handle.


Figure 2. Improper welding

(2) The HEAT lamp is blinking.

- The screws are loose. Turn off the POWER switch and remove the cord from the receptacle. Then retighten the screws securely and turn on the POWER switch. If the problem still persists, contact an authorized sales or service representative.

(3) Although the HEAT lamp lights up, the seal tube is not welded.

- When you pull the handle up and down time after time, although the HEAT lamp lights up, the seal tube may not be welded. If the seal tube is not welded, keep pulling the handle up for 40 seconds and then repeat sealing the tube (see “6. Operation (10)”).

(4) Sample leakage occurs during centrifugation even though the tube is sealed properly.

Check as follows:

① Rotor

- Is the rotor cap tightened securely?
- Is the cap packing (O-ring) deteriorated or scratched?

② Space cap and Float spacer

- Is the specified space cap (Refer to the “ROTORS, TUBES, BOTTLES AND CAPS Instruction Manual (Part No.S999204)”.) used?
- Is the specified float spacer (Refer to the “ROTORS, TUBES, BOTTLES AND CAPS Instruction Manual (Part No.S999204)”.) used?
- Is the space cap deformed or damaged?
- Is the float spacer deformed or damaged?

③ Welding

- Is the sample amount enough?
- Is any liquid adhered on the surface of the tube inlet?
- Recheck the height of the welded part according to “4. Operation (13)”.
- If the welded part is not smooth, the heater may be damaged. Contact an authorized sales or service representative.

⚠ CAUTION:

- If the sample adheres to heater unit, the heater may be corroded and a poor contact may occur. If any sample adheres, immediately remove the cord from the receptacle. Clean the heater using cotton swabs moistened with water and dry it well.

(5) The POWER lamp does not light up when the POWER switch is turned on.

The fuse is blown.

Turn off the POWER switch and remove the power cord from the receptacle.


Insert the screwdriver etc. in the cover of the Inlet and remove the cover of the Inlet.

The fuse is mounted on the cover of the Inlet.



Replace the fuse with the new one (accessory) and put the cover of the Inlet to the original position.

⚠ WARNING:

- Turn off the POWER switch and disconnect the power cord from the receptacle before replacing the fuse with the new one. Be sure to use the fuse that is designated for this equipment by us.

10. Packing List

| Item name | Part No. | Quantity | Item drawing | Note |
|--------------------|-----------------------|----------|--|--|
| Instruction manual | S998379 | 1 |  | |
| Power cord* | S204816 (200-240V) | 1 |  | |
| | S204746 (110-120V) | 1 |  | |
| Fuse** | S413093 | 2 |  | FSL 250V 5A Rating: 250V 5A Melting Time: 210% 2min. Max. |

*Use the provided power cord that is designated for this equipment by us.
Do not use the provided power cord for a machine that is not designated by us.

**Because a current exceeding 2A flows temporarily through this equipment, 5A fuse is used in this equipment, though the rated current of this equipment is 2A.

11. Decontamination

 WARNING:

- If the centrifuge, rotor or an accessory (this equipment etc.) is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory (this equipment etc.) is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory (this equipment etc.) properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory (this equipment etc.) unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, the accessory (this equipment etc.), or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the centrifuge, rotor, the accessory (this equipment etc.), or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination.
Note that we cannot repair or inspect the centrifuge, rotor or the accessory (this equipment etc.) unless sterilization or decontamination is completed.

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Copyright © 2022 Eppendorf Himaс Technologies Co., Ltd.
Copyright © 2022 Eppendorf Himac Technologies Co., Ltd. All rights reserved.