

Stainless steel rotor lid for the Rotor F-45-12-11 – Instructions for use

These instructions for use for the stainless steel rotor lid for the rotor F-45-12-11 are a supplement to the operating manual for the centrifuge and are not a substitute, so also read the operating manual before using the rotor lid for the first time. The current version of the operating manual can be found on the Internet at www.eppendorf.com.

1. Introduction

The Eppendorf stainless steel rotor lid has been designed for the Rotor F-45-12-11 of the MiniSpin[®], MiniSpin[®] plus and MiniSpin[®] SPACE.



Please read this supplement sheet prior to using the rotor lid for the first time.

2. Safety instructions and applicational limitations

For your own safety please be sure to observe the following instructions:

- The stainless steel rotor lid may only be used in the centrifuges MiniSpin[®] / MiniSpin[®] plus and MiniSpin[®] SPACE for the intended use specified in the operating manual. Only original replacement parts, rotors and rotor lids from Eppendorf should be used.
- Rotors must always be fixed in position according to the instructions. Centrifugation may only be carried out with a firmly tightened rotor.
- Be sure to protect the rotor lid from mechanical damage. Even slight scratches and tears can lead to serious internal material damage. Avoid damage caused by aggressive chemicals, including among others: Strong and weak alkali, strong acids, solutions with mercury, copper and other heavy metals, chlorinated hydrocarbons, concentrated saline solutions and phenol. Where contamination from aggressive media occurs, clean the rotor immediately with a neutral detergent.
- Rotor lids showing visible signs of corrosion or mechanical damage should not be used. Inspect the rotor lid regularly.
- Be sure to close the tube lids tightly prior to centrifugation. Open tube lids can be torn off during centrifugation and damage the rotor lid or centrifuge.
- When using poisonous, radioactive liquids or pathogenic germs of Risk Group II (see World Health Organization: Laboratory Biosafety Manual), the appropriate national regulations are to be observed. Careful, proper cleaning of the centrifuge must be carried out following accidental spillage of such liquids into the rotor or the rotor area. Before using any other cleaning or decontamination method than that given in point 5 "Cleaning", please be sure to ascertain from Eppendorf that the intended method will not damage the lid.

3. Delivery package

1 stainless steel rotor lid incl. special rotor nut (see also point 4: Operation).

4. Operation

The accompanying special rotor nut is required for using the stainless steel rotor lid. Please replace your standard rotor nut with this special rotor nut. Tighten the special rotor nut clockwise to lock the rotor.



When you use a MiniSpin[®], MiniSpin[®] plus or MiniSpin[®] SPACE manufactured prior to 10/2003 and delivered with a polypropylene rotor lid, the stainless-steel rotor lid can only be used when the original rotor nut is exchanged for the new special rotor nut. Please only use the special rotor nut together with the stainless-steel rotor and **under no circumstances** together with the polypropylene rotor lid.

5. Cleaning

Use a mild cleaning agent for cleaning.

Select a disinfection method which complies with the legal regulations and guidelines for your range of application. For example, use alcohol (ethanol, isopropanol) or other disinfectants containing alcohol.

6. Technical data

Rotor lid material: Stainless steel V2A
Rotor nut material: Stainless steel V2A

eppendorf

Your local distributor: www.eppendorf.com/worldwide

Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: eppendorf@eppendorf.com
Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA
Tel: +1 516 334 7500 · Toll free phone: +1 800-645-3050 · Fax: +1 516 334 7506 · E-mail: info@eppendorf.com

Application Support Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)
support@eppendorf.com
North America: Tel: +1 800 645 3050 · E-mail: techserv@eppendorf.com
Asia Pacific: Tel: +60 3 8023 6869 · E-mail: support_asiapacific@eppendorf.com