

## **Instructions for Use** Thermoblock DWP 2000 English (EN)

### **1 About this manual**

#### **1.1 Notes about this manual**

This document supplements the operating manual of the product. This document does not replace the operating manual.

1. Read these instructions for use before using this product.
2. In addition, read the operating manual of the device with which you are using the accessories.

#### **1.2 Other applicable documents**

The following documents supplement this manual:

- Hardware manual of the device
- Software manual epBlue with MultiCon PC

### **2 Safety**

#### **2.1 Residual risks when used as intended**

If the product is not used as intended, the installed safety devices may not function correctly. To reduce the risk of personal injury and material damage and to avoid dangerous situations, please observe the general safety instructions.

##### **2.1.1 Risks of burns**

The thermoblock can reach temperatures of up to 96 °C. Persons may burn themselves on the thermoblock.

- Do not touch the hot thermoblock.
- Before working on the device, check that the thermoblock has cooled down.

### **3 Product description**

#### **3.1 Features**

The following labware can be inserted into the Thermoblock DWP 2000:

- Eppendorf Deepwell Plates 96/2000 µL

In epBlue, the Thermoblock DWP 2000 is not permanently combined with a labware. A separate labware definition is available for the Thermoblock DWP 2000.

### 3.2 Product overview




Fig. 3-1: Thermoblock DWP 2000

## 4 Maintenance

### 4.1 Cleaning

Information on cleaning can be found in the hardware manual of the device.

## 5 Disposal

### 5.1 Preparing for disposal

#### Preparing disposal in accordance with legal regulations



For information on the legal regulations that apply in your country, please contact your local competent authority or your Eppendorf partner.



Dispose of non-decontaminable devices as hazardous waste.

1. Check which legal regulations apply to disposal in your country.
2. Choose a certified disposal company or contact your Eppendorf partner.

#### Creating a decontamination certificate

Prerequisites:

- The device has been decontaminated.
1. Download a decontamination certificate from our webpage [www.eppendorf.com](http://www.eppendorf.com).
  2. Complete the decontamination certificate.

## 5.2 Handing over the device to the disposal company

1. Inform the disposal company of any hazards posed by the device, e.g., locking devices, flammable substances.
2. Hand over the device and the decontamination certificate to the certified disposal company.

## 6 Ordering information

| Description   | Order no.    |
|---|--------------|
| <b>Thermoblock DWP 2000</b><br>for 2 mL Eppendorf Deepwell Plates 96/2000 | 5075 751 330 |