

# Automating Cell-based Apoptosis Assays with the Eppendorf epMotion® 5075t\* increases the reproducibility

Eric Gancarek, Sandrine Hamels, Muriel Art

Eppendorf Application Technologies S.A., Namur, Belgium

## Abstract

Cell-based assays have become an established and widely used tool in modern life sciences, whether basic research or routine drug screening. Performing cell-based assays can be laborious and time consuming when multiple parameters such as different cell types, compounds, concentrations etc. have to be assessed at the same time. Medium scale benchtop automation can significantly facilitate this work, especially when multiple 96 or 384 well plates have to be processed, while simultaneously increasing quality and results consistency, making different sets of experiments more comparable.

In this study HeLa and Jurkat cell lines were cultivated

and treated with different concentrations of staurosporine, a compound well characterized to induce apoptosis. The level of apoptosis was determined via the Apo-ONE® Homogeneous Caspase-3/7 Assay (Promega®) that generates a fluorescent signal. A comparison of a manual and automated approach shows that automating the entire workflow of an apoptosis assay including cell seeding, compound dilution and addition, and assay setup with the Eppendorf epMotion 5075t automated liquid handling system generates reliable and consistent results without the risk of contaminating the cultured cells and also increases the reproducibility.

## Introduction

Originally, drug discovery researches were conducted on animal models. Over the past few years, the emergence of cell-based assays specifically developed for high or medium-throughput screening allowed to limit this practice to an absolute minimum, especially for ethical and economic reasons. Compared to biochemical assays, also commonly used for therapeutic agent screening, *in vitro* information provided by cell-based assays is considered as more pertinent as the drug response is evaluated in a cellular context which is more representative of the real-life [1]. Among all cell responses which can be evaluated, apoptosis is a subject being extensively studied since the early 1990s.

This natural physiological process, also called programmed cell death, contributes to regulation of cell populations in the human body by balancing cell division. Apoptosis plays an essential role during nervous and immune system development, takes part in wound healing and is also needed to destroy cells representing a threat such as auto-aggressive immune cells or cells infected with virus. When this normal process is halted or impaired, it can also be involved in disease formation. Some cancers can be induced by cells not correctly removed by apoptosis, which are subsequently maintained and finally becoming immortal. Excessive apoptosis can also inflict tissue damages as

observed in Alzheimer's, Huntington's and Parkinson's diseases [2]. Mechanisms of apoptosis are complex and characterized by specific morphological and biochemical features, including activation of caspases. Just like a large variety of assays developed for detecting apoptosis, the Apo-ONE Homogeneous Caspase-3/7 Assay provided by Promega and used in this study is based on caspase activity detection via cleaved fluorescent substrates [3]. If cell-based assays development is crucial for the drug screening process,

automation has also to be considered. Compatibility between liquid handling workstation and cell-based assay has to be demonstrated. All instrument settings have to be defined to ensure correct solution dispensing without damaging mono-layers of adherent cells by mechanical forces. Finally, assay reproducibility is an essential factor which has to be evaluated [4]. This study shows the possibility to implement a widely used cell-based apoptosis assay on the Eppendorf epMotion 5075t

## Material and Methods

### Methods for automated Cell-based assay

The complete workflow is programmed to process three 96-well plates in parallel and is divided into three epMotion methods. For the first two methods, the epMotion 5075t surfaces and tools are cleaned using a disinfection solution. Both UV-lights and air filters are started 15 minutes before using the system, whereas UV stopped automatically after this time span. At the end of each method, a user intervention is requested to handle the plates to downstream steps (incubation, reading).

#### Method 1:

Cell seeding is performed during this first protocol step. Adherent HeLa cells are seeded into 96-well cell culture plates with transparent bottom at 30,000 cells per well. With suspension Jurkat cells, a density of 10,000 cells per well is dispensed. Before starting the epMotion method 1, a cell solution at the appropriate concentration (10,000 cells per 90  $\mu$ L for Jurkat cells and 30,000 cells per 90  $\mu$ L for HeLa cells) is prepared in a sterile tube and transferred to a sterile autoclaved epMotion reservoir of 100 mL on the epMotion 5075t. Cell culture media, used as blank, is transferred to a sterile epMotion reservoir of 30 mL size. These steps need to be performed under sterile cell culture conditions. The method 1 generates the dispensing of 90  $\mu$ L of cells per well in column 1 to 11 in three 96-well microplates and 90  $\mu$ L of culture medium without cells in column 12 in these three plates. At the end of the program, the lid is manually replaced on the plate and the plates are placed into the CO<sub>2</sub> incubator at 37 °C for 24 hours.

#### Method 2:

The goal of the second method is to generate the concentration curve of the apoptosis inducer. Compound concentration range is between 1 nM and 20  $\mu$ M for HeLa cells and between 3 nM and 100  $\mu$ M for Jurkat cells.

16 increasing staurosporine (STS) concentrations are used for inducing apoptosis. In order to generate a concentration curve with the same amount of vehicle solution in each well, the dilution curve is produced by using two different dilution steps. The first dilution is performed into DMSO 99.9 % vehicle solution in a Deepwell plate for generating 16 STS concentrations 10 times more concentrated than desired. Culture medium without STS and DMSO being used as negative control is also added into the Deepwell Plate. At the end of this first part, the 96-well plates seeded with the cells are placed on the worktable. The second dilution step is performed into the 96-well plates and allows to obtain the 16 final STS concentrations. 10  $\mu$ L of each STS concentration is dispensed from the Deepwell plate to the 96-well plates. 10  $\mu$ L of culture medium without STS and DMSO is also added as control. The final vehicle (DMSO) concentration did not exceed 1 % and was equivalent for all staurosporine concentrations tested. The plates are transferred on the TMX of the epMotion 5075t and mixed at 500 rpm for 30 sec. At the end of the method, the lid is replaced on the plates manually. The optimal apoptosis induction time is cell-type dependent, plates seeded with HeLa cells were incubated 24 hours at 37 °C (5 % CO<sub>2</sub>) while 5 hours were optimal for stimulating Jurkat cells.

#### Method 3:

After incubation, the plates are removed from the CO<sub>2</sub> incubator. This method includes the dispensing of 100  $\mu$ L freshly prepared Apo-ONE Caspase-3/7 reagent to each well of three 96-well plates. The solution is first dispensed into the blank and negative controls. The plates are transferred to the ThermoMixer® and mixed at 500 rpm for 30 sec. At the end of this method, the lid is replaced on the plates manually and the plates are incubated for 1 hour at room temperature protected from light. Fluorescence generated by the caspase substrate cleavage was read in each well at two wavelengths (excitation at 485 nm and emission at 535 nm).

Those methods allow the generation of three 96-well plates, but can be easily adapted for a different number of plates. In each plate, 16 apoptosis inducer concentrations are

evaluated in five replicates. Blank and negative controls are also included. The final 96-well plate layout is presented below:

	1	2	3	4	5	6	7	8	9	10	11	12
A	CC 1	CC 1	CC 1	CC 1	CC 1	CC 9	CC 9	CC 9	CC 9	CC 9	Negative control	Blank
B	CC 2	CC 2	CC 2	CC 2	CC 2	CC 10	CC 10	CC 10	CC 10	CC 10		
C	CC 3	CC 3	CC 3	CC 3	CC 3	CC 11	CC 11	CC 11	CC 11	CC 11		
D	CC 4	CC 4	CC 4	CC 4	CC 4	CC 12	CC 12	CC 12	CC 12	CC 12		
E	CC 5	CC 5	CC 5	CC 5	CC 5	CC 13	CC 13	CC 13	CC 13	CC 13		
F	CC 6	CC 6	CC 6	CC 6	CC 6	CC 14	CC 14	CC 14	CC 14	CC 14		
G	CC 7	CC 7	CC 7	CC 7	CC 7	CC 15	CC 15	CC 15	CC 15	CC 15		
H	CC 8	CC 8	CC 8	CC 8	CC 8	CC 16	CC 16	CC 16	CC 16	CC 16		

CC: Compound Concentration

Figure 1: Final plate layout.

## Results and Discussion


### Cell-based Assay Automation

Cell-based assays commonly contain three successive stages: sterile cell seeding, specific cell treatment including compound dilution and assay reagent addition. Essentially based on liquid handling, implementation of cell-based assays on automation platforms has been demonstrated as feasible. Depending on integrated equipment (shaker, incubator, plate reader), the assay procedure can be partly [5][6], or completely automated. Through the CleanCap configuration available for the Eppendorf epMotion 5075t, the workstation enclosure can be used to maintain a clean environment and the entire cell-based assay protocol can be carried out with the automation platform. The complete workflow has been divided into three methods. For steps requiring sterile environment, epMotion 5075t surface and tools were disinfected beforehand (steps 1 and 3). Each method corresponds to a protocol part which is independently processed by the epMotion 5075t. Once a method is running, all manipulations are performed by liquid handling workstation. Only plate incubation at 37 °C and data reading were performed offline.

### Comparison of manual versus automated cell-based assay performances.

The transferability of the Apo-ONE Homogeneous Caspase-3/7 Assay to the epMotion 5075t was assessed with HeLa and Jurkat cell lines, both treated with a staurosporine dose-response curve. In parallel, the assay was performed manually using the same cellular models. Five replicates of each staurosporine concentration were tested and three plates were processed in parallel. The data obtained confirmed the possibility to automate the Apo-ONE Homogeneous Caspase-3/7 Assay on the epMotion 5075t. As expected, apoptosis is induced by staurosporine in both cell lines. Automated assays provided caspase 3/7 activity profiles comparable to those achieved when assays were manually performed (Figures 2 and 3).

A



B





Figure 2: Dose-dependent apoptosis induction of HeLa cells. Cell-based assay performed manually (A) or on epMotion 5075t (B)

A



B




Figure 3: Dose-dependent apoptosis induction of Jurkat cells. Cell-based assay performed manually (A) or on epMotion 5075t (B)

The Z'-factor, a statistical value comparing the assay's dynamic range to its data variation, was used to estimate and compare the assay quality (Table 1). Z'-factor is commonly used as an indicator of assay robustness [7]. A value above 0.5 is the sign of an excellent assay quality. All assays had a Z'-factor value equal or higher than 0.75 indicating that automated assays are as robust as assays performed manually

Table 1: Z'-Factor for manual and automated Apo-ONE Homogeneous Caspase-3/7 assays

Z'-Factor Analysis		epMotion 5075t	Manual
HeLa cells	Plate 1	0.82	0.75
	Plate 2	0.8	0.89
	Plate 3	0.9	0.81
Jurkat cells	Plate 1	0.86	0.85
	Plate 2	0.88	0.81
	Plate 3	0.8	0.77

Assay automation does not only guarantee robustness, performing cell-based assays with the epMotion 5075t also permits to reduce assay variability. Comparison of intra- and inter-plate precision calculated for manual and automated assay clearly proves that using the epMotion results in a significantly higher precision (Table 2). By reducing human intervention, automation of cell-based assays allows to eliminate human error as one of the major sources of variability. Errors can also be efficiently reduced especially

when a broad concentration range of apoptotic compound has to be prepared.

Through the flexible and intuitive epBlue™ software, this error prone protocol step can completely be handled by the workstation ensuring reproducible dilutions and correct dispensing. Depending on the cellular models used for this study, the cell-based assay variability can be decreased by 30 % to 73 % when automated.


**Table 2:** Global intra- and inter-plate coefficient of variation (CV) calculated for manual and automated Apo-ONE Homogeneous Caspase-3/7 assays

		epMotion 5075t		Manual	
		Global intra-plate CV	Global inter-plate CV	Global intra-plate CV	Global inter-plate CV
HeLa cells	Plate 1	6.11 %		11.90 %	
	Plate 2	7.15 %		7.69 %	
	Plate 3	5.80 %		15.84 %	
	Plates 1-2-3		5.53 %		19.87 %
Jurkat cells	Plate 1	9.43 %		12.71 %	
	Plate 2	6.15 %		12.71 %	
	Plate 3	7.83 %		13.50 %	
	Plates 1-2-3		8.47 %		11.81 %

### Automation Benefit

Cell-based assays are often used to screen chemical compound libraries. To rapidly identify active molecules, scientists want to perform large series of assays reliably and quickly. Assay variability is one of the parameters evaluated to determine reliability. A simple approach for estimating assay variability is to perform a large number of replicates (at least 20) in a single run on a single day. Automated Apo-ONE Homogeneous Caspase-3/7 assay variability was

determined by seeding 30,000 HeLa cells per well and treating those cells with two staurosporine concentrations. The first concentration belonged to the exponential phase (1  $\mu$ M) and the second was at the beginning of the plateau phase (10  $\mu$ M) as determined in previous experiments. For each concentration, 22 replicates were tested and three plates were included (Figure 4).



**Figure 4:** Caspase 3/7 activity induced by 1  $\mu$ M (A) or 10  $\mu$ M (B) of staurosporine (STS), determined in three plates with 22 replicates each.

**Table 3:** Apo-ONE Homogeneous Caspase-3/7 assay variability evaluation on the epMotion 5075t

	Intra-plate 1 CV	Intra-plate 2 CV	Intra-plate 3 CV	Mean intra-plate CV	Inter-plates CV
1 $\mu$ M STS	6.7 %	7.5 %	10.0 %	8.0 %	8.5 %
10 $\mu$ M STS	4.3 %	3.5 %	3.7 %	3.8 %	1.5 %

The variability when the Apo-ONE Homogeneous Caspase-3/7 assay is implemented on the epMotion 5075t is excellent. For the lowest staurosporine concentration tested (1  $\mu$ M), the mean variability observed within a plate is 8.0 %, while the inter-plate precision is 8.5 %. When a higher staurosporine concentration (10  $\mu$ M) is used, variability is reduced to maximum 3.8 % and 1.5 % respectively for the intra- and inter-plate precision (Table 3). With all Z'-factors

being above 0.5, the assay robustness is also guaranteed when a large number of replicates is handled within an assay run (data not shown). Besides the high precision offered to researchers by automation, assay implementation on a platform as the epMotion 5075t also allows to reduce time spent for manipulation. As a result, 75 % of the time needed to perform the assay does not require the presence of an operator.

## Conclusion

In the present Application Note, we demonstrate the possibility to automate the Apo-ONE Homogeneous Caspase-3/7 assay on the Eppendorf epMotion 5075t. Because this workstation can be equipped with the CleanCap configuration, a clean environment can be maintained and the complete cell-based assay protocol can be carried out on the automation platform. Assay transferability has been evaluated with two cell lines and all assays performed had a Z'-factor value equal or higher than 0.75, indicating an excellent robustness.

Comparison to the manual method showed that automation with the epMotion 5075t significantly increases the assay reproducibility. Depending on the cellular models used for this study, the cell-based assay variability can be decreased by 30 % to 73 % when automated. Excellent intra-plate and inter-plate precision is preserved when the number of samples handled in parallel increases. This clearly shows that automation of cell-based assays on the epMotion 5075t is an excellent solution for scientists interested in a low to a medium-throughput screening.

## Literature

- [1] Zang R, Li D, Tang I-C, Wang J, Yang S-T. Cell-based Assays in High-Throughput Screening for Drug Discovery. *International journal of Biotechnology for Industries* 2012; 1:31–51.
- [2] Elmore S,. Apoptosis: A Review of Programmed Cell Death. *Toxicologic Pathology* 2007; 35:495–519.
- [3] Apo-ONE® Homogeneous Caspase-3/7 Assay *Technical Bulletin* #TB295, Promega Corporation, 2009
- [4] An W.F, Tolliday N. Cell-based Assays for High-Throughput Screening. *Mol Biotechnol* 2010; 45:180–186.
- [5] Worzella T, Larson B. Automating Promega Cell-Based Assays in Multiwell Formats. *Promega Notes*, number 85, 2005.
- [6] Bracht D, Baranowski S. Apoptosis assay with the Promega Apo-ONE® Homogeneous Caspase-3/7 Assay on the epMotion® 5075 TMX. Eppendorf Application Note 209, [www.eppendorf.com](http://www.eppendorf.com)
- [7] Zhang J, Chung T, Oldenburg K. A Simple Statistical Parameter for Use in Evaluation and Validation of High Throughput Screening Assays. *J. Biomol. Screening* 2009; 4:67–73.

### Ordering information

Description	Order no. international	Order no. North America
epMotion® 5075t	5075 000 042	5075000042
epMotion® 5075t CleanCap upgrade set	5075 001 888	5075001888
Gripper	5282 000 018	960002270
TS50 pipetting tool	5280 000 010	960001010
TM50-8 pipetting tool	5280 000 215	960001044
TM1000-8 pipetting tool	5280 000 258	960001061
ReservoirRack for epMotion®	5075 754 002	960002148
CellXpert® C170i (Cell Culture Incubator)	<a href="#">see online</a>	<a href="#">see online</a>
epT.I.P.S.® Motion 50 µL Filter	0030 015 215	0030015215
epT.I.P.S.® Motion 1000 µL Filter	0030 015 258	0030015258
epMotion® Reservoir	0030 126 505	960051009
30 mL	0030 126 513	960051017
100 mL		
Reservoir Rack Module TC, for use in epMotion® Reservoir Racks, temperable, 4 x Safe-Lock tubes 0.5/1.5/2.0 mL	5075 799 081	960002620
Eppendorf Safe-Lock Tubes, 1.5 mL	0030 120 086	022363204
Eppendorf Deepwell Plate 96/1000 µL	0030 502 205	951032701
Waste bag	5075 752 034	5075752034

\*Developed on a predecessor model, but thanks to the migration feature, this method can easily be transferred to the newest generation of epMotion®.

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf SE · 22331 Hamburg · Germany

[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

Sigma-Aldrich® is a registered trademark of Sigma-Aldrich Co., LLC Lim., USA. Promega® and Apo-ONE® are registered trademarks of Promega Corporation, USA. Umonium38 ® is a registered trademark of Etienne Huckert, Belgium. Eppendorf®, the Eppendorf brand design, epMotion®, epT.I.P.S.® and Eppendorf ThermoMixer® are registered trademarks of Eppendorf SE, Germany. U.S. Design Patents are listed on [www.eppendorf.com/ip](http://www.eppendorf.com/ip). All rights reserved, including graphics and images. Copyright © 2024 by Eppendorf SE, Germany.

Methods are intended for molecular research applications. They are not intended, verified or validated for use in the diagnosis of disease or other human health conditions.