Keep this rotor instruction manual and the centrifuge manual in the file.

## Swing Rotor for Refrigerated Centrifuge

# R5S4

Before using these rotors, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

### Contents

1. Specifications	_____	1
2. Microplates	_____	3
3. How to use	_____	4
4. Maintenance	_____	7
5. Packing list	_____	9
6. Spare parts	_____	9
7. Decontamination	_____	10
8. Rotor retirement	_____	10
(Appendix) K Factors	_____	11

## Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.  
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.  
These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

### **WARNING**


- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor body).  
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density, kinds of microplates, etc.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

### **CAUTION**

- Balance the microplates and the samples including the MP adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.  
Contact an authorized sales or service representative.

# 1. Specifications

## 1.1 Construction

## 1.2 Specifications

### Applicable centrifuge\*1

(Some centrifuges are on sale in only Japan.  
Visit to our web site.)

Max. speed: 4,700 rpm

Max. RCF: 3,100 × g

Allowable mean density of sample: 1.2 g/mL

Max. dimensions of microplate:  
128(W)X86(D)X90(H) mm

Number of plates:  
Number of microplate (height: 1.4 cm) : 6 pcs./bucket  
Number of deep well plate (height: 4.5 cm) : 2 pcs./bucket

Rotor Size:  $\phi 422 \times 208$  mm

Weight: 5.2 kg (including buckets)


Material of rotor: Aluminum alloy

Material of bucket: Aluminum alloy

Rotor No.: 78

\*1: For applicable centrifuges, refer to the "Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))".

## 1.3 Cross-sectional view of rotor

\*: Rmin. is 4.50 cm when you use two deep well plates.

Figure 1-2 Cross-sectional view of R5S4 rotor

## 1.4 Characteristic of rotor

Table 1-1 Characteristics

Rotor speed (rpm)	RCF (× g)		K factor ※
	R = Rmin. (4.5 cm)	R = Rmax. (12.55 cm)	
1,000	50	140	259,798
1,500	113	316	115,466
2,000	201	561	64,950
2,500	314	877	41,568
3,000	453	1,263	28,866
3,500	616	1,719	21,208
4,000	805	2,245	16,237
4,700	1,111	3,100	11,761

※For details about K factor, refer to "(Appendix) K Factors"

### Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: speed (rpm)

### 1.5 Allowable speed of the rotor

 **WARNING :**  
 Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic, kinds of microplates, and centrifuge model.  
 Do not exceed the allowable speed of the rotor.

The speed must be limited depending on the total weight of contents (MP adapter, microplate and sample) in the bucket. When the total weight is over 825g, calculate the allowable speed according to the following equation.

(1) Allowable speed depending on the total weight of contents in the bucket.

 **WARNING :**  
 the total weight of contents : MP adapter, microplate, and sample

$$\text{Allowable speed (rpm)} = 4,700 \text{ (rpm)} \times \sqrt{\frac{825 \text{ (g)}}{\text{Total weight of contents in the bucket (g)}}}$$

### 1.6 Rotor registration data for the CR-G II series centrifuges

If this rotor is not registered in the CR-G II series centrifuges, register this rotor by entering the following data according to “3-5-5 Entry of New Rotor” on the instruction manual of the CR-G II series centrifuges.

Item/No.	DATA	Item/No.	DATA	Item/No.	DATA	Item/No.	DATA
ROTOR No.	78	6	1194	12	83052	18	0
1	24000	7	3132	13	32032	19	0
2	18096	8	100	14	32032	20	0
3	35142	9	18000	15	18092	21	5000
4	4231	10	11000	16	2154	22	6255
5	3085	11	82053	17	0	23	128064

## 2. Microplates

### 2.1 Applicable microplates

Microplates are not included in this rotor. Please purchase applicable ones on the market.

#### CAUTION:

- Use the microplates at a temperature between 4°C and 25°C. Otherwise, microplates may be broken or deformed during operation.
- Before using microplates on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- The maximum height of a microplate that can be held in a bucket is 90 mm.  
Do not use any microplate that is higher than 90 mm. Otherwise, the shield cap can not be mounted or improper swinging may occur. And it may result in damage to the rotor or the centrifuge.
- If it is necessary to set two or more microplates in one bucket, test them by filling with water instead of sample and run them at the intended speed to ensure that there is no break before actual operation.
- When setting two or more microplates in one bucket, use the microplates with rubber lids for sealing. Otherwise, the microplates may be damaged during operation.
- Use the microplates whose well-bottom can touch the MP adapter when the microplate is set into the MP adapter. Otherwise, the microplates may be damaged during operation.

Table2-1 Applicable microplates

Part No.	Name	Material * 2	Size (mm)	Actual Capacity	Max. Speed (rpm)	Max. RCF (×g)
Marketed	Microplate	PP	128×86×14	—	4,700	3,100
Marketed	Deep well plate	PP	128×86×43	—	4,700	3,100
Marketed	Deep well filter plate	PP	128×86×74	—	4,700	3,100

\* 1: When using microplates on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer. Otherwise the microplates may be broken during operation.

\* 2: PP: Polypropylene

The rotor speed is limited by the factors (the total weight of contents in the bucket etc.) (see Section1.5).

- Maximum speeds listed are guidelines only. Because of variances in user methodologies, no guarantee of performance is expressed or implied.

### 2.2 Cleaning and sterilization of microplates

Clean and sterilize microplates on the market according to the instructions of the manufacturer.

### 2.3 Life expectancy of microplates

In spite of the sterilization condition, the characteristics of the samples, speed of the rotor used, temperature, and the run time, the life expectancy of the microplates on the market is 1 operation.

## 3. How to use

---

### 3.1 Preparation of the rotor

- (1) Check that there is no crack, corrosion or deformation on the rotor and the buckets.
- Check that the drive pin of the drive hole (crown hole) is not deformed.
  - Check that there is no corrosion, foreign matters, dirt, and scratches in the inside bottom of the bucket and in the drive hole (crown hole).

**CAUTION:**

If there is any abnormality such as corrosion or cracks, stop use of the rotor and the buckets and contact an authorized sales or service representative.

- (2) Check that there is no foreign matter in the buckets.

### 3.2 Preparation of the sample and the microplates

Inject the sample into the microplates.

**WARNING :**

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

**CAUTION:**

- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor (including buckets), MP adapters, shield cap (including packing and sliders), and microplates. Using such a sample could corrode or deteriorate them.
- Chemical resistance and the strength of microplates varies with speed, temperature, and so on. Before using sample, fill microplates with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise the sample leakage may occur during operation.
- Shield caps (including packing and sliders) are consumables. Replace them with new ones if damage or deterioration is found. It is recommended to periodically replace them with new ones once a year.
- Use the microplates that are same products. Otherwise, the microplates may be damaged due to their different form or the imbalance operation may occur, and it may result in the damage to the centrifuge and the rotor.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, buckets, microplates, and samples are warmed or cooled to the desired temperature before operation.

### 3.3 Balancing the microplates

The allowable imbalance of this rotor is 5g (Allowable imbalance: Difference in weight of contents (MP adapters, microplates, and samples) in the symmetrically arranged buckets).

**CAUTION:**

- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise, the precipitation levels may be different by centrifugation and there may be variations in position of center of gravity, and it may result in damage to the centrifuge and the rotor due to the imbalance operation.
- Balance between symmetrically arranged microplates filled with the same-density sample. Variations in sample density cause imbalance operation, and it may result in the damage to the centrifuge and the rotor.


Use of shield caps is recommended to prevent contamination of the sample with dust or water.

Use of microplates with films for sealing (on the market) (or rubber lids (on the market) ) is recommended to prevent contaminating the sample of the well with the sample of another well.

### 3.4 Loading the microplates into the buckets


Place the balanced microplates into the MP adapters, and load them into the buckets.

**CAUTION:**  
When setting two or more microplates in one bucket, check that the microplates are piled with their corresponding wells aligned.

### 3.5 Mounting the shield cap

Mount the shield cap to the bucket, and then align the mark (triangle) of the slider with the mark (triangle) of the shield cap by sliding the slider to the “LOCK” position from the right side (the state of “UNLOCK”) to the left side (see the below). (There are four “LOCK” positions .)

**CAUTION:**

- While the bucket is set to the rotor, do not slide the slider. Otherwise it may result in the damage to the drive shaft.
- Use the microplates with the attached MP adapters (material: stainless steel). Otherwise the microplates may be deformed or cracked during operation.
- Do not hold the attached MP adapters (material: stainless steel) strongly. Otherwise the MP adapters may be deformed and the microplates can not be placed into them.
- Be sure to set all the shield caps to the buckets. Operation with only one shield cap is very dangerous because it may cause a severe imbalance state.
- Align the marks (triangle) of all four sliders with the marks (triangle) of the shield caps by sliding the sliders to all the “LOCK” positions. (There are four “LOCK” positions.)  
Otherwise the shield caps will disengage from the buckets while the rotor is running. And the rotor and the centrifuge may be damaged.


### 3.6 Mounting the rotor and the buckets

Mount the balanced buckets with respect to the center of the rotor, and then mount the rotor to the drive shaft.

#### WARNING:

- Be sure to set all the buckets of the same type. Operation with only one bucket is very dangerous because it may cause a severe imbalance state. Never use any bucket or adapter that is not designated for this rotor by us.
- Be sure to set all the buckets of the same serial number (The serial number is shown on the side of the bucket). Otherwise, it may cause the imbalance operation and it may result in the damage to the centrifuge and the rotor.


- Securely hook the buckets on the rotor pins. Otherwise the buckets may come off while the rotor is running.

#### CAUTION:

- Set the microplates so that the samples in the microplates are arranged symmetrically with respect to the center of rotation. Otherwise, imbalance operation may be caused and the centrifuge may be damaged.

Check that the buckets swing smoothly. If the swinging is not smooth, clean both the rotor and the pin holes of the buckets.

### 3.7 Operation

Run the rotor according to the procedures specified in the centrifuge instruction manual.

### 3.8 Taking out samples

- (1) After the operation, carefully remove the rotor from the centrifuge.
- (2) Remove the shield caps from the buckets by sliding the slider and take out MP adapters and microplates. Then take out the separated sample. If it is not easy to remove the shield cap from the bucket, deform the shield cap by pressing the slider horizontally while keeping the slider placed the state of "UNLOCK".

## 4. Maintenance

---

### 4.1 Corrosion resistance of the rotor, buckets, shield caps, and MP adapters

**WARNING:**

Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor, buckets, shield caps, and MP adapters.

The material of the shield caps is polypropylene. The material of the MP adapters is stainless steel and chloroprene rubber.

The rotor and the buckets are made of aluminum alloy. They have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode them and decrease their strength. Use samples that will not affect the rotor referring to the attached chemical resistance chart.

### 4.2 Maintenance of the rotor and the buckets

After the operation, properly maintain the rotor and the buckets to prevent corrosion that can cause them breakage.

- Normal maintenance

**CAUTION:**

Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor, buckets, shield caps, and MP adapters with tap water or a dilute solution of neutral detergent and rinse them out with distilled water. Wipe them with a soft cloth. Dry the rotor which the buckets, shield caps, and MP adapters are removed from. Check that the rotor is completely dried and then put a light coat of silicone grease (vacuum grease)(standard accessory of the centrifuge) on the rotor. Store the rotor in a dry place that is not exposed to the direct sunlight.

- Maintenance after use of a corrosive sample

Wash the rotor, buckets, shield caps, and MP adapters with tap water immediately after the operation. Then perform the normal maintenance.

- Maintenance when foreign substances are adhered to the rotor and buckets

**CAUTION:**

Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and buckets can be discolored or corroded.

Soak the rotor and the buckets in warm water for one or two hours and wash the inside of the buckets and the crown hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the crown hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

### 4.3 Sterilizing

**WARNING:**

Never sterilize rotor, buckets, shield caps, MP adapters, packing, and sliders by boiling, otherwise they will be brittle.

(1) Sterilizing the rotor

Sterilize the rotor in any of the following methods; gas sterilization (ethylene oxide or formaldehyde), chemicals sterilization (ethanol (70%), hydrogen peroxide (3%) or formalin (3%)), or ultraviolet rays sterilization(200-300nm).

**CAUTION:**

- Never sterilize the rotor by autoclaving (121°C, 20 minutes). Otherwise the strength of the rotor may be significantly decreased.
- Do not dip the rotor in the formalin (3%) solution more than 2 hours.

(2) Sterilizing the buckets

Sterilize the buckets in any of the following methods; autoclaving (121°C, 20 minutes), gas sterilization (ethylene oxide or formaldehyde), chemicals sterilization (ethanol (70%), hydrogen peroxide (3%) or formalin (3%)), or ultraviolet rays sterilization(200-300nm). Turn the buckets upside down to dry well.

**CAUTION:**

- Do not dip the rotor in the formalin (3%) solution more than 2 hours.
- Sterilize the buckets after removing the shield caps from the buckets.
- Sterilize the buckets after removing the shield caps from the buckets. Take out the buckets carefully from the autoclaving chamber as they are wet and slippery.

(3) Sterilizing the shield caps (including the packing and the sliders)

Sterilize the shield caps in any of the following methods; autoclaving (121°C, 20 minutes), gas sterilization (ethylene oxide or formaldehyde), or chemicals sterilization (hydrogen peroxide (3%) or formalin (3%)).

**CAUTION:**

- Never sterilize the shield caps by ultraviolet radiation. Otherwise they may be significantly deteriorated and they may be damaged.
- Do not use the shield caps more than ten times when sterilizing them by autoclaving (121°C, 20 minutes).
- Inspect the packing of the shield cap. If the packing is deteriorated, it may cause the damaged packing and bucket seal can not work properly when the sample leaks from the buckets.

(4) Sterilizing the MP adapters

Sterilize the MP adapters in any of the following method; chemicals sterilization (ethanol (70%) or hydrogen peroxide (3%)). MP adapters can not be sterilized by autoclaving due to their specifications.

### 4.4 Inspecting rotor

Periodically check every 100 hours or once a month whether the rotor and buckets are not corroded. Especially be careful of the inside bottom of the bucket and the crown hole of the rotor bottom, because the rotor and the buckets can be brittle if these portions are corroded. If their surfaces are discolored, dented or cracked, they are corroded. Do not use such a corroded rotor or a bucket. Contact an authorized sales or service representative immediately for further inspection.

**CAUTION:**


If there is any abnormality such as corrosion or cracks, stop use of the rotor and the buckets and contact an authorized sales or service representative.

## 5. Packing list

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	Remark
①	Rotor	1	—
②	Bucket	2	—
③	MP adapter	2	—
④	Shield cap	2*	—
⑤	Packing (for the shield cap)	2*	—
⑥	Rotor instruction manual	1	Part No.S999924

\*: The quantities of the parts that are mounted to the rotor.

## 6. Spare parts

Figure	PART No.	Name	Qty
	S411670A	Packing (for the shield cap)	2
	S308646A	Shield Cap Ass'y (This includes the packing [for the shield cap].)	2
	483719	Silicone grease (Vacuum grease)	1
	S999924	Rotor instruction manual	1

## 7. Decontamination

---

### WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

## 8. Rotor retirement

---

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

## (Appendix) K Factors

---

### • Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit (S).

Namely, the time taken for particles of 1 S to settle from Rmin. to Rmax. is K factor.

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max.}) - \ln(R_{\min.})}{\omega^2} \cdot \frac{10^{13}}{3,600}$$

$\omega$ : Angular velocity  $\omega = 2\pi/60 \times N$

N: Speed (rpm)

Rmin.: Distance from axis of rotation to top of solution (cm)

Rmax.: Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t (in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 4,700 rpm of Model R5S4 rotor is 11,761. The time taken for particles of 1,000 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{11,761}{1,000} = 11.8 \text{ hours}$$

The time taken for settling of particles of 1,000 S by the same rotor at a speed of 3,000 rpm is:

$$K = \left[ \frac{4,700}{3,000} \right]^2 \times 11,761 = 28,867 \quad t = \frac{28,867}{1,000} = 28.9 \text{ hours}$$

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Copyright © 2022 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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