Keep this rotor instruction manual and the centrifuge manual in the file.

## Angle Rotor for Refrigerated Centrifuge

# R9A

(Serial No. 4001 or later)

Before using this rotor, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.  
Keep this instruction manual for your reference and refer to it as required.

### Contents

1.	Specifications	_____	1
2.	Bottles	_____	3
3.	How to Use the Rotor	_____	4
4.	How to Use the Bottle Vise	_____	7
5.	Maintenance	_____	8
6.	Packing List	_____	10
7.	Spare Parts of the Rotor	_____	10
8.	Decontamination	_____	11
9.	Rotor Retirement	_____	11
	(Appendix) K Factors	_____	12

## Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.  
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.  
These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

### **WARNING**

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the surface of the cover-locking knob).  
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density or kinds of tubes/bottles.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

### **CAUTION**

- Do not run the rotor without the rotor cover in position.  
Tighten the cover-locking knob securely.
- Balance the tubes/bottles and the samples including the caps and the adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the rotor tubes and bottles within their actual capacities.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge. The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature.  
Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/ bottles if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.  
Contact an authorized sales or service representative.

# 1. Specifications

## 1.1 Construction


Fig.1-1 Construction of the R9A Rotor

## 1.2 Specifications

Applicable centrifuge\*1

(Some centrifuges are on sale in only Japan.  
Visit our web site.)

Max. Speed\*2: 9,000 rpm

Max. RCF: 15,300 × g

Nominal Rotor Capacity:

1,000 mL × 4 bottles = 4,000 mL

Rotor Size: φ 356 × 207 mm

Weight: 13.7 kg

Material: Aluminum alloy

Rotor No.: 76

\*1: For applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges) (Part No. S998611)”.

\*2: The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges) (Part No. S998611)”.

## 1.3 Cross-Sectional View of the Rotor


Fig.1-2 Cross-Sectional View of the R9A Rotor

## 1.4 Characteristics of the Rotor

Table1-1 Characteristics

Rotor speed (rpm)	RCF (×g)		K factor*
	Rmin. (6.8cm)	Rmax. (16.9cm)	
1,000	76	189	230,605
2,000	304	756	57,651
3,000	684	1,700	25,623
4,000	1,220	3,020	14,413
5,000	1,900	4,720	9,224
6,000	2,740	6,800	6,406
7,000	3,730	9,260	4,706
8,000	4,870	12,100	3,603
8,500	5,490	13,700	3,192
9,000	6,160	15,300	2,847

\*For details about K factor, refer to “(Appendix) K Factors”

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)  
N: speed (rpm)

### 1.5 Allowable speed of the rotor

**⚠ WARNING :**  
 Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic and kinds of tubes/bottles.  
 Do not exceed the allowable speed of the rotor.

The maximum speed marked on the surface of the cover-locking knob is permitted under the conditions as follows.

-Average sample density : less than 1.2g/mL

Depending on sample density, kinds of bottles and centrifuge model, the rotor speed may be limited as follows. Do not exceed the allowable speed of the rotor.

#### (1) Allowable speed for sample density

**⚠ WARNING :**  
 When sample density exceed 1.2g/mL, calculate the allowable speed according to the following equation.

$$\text{Allowable speed (rpm)} = 9,000(\text{rpm}) \times \sqrt{\frac{1.2(\text{g/mL})}{\text{Average density of sample}(\text{g/mL})}}$$

### 1.6 Rotor registration data for the CR-G II series centrifuges

If this rotor is not registered in the CR-G II series centrifuges, register this rotor by entering the following data according to “3-5-5 Entry of New Rotor” on the instruction manual of the CR-G II series centrifuges.

Item/No.	DATA	Item/No.	DATA	Item/No.	DATA	Item/No.	DATA
ROTOR No.	76	12	65032	30	76	42	238
1	76000	13	32032	31	7208	43	0
2	35044	14	32032	32	0	44	0
3	67142	15	35040	33	30	45	0
4	6157	16	9057	34	18152	46	0
5	6137	17	0	35	4044	47	0
6	2168	18	0	36	62	48	0
7	5220	19	3186	37	30130	49	35040
8	29	20	0	38	6075	50	7208
9	16000	21	7002	39	102	51	10
10	1001	22	10255	40	255255	52	132225
11	82057	23	16125	41	10193	—	—

## 2. Bottles

### 2.1 Applicable bottles

Table 2-1 Applicable bottles

Part No.	Name	Size (mm)	Actual capacity	Max. speed	Max. RCF
S307700A	1000 PP bottle (WM)	φ98 × 170	900mL	9,000rpm	15,300 × g
S308357A	1000 PC bottle (WM)	φ98 × 170	900mL	9,000rpm	15,300 × g

Refer to the attached Chemical resistance chart (Part No. S999313) for chemical resistance of each material.

### 2.2 Cleaning and sterilization of bottles and caps of bottles

Clean and sterilize bottles according to Table 2-2.

#### ● Autoclaving PP and PC bottles

Wash the bottles well. Put them on a heat-resistant bottle stand facing their mouths upward. Remove the cap of the bottle to prevent deformation or breakage during autoclaving. After autoclaving, wait until the temperature in the autoclaving chamber reduces to the room temperature, then take out the bottles.

**⚠ CAUTION:**  
Do not place the cap on the mouth of the bottle. Otherwise the bottle may be deformed or damaged during autoclaving.

Table 2-2 Washing and sterilization of bottles

○: Applicable    x: Inapplicable

Conditions		Material	PC	PP	
Washing	Washing solution	pH 5 or lower	x	x	
		pH 5 to 9	○	○	
		pH 9 or higher	x	○	
		Warm water (50°C or lower)	○	○	
	Ultrasonic washing	Neutral detergent (pH 7)	○	○	
Sterilization	High-pressure steam sterilization(Autoclaving)	115°C(0.7kg/cm <sup>2</sup> ), 30 minutes	○	○	
		121°C (1.0kg/cm <sup>2</sup> ), 20 minutes	○	○	
		126°C (1.4kg/cm <sup>2</sup> ), 15 minutes	x	x	
		Boiling	For 15 to 30 minutes	○	○
		Ultraviolet ray	200-300nm	x	x
		Gas	Ethylene oxide	x	○
	Formaldehyde		○	○	
		Chemical	Ethanol (70%)	x	○
	Hydrogen peroxide (3%)		○	○	
	Formalin (3%)		○	○	

### 2.3 Life expectancy of plastic bottle and cap of bottle

The life of plastic bottle and cap of bottle depends upon the sample, speed of rotor, temperature, etc. The standard life of plastic bottle when used with aqueous samples (of pH5 to pH9) is specified as follows. (in case they are used repeatedly for 1 hour at the maximum speed of rotor):

- Bottle (PP,PC) ..... 50 times

However, the standard life of the PC bottle is specified as follows depending upon the pre-treating conditions such as washing, sterilization, etc.:

Table 2-3 Standard Life of tube and the PC bottle

Sample	Washing and sterilization	Sterilization in gas or washing in hot water	Autoclaving at 121°C x 20min
Neutral(pH7)		50 times	10 times
Weak alkali(pH7 to pH9)		30 times	5 times

The standard life of plastic bottles as the above are an approximate guide when they are used repeatedly for 1 hour at the maximum speed of rotor. We do not warrant the life expectancies of bottles.

#### ● Washing :

Wash the bottles with tap water or a dilute solution of neutral detergent. If the bottles stain heavily, soak them in a dilute solution of neutral detergent or wash them with a soft brush being careful not to scratch them. Then rinse the tubes/bottles in tap water and finally rinse them with distilled water. Dry the bottles in the air. Use a detergent of pH 5 to 9 because the PC bottles are susceptible to alkaline solutions. Avoid soaking the bottles in a dilute solution of neutral detergent for a long time. Otherwise, the bottles may be degraded in strength. For the caps of the bottles, remove them from the bottles. Then disassemble and wash them. And dry and store the caps in the same manner as the bottles.

## 3. How to Use the Rotor

### 3.1 Preparation of bottles

Check that the bottles are free from crack or deformation. Do not use remarkably deformed or cracked bottles.

**⚠ WARNING :**

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

**⚠ CAUTION:**

- Do not use any sample inapplicable to the bottles and the caps in use referring to the attached Chemical resistance chart. Otherwise, the bottles and the caps may be degraded.
- Chemical resistance and the strength of bottles varies with speed, temperature, and so on. Before using sample, fill bottles with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the specified sample capacity of the bottles, otherwise sample leakage or deformation of bottles may occur.
- Balance the sample within the allowable imbalance including the weight of the bottles and caps. Excessive imbalance can cause damage to the centrifuge and the rotor.
- Balance between symmetrically arranged bottles filled with the same-density sample. Variations in sample density cause imbalance operation.
- Run the rotor at the proper temperature (from 4 °C to 25°C). Otherwise the bottles may be damaged or deformed during operation. When running the rotor beyond this limit (under 4°C or above 25°C), test the bottles to ensure that there is no abnormality before actual operation.
- Tighten the cap of the bottle securely, otherwise the sample leakage during operation can cause the excessive imbalance, and it may result in damage to the centrifuge and the rotor.

(1) Inject the sample into the bottles and tighten the cap of the bottle securely. Do not exceed the specified sample capacity of the bottles, otherwise sample leakage may occur during operation.

If you can not tighten the cap properly because of using laboratory gloves etc., we recommend that you use the bottle vise (option) (see "4. How to use the bottle vise").

(2) Balance the bottles that are symmetrically arranged. The allowable imbalance of this rotor is within 6g/bottles.

### 3.2 Preparation of the rotor

(1) Inspection of the rotor and the tube holder

- Check that the bottom of the tube holder and the drive hole of the rotor are not cracked and are not stained. And check that there is no foreign matter in the bottom of the tube holder and the drive hole of the rotor.
- Check that the drive pins of the drive hole are not deformed.
- Check that the screw portion of the rotor, the cover-locking knob and the rotor are not worn.

(2) Apply the lubricant for screw (Part No.84810601) to the screw portion of the cover-locking knob and the rotor slightly.


### 3.3 Setting of bottles into rotor

**⚠ CAUTION:**


- Do not load only the one bottle or three bottles to this rotor. Otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- Before using the tube holder, be sure to wipe off the water and the sample which remained into the tube holder. Otherwise it may be hard to remove the bottle from the tube holder.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, bottles, and sample are warmed or cooled to the desired temperature before operation.

Check that the tube cavities of the tube holder are free from foreign substances. Set the balanced bottles symmetrically with respect to the drive shaft as shown in the figure below. Load the two or four bottles.

Center of the drive shaft

Setting of two bottles

Setting of four bottles

Fig. 3-1 Setting of bottles


### 3.4 Setting the cap to the tube holder

#### CAUTION:

- Use the caps, whose color is silver, and the rotor, whose serial number is 4001 or later (serial number of the rotor is mentioned on the upper surface of the rotor), otherwise the caps or packing may come off during operation and it may cause the centrifuge and the rotor to be damaged.
- Be careful not to get your fingers caught between the caps when setting the cap.
- Do not set only the one cap or three caps to the tube holder. Otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- Be sure to set all caps and tighten them securely, otherwise samples may evaporate during operation.

- (1) Apply the lubricant for screw (Part No.84810601) to the screw portion of the caps.
- (2) Apply the silicone grease (vacuum grease) (the standard accessory of the centrifuge) to the packing and then fit it into the groove of the cap. Replace the packing with the new one if damaged or deteriorated.
- (3) Tighten the cap to the tube holder securely by the cap opener (Part No.S204715)(the standard accessory of this rotor). Tighten the cap to the tube holder until the caps touches the tube holder.

Fit the packing into the groove of the cap so that the protrusion of the packing faces downward

Fig. 3-2 Packing and cap

### 3.5 Mounting the rotor cover

#### CAUTION:

Be sure to mount the rotor cover to the rotor body and tighten the cover-locking knob securely before operation, otherwise the rotor cover or the rotor cover may be removed during operation and result in damage to the centrifuge and the rotor.

- (1) Apply the lubricant for screw (Part No.84810601) to the screw portion of the cover-locking knob.
- (2) Mount the rotor cover to the rotor body and tighten the cover-locking knob securely.

### 3.6 Operation

- (1) Mount the rotor onto the drive shaft of the centrifuge gently and securely. Check that the rotor cover is mounted and the cover-locking knob is securely tightened.
- (2) Perform operation according to the instructions of the centrifuge instruction manual.

### 3.7 Taking out samples

After the operation, gently take out the rotor from the centrifuge and remove the rotor cover by turning the cover locking knob. Take out the bottles/tubes from the rotor and then take out the sample in a proper manner.

#### CAUTION:

- Be careful not to get your fingers caught between the caps when removing the cap.
- Remove the bottle holding the rotor. If there is water or the sample in the tube holder, it may be hard to remove the bottle from the tube holder.

- (1) After the operation, gently take out the rotor from the centrifuge and remove the rotor cover by turning the cover-locking knob.
- (2) Remove the cap from the tube holder by the cap opener (Part No.S204715)(the standard accessory of this rotor), and then remove the bottles. Take out the sample in a proper manner. Mount small pin of the cap opener in the notch of the cap and then turn the cap opener counterclockwise (See Fig. 3-3).(When opening the cap of the bottle, use the large pins of the cap opener .)



Fig. 3-3 How to use the cap opener

## 4. How to Use the Bottle Vise

- The bottle vise is not included in this rotor. If necessary, please order the bottle vise (Part No. S205037A).

### 4.1 Mounting the bottle vise

- (1) Check the thickness of the desk. If the distance between the bottom of the bottle vise and the handle (B) is shorter than the thickness of the desk, turn the handle (B) counterclockwise to make the distance longer between the bottom of the bottle vise and the handle (B).
- (2) Position the bottle vise on the desk and turn the handle (B) clockwise to mount the bottle vise on the desk.

### 4.2 Tightening the bottle

#### CAUTION:

- Do not tighten the handle (A) too tightly. Failure to do so might cause the bottle to be deformed.

- (1) Set the bottle to the middle of the bottle vise.
- (2) If there is not enough space for the bottle, turn handle (A) counterclockwise to open the bottle vise and give the bottle more space. (Fig. 4-1)
- (3) Push the bottle against the fixing portion of the bottle vise and turn the handle (A) clockwise to move the fixing plate until the fixing plate touches the bottle lightly. (Fig. 4-2)
- (4) As soon as the fixing plate touches the bottle, turn another 90° to have the bottle fixed tightly in place. (Fig. 4-3)


Fig. 4-1


Fig. 4-2



Fig. 4-3

---

### 4.3 Tightening the cap

Tighten the cap with the provided cap opener.

Mount the two pins of the cap opener in the holes of the cap and then turn the cap opener clockwise.

## 5. Maintenance

---

### 5.1 Corrosion resistance of rotor and cap

**⚠ WARNING:**  
Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor and the caps.

This rotor body and caps is made of aluminum alloy. The rotors made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode the rotor body and caps and decrease the strength of them. Use samples that will not affect the rotor and caps referring to the attached chemical resistance chart.

---

## 5.2 Maintenance of rotor and cap

After the operation, properly maintain the rotor and cap to prevent corrosion that can cause the rotor and caps breakage.

- Normal maintenance

 **CAUTION:**  
Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor and caps with tap water or a dilute solution of neutral detergent and rinse them out with distilled water. Wipe the rotor and caps with a soft cloth. Turn the rotor upside down with the rotor cover and caps removed to dry it well. Check that the rotor and caps is completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor and caps. Store the rotor and caps in a dry place.

- Maintenance after use of a corrosive sample

Wash the rotor and caps with tap water immediately after the operation. Then perform the normal maintenance.

- Maintenance when foreign substances are adhered to the rotor and caps

 **CAUTION:**  
Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and caps can be discolored or corroded.


Soak the rotor and caps in warm water for one or two hours and wash the inside of the tube cavities and the drive hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the drive hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

## 5.3 Sterilizing rotor and cap

Sterilize this rotor and caps in any following methods; autoclaving(120 °C , 20minutes), gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).


 **WARNING:** Be sure to remove the rotor cover and the caps from the rotor before autoclaving, otherwise the rotor, the rotor cover or the caps may be deformed.

 **WARNING:** Never sterilize the rotor, the rotor cover and the cap by boiling, otherwise they will be brittle.

 **CAUTION :** After autoclaving, wait until the temperature in the autoclaving chamber reduces to the room temperature, then take out the rotor, the rotor cover, and the caps.

 **CAUTION :** Take out the rotor , the rotor cover and the caps carefully from the autoclaving chamber as they are wet and slippery.

 **CAUTION :** Inspect the packing. If the packing is deteriorated, it may cause the damaged packing and rotor seal can not work properly when the sample leaks from the bottle.

 **CAUTION :** Do not dip the rotor in the formalin (3%) solution more than 2 hours.

## 5.4 Inspecting rotor


Periodically check that the rotor is free from corrosion every 100 hours use. Especially be careful of the tube cavities of the rotor and the drive hole at the rotor bottom because the rotor can be brittle if these portions are corroded. If the rotor surface is discolored, dented or cracked, the rotor is corroded. Do not use such a corroded rotor and contact an authorized sales or service representative immediately for further inspection.

## 6. Packing List

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	Part No. of spare parts (Qty)
①	Rotor	1	—
②	Cap	4*	S414217A (2)
③	Packing	4*	S409741A (2)
④	Cap opener	1	S204715 (1)
⑤	1000 cap (WM) ass'y	4	S410870A (2)
⑥	1000 bottle (WM) ass'y	4	S307700A (6)
⑦	Bottle setter	1	84550301 (1)
⑧	Lubricant for screw	1	84810601 (1)
⑨	Rotor stand	1	S205816 (1) : For the rotor of which the red adapter is equipped at the bottom
⑩	Rotor instruction manual	1	S999828

\*: the quantities of the parts that are mounted to the rotor.

## 7. Spare Parts of the Rotor

Figure	Part No.	Name	Qty
	S414217A	1000AL cap (WM) ass'y2	2
	S409741A	Packing ass'y	2
	84810601	Lubricant for screw	1
	483719	Silicone grease (vacuum grease)	1
	S204715	Cap opener	1
	84550301	Bottle setter	1
	S205816 (*1)	Rotor stand	1
	S999828	Rotor Instruction manual	1

(\*1): For the rotor of which the red adapter is equipped at the bottom

## 8. Decontamination

---

 **WARNING:**

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

## 9. Rotor Retirement

---

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

## (Appendix) K Factors

---

- Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

Namely, the time taken for particles of 1S to settle from Rmin. to Rmax. is K factor.

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max.}) - \ln(R_{\min.})}{\omega^2} \cdot \frac{10^{13}}{3600}$$

$\omega$  : Angular velocity  $\omega = 2\pi/60 \times N$

N : Speed (rpm)

Rmin. : Distance from axis of rotation to top of solution (cm )

Rmax. : Distance from axis of rotation to bottom of tube (cm )

Using the K factor, the settling time t(in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 9,000rpm of Model R9A rotor is 2,847. The time taken for particles of 1000 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{2,847}{1,000} = 2.8 \text{ hours}$$

The time taken for settling of particles of 1000 S by the same rotor at a speed of 7,000 rpm is:

$$K = \left[ \frac{9,000}{7,000} \right]^2 \times 2,847 = 4,706 \quad t = \frac{4,706}{1,000} = 4.7 \text{ hours}$$

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Copyright © 2022 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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