Keep this rotor instruction manual and the centrifuge manual in the file.

Horizontal Rotor for Refrigerated Centrifuge

R10H

Before using these rotors, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Plates	_____	3
3. How to use	_____	5
4. Maintenance	_____	8
5. Packing list	_____	10
6. Spare parts	_____	10
7. Decontamination	_____	11
8. Rotor retirement	_____	11

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.

These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density, kinds of microplates, etc.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Balance the plates and the samples including the plate adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.
Contact an authorized sales or service representative.

1.Specifications

1.1 Construction


Figure1-1 Construction of R10H rotor

1.2 Specifications

Applicable centrifuge*1

(Some centrifuges are on sale in only Japan. Visit to our web site.)

Max. speed: 10,000 rpm

Max. RCF: 13,000 × g

Tube hole angle: 90°

Number of the microplate for one compartment of the rotor: 1 microplate (maximum 3 microplates 1.5cm in thickness)

Total number of the microplate of the rotor: 4 microplates (maximum 12 microplates 1.5cm in thickness)

Dimensions of the microplate: 12.8(W)X8.6(D)X1.5-4.5(H) cm

Number of the PCR plate for one compartment of the rotor: 1 PCR plate (The PCR tube rack is required.)

Total number of the PCR plate of the rotor: 4 PCR plates

Dimensions of the PCR plate: 12.8(W)X8.6(D)X2.3(H) cm

Allowable mean density of sample: 1.2 g/mL

Weight: 10.6 kg

Material of the rotor body: Aluminum alloy

Material of the rotor cover: Aluminum alloy

Rotor No.: 53

*1: For applicable centrifuges, refer to the "Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))".

⚠ CAUTION :

Only one microplate can be loaded in one compartment of the rotor when running the centrifuge at 2,700 rpm or higher. Use of two or more microplates in one compartment of the rotor can cause damage to the microplates and the rotor.

1.3 Characteristic of rotor

Rotor Speed (rpm)	RCF (xg)	K factor
	Rmax.	
2,000	516	7,527
4,000	2,080	1,882
6,000	4,670	836
7,000	6,350	614
8,000	8,300	470
10,000	13,000	301

Calculation formula

$$\text{RCF} = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: speed (rpm)

1.4 Allowable speed of the rotor

WARNING :

Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on the weight of contents of one compartment of the rotor, kinds of microplates, and centrifuge model. Do not exceed the allowable speed of the rotor.

The speed must be limited depending on the total weight of contents in one compartment of the rotor. When the total weight is over 340g, calculate the allowable speed according to the following equation.

(1) Allowable speed depending on the total weight of contents in one compartment of the rotor.

WARNING :

- Microplates

Total weight of contents in one compartment of the rotor: MP plate adapter, microplate, and sample

$$\text{Allowable speed (rpm)} = \text{Maximum speed of the rotor (rpm)} \times \sqrt{\frac{340 \text{ (g)}}{\text{Total weight of contents in one compartment of the rotor (g)}}}$$

- PCR plates

Total weight of contents in one compartment of the rotor: PCR plate adapter, PCR plate, PCR tube rack, and sample

$$\text{Allowable speed (rpm)} = \text{Maximum speed of the rotor (rpm)} \times \sqrt{\frac{340 \text{ (g)}}{\text{Total weight of contents in one compartment of the rotor (g)}}}$$

2. Plates

2.1 Applicable plates

Microplates and PCR plates are not included in this rotor. Please purchase applicable ones on the market.

CAUTION:






- Be sure to use the microplates together with a hermetic seal or mat.
- Use the microplates at a temperature between 4°C and 25°C. Otherwise, microplates may be broken or deformed during operation.
- Before using microplates on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- The maximum height of a microplate that can be loaded in one compartment of the rotor is 4.5 cm. Do not use any microplate that is higher than 4.5 cm.
- If it is necessary to set two or more microplates in one compartment of the rotor, test them by filling with water instead of sample and run them at the intended speed to ensure that there is no break before actual operation.
- Use the microplates whose well-bottom can touch the M plate adapter. Otherwise, the microplates may be damaged during operation.

CAUTION:

- Be sure to use the PCR plate together with a hermetic seal or a hermetic mat.
- Be sure to use the PCR tube rack when using the PCR plate. Only one PCR plates and one PCR tube rack can be loaded in one compartment of the rotor.
- Use the PCR plates at a temperature between 4°C and 25°C. Otherwise, PCR plates may be broken or deformed during operation.
- Before using PCR plates on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- The maximum total height of one PCR plate including one PCR tube rack that can be loaded in one compartment of the rotor is 4.5 cm. Do not exceed the total height that is higher than 4.5 cm.
- Use the PCR plates whose bottom can fit the PCR tube rack. Otherwise, PCR plates may be damaged during operation.

Table2-1 Applicable plates

Plates			Racks		Adapters			Max. Speed (rpm)	Max. RCF (xg)
Part No. (Qty)	Name	Size (cm)	Part No. (Qty)	Name	Part No. (Qty)	Name	Material		
Marketed	Microplate 	12.8X8.6X1.5-4.5 (W) X (D) X (H)	---	---	S410604A (4 pcs.)	M plate adapter 	* 2	10,000	13,000
Marketed	PCR plate 	12.6 X 8.6 X 2.3 (W) X (D) X (H)	Marketed	PCR tube rack 	S410605A (4 pcs.)	PCR plate adapter 	* 3	10,000	13,000

* 1: When using microplates and PCR plates on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer. Otherwise they may be broken during operation.

* 2: ABS

* 3: Polypropylene

The rotor speed is limited by the factors (the total weight of contents in the microplate or the PCR plate etc.) (see Section1.5).

- Maximum speeds listed are guidelines only. Because of variances in user methodologies, no guarantee of performance is expressed or implied.

2.2 Cleaning and sterilization of microplates and PCR plates

Clean and sterilize microplates and PCR plates on the market according to the instructions of the manufacturer.

2.3 Life expectancy of plates

In spite of the sterilization condition, the characteristics of the samples, speed of the rotor used, temperature, and the run time, the life expectancy of the microplates and PCR plates on the market is 1 operation.

3. How to use

3.1 Preparation of the rotor

(1) Check that there is no crack, corrosion or deformation on the rotor.

- Check that the drive pin of the drive hole (crown hole) is not deformed.
- Check that there is no corrosion, foreign matters, dirt, and scratches in the drive hole (crown hole).

CAUTION:

If there is any abnormality such as corrosion or cracks, stop use of the rotor and contact an authorized sales or service representative.

(2) Check that there is no foreign matter in the rotor.

3.2 Preparation of the sample and the microplates/PCR plates

Check that the microplates/PCR plates are free from cracks or deformation. Do not use cracked or deformed microplates/PCR plates.

Inject the sample into the microplates or the PCR plates.

WARNING :

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

CAUTION:

- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor, M plate adapters/PCR plate adapters, microplates/PCR plates, and PCR tuber racks. Using such a sample could corrode or deteriorate them.
- Chemical resistance and the strength of microplates/PCR plates varies with speed, temperature, and so on. Before using sample, fill microplates/PCR plates with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise the sample leakage may occur during operation.
- Use the microplates/PCR plates that are same products. Otherwise, the microplates/PCR plates may be damaged due to their different form or the imbalance operation may occur, and it may result in the damage to the centrifuge and the rotor.

3.3 Balancing the microplates

The allowable imbalance of this rotor is 4g

- Microplates

Allowable imbalance: Difference in weight of contents (M plate adapters, microplates, and samples) in the symmetrically arranged compartment of the rotor

- PCR plates

Allowable imbalance: Difference in weight of contents (PCR plate adapters, PCR plates, PCR tube racks, and samples) in the symmetrically arranged compartment of the rotor

CAUTION:

- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise, the precipitation levels may be different by centrifugation and there may be variations in position of center of gravity, and it may result in damage to the centrifuge and the rotor due to the imbalance operation.
- Balance between symmetrically arranged microplates/PCR plates filled with the same-density sample. Variations in sample density cause imbalance operation, and it may result in the damage to the centrifuge and the rotor.

3.4 Loading the plates in the rotor

Check that the rotor is free from foreign substances. Adjust the position of the fix plate according to the size of the plates. Load the balanced plates and plate adapters symmetrically with respect to the drive shaft as shown in Fig. 3-1.


Fig. 3-1 Loading plates

CAUTION:

- Do not load only one plate nor arrange the plates asymmetrically. Otherwise, imbalance operation resulting in a failure of the centrifuge or rotor can be caused.
- Mount the fix plates symmetrically with respect to the drive shaft and securely tighten them with the mounting screws. If the mounting screws are loose, they may be removed during operation and it may cause imbalance operation resulting in damage to the centrifuge and the rotor.
- Be sure to mount all of the M plate adapters or the PCR plate adapters regardless of absence or presence of samples.
- Only one microplate can be loaded in one compartment of the rotor when running the centrifuge at 2,700 rpm or higher. Use of two or more plates in one compartment can cause damage to the plates and the rotor.
- Do not use the PCR tube rack with the M plate adapter. Otherwise, the PCR tube rack may deform or break due to mismatch in shape.

Be sure to use the PCR tube rack when using the PCR plate. (Remove the lid of the PCR tube rack before centrifugation.)

3.5 Mounting the rotor cover


CAUTION:

- Be sure to mount the rotor cover to the rotor body before operation, otherwise the rotor may be removed during operation and result in damage to the centrifuge and the rotor.

(1) Apply silicone grease (vacuum grease) (standard accessory of the centrifuge) to the packing (for the rotor) and then fit the packing (for the rotor) in the groove of the rotor body. Replace the packing (for the rotor) with new one if damaged or deteriorated.

(2) Mount the rotor cover on the rotor body.

3.6 Operation

Mount the rotor onto the drive shaft of the centrifuge gently and securely. Perform operation according to the instructions of the centrifuge instruction manual.

3.7 Taking out samples

After the operation, confirm that the rotor is completely stopped. Gently take out the rotor from the centrifuge and remove the rotor cover. Take out the plates from the rotor and then take out the sample in a proper manner.

4. Maintenance

4.1 Corrosion resistance of rotor

 **WARNING:**
Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.

This rotor is made of aluminum alloy. The rotors made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode these rotors and decrease the strength. Use samples that will not affect the rotor referring to the attached chemical resistance chart.

4.2 Maintenance of rotor

After the operation, properly maintain the rotor to prevent corrosion that can cause the rotor breakage.

- Normal maintenance

 **CAUTION:**
Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor with tap water or a dilute solution of neutral detergent and rinse it out with distilled water. Wipe the rotor with a soft cloth. Turn the rotor upside down with the rotor cover removed to dry it well. Check that the rotor is completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor. Store the rotor in a dry place.

- Maintenance after use of a corrosive sample

Wash the rotor with tap water immediately after the operation. Then perform the normal maintenance.

- Maintenance when foreign substances are adhered to the rotor

 **WARNING :**
Do not allow the rotor temperature to rise over 100°C. Otherwise, the rotor can be brittle.

 **CAUTION:**
Use a neutral detergent having a pH between 5 and 9, otherwise the rotor can be discolored or corroded.

Soak the rotor in warm water for one or two hours and wash the inside of the rotor and the crown hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the crown hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

4.3 Sterilizing rotor

- Sterilize the aluminum alloy rotors according to either gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).

WARNING :

- Never sterilize the rotor by autoclaving or boiling. Otherwise, the rotor will be brittle. It is very dangerous.

CAUTION:

- Do not dip the rotor in the formalin (3%) solution more than 2 hours.

4.4 Inspecting rotor


Periodically check that the rotor is free from corrosion every 100 hours use. Especially be careful of the inside of the rotor and the crown hole at the rotor bottom because the rotor can be brittle if these portions are corroded. If the rotor surface is discolored, dented or cracked, the rotor is corroded. Do not use such a corroded rotor and contact an authorized sales or service representative immediately for further inspection.

5. Packing list

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	Remark
①	Rotor	1	—
②	Packing (O-ring)(for the rotor)	1*	—
③	Fix plate	8*	—
④	Hex. bar wrench (for tightening the mounting screws of the fix plate)	1	—
⑤	M plate adapter	4	—
⑥	PCR plate adapter	4	—
⑦	Rotor instruction manual	1	Part No.S999429

*: The quantities of the parts that are mounted to the rotor.

6. Spare parts

Figure	Part No.	Name	Qty
	S401867	Packing (O-ring)(for the rotor)	1
	8046005	Hex. bar wrench (for tightening the mounting screws of the fix plate)	1
	S410604A	M plate adapter	4
	S410605A	PCR plate adapter	4
	483719	Silicone grease (Vacuum grease)	1
	S999429	Rotor instruction manual	1

7. Decontamination

WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

8. Rotor retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

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Copyright © 2022 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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