



# Remote and quantitative cell culture monitoring: Vero cells

## Performance of Olympus Provi™ CM20 incubation monitoring system in the CellXpert® C170i CO<sub>2</sub> incubator

**Introduction & Aim:** Systems for remote monitoring of confluency and cell count can significantly reduce efforts in cell culture. As CO<sub>2</sub> incubators are thermally insulated and typically do not have active cooling, additional heat introduced by monitoring devices may interfere with the CO<sub>2</sub> incubators temperature control. Furthermore, the size of the device may interfere with the airflow. Both factors can lead to a higher risk of condensation and temperature control failure. Therefore, it is recommended to test the combination of systems to ensure optimal performance.




**Procedure:** Vero cells were seeded in appropriate culture medium (RPMI 1640, 10 % FBS). Two cell culture vessels were placed each on one Provi™ CM20 incubation monitor head, installed in the CellXpert C170i CO<sub>2</sub> incubator from Eppendorf in a room at a temperature of 21.2 °C. Cells were incubated at 37 °C and 5 % CO<sub>2</sub>. Cell growth was monitored for approximately 60 hours after seeding. CO<sub>2</sub> concentration and temperature homogeneity were simultaneously documented. The temperature homogeneity in the CO<sub>2</sub> incubator was measured by external probes at 27 measuring points distributed on the top, middle, and bottom shelf according to DIN 12880:2007-05 standard.



Representative microscopic images of Vero cells (A) scanned regularly with the CM20 system during the incubation period to determine confluency and number of cells (B).

**Results:** The two monitoring heads did not reduce the specified temperature homogeneity of  $\pm 0.3\text{ }^{\circ}\text{C}$  (at  $37\text{ }^{\circ}\text{C}$ ) inside the CellXpert C170i incubator. Furthermore, neither interference with  $\text{CO}_2$  control nor condensation due to airflow interference was detected.

**Conclusion:** Here, we have shown the compatibility of two Olympus Provi™ CM20 incubation monitoring systems operating inside the CellXpert C170i  $\text{CO}_2$  incubator from Eppendorf. Thus, this setup ensures equal growth conditions for cells at different locations inside the  $\text{CO}_2$  incubator, even on the top shelf (fanless design). It also proves comparability of growth conditions between experiments performed with the CellXpert C170i with vs. without the use of the Olympus Provi™ CM20 system.



The set  $\text{CO}_2$  as well as temperature parameters and their homogeneity are maintained throughout the experimental incubation period.

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)