



Visionize[®] Lab Suite Product Guide

You and Your Lab – Connected.

Table of Contents

Introducing VisioNize Lab Suite

Introduction	3
> Connectable devices	4
> Video: Experience the VisioNize Lab Suite	5

Benefits of using VisioNize Lab Suite

> Features of VisioNize Lab Suite

> Benefits of connecting Your	
> ULT freezers	11
> Incubators	12
> Shakers.....	13
> Cyclers, mixers, and centrifuges.....	14
> Enabling connectivity with VisioNize box	15
> Enabling connectivity with VisioNize sense	16

Getting started

> Getting you and your lab connected	17
> Sign up: Choose the plan that meets your and your labs needs	18
> Connect your laboratory devices.....	19

Integrations - Leveraging device monitoring data

> eLabNext: Improve your lab workflows	21
> Clustermarket: Improve your equipment planning	22

Reference cases	24
------------------------------	-----------

List of helpful documents	25
--	-----------



Monitor Your Laboratory Equipment in Real Time and Receive Notifications Straight to Your Phone when Issues Occur

VisioNize® Lab Suite ensures you access to critical parameters of your connected lab devices to ensure the safety of your samples and the reliability of your research. The user-friendly platform provides easy, quick access to all relevant data so you can react rapidly to avoid severe problems in the lab.



Remote monitoring

- > Monitor connected lab devices and track status, parameters, and performance.



Real-time-notification

- > Get alert and event notifications to address issues quickly.



Document management & device data export

- > Quickly retrieve relevant device documentation, e.g. operating manuals and other SOPs. Export device parameter, event data, and any user interactions.



Maintenance management

- > Plan your upcoming maintenance tasks in one tool and receive reminders when service is due.



Connectable Devices

Do you already have a new Eppendorf device with a touch screen?:

Then you're ready to start VisioNize® Lab Suite. Out-of-the-box connectivity with VisioNize onboard devices:

Eppendorf ULT Freezers

CryoCube® F740i
CryoCube® F740hi
CryoCube® F740iW
CryoCube® F740hiw

Eppendorf CO₂ Incubators

CellXpert® C170i

Eppendorf Shakers

CellXpert® CS220
Innova® S44i*
**only compatible after serial number S44IJH302546*

Eppendorf Centrifuge

Centrifuge 5910 Ri

Eppendorf Cyclers

Mastercycler® X50 series*
Mastercycler® X40

*Mastercycler X50 eco not supported without Mastercycler X50 master device

No VisioNize® onboard devices?

Don't worry, easily connect with our [connect ware](#) VisioNize box or VisioNize sense*:

VisioNize box enabled

The [VisioNize box](#) acts as a hardware hub and works as a translation system enabling data collection of the device data in VisioNize Lab Suite.

VisioNize sense

Connect your non-Eppendorf devices to the VisioNize® Lab Suite. Sensors, [monitor](#), e.g., ambient temperature, humidity, pressure & light, and temperature of any device (including ULT freezers) within the range of +60 to -90 °C.

Note: VisioNize sense monitoring kit not available worldwide. Please contact us for local availability.

VisioNize box enabled devices

Eppendorf Centrifuges

Centrifuge 5910 R
Centrifuge 5424 R
Centrifuge 5425
Centrifuge 5425 R
Centrifuge 5427R

Eppendorf CO₂ Incubators

New Brunswick S41i
CellXpert® C170
Galaxy 48R

Eppendorf Shakers

Innova® 40/40R
Innova® 42/42R
Innova® 43/43R
Innova® 44/44R

Eppendorf Thermomixer

Eppendorf ThermoMixer C,
F0.5, F1.5, F2.0, Fp
Eppendorf Thermostat C

Eppendorf Cyclers

Mastercycler® Nexus
Mastercycler® Nexus Gradient
Mastercycler® Nexus flat
Mastercycler® Nexus X2
Mastercycler® nexus GX2
Mastercycler® Nexus SX1
Mastercycler® Nexus GSX1

Eppendorf ULT Freezer

CryoCube® F440n
CryoCube® F440h
CryoCube® F570

Eppendorf ULT Freezer

CryoCube® F570h
CryoCube® F570hw
CryoCube® FC660
CryoCube® FC660h
CryoCube® F740
CryoCube® F740h
Innova® U101
Innova® U535
Innova® U570
Innova® U725
Innova® U725-G
Innova® C760
Innova® C585
Premium U410
Premium U570
HEF U410



Experience the Visionize Lab Suite



Learn more about the benefits of switching to a connected digital lab. Our specialists will give you insights into each service of the Visionize Lab Suite.

Join product tour:



**Eppendorf digital platform –
Demo of the Visionize Lab Suite**



Features of VisioNize Lab Suite



Real-Time Monitoring

Monitor your devices remotely and track device status, parameters (e.g. temperature), and history from anywhere



Alert Management

Track all alert acknowledgements with time-stamped data logs



Task Management

Plan and assign device-related maintenance tasks to ensure optimal device performance



Instant Notifications

Receive alert and event notifications and quickly address any issues



Device Connectivity

Connect and monitor all your temperature-controlled devices including 3rd party with VisioNize sense monitoring kit



Response-Based Alarm Escalation

Define escalation schemes to individual or multiple recipient groups based on presence in the lab



Easy Data Access

Export device and event data, as well as logged user interactions




Protecting your samples
Ensure the safety of your samples 24/7



More Peace of Mind with Real-Time Device Monitoring & Notifications

Overview of all connected devices, e.g. freezers and incubators.


Monitoring dashboard of a CryoCube® F740hi freezer showing:
> Device status
> Parameters
> History chart





Meet Compliance Needs with Device Data Export

Digitalize your audit preparations with VisioNize Lab Suite and eliminate paper-based documentation.





Maintenance Management for an Effective Lab

VisioNize Lab Suite supports you in managing your laboratory equipment to ensure your lab processes keep going. Oversee the maintenance status and plan and assign upcoming tasks to maintain your equipment.

The screenshot displays the VisioNize task management interface. On the left, a sidebar shows a list of devices under 'Overview Upcoming Tasks'. The 'HH S44i Innova S44i' device is highlighted in blue. The main panel shows the details for this device, including its maintenance status ('Overdue'), a task list ('Shaker yearly maintenance' due 3 days overdue), and device information (Serial Number S44iIM102404, Registration Number -, Location Building1 Cambridge). A red arrow points from a text box to the 'Maintenance status: Overdue' indicator. At the bottom, an 'Add Task' button is highlighted with a blue arrow pointing to another text box.

VisioNize | task management

Device | Task

Overview Upcoming Tasks

Search

Sort by Maintenance Status

Overdue

HH Innova 40 Innova 4x Shaker

HH S44i Innova S44i

Ready to work

1006GI300054 VisioNize box

1006GM100159 VisioNize box

1006HJ000223 VisioNize box

HH S44i Innova S44i

Maintenance status: Overdue

Tasks

Shaker yearly maintenance 3 days overdue

Due Date 2021-01-08

Description -

Responsible Person visionize.user@eppendorf.com

Edit Task Perform Task

Device Information

Serial Number S44iIM102404

Registration Number -

Location Building1 Cambridge

Edit Device

Related Documents

User Manuals (3)

Certificates (3)

Add Task

Add customizable tasks and assign them to lab members to effectively manage your lab.

Check the maintenance status of your device and if you can use the device for your experiment.



Manage Your Activities

Plan your **upcoming tasks in one tool** and receive reminders when service is due. **Never lose sight** of your device maintenance needs with the VisioNize Lab Suite. All new software updates for VisioNize touch devices, VisioNize box and VisioNize sense are available immediately through VisioNize Lab Suite. Schedule the update based on your convenience and your lab's schedule!

The screenshot displays the 'Pipette calibration' task planning interface. It is divided into three main sections: Task Description, Task Planning, and Devices.

- Task Description:** Contains a text box with the description 'Send Pipettes to Eppendorf calibration service.'
- Task Planning:** Includes fields for 'Recurring Task' (set to '6 Months(s)'), 'Remind' (set to '2 Week(s) before'), 'Responsible' (set to 'VisioNize User'), 'Email-notification' (checked), and 'SMS-notification' (unchecked).
- Devices:** Lists two devices: 'M45090I Research plus 3120' and 'N47006I Research plus 3122', both marked as 'due in 182 days'. Below the list is a button labeled 'Add / Remove Device'.

Annotations with blue callout boxes provide additional context:

- A callout points to the 'Task Description' section: 'Plan one-time or recurring maintenance tasks for all of your devices.'
- A callout points to the 'Email-notification' and 'SMS-notification' options: 'Get a reminder for a task via email and/or SMS text message.'
- A callout points to the 'Add / Remove Device' button: 'Get notified when a maintenance task is due.'

At the bottom of the interface, there are buttons for 'Perform Task', 'Delete task', 'Cancel', and 'Save'.



Benefits of Connecting Your Freezers to VisioNize Lab Suite

- > **Monitor your entire freezer fleet** and **document** all related **device parameters** like temperature and ambient temperature as well as all device interactions, e.g., set-point changes
- > Export parameter and event data
- > Receive **real-time notifications** via email and SMS:
 - "Temperature above alert limit -60 °C"
 - "Door open longer than 2:00 min."
 - "Main power failure. Check power supply to the device."





Benefits of Connecting Your Incubators to VisioNize Lab Suite

- > **Monitor your incubators** and **track all related device parameters** such as **temperature, O₂ and CO₂** as well as all device interactions, e.g., set-point changes
- > Export parameter and event data
- > Receive **real-time notifications** via email and SMS:
 - "Temperature below alert limit 30 °C"
 - "CO₂ pressure below warning level 0.4 bar. Check gas supply."
 - "Long door opening: Close the door."





Benefits of Connecting Your Shakers to VisioNize Lab Suite

- > **Monitor your shaker runs** and **track all related device parameters** such as temperature, **agitation**, and **remaining runtime** as well as all device interactions like set-point changes
- > Export parameter and event data
- > Receive **real-time notifications** via email and SMS:
 - "Program 'pre-culture' has finished."
 - "Temperature below alarm limit 30°C"
 - "Door open longer than 5:00 min."





Benefits of Connecting Your Cyclers, Mixers and Centrifuges to VisionNize Lab Suite

- > **Monitor** all **device-related parameters** like temperature and agitation as well as all device interactions like set-point changes
- > Check the current status of your experiment
- > Export parameter and event data
- > Receive **real-time notifications** via email and SMS:
 - "Execution of program ended"
 - "Imbalance: imbalance during run"





Enabling Connectivity with VisioNize box

The easy-to-install hardware hub connects your VisioNize box enabled devices to the VisioNize Lab Suite.

The hub works as a translation system, **enabling data collection** of the device data in the platform. The box offers connections for **up to four devices** depending on the lab infrastructure/device location.

- > Connect VisioNize box **enabled (Eppendorf) devices** like freezers, incubators, shakers, centrifuges and thermomixers
- > **Monitor** all device-related parameters and run statuses
- > Have all parameter **data ready for export**

- > **Receive notifications** for connected devices covering:
 - > Significant parameter changes
 - > Errors
 - > Device status changes

Order VisioNize box and additional accessories such as cables via the Eppendorf e-shop or contact your local Eppendorf sales representative.





Enabling Connectivity with VisioNize sense

VisioNize sense can capture, record and transfer even the smallest changes in your environmental conditions to VisioNize Lab Suite. Use VisioNize sense to monitor your lab rooms, CO₂ incubators, fridges or freezers from any supplier.

- > Sensors can be used for any kind of **temperature-controlled device monitoring** (from -90°C to +60°C)
- > Configure **temperature limits** in VisioNize Lab Suite and receive notifications via SMS or email.
- > All temperature data will be documented for audits
- > Use VisioNize sense as an independent (redundant) temperature probe for a fridge or freezer
- > The environmental sensor tracks additional parameters, e.g., **ambient temperature** (-10°C to +59°C), **humidity** (0 % to 95 %, non-condensing), pressure (700 mbar to 1.200 mbar) and **brightness** (0 % to 100 %).

Order VisioNize sense and additional accessories such as cables via the Eppendorf e-shop or contact your local Eppendorf sales representative.



VisioNize Lab Suite Setup Guide



Getting You and Your Lab Connected

Make the **right preparations** for your digital lab. We recommend **checking** the following **before subscribing** to the VisioNize Lab Suite and its services:

1. Check the **IT prerequisites** for VisioNize Lab Suite with the [VisioNize Lab Suite Setup Guide](#)*
2. Check the **compatibility of your device** by having a look in the [VisioNize Lab Suite Device Overview](#)*
3. Check **what else is required** for connecting your lab devices, e.g., VisioNize box or VisioNize sense

For more guidance and precise planning, contact your local Eppendorf sales representative.



Find helpful tips and tricks to help you navigate your lab through digital transformation



*A [list of helpful documents](#) can be found in the back including the links to the documents online.



Choose the Plan That Meets Your and Your Labs Needs



You can easily subscribe online with your myEppendorf-Account to receive access to the VisioNize Lab Suite: www.eppendorf.com/visionize-subscription

Sign up for free!

For starters:

VisioNize Lab Suite Essential FREE (for first 3 devices)

Gain free, hands-on experience with VisioNize Lab Suite and enjoy a glimpse into the digital, connected lab. The Essentials plan covers all necessary aspects to get started to digitally manage the lab for increased sample safety.

- > Real-Time Monitoring
- > Email Notification
- > Data Access for 7 days
- > Device Data Export
- > Document Management
- > Task Management

Please note:

This plan already includes 3 licenses for connected devices right from the start. Each additional device requires an additional license.

Hardware and software upgrades are not included and may require an additional fee.

For more peace of mind:

VisioNize Lab Suite Sample Safety+ Price / device / year

Be confident about your sample safety and guarantee smooth operation around the clock with the Sample Safety+ plan.

Enjoy the enriched feature set in VisioNize Lab Suite for more convenience!

- > All features of Essentials plan
- > SMS Notifications
- > Data Access for 90 days
- > Response-Based Alarm Escalation
- > Audit Trail

Please note:

The price per device license depends on total number of connected devices.



Connect Your Laboratory Devices

Connecting VisionNize touch enabled devices

1. Connect your device:


Connect the Ethernet interface at the rear side of the device to a patched wall port with a standard Ethernet cable.

2. Add the device to your device list:

Navigate to the “device management” in the VisionNize Lab Suite, click “add device”, select “A VisionNize touch enabled device” and enter the serial number.

3. Accept your device:

Navigate to the “Requests” tab in the device management app, select the device you want to add and click “Accept devices”.



VisionNize Lab Suite Setup Guide



Connecting VisioNize box enabled devices

1. Connect the Eppendorf devices to the VisioNize box:

Use an appropriate adapter cable. One VisioNize box can manage up to 4 devices.

2. Update the VisioNize box:

Download the latest software version here, unpack the ZIP file on an empty USB stick and connect it to the VisioNize box.

3. Add your VisioNize box to the network:

Connect the Ethernet interface at the rear side of the device to a patched wall port with a standard Ethernet cable.

4. Add the device to your device list:

Navigate to the “device management” in the VisioNize Lab Suite, click “add device”, select “VisioNize box or VisioNize sense gateway router” and enter the serial number.

5. Accept your device:

Navigate to the “Requests” tab in the device management app, select the device you want to add and click “Accept devices”.



VisioNize Box Manual

Connecting VisioNize sense

1. Connect the VisioNize sense gateway:

Connect the Ethernet interface at the rear side of the device to a patched wall port with a standard Ethernet cable.

2. Add the device to your device list:

Navigate to the “device management” in the VisioNize Lab Suite, click “add device”, select “VisioNize box or VisioNize sense gateway router” and enter the serial number.

3. Accept your device:

Navigate to the “Requests” tab in the device management app, select the device you want to add and click “Accept devices”.

VisioNize Sense Manual



Step-by-step Setup Tutorial





eLabJournal

eLabJournal is an intuitive and flexible software solution to document research and track research data. The Electronic Lab Notebook transfers the proven paper-based documentation to the more safe, collaborative, and efficient digital world.




[Show more product information](#)

eLabInventory

With eLabInventory software, you can efficiently track any sample in the lab. Whether you want to track cell lines, tissue specimens, bacterial strains, DNA, RNA, proteins, or any biobank, you can set up templates to track any information about any item in your lab inventory.



[Show more product information](#)



Are you an eLabInventory user? Try out the VisioNize sense add-on!

The ability to monitor the viability and condition of samples in an inventory management system is imperative for laboratories, especially in GxP-compliant labs.

The VisioNize sense add-on for eLabNext eLabInventory lets users enjoy a seamless experience of their inventory management system while using environmental condition monitoring sensors from Eppendorf.

The add-on is available for free in the eLabNext Marketplace.



[Learn more](#)



Clustermarket: Improve Your Equipment Planning

Improve your lab workflows and planning of equipment

Clustermarket provides the easiest way to manage shared lab device scheduling and get all your team members' schedules aligned. Lab managers add the equipment and define booking rules, allowing team members to book the equipment by simply dragging and dropping.



[Learn more about Clustermarket](#)



[See how to activate your integration in Clustermarket.](#)





Lab Innovation Drivers & Proud VisioNize Lab Suite Users



Making data simple and accessible has become much easier thanks to increasing digitalization. Laboratories like GeneQuine Biotherapeutics are striving not only for sample security but also seek solutions allowing data access in a flexible way and moreover supporting them to keep their lab running.

Meet Kilian, Co-Founder and CEO sharing his experience using VisioNize Lab Suite.



Reference Case: Data management and performance monitoring for future-proofing the laboratory



'VisioNize Lab Suite adds value to the lab. We have a lot of sensitive samples and it is really important to not waste time on things that should not be problems. VNLS has been great in monitoring our devices enabling users to regulate things from afar.'

**Ashmita Jha, Research Assistant, MIT,
Boston, United States**



'VisioNize Lab Suite helps us automatically log device data and use it for further purposes. This makes our evaluations extremely robust and our protocolling complete.'

**Process Development Drug Substance,
CureVac AG, Tübingen, Germany**

'An unreliable freezer that I could not monitor from outside has been one of my negative experiences in the past. Thanks VisioNize® Lab Suite, I am confident that our samples are safe as it allows monitoring from anywhere and warns irregularities immediately.'

**Prof. Dr. Sabine Ellinger,
Institute of Nutrition and Food Sciences,
Human Nutrition, Uni Bonn**





'We expect lab instruments to be connectable which is why we chose Eppendorf for our -80°C freezers and cell culture incubators.'


**Kilian Guse, PhD, CEO and Co-Founder,
GeneQuine Biotherapeutics,**





Everything at a Glance


 Get ready to setup the VisioNize Lab Suite:
www.eppendorf.com/visionize-setup-guide


 Learn what is required to connect your device:
www.eppendorf.com/visionize-device-overview

 Find useful tips & tricks in our knowledge base:
www.eppendorf.com/visionize-lab-suite-knowledge-base

 Request your personal demo here:
www.eppendorf.com/visionize-demo

 Register and manage your subscription online:
www.eppendorf.com/visionize-subscription

 Explore all Eppendorf digital lab solutions here:
www.eppendorf.com/digitalabsolutions

 Explore our ecosystem partners:
www.eppendorf.com/vnls-partners



Your local distributor: www.eppendorf.com/contact

Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com/visionize

iPad, iPhone and Apple Watch are registered trademarks of Apple Inc., USA. Microsoft and Azure are registered trademarks of Microsoft Corporation, USA.

Eppendorf®, the Eppendorf Brand Design, VisioNize®, the VisioNize logo, CellXpert®, Mastercycler®, CryoCube®, Eppendorf Research® plus, Eppendorf Xplorer®/ Eppendorf Xplorer® plus are registered trademarks of Eppendorf SE, Germany. Innova®, Galaxy® are registered trademarks of Eppendorf Inc., USA. All rights reserved, including graphics and images · Copyright © 2024 by Eppendorf SE.