

High Throughput Purification of Plasmid DNA and *in vitro* Transcribed mRNA Using the High-Speed Centrifuge CR22N

Authors

Pascal Rowart¹, Vincent Dufey¹, Jan Knop², Françoise De Longueville¹

¹Eppendorf Application Technologies S.A., Namur, Belgium, ²Eppendorf SE, Hamburg, Germany

* Corresponding author: rowart.p@eppendorf-eat.be

Abstract

The generation of large amounts of plasmid DNA is important for various applications including the downstream processing of *in vitro* transcribed mRNA. This workflow is central to instant vaccine research, where plasmids are stored in large DNA library and the inserted target gene on the plasmid can be used as a template for, for example, *in vitro* transcription, to generate mRNA. To obtain consistent results, an optimized process is important. In this Application Note, a high throughput method is shown to harvest bacteria with the unique

4 x 1.5 L capacity of Rotor R9A2 in the Centrifuge CR22N. Also using Centrifuge CR22N, plasmid DNA is purified using the 10 x 50 mL and 10 x 15 mL capacity of Rotor R15A and the 30 x 2 mL capacity of Rotor R22A4. This results in the generation of a high quality and consistent DNA library that can be used for the production of mRNA by *in vitro* transcription, which can then be purified on the same device.

Introduction

Optimal plasmid production requires efficient processes to produce high-yield and high-purity products. This process has recently gained much attention because of its central role in the biotechnology and pharmaceutical industries. Here, high throughput and high-quality plasmid production are important for the development of mRNA therapeutics and vaccines for the pharmaceutical industry. In this workflow, it is important to grow large numbers of bacteria containing plasmid to generate a unique, consistent, and pure library of purified DNA for subsequent steps, to the final product. Previously, we have shown that the Centrifuge CR22N with its unique 1.5 L bottle can save up to 32% handling time when processing large bacterial batches¹⁻⁴.

In this Application Note, we show how the Centrifuge CR22N can be further used to purify plasmid DNA and PCR products, which serves as a first step in the mRNA production workflow of *in vitro* transcription. To this end, we harvested a large bacterial culture containing the pluri-potency-associated transcript target gene, LIN28A, on a plasmid DNA, grown in the Shaker Innova[®] S44i. Rotor R9A2 was used to pellet 3 L (2 x 1.5 L) of bacterial culture. The DNA purification of the entire bacterial pellet obtained from one 1.5 L bottle (1500PP bottle) was performed using a combination of the Rotor R15A, which can hold up to 10 x 50 mL and 10 x 15 mL, and the Rotor R22A4 which can hold up to 30 x 2 mL. This combination allows a reduced number of spins due to their high capacity. Finally, we show that high-quality transcription processes can produce the mRNA by *in vitro* transcription (IVT)^{5,6}.

Focus: High Throughput Production of Plasmid DNA and Purification Using High-speed Centrifuge CR22N




Fig. 1: Workflow from the bacterial culture to *in vitro* transcription using the Centrifuge CR22N and a combination of Rotor R9A2, Rotor R15A, and Rotor R22A4.

Materials and Methods

The entire procedure followed the protocol and recommendations of the research papers by Rossi et al ^{7,8}.

Plasmid

A plasmid encoding human LIN28A was purchased from the Addgene plasmid repository (Addgene plasmid #*26819; <http://n2t.net/addgene:26819>; RRID: Addgene_26819). The plasmid construct was designed by Derrick Rossi in the DH5alpha growth strain. It contained the LIN28A sequence (813 bp), the 5' sequencing primer "CMV Forward" and the 3' sequencing primer "TK polyA reverse" (5'UTR-LIN28A-3'UTR). The total plasmid backbone size with the insert was 6213 bp.

* # is the symbol for „order number“

Culturing of bacteria

Amplification in dish and cryopreservation in glycerol

Bacteria were inoculated on agar (Invitrogen™, #11518916) containing ampicillin at 100 µg/mL (Sigma-Aldrich, #A5354) in a 10 mm dish and grown for 24 h at 37 °C in an incubator. A single colony was picked and transferred to a 15 mL conical tube filled with 10 mL of LB Broth (Sigma Aldrich, #L2897) containing ampicillin at 100 µg/mL and amplified overnight at 37 °C and 200 rpm in an Innova® S44i shaker (Eppendorf, #S44I311001). Bacteria were spun at 1,400 × g for 5 min and the pellet was resuspended in 10 mL of LB Broth and glycerol 100% (v/v). Aliquots of 500 µL were frozen in cryovials at -80 °C for further experiments.

The master mix was aliquoted into 8 x 25 μ L PCR strips and run on a MasterCycler[®] X50 Aluminum (Eppendorf, #6313000018) with the following conditions: 1 cycle (Denaturation: 95 °C, 3 min) – 30 cycles (Denaturation: 98 °C, 20 s / Annealing: 60 °C, 15 s / Extension: 72 °C, 60 s) – 1 cycle (Denaturation: 72 °C, 3 min). PCR products were purified using the QIAquick PCR Purification Kit (Qiagen, #28106) as described above. Tailed PCR products were measured using the Nanodrop[™] (Thermo Fisher, #ND-2000) and stored at -20 °C for further analysis.

***In vitro* Transcription**

Before starting the synthesis of RNA, the working area and pipettes were cleaned with RNaseZap[™] RNase Decontamination Solution (Invitrogen[™], #10708345). RNA was synthesized with the MEGAscript[™] T7 Transcription Kit (Invitrogen[™], #10065754), with 40 ng/ μ L (final concentration) of purified tailed PCR product to template each 40 μ L reaction. A custom ribonucleoside blend was added comprising 4 μ L of 3'-O-Me-m7G(5')ppp(5')G RNA Cap Structure (New England Biolabs, #S1411S), 4 μ L of adenosine triphosphate, and 0.8 μ L of guanosine triphosphate (from the MEGAscript T7 kit), 3 μ L of each 5-methylcytidine triphosphate and pseudouridine triphosphate (TriLink Biotechnologies, #N-1014 and #N-1019). The custom ribonucleoside blend was vortexed and spun briefly before adding the tailed PCR product, 4 μ L of each the T7 Buffer and T7 enzyme (from the MEGAscript T7 kit). Water was added to reach 40 μ L final volume. The solution was incubated for 4 h at 37 °C in a MasterCycler[®] X50 Aluminum (Eppendorf, #6313000018). After incubation, 2 μ L of Turbo DNase (from the MEGAscript T7 kit) was added to the sample and incubated for 15 min at 37 °C. RNA was purified using the MEGAclean[™] Transcription Clean-Up Kit (Invitrogen[™], #AM1908) following the manufacturer's instructions. Briefly, samples were mixed with 60 μ L of Elution buffer, 350 μ L of Binding solution, and 250 μ L of 100 % ethanol, transferred to a cartridge, and centrifuged at 10,000 $\times g$ for 1 min at 4 °C using the Centrifuge 5702R (Eppendorf, #5703000010). Cartridges were washed in 2 steps consisting of 500 μ L of Wash solution and centrifugation at 10,000 $\times g$ for 1 min at 4 °C. RNA was collected in 2 steps consisting of 50 μ L of pre-heated (95 °C) Elution solution and centrifugation at 10,000 $\times g$ for 1 min at 4 °C. Purified RNA was treated with 11 μ L of 10 \times Antarctic phosphatase buffer and 2 μ L of Antarctic phosphatase (New England Biolabs, #M0289S) for 1 h at 37 °C using a Thermo-Mixer[®] C (Eppendorf, #5382000015). Phosphatase-treated RNA was then purified again using the MEGAclean[™] Tran-

scription Clean-Up Kit (Invitrogen[™], #AM1908) as described above. RNA was measured with the Nanodrop[™] (Thermo Fisher, #ND-2000) and stored at -80 °C for further analysis.

Electrophoresis

The analysis of linearized plasmid, PCR products, and RNA was performed using the automated platform TapeStation 4150 system (Agilent, #G2992AA) according to the manufacturer's instructions. Digested plasmid and PCR products were run through the D5000 ScreenTape System (Agilent, #5067-5588, #5067-5589). The RNA ScreenTape System (Agilent, #5067-5576, #5067-5577, #5067-5578) was used to visualize the RNA.

Result and Discussion

Bacterial culture in a volume of up to 6 liters

To increase the culture volume and quantity of bacteria in order to generate a DNA library, we used Ultra Yield 2.5 L flasks and Rotor R9A2 in the Centrifuge CR22N. In this Application Note, six Ultra Yield 2.5 L flasks were used to collect 3 liters of media to fill two bottles. Of note, users can increase the bacterial culture to up to twelve Ultra Yield 2.5 L flasks as the rotor can accommodate up to four 1.5 L bottles, which represents 6 liters of media in one run. The 1.5 L bottles have a wide mouth which makes it easy to transfer the media from Ultra Yield 2.5 L flask to the bottle. Bottles were closed using the included tools such as bottle vise and cap opener. The bottles were loaded with the bottle setter. After the centrifugation, the pellet was easy to resuspend and collect due to the unique triangular shape and the wide mouth of the bottle (Figure 1). In this Application Note, two bottles were used to harvest bacteria. Nevertheless, only one bottle was used to generate data for DNA purification, PCR, and IVT.

Plasmid DNA purification

A pellet from one bottle was used to generate a library of DNA, enough to run multiple PCRs and IVTs. The pellet was resuspended, mixed with a buffer from the CompactPrep Plasmid kit from Qiagen, and distributed to 5 x 50TC tubes. The final volume of 15 mL was centrifuged at $20,000 \times g$ for 30 min at 4 °C using the Rotor R15A. Each supernatant was then distributed to 3 extenders attached to the CompactPrep column and to the vacuum. Each column was then centrifuged at $20,000 \times g$ for 1 min at 4 °C using Rotor R22A4. The average DNA concentration at the end of the purification in columns was $1.5 \mu\text{g}/\mu\text{L}$. Taking into consideration the total elution volume and the number of columns, an estimated amount of 7.23 mg of DNA can be purified from a single bottle.

Moreover, the unique designs of these two fixed-angle rotors (Rotor R15A: 10 x 15TC and 10 x 50TC conical tubes - Rotor R22S4: 30 x 2 mL) would allow the processing (the first steps) of four bottles in a single run using Rotor R15A, with the final steps in two runs using Rotor R22S4.




Fig. 2: Representative workflow from bacterial culture to bacterial pellet. Media were transferred to the 1.5L bottle and spun at $3,340 \times g$ for 30 min at 4 °C using the Rotor R9A2 in the Centrifuge CR22N. Media was discarded, and the pellet was resuspended and collected for DNA purification

The bacterial pellet weights for both bottles were 8.92 g and 8.55 g. Considering the capacity of Rotor R9A2 and its 4 x 1.5 L bottles, users can harvest around 35 g of total bacteria to proceed to DNA purification and later IVT.

This drastically reduces operation time, and it also reduces any differences that might appear when samples are in line for the next centrifugation run when there is not enough space in the rotor. From the entire 6 L bacterial culture, users would be able to generate and purify around 29 mg of DNA (Figure 2).




Fig. 3: Representative workflow of the maximum capacity of all rotors used in the Centrifuge CR22N to isolate bacterial plasmid DNA. Using the Rotor R9A2, four pellets can be obtained from 6 L of bacterial culture. Each pellet is resuspended and distributed into 5 different conical tubes (either 15 or 50 mL). The 20 tubes can be spun all at once using Rotor R15A. Each supernatant from the tubes is divided into 3 extenders and columns. The CompactPrep Plasmid kit (Qiagen, #12863) was used as it is designed for large volumes thanks to its extenders. DNA is eluted from each column using Rotor R22S4, in two spins, generating around 29 mg of DNA. Nevertheless, in this Application Note, two bottles were spun and only one was used to purify the DNA.

Plasmid linearization and purification

Linearization of the plasmid was performed using the *SpeI* restriction enzyme. This enzyme was used as it cut the plasmid in two places creating a short (~684 bp) and a long fragment containing the inserted gene (~5,529 bp). One cut was situated a few bp downstream of the inserted gene. After purification using the R22A4 fixed-angle rotor, the digestion was verified using the Agilent D5000 ScreenTape in the TapeStation. The band at the expected size confirmed the proper digestion by the enzyme (Figure 4-A). The mean DNA concentration was 70 ng/μL.

Amplification of the target gene by PCR

To generate the template for the IVT, the addition of a poly-(A) tail to the template DNA (target gene) by tail-PCR was performed using primers Xu-F1 and Xu-T120 in the MasterCycler® X50 Aluminum. PCR products were purified using Rotor R22A4. The expected length of PCR products, based on data generated previously by Mandal et al.⁷, was confirmed using Agilent D5000 ScreenTape in the TapeStation (Figure 4-B). The mean concentration of the PCR products was 686 ng/μL.

IVT and purification of the mRNA

The final steps, which consist of the production of the mRNA, were carried out using the MasterCycler® X50 Aluminum. After DNase and phosphatase treatment, the mRNA was purified and then analyzed using the Agilent RNA ScreenTape. A single band below 1,000 bp was observed which confirmed the IVT process and the single amplification of the targeted gene. The mean concentration of the mRNA was 1,023 μg/μL. The RNA integrity number (RIN) assigns a value of 1 to 10, with 1 being strongly degraded and 10 being the least degraded (intact). The RIN score here was shown to be 10, meaning pure and intact mRNA was produced.




Fig. 4: A) Digested plasmid, B) PCR amplifying the target gene, C) IVT producing single and intact mRNA.

Conclusion

Here, we have shown that the Centrifuge CR22N can be used to consistently produce high-quality DNA from large bacterial cultures. Using different fixed-angle rotors, Centrifuge CR22N allows harvesting from 6 L of bacterial culture and obtaining 35 g of bacteria in a single centrifuge run with Rotor R9A2 and the 4 x 1.5 L triangular bottles. This rotor reduces the number of vessels to be handled, therefore greatly improving process efficiency without impacting pelleting efficacy. The bottles are designed for specific harvesting applications with many advantages described and explained in the Eppendorf white paper N°64⁴. After resuspension, the pellet from four bottles can be redistributed into 50TC and 15TC mL tubes to start the purification steps using the Rotor R15A. This rotor allows the use of both 10 x 15TC and 10 x 50TC conical tubes, holding up to 650 mL in a single run. The first steps of the DNA purifica-

tion of the four pellets can therefore be achieved using only one run. To finish the purification in 1.5/2 mL tubes using a column, the Rotor R22A4 was used as it can hold up to 30 samples. The subsequent PCR product was also purified using the Rotor R22A4.

This combination of rotors, all compatible with the Centrifuge CR22N, reduces drastically the time of sample handling and processing by one-third⁴. It simultaneously enables fast and efficient DNA purification from a single bacterial culture on a single, versatile device, under consistent conditions, thus reducing variability and inconsistencies between DNA libraries and allowing the production of high quantity and quality mRNA.

Literature

- [1] Helinski DR. A Brief History of Plasmids. *EcoSal Plus*. 2022;10(1). doi:10.1128/ECOSALPLUS.ESP-0028-2021
- [2] Liu MA. A Comparison of Plasmid DNA and mRNA as Vaccine Technologies. *Vaccines (Basel)*. 2019;7(2). doi:10.3390/VACCINES7020037
- [3] Lederberg J. Cell genetics and hereditary symbiosis. *Physiol Rev*. 1952;32(4):403-430. doi:10.1152/PHYSREV.1952.32.4.403
- [4] Eppendorf AG. Unique 4 x 1.5 L Capacity Rotor for High- Speed Centrifuges CR22N and CR30NX.
- [5] Beckert B, Masquida B. Synthesis of RNA by *In Vitro* Transcription. In: *Methods in Molecular Biology*. Vol 703. Humana Press Inc.; 2011:29-41. doi:10.1007/978-1-59745-248-9_3
- [6] Karikó K. *In vitro*-Transcribed mRNA Therapeutics: Out of the Shadows and Into the Spotlight. *Molecular Therapy*. 2019;27(4):691-692. doi:10.1016/j.ymthe.2019.03.009
- [7] Mandal PK, Rossi DJ. Reprogramming human fibroblasts to pluripotency using modified mRNA. *Nat Protoc*. 2013;8(3):568-582. doi:10.1038/nprot.2013.019
- [8] Warren L, Manos PD, Ahfeldt T, et al. Highly efficient reprogramming to pluripotency and directed differentiation of human cells with synthetic modified mRNA. *Cell Stem Cell*. 2010;7(5):618-630. doi:10.1016/J.STEM.2010.08.012

Ordering information

Description	Order No.
Centrifuge CR22N	Inquire*
Rotor R15A	5721 221 007
Rotor R22A4	5721 221 003
Rotor R9A2	5721 221 014
1500PP bottle	5721 411 030
Innova® S44i	S44I 300 001
Mastercycler® X50	6313 000 018
50TC tubes	5721 411 170
15TC tubes	5721 411 049

*Inquire the part number for your country

Your local distributor: www.eppendorf.com/contact

Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Nanodrop™, RNAzap™, MEGAclean™, MEGAscript™ and Invitrogen™ are trademarks of Thermo Fisher Scientific.

Himac® is a registered trademark of Eppendorf Himac Technologies Co., Ltd., Japan.

Eppendorf® and the Eppendorf Brand Design, ThermoMixer® C, Innova® and Mastercycler® are registered trademarks of Eppendorf SE, Germany. Eppendorf SE reserves the right to modify its products and services at any time.

This application note is subject to change without notice. Although prepared to ensure accuracy, Eppendorf SE assumes no liability for errors, or for any damages resulting from the application or use of this information. Viewing application notes alone cannot as such provide for or replace reading and respecting the current version of the operating manual.

All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE, Hamburg, Germany.