Keep this rotor instruction manual and the centrifuge manual in the file.

Angle Rotor for Refrigerated Centrifuge

R19A2

Before using these rotors, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Tubes	_____	3
3. How to use	_____	4
4. Maintenance	_____	6
5. Packing list	_____	8
6. Spare parts	_____	8
7. Decontamination	_____	9
8. Rotor retirement	_____	9
(Appendix) K Factors	_____	10

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.
These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density or kinds of tubes.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Do not run the rotor without the rotor cover in position.
Tighten the cover-locking knob securely.
- Balance the tubes and the samples including the caps and the adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the rotor tubes within their actual capacities.
- Do not use tubes that have exceeded their life expectancy. Failure to do so could result in damage of tubes and the rotor and the centrifuge. The life expectancy of tubes depends on factors such as the characteristics of samples, speed of the rotor used, and temperature.
Always check for deterioration and damage (cracks, deformation, and so on) on tubes before using them. Do not use the tubes if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.
Contact an authorized sales or service representative.

1. Specifications

1.1 Construction


Figure1-1 Construction of R19A2 rotor

1.2 Specifications

Applicable centrifuge*1

(Some centrifuges are on sale in only Japan.
Visit our web site.)

Max. speed*2:	19,000 rpm
Max. RCF:	50,000 × g
Nominal rotor capacity:	50 mL × 8 tubes = 400 mL
Rotor Size:	φ 255 × 180.5 mm
Weight:	9.2 kg
Material:	Aluminum alloy
Rotor No.:	63

*1: For applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

*2: The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

1.3 Cross-sectional view of rotor


Figure1-2 Cross-sectional view of R19A2 rotor

1.4 Characteristic of rotor

Table1-1 Characteristics (when himac 50TC tubes are used)

Rotor Speed (rpm)	RCF (×g)		K factor*
	Rmin. (6.67cm)	Rmax. (12.4cm)	
4000	1190	2220	9817
6000	2680	4990	4363
8000	4770	8870	2454
10000	7460	13900	1571
12000	10700	20000	1091
14000	14600	27200	801
16000	19100	35500	614
18000	24200	44900	485
19000	26900	50000	435

*For details about K factor, refer to “(Appendix) K Factors”

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)
N: Speed(rpm)

1.5 Allowable speed of the rotor

WARNING :

Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic, kinds of tubes and centrifuge model.

Do not exceed the allowable speed of the rotor.

The maximum speed marked on the surface of the rotor cover is permitted under the conditions as follows.

-Average sample density : less than 1.2g/mL

-Standard tubes are used.

Depending on sample density, kinds of tubes and centrifuge model, the rotor speed may be limited as follows. Do not exceed the allowable speed of the rotor.

(1) Allowable speed for sample density

WARNING :

When sample density exceeds 1.2g/mL, calculate the allowable speed according to the following equation.

$$\text{Allowable speed (rpm)} = 19,000(\text{rpm}) \times \sqrt{\frac{1.2(\text{g/mL})}{\text{Average density of sample}(\text{g/mL})}}$$

(2) Allowable speed for combination of tube, cap, adapter, etc.

The rotor speed should be limited by combination of tube, cap, adapter, etc.

For details, refer to “2.1 Applicable tubes and adapters”.

Perform operation under the allowable speed, otherwise tubes may be broken during operation.

When using tubes on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer. Otherwise tubes may be broken during the operation.

2. Tubes

2.1 Applicable tubes and adapters

CAUTION:

The tube cavities of this rotor and 15TA50T adapter (C) are specifically designed for conical culture tubes. Do not use round-bottom or self-standing culture tubes.

Use tubes which match the tube cavities of this rotor. Failure to do so, tube might be damaged and it might be hard to take out the tubes from the rotor.

Use tubes which match the adapters. Failure to do so, tube might be damaged and it might be hard to take out the tubes from the adapters.

- When using tubes on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tubes whose allowable speed or RCF is not specified, otherwise the tubes may be broken during operation.
We will check the applicability of tubes upon request. Contact our service representative.
- Clean and sterilize tubes on the market according to the instructions of the manufacturer.
- Tubes on the market, himac 50TC tubes, and himac 15TC tubes should be discarded after one use.

Table2-1 Applicable tubes and adapters

Adapter			Tubes				Max. speed (Max. RCF)
Part No.	Name	Capacity × Qty	Part No.	Name	Actual capacity	Size (mm)	
—	—	—	S311187A	himac 50TC tube(*2) (100 pcs./box)	45 mL	φ 30 (*3) × 115	19,000 rpm (50,000 × g)
—	—	—	S307904A	himac 50TC tube(*2) (300 pcs./box)	45 mL	φ 30 (*3) × 115	19,000 rpm (50,000 × g)
—	—	—	(Marketed)	50mL PP culture tube (conical bottom)	45 mL	φ 30 × 115	11,600 rpm (18,650 × g)
S411576A	15TA50T adapter (C) (*1) (2 pcs.)	15mL × 1	S311184A	himac 15TC tube(*2) (100 pcs./box)	15 mL	φ 17 (*3) × 120	9,650 rpm (12,610 × g)
S411576A	15TA50T adapter (C) (*1) (2 pcs.)	15mL × 1	(Marketed)	15mL PP culture tube (conical bottom)	13 mL	φ 17 × 120	9,650 rpm (12,610 × g)

PP: Polypropylene

(*1): Material of the adapter: Polyacetal

(*2): Material of the tube: Polypropylene Material of the cap: Polyethylene

(*3): Excluding the cap

The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”.

The rotor speed is limited by the sample density and so on (see Section1.5).

Maximum speeds listed are guidelines only. Because of variances in user methodologies, no guarantee of performance is expressed or implied.

3. How to use

3.1 Preparation of tubes

Check that the tubes are free from crack or deformation. Do not use remarkably deformed or cracked tubes.

WARNING :

- This rotor and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

CAUTION:

- Do not use any sample inapplicable to the tubes and the caps and so on in use referring to the attached "Chemical resistance chart". Otherwise, the tubes and the caps may be degraded.
- Chemical resistance and the strength of tubes varies with speed, temperature, and so on. Before using sample, fill tubes with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the specified sample capacity of the tubes, otherwise sample leakage or deformation of tubes may occur.
- Balance the sample within the allowable imbalance including the weight of the tubes, caps and adapters. Excessive imbalance can cause damage to the centrifuge and the rotor.
- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise the precipitation levels may be different by centrifugation and such operation may increase the level of imbalance, and it may result in damage to the centrifuge and the rotor.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise the sample leakage may occur during operation.
- Use the tubes that are same products. Otherwise, the tubes may be damaged due to their different form or the imbalance operation may occur, and it may result in the damage to the centrifuge and the rotor.
- Use the tubes at a temperature between 4°C and 25°C. Otherwise, tubes may be broken or deformed during operation.
- Before using tubes on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.

Inject the samples into the tubes and cap the tubes properly. Do not exceed the specified sample capacity of the tubes, otherwise sample leakage or deformation of tubes may occur.

Balance the tubes that are symmetrically arranged. The allowable imbalance of this rotor is within 2 g.

3.2 Setting of tubes into rotor


CAUTION:

Do not load only one tube nor arrange the tubes asymmetrically with respect to the drive shaft, otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.

Check that the tube cavities of the rotor are free from foreign substances. Set the balanced tubes symmetrically with respect to the drive shaft.

Pre-cool the rotor when the operating temperature of the rotor is lower than the room temperature or pre-heat the rotor when the operating temperature of the rotor is higher than the room temperature before setting the tubes in the rotor.

Examples of incorrect tube setting are shown below.

3.3 Mounting the rotor cover


⚠ CAUTION:

- Be sure to mount the rotor cover to the rotor body and tighten the cover locking knob securely before operation, otherwise the rotor or the rotor cover may be removed during operation and result in damage to the centrifuge and the rotor.
- Do not mount the rotor cover into the rotor while the rotor is placed in the rotor chamber. Doing so might cause the rotor cover to be tightened insufficiently.
Do not remove the rotor cover from the rotor while the rotor is placed in the rotor chamber. Doing so might cause the parts of the rotor cover to be damaged.

- (1) Apply silicone grease (vacuum grease) (standard accessory of the centrifuge) to the packing (for the rotor) and then fit the packing (for the rotor) in the groove of the rotor body. Replace the packing (for the rotor) with new one if damaged or deteriorated.

If you neglect to apply silicone grease (vacuum grease), it might be hard to mount or remove the rotor cover. When it is hard to mount or remove the rotor cover, apply silicone grease (vacuum grease) to the packing (for the rotor).

- (2) Apply lubricant for screw (Part No. 84810601) to the threaded portion of the cover locking knob. Put the rotor cover on the rotor body. Turn the cover locking knob to screw in the rotor body securely.

3.4 Operation

Mount the rotor onto the drive shaft of the centrifuge gently and securely. Check that the rotor cover is mounted and the cover-locking knob is securely tightened. Perform operation according to the instructions of the centrifuge instruction manual.

3.5 Taking out samples

After the operation, gently take out the rotor from the centrifuge and remove the rotor cover by turning the cover locking knob. Take out the tubes from the rotor and then take out the sample in a proper manner.

4. Maintenance

4.1 Corrosion resistance of rotor

WARNING:

Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.

This rotor is made of aluminum alloy. The rotors made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode these rotors and decrease the strength. Use samples that will not affect the rotor referring to the attached chemical resistance chart.

4.2 Maintenance of rotor

After the operation, properly maintain the rotor to prevent corrosion that can cause the rotor breakage.

- Normal maintenance

CAUTION:

Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor with tap water or a dilute solution of neutral detergent and rinse it out with distilled water. Wipe the rotor with a soft cloth. Turn the rotor upside down with the rotor cover removed to dry it well. Check that the rotor is completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor. Store the rotor in a dry place.

- Maintenance after use of a corrosive sample

Wash the rotor with tap water immediately after the operation. Then perform the normal maintenance.

- Maintenance when foreign substances are adhered to the rotor

WARNING :

Do not allow the rotor temperature to rise over 100°C. Otherwise, the rotor can be brittle.

CAUTION:

Use a neutral detergent having a pH between 5 and 9, otherwise the rotor can be discolored or corroded.

Soak the rotor in warm water for one or two hours and wash the inside of the tube cavities and the crown hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the crown hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

4.3 Sterilizing rotor

- Sterilize this rotor according to either gas sterilization method (ethylene oxide or formaldehyde), chemical sterilization method (70% ethanol, 3% hydrogen peroxide, or 3% formalin), autoclaving (121°C, 20 minutes), or ultraviolet rays sterilization.
- Sterilize the plastic or rubber adapters according to either gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, or 3% formalin).

- | |
|--|
| <p> WARNING: Be sure to remove the rotor cover from the rotor before autoclaving, otherwise the rotor or the rotor cover may be deformed.</p> <p> WARNING: Never sterilize the rotor by boiling, otherwise the rotor will be brittle.</p> <p> CAUTION : Never sterilize the adapters by autoclaving or boiling, otherwise the adapters may be deformed. Such adapters cannot be used.</p> <p> CAUTION : After autoclaving, wait until the temperature in the autoclaving chamber reduces to the room temperature, then take out the rotor and the rotor cover.</p> <p> CAUTION : Take out the rotor and the rotor cover carefully from the autoclaving chamber as they are wet and slippery.</p> <p> CAUTION : Check the packing (for the rotor) after sterilizing. The packing (for the rotor) may be damaged during operation due to deterioration. Replace the new one if a damage or deterioration is observed.</p> <p> CAUTION : Do not dip the rotor and the adapters in the formalin (3%) solution more than 2 hours.</p> |
|--|

4.4 Inspecting rotor


Periodically check that the rotor is free from corrosion every 100 hours use. Especially be careful of the tube cavities of the rotor and the crown hole at the rotor bottom because the rotor can be brittle if these portions are corroded. If the rotor surface is discolored, dented or cracked, the rotor is corroded. Do not use such a corroded rotor and contact an authorized sales or service representative immediately for further inspection.

5. Packing list


Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	Remark
①	Rotor	1	—
②	Packing (for the handle) (O-Ring)	1*	—
③			
④	himac 50TC tube	1	50TC Tube: 100 pcs. 50TC Cap: 100 pcs.
⑤	Rotor instruction manual	1	Part No. S998652

*: the quantities of the parts that are mounted to the rotor

6. Spare parts

Figure	Part No.	Name	Qty
	84521017	Packing (for the handle) (O-Ring)	3
	S409676A	Packing (for the rotor)	3
	S311187A	himac 50TC tube	100
	84810601	Lubricant for screw	1
	483719	Silicone grease (vacuum grease)	1
	S998652	Rotor instruction manual	1

7. Decontamination

WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

8. Rotor retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

(Appendix) K Factors

- Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

Namely, the time taken for particles of 1S to settle from Rmin. to Rmax. is K factor.
The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R \text{ max. }) - \ln(R \text{ min. })}{\omega^2} \cdot \frac{10^{13}}{3600}$$

ω : Angular velocity $\omega = 2\pi/60 \times N$

N : Speed (rpm)

Rmin. : Distance from axis of rotation to top of solution (cm)

Rmax. : Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t (in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 19,000rpm of Model R19A2 rotor is 435.
The time taken for particles of 100 S to be settled from Rmin. to Rmax. using this rotor is as follows:

$$t = \frac{435}{100} = 4.4 \text{ hours}$$

The time taken for settling of particles of 100 S by the same rotor at a speed of 15,000 rpm is:

$$K = \left(\frac{19,000}{15,000} \right)^2 \times 435 = 698 \quad t = \frac{698}{100} = 7.0 \text{ hours}$$

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Copyright © 2021 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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