Keep this rotor instruction manual and the centrifuge manual in the file.

Swing Rotor for Refrigerated Centrifuge

R25ST

Before using this rotor, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Tubes/Bottles	_____	2
3. How to Use	_____	3
4. Maintenance	_____	6
5. Using Cesium Chloride	_____	7
6. Using Density Gradient Centrifugation	_____	7
7. Packing List	_____	8
8. Rotor Life	_____	9
9. Decontamination	_____	10
10. Rotor Retirement	_____	10
(Appendix) K Factors	_____	11


2021.08
S99865509

Original instruction manual

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.
These notes are defined as follows.

-  **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.
-  **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING

- Never use any material capable of producing flammable or explosive vapors. Your centrifuge is not explosion proof centrifuge and your rotor is not explosion proof rotor.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density or kinds of tubes.
- Check the chemical resistance chart attached the rotor, and do not use any sample inapplicable to the rotor (including buckets). Using such a sample could corrode the rotor (including buckets).
- Do not allow the rotor temperature to rise 100°C.: This would cause the material brittle.
- Never sterilize the rotor by autoclaving or boiling. Otherwise the strength of rotor may be significantly decreased.
- When using swing rotor, be sure to set all buckets of the same type whether or not samples are put in buckets: Failure to do so could not only cause the rotor to vibrate, but could result in the rotor being deformed and the buckets becoming detached, which is very dangerous.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Check the chemical resistance chart attached to the rotor, and do not use any sample inapplicable to the tubes. Using such a sample could deteriorate them.
- The allowable speed may be lower depending on the density of sample, use of salt-such as cesium chloride-or the combination of tubes and tube adapters. Refer to this instruction manual.
- Always balance the tubes and samples within the allowable imbalance of rotor. Do not exceed the allowable imbalance. Otherwise the rotor and the centrifuge may be damaged.
- Do not load only one tube or load tubes asymmetrically: Asymmetrically loading may cause imbalance operation and damage the centrifuge and rotor.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the rotor tubes and bottles within their actual capacities.
- Do not use tubes that have exceeded their life expectancy. Failure to do so could result in damage of tubes and the rotor and the centrifuge. The life expectancy of tubes depends on factors such as the characteristics of samples, speed of the rotor used, and temperature. Always check for deterioration and damage (cracks, deformation, and so on) on tubes before using them. Do not use the tubes if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.
Contact an authorized sales or service representative.

1. Specifications

1.1 Construction


Fig. 1-1. Construction of R25ST rotor

1.2 Specifications

- Applicable centrifuge*1

(Some centrifuges are on sale in only Japan. Visit to our web site.)

- Max. speed*2: 25,000 rpm
- Max. RCF: 110,000 × g
- Allowable density at maximum speed of the rotor: 1.2g/mL
- K Factor at maximum speed of the rotor: 325
- Max. rotor capacity(nominal): 40mL × 6tubes = 240mL
- Material of rotor: Ti alloy
- Tube size: $\phi 2.6 \times 9.0$ cm
- Weight(includes bucket and cap): 7.7 kg
- (without bucket and cap) : 6.9 kg
- Max. diameter: 21.7(31.8)cm
- Max. height: 17.6 cm
- Rotor No.: 72

(The numerical value in parentheses is maximum diameter measured with the bucket swinging.)

*1: For applicable centrifuges, refer to the "Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))."

*2: The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the "Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))."

1.3 Cross-sectional view of rotor


Fig. 1-2. Cross-sectional view of R25ST rotor

⚠ CAUTION:

Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

1.4 Characteristic of rotor

Table 1-1. Characteristic of R25ST

Rotor Speed (rpm)	RCF (×g)		K factor*
	Rmin. (7.05cm)	Rmax. (15.72cm)	
5000	1970	4390	8125
10000	7880	17600	2031
15000	17700	39500	903
20000	31500	70300	508
21000	34800	77500	461
22000	38100	85100	420
23000	41700	93000	384
24000	45400	101000	353
25000	49300	110000	325

*For details about K factor, refer to "(Appendix) K Factors"

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: Speed (rpm)

2. Tubes/Bottles

2.1 Applicable tubes/bottles

Table 2-1. Applicable tubes/bottles

Nominal capacity (mL)	Tubes/Bottles				Tube adapter		Cap		Max. speed (rpm)	Max. RCF (xg)
	Part No.	Name	Size (ϕ X L cm)	Actual capacity (mL/tube)	Part No.	Name	Part No.	Name		
40	329607A	40PA tube	2.6 X 9.0	35.6	---	---	---	---	25,000	110,000
	325754A	40PC tube	2.6 X 9.0	38.7	---	---	---	---	25,000	110,000
	S303279A	40PET tube	2.6 X 9.0	38.6	---	---	---	---	25,000	110,000
	S412571A	40SST tube (B)	2.6 X 9.0	38.4	---	---	---	---	15,800	43,900
33	S311544A	33PA seal tube	2.6 X 8.8	34.0	---	---	S413983	E4-space cap	25,000	110,000

1: When using 33PA seal tubes, we recommend the purchase of the tube kit, what is the set of tubes and necessary tools.

Refer to Chapter 7 for the tube kit.

2: Material of the space cap: Modified polyphenylene ether


CAUTION :

Chemical resistance and the strength of tubes varies with speed, temperature, and so on. Before using sample, fill tubes with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.

3. How to Use

The brief procedure of tubes/bottles handling is described as follows.

=For handling=

● Allowable speed of rotors

1. To centrifuge a sample of an average density more than 1.2g/mL.
Determine the allowable speed from the following equation.

$$\text{Allowable speed (rpm)} = 25,000 \text{ (rpm)} \times \sqrt{\frac{1.2 \text{ (g/mL)}}{\text{Mean density of sample (g/mL)}}}$$

2. The allowable speed may be lower depending on the density of a sample when using of salt such as cesium chloride. Refer to "5. Using Cesium Chloride" for the detail.
3. The allowable speed may be lower depending on the tubes, etc. Refer to "2. Tubes/Bottles"

● Tube

Fill the tube up to within 3mm from the top of the tube with sample.

● Seal tube

These tubes must be filled up the sample. Weld the seal tube completely.

● Handling the rotor

WARNING :

- Never exceed the maximum speed of the rotor (mentioned on the rotor). Otherwise, rotor may be damaged and result in damage to the centrifuge.
- Check the chemical resistance chart attached the rotor (including buckets), and do not use any sample inapplicable to the rotor. Using such a sample could corrode the rotor (including buckets).

CAUTION :

- Inspect the rotor before use (see the "4. Maintenance · Inspecting rotor").
- Do not load only one tube or load tubes asymmetrically : Asymmetrically loading may cause imbalance operation and damage the centrifuge and rotor.
- Balance the tubes and the samples within the allowable imbalance of the rotor.
Do not exceed the allowable imbalance. Otherwise, the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- If foreign matter adheres to the rotor (including buckets), wipe it off promptly. Otherwise, the rotor (including buckets) may be corroded.

Balance the tubes that are symmetrically arranged. The allowable imbalance of this rotor is within 1 g/tube.


Set the balanced tubes to the rotor symmetrically.

<Tightening the cap to the bucket>

CAUTION :

Imperfect tightening of the cap may cause imbalance operation and it may happen the serious accident.

Tighten the numbered cap to the same numbered bucket to match each mark of the bucket and the cap.


Tighten the cap until the bottom of the cap bumps to the bucket and the cap must be adjusted to the match mark of the bucket.

<Setting the bucket to the rotor>

CAUTION :

Check that "R25ST" is marked on the buckets, caps, and the rotor before use.


Insert a bar of the bucket to vertical groove of the rotor.

Match the numbered buckets with number of rotor and set the buckets to the rotor.

⚠ WARNING :

- Be sure to set all buckets of the same type whether or not samples are put in buckets: Failure to do so could not only cause the rotor to vibrate, but could result in the rotor being deformed and the buckets becoming detached, which is very dangerous. Never use buckets made by other companies, or any other type of bucket that is not exclusively made for the rotor even if it is made by us.
- Use the buckets (the bucket cap is red) which are exclusively for the R25ST rotor, otherwise the rotor and buckets may be damaged.

<Examples of incorrect bucket setting>

Tighten the cap to the bucket securely. Set the bucket to the rotor securely.

⚠ WARNING :

Before setting a rotor in the centrifuge, make sure that buckets have been securely installed in the rotor: Incorrect installation of buckets in rotor could damage the rotor or detach the buckets, which is very dangerous.

⚠ CAUTION :

Inspect and maintain the rotor after use. If abnormality is observed, do not use it. Contact an authorized sales or service representative.

● Sample limitation

⚠ WARNING :


- Never use any material capable of producing flammable or explosive vapors. Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage.
- When using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.

4. Maintenance

● Maintenance of rotor

⚠ WARNING :
Do not allow the rotor temperature to rise over 100°C.
Otherwise, the rotor can be brittle.

⚠ CAUTION :
Use a neutral detergent having a pH between 5 and 9.
Otherwise, the rotor can be corroded.

Wash the bucket with tap water or a dilute solution of neutral detergent and rinse it out with distilled water. Turn the bucket upside down with the cap removed to dry it well. Check the bucket is completely dried and then put the light coat of silicone grease on the surface of the bucket, the light coat of lubricant for screw on the thread portion, the light coat of silicone grease (vacuum grease) on the O-ring and the gasket. Store them in a dry place. If foreign matter adheres to the bucket, soak the bucket in warm water in one or two hour and remove foreign matter with a soft brush etc. Use a neutral detergent having a pH between 5 and 9. If you can not remove it, contact an authorized sales or service representative.

Wipe the stain off the rotor body and then put the light coat of silicone grease (vacuum grease) on the surface of the rotor body. Store it in a dry place.

Put the light coat of silicone grease inside the drive hole (crown hole) and then wipe it out with a soft cloth.

● Sterilizing rotor

⚠ WARNING :
Never sterilize the rotor(including buckets) by autoclaving or boiling. Otherwise the strength of rotor(including buckets) may be significantly decreased.

Sterilize the rotor according to either gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).

* Do not dip the rotor in the formalin (3%) solution more than 2 hours.

● Inspecting rotor

Check the rotor after every run as follows.

- The thread portion is worn. Contact an authorized sales or service representative.
- The tube cavity is corroded. Contact an authorized sales or service representative.
- The O-ring, the gasket, etc. are deteriorated. Replace it with the new one.
- The adapter is stained. Wipe it out with a soft cloth.

Check the rotor every 100 hours as follows.

- The rotor pin is bent. Contact an authorized sales or service representative.

Check the rotor once a month as follows..

- The inside of drive hole (crown hole) is stained. Wipe it out with a soft cloth dampened with water.

⚠ CAUTION :
Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

5. Using Cesium Chloride

The cesium chloride (CsCl) solution is frequently used as a density gradient medium, but if the CsCl solution with a high density is used, it may saturate during the rotation of the rotor depending on the rotation condition which may cause CsCl crystals to be educed. Educated crystals have a high density (approx. 4g/mL) and apply an excessive load to the rotor which is very dangerous.

The crystallization will also greatly affect the density gradient of the CsCl solution and separation status of the sample. Therefore, always use a solution which will not crystallize the allowable speed.

Fig.5-1 shows the relationship between the maximum allowable density and the speed at which cesium chloride does not crystallize in the condition that standard tubes is filled with the cesium chloride solution. In this diagram, crystallization will occur above the curves.

Be sure to select the density of the cesium chloride solution and the rotor speed at which cesium chloride does not crystallize. The curve in this diagram vary in temperature etc. Before actual operation, centrifuge the cesium chloride solution with the same density to ensure that crystallization will not occur.


Fig.5-1. Relationship between density of cesium chloride and speed

6. Using Density Gradient Centrifugation

● Notes for density gradient centrifugation

Depend on a centrifugal condition, the desired density gradient curve might not be formed.

Before using sample, it is recommended that you check the desired density gradient curve is formed at the intended centrifugal condition.

If the desired density gradient curve is not stably formed, it may be improved by appropriately adding a high density liquid having a higher density than the density gradient medium to the bottom of the tube.

However, adding high density liquid might lower the allowable speed of this rotor (see “● Allowable speed of rotors” in this manual).

7. Packing List

Check the parts and accessories of rotor when receiving the rotor.
Contact an authorized sales or service representative when the need arises.

No.	Name	Qty	Remark	Part No. (Qty)	No.	Name	Qty	Remark	Part No. (Qty)
①	R25ST Rotor	1			⑧				
②	Rotor body	1			⑨	Tool box(C)	1		
③	Bucket body	6*			⑩	Tweezers	1		452701 (1)
④	NBR O-Ring	6*		S401801A (12)	⑪	Bucket stand	1		209641A (1)
⑤	Cap	6*			⑫	Lubricant for screw	1		84810601 (1)
⑥	NBR O-Ring	18	For spare		⑬	Silicone grease (vacuum grease)	1		483719 (1)
⑦	40PA tube	1	50 pcs./set	329607A (50)	⑭	Instruction manual	1		S998655 (1)


*: the quantities of the parts that are mounted to the rotor

When using 33PA Seal Tubes, we recommend the purchase of under-mentioned tube kit, what is the set of tubes and necessary tools.

● 33PA Seal Tube Kit (Part No.S311546A)
(component)

S311544A 33PA Seal Tube : 2 (50 pcs./pack)
S413983 E4-Space Cap : 6
S407157 Tube Setter : 1

※When using 33PA Seal Tubes, Tube Sealer (STF3, STF2) and Tube Rack (E2) (Part No.S206017D) are needed.

8. Rotor Life

While using rotor repeatedly, its strength decreases gradually due to fatigue and creep of material, by which the rotor life is determined.

Sum up the number of runs. if the value reaches the figures shown in Table 8-1, don't use the rotor any more.

Table 8-1. Rotor life

Number of runs	10,000 operations
----------------	-------------------

WARNING:

Care must be taken not to exceed the specified life of a rotor since it could lead to serious accident.

WARNING:

If you found corrosion and some scratches in the rotor and buckets, do not use them. Failure to do so might cause them to be damaged during centrifugation because of decrease in strength.

9. Decontamination

⚠ WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

10. Rotor Retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Retire After Years
R25ST rotor	10

(Appendix) K Factors

• Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

Namely, the time taken for particles of 1S to settle from Rmin. to Rmax. is K factor.

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max.}) - \ln(R_{\min.})}{\omega^2} \cdot \frac{10^{13}}{3600}$$

ω : Angular velocity $\omega = 2\pi/60 \times N$

N : Speed (rpm)

Rmin. : Distance from axis of rotation to top of solution (cm)

Rmax. : Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t(in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 25,000rpm of Model R25ST rotor is 325. The time taken for particles of 100 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{325}{100} = 3.3 \text{ hours}$$

The time taken for settling of particles of 100 S by the same rotor at a speed of 15,000 rpm is:

$$K = \left[\frac{25,000}{15,000} \right]^2 \times 325 = 903 \quad t = \frac{903}{100} = 9.0 \text{ hours}$$

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Copyright © 2021 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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