Keep this rotor instruction manual and the centrifuge manual in the file.

## Angle Rotor for Refrigerated Centrifuge

# R9A2

(Serial No. 4001 or later)

Before using this rotor, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.  
Keep this instruction manual for your reference and refer to it as required.

### Contents

1.	Specifications	_____	1
2.	Bottles	_____	3
3.	How to Use the Rotor	_____	5
4.	How to Use the Bottle Vise	_____	9
5.	Maintenance	_____	9
6.	Packing List	_____	11
7.	Spare Parts of the Rotor	_____	12
8.	Decontamination	_____	14
9.	Rotor Retirement	_____	14
	(Appendix) K Factors	_____	15

## Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.  
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.  
These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

### **WARNING**

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor).  
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density or kinds of tubes/bottles.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

### **CAUTION**

- Do not run the rotor without the rotor cover in position.  
Tighten the cover-locking knob securely.
- Balance the tubes/bottles and the samples including the caps and the adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the rotor tubes and bottles within their actual capacities.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge. The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature.  
Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/bottles if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.  
Contact an authorized sales or service representative.

# 1. Specifications

## 1.1 Construction


Fig.1-1 Construction of the R9A2 Rotor

## 1.2 Specifications

### Applicable centrifuge\*1

(Some centrifuges are on sale in only Japan.  
Visit our web site.)

Max. Speed\*2: 8,500 rpm

Max. RCF: 15,100 × g

Nominal Rotor Capacity:

1,500 mL × 4 bottles = 6,000 mL

Rotor Size: φ 400 × 247 mm

Weight: 17.1 kg

Material: Aluminum alloy

Rotor No.: 88

\*1: For applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges) (Part No. S998611)”.

\*2: The maximum allowable speed varies with centrifuge models and so on. For further details, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges) (Part No. S998611)”.

## 1.3 Cross-Sectional View of the Rotor


Fig.1-2 Cross-Sectional View of the R9A2 Rotor

## 1.4 Characteristics of the Rotor

Table1-1 Characteristics

(when 1500PP bottles (WM) are used)

Rotor Speed (rpm)	RCF (×g)		K Factor*
	Rmin. (5.4cm)	Rmax. (18.6cm)	
1,000	60	210	311,621
2,000	240	830	77,905
3,000	550	1,880	34,625
4,000	970	3,340	19,476
5,000	1,520	5,210	12,465
6,000	2,190	7,510	8,656
7,000	2,990	10,200	6,360
8,000	3,900	13,300	4,869
8,500	4,400	15,100	4,313

\*For details about K factor, refer to “(Appendix) K Factors”

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: speed (rpm)

---

## 1.5 Allowable Speed of the Rotor

**WARNING:**

Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic, kinds of tubes/bottles centrifuge model.  
Do not exceed the allowable speed of the rotor.

The maximum speed marked on the surface of the rotor cover is permitted under the conditions as follows.

-Average sample density: less than 1.2 g/mL

-Standard tubes/ bottles are used.

Depending on sample density, kinds of tubes/bottles and centrifuge models, the rotor speed may be limited as follows. Do not exceed the allowable speed of the rotor.

### (1) Allowable speed for sample density

**WARNING:**

When sample density exceeds 1.2 g/mL, calculate the allowable speed according to the following equation.

$$\text{Allowable speed (rpm)} = 8,500 \text{ (rpm)} \times \sqrt{\frac{1.2 \text{ (g/mL)}}{\text{Average density of sample (g/mL)}}}$$

## 2. Bottles

### 2.1 Applicable Bottles

Table 2-1 Applicable Bottles

Bottles					Cap			Max. speed (rpm)	Max. RCF (xg)
Part No. (Qty)	Name	Size (mm)	Material	Actual capacity (mL/tube)	Part No. (Qty)	Name	Material		
S310087A (4 pcs.)	1500PP Bottle (WM)	126(W) X 175(H)	Polypropylene	1,500	S310071B (2 pcs.)	1500AL Cap	Aluminum alloy	8,500	15,100
					S412932A (4 pcs.)	Cap Support	Polycarbonate including glass fibers		
					S412933A (4 pcs.)	1000 Inner Cap (WM)2	Polypropylene including glass fibers		
					S412935A (10 pcs.)	O-Ring	Silicone rubber		

For further details of the cap, see the section 2.2.

Refer to the provided "Chemical resistance chart (Part No. S999313)", for chemical resistance of each material.

\* Although there are wrinkle patterns on the inside of the 1500PP Bottle (WM), there is no problem. These wrinkle patterns grow during the molding process.


Fig. 2-1 Wrinkle patterns on the inside of the 1500PP Bottle (WM)

### 2.2 Construction of the cap

The construction of the cap is illustrated below.


Fig. 2-2 Construction of the cap

## 2.3 Sterilization and Washing of Bottles and Caps

### ● Sterilization

Sterilize bottles and caps in any following methods; autoclaving(121°C, 20minutes), gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).

Put bottles, whose mouths face upward, on a heat-resistance tube stand when autoclaving bottles.

#### ⚠ CAUTION:

- Do not sterilize the cap support (material: polycarbonate) with organic solvent and alkaline solution. Do not wash the cap support (material: polycarbonate) with organic solvent and alkaline solution.
- Be sure to remove the caps from the bottles before autoclaving. Failure to do so can cause the bottles or the caps to become deformed and damaged.
- Do not place the caps on the mouths of the bottles during autoclaving. Failure to do so may cause the bottles or the caps to become deformed and damaged.
- After autoclaving, wait until the temperature in the autoclaving chamber is the same as the room temperature and then take out the rotor, the rotor cover, and the caps.
- Do not dip the rotor in the formalin (3%) solution more than 2 hours.
- Inspect the O-ring of the cap after sterilization. If the O-ring is deteriorated, it may cause the damaged O-ring during centrifugation. If you observe that the O-ring is damaged and is deteriorated, replace it with the new one.

### ● Washing

Wash the bottles with tap water or a dilute solution of neutral detergent. If the bottles stain heavily, soak them in a dilute solution of neutral detergent or wash them with a soft brush being careful not to scratch them. Then rinse the bottles in tap water and finally rinse them with distilled water. Dry the bottles in the air. Avoid soaking the bottles in a dilute solution of neutral detergent for a long time. Failure to do so may cause the bottles to be degraded in strength. For the caps of the bottles, disassemble caps and then wash, dry and store the caps in the same manner as the bottles.

## 2.4 Inspecting the Bottles and Caps

The bottle, cap support, inner cap, and O-ring are consumable parts. Check these consumable parts every time they are used. If any abnormality such as deformation is found, replace the abnormal part with the new one. Periodical replacement is recommended.


The life expectancy of plastic products, such as the bottle, cap support, inner cap, and O-ring, depends upon the sample, speed of the rotor, temperature, etc. The life expectancy of plastic products when used with aqueous samples (of pH5 to pH9) is specified as follows. (In case they are used repeatedly for 1 hour at the maximum speed of rotor):

- Bottle, cap support, inner cap, O-ring ..... 50 times

However, the life expectancy of plastic products is an approximate guide (see above). We do not warrant the life expectancy of plastic products, such as the bottle, cap support, inner cap, and O-ring.

#### ⚠ CAUTION:

If the upper portion of the 1500PP Bottle (WM) is cracked or deformed, do not use the bottle. (Although the upper portion of the 1500PP Bottle (WM) might whiten depending on sample, speed of the rotor, temperature, etc., there is no problem.)

Check the 1500PP Bottle (WM) before use.  
If the upper portion of the 1500PP Bottle (WM) is deformed excessively, do not use the bottle.

### 3. How to Use the Rotor

#### 3.1 Preparation of the Caps

Push the outer cap on the cap support until you can hear the click sound (see fig. 3-1).

Mount the O-ring in the groove of the inner cap (see fig.3-2)


Fig. 3-1 Preparation of the 1500AL Cap and Cap Support




Fig. 3-2 Preparation of the Inner Cap

#### ⚠ CAUTION:

- Push the outer cap on the cap support so that the flat surface of the cap support faces the outer cap (see below).

- Check that the O-ring is put in the groove of the inner cap before tightening the cap (see below). When tightening the cap, if the O-ring is put in its improper position, the cap can not be tightened securely.

### 3.2 Preparation of the Bottles

(1) Check that the bottles are free from a crack or deformation. Do not use remarkably deformed or cracked bottles.

**⚠ WARNING:**

- This rotor and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

**⚠ CAUTION:**

- Do not use any sample inapplicable to the bottles and the caps in use referring to the provided chemical resistance chart. Failure to do so can cause the bottles and caps to be deteriorated.
- Chemical resistance and the strength of bottles varies with speed, temperature, and so on. Before using sample, fill bottles with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the specified sample capacity of the bottles. Failure to do so can cause sample to leak or can cause the bottles to be deformed.
- Balance the sample within the allowable imbalance including the weight of the tubes/bottles and caps. Excessive imbalance can cause damage to the centrifuge and the rotor.
- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise the precipitation levels may be different by centrifugation and such operation may increase the level of imbalance, and it may result in damage to the centrifuge and the rotor.
- Run the rotor at the proper temperature (from 4 °C to 25°C). Failure to do so may cause the bottles to be damaged or deformed during operation. When running the rotor beyond this limit (under 4°C or above 25°C), before actual operation, test the bottles to ensure that there is no abnormality.
- Tighten the cap of the bottle securely. Failure to do so can cause the sample to leak during operation and can cause the excessive imbalance. As a result, it can cause the centrifuge and the rotor to be damaged.

(2) Fill the bottles with the sample and then balance the bottles that are symmetrically arranged. Do not exceed the specified sample capacity of the bottles. Failure to do so can cause sample to leak or can cause the bottles to be deformed.

The allowable imbalance of this rotor is within 10 g.

### 3.3 Tightening the cap

Tighten the cap with the provided bottle vise and the provided cap opener.

Mount the bottle vise on the desk etc. before use (see Chapter 4).

- (1) Fill the bottle with the sample and then place the inner cap, in which the O-ring is mounted, on the mouth of the bottle gently.
- (2) Screw the 1500AL cap, which is mounted on the cap support, on the bottle lightly.
- (3) Mount the two pins of the cap opener in the holes of the cap and then turn the cap opener clockwise to tighten the cap properly.


Fig. 3-3 Operating Procedures for Tightening the Cap

**CAUTION:**

Tighten the cap securely with the provided cap opener. Failure to do so can cause the sample to leak and can cause the bottle to become deformed.

### 3.4 Preparation of the Rotor

#### (1) Inspection of the rotor

- Check that the drive pins of the drive hole are not deformed.
- Check that the screw portion of the cover-locking knob and the rotor are not worn.

#### (2) Apply the lubricant for screw (Part No.84810601) to the screw portion of the cover-locking knob and the rotor slightly.

### 3.5 Setting the Bottles into the Rotor

Load the bottle into the tube cavity gently. The provided bottle setter is useful to load the bottle into the tube cavity if the sample is heavy.


Fig. 3-4 How to use the bottle setter

**CAUTION:**

- Be sure to support the bottle by your hand when using the bottle setter.
- Do not load only the one bottle or three bottles to this rotor. Otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- Be sure to wipe off the water and the sample which remained into the tube cavities of the rotor before use. Otherwise it may be hard to load the bottle into the tube cavity.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, bottles, and sample are warmed or cooled to the desired temperature before operation.


Fig.3-5 Correct Setting of the Bottles

### 3.6 Installing the Rotor in the Centrifuge

Secure the provided rotor set handle to the rotor center hole.

Hold the rotor set handle by both hands and then mount the rotor onto the drive shaft of the centrifuge gently and securely.

After installing the rotor in the centrifuge, remove the rotor set handle from the rotor.

### 3.7 Mounting the Rotor Cover

**⚠ CAUTION:**

- Before mounting the rotor cover on the rotor, check that the flat surface of the cap support faces the outer cap. If the flat surface of the cap support does not face the outer cap, the rotor and the bottle may be damaged during centrifugation. Before mounting the rotor cover on the rotor, check that the rotor set handle is removed from the rotor.
- Be sure to mount the rotor cover to the rotor body and tighten the cover-locking knob securely before operation, otherwise the rotor or the rotor cover may be removed during operation and result in damage to the centrifuge and the rotor.

(1) Apply silicone grease (vacuum grease) (standard accessory of the centrifuge) to the packing (for the cover) and then fit the packing (for the cover) in the groove of the cover. Replace the packing (for the cover) with new one if the packing (for the cover) is damaged or deteriorates.

If you neglect to apply silicone grease (vacuum grease) to it, it might be hard to mount or remove the rotor cover. When it is hard to mount or remove the rotor cover, apply silicone grease (vacuum grease) to the packing (for the cover).

(2) Apply lubricant for screw (Part No. 84810601) to the threaded portion of the cover-locking knob. Put the rotor cover on the rotor body. Turn the cover-locking knob to screw in the rotor body securely.


Fig.3-6 Construction of the Rotor Cover

### 3.8 Operation

Before the operation, check that the rotor cover is mounted and the cover-locking knob is securely tightened. Perform operation according to the instructions of the centrifuge instruction manual.

### 3.9 Taking out Samples

After the operation, remove the rotor cover by turning the cover-locking knob and then remove the bottles. If it is hard to take out the bottles from the rotor, use the provided bottle setter.

Open the cap with the provided bottle vise and the provided cap opener and then take out the sample in a proper manner.

After removing the bottles from the rotor, secure the provided rotor set handle to the rotor center hole.

Hold the rotor set handle by both hands and then remove the rotor from the drive shaft of the centrifuge gently. After removing the rotor from the centrifuge, remove the rotor set handle from the rotor.

**⚠ CAUTION:**

Be sure to support the bottle by your hand when using the bottle setter. Be careful not to drop the bottle when the surface of the bottle is slippery, because there might be condensation on the surface of the bottle when the sample temperature is lower than the room temperature.

## 4. How to Use the Bottle Vise

- (1) Check the thickness of the desk. If the distance between the bottom of the bottle vise and the handle is shorter than the thickness of the desk, turn the handle counterclockwise (if viewing the handle from the lower side) to make the distance longer between the bottom of the bottle vise and the handle.
- (2) Position the bottle vise on the desk and turn the handle clockwise (if viewing the handle from the lower side) to mount the bottle vise on the desk.

### CAUTION:

- Before use, check that the bottle vise is mounted on the desk etc. securely.
- Be careful not to injure yourself with the edge of the bottle vise. If you do not use the bottle vise for a long time, remove the bottle vise from the desk etc.


Fig. 4-1 Mounting the bottle vise

## 5. Maintenance

### 5.1 Corrosion Resistance of Rotor

### WARNING:

Check the provided chemical resistance chart, and do not use any sample inapplicable to the rotor.

This rotor body is made of aluminum alloy. The rotors made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode the rotor body and decrease the strength of the rotor. Use samples that will not affect the rotor (see the provided chemical resistance).

### 5.2 Maintenance of Rotor and Cap

After the operation, properly maintain the rotor to prevent corrosion that can cause the rotor breakage.

- Normal maintenance

### CAUTION:

- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Do not soak the rotor in a dilute solution of neutral detergent or tap water. Failure to do so can cause the rotor to be damaged due to the remaining liquid inside of the rotor.

After operation, wipe the rotor and the rotor cover with a soft cloth which is dampened with a dilute solution of neutral detergent or tap water and wipe the rotor and the rotor cover with a dry cloth. Dry the rotor and the rotor cover separately.

Check that the rotor and the rotor cover are completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor and the rotor cover. Store the rotor and the rotor cover in a dry place.

- Maintenance after use of a corrosive sample

Wipe the rotor and the rotor cover with a soft cloth which is dampened with tap water soon after operation. Then perform the normal maintenance.

- 
- Maintenance when foreign substances are adhered to the rotor and the rotor cover

 **CAUTION:**  
Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and caps can be discolored or corroded.

Wipe the rotor and the rotor cover with a soft cloth which is dampened with tap water to remove foreign substances.

Be sure to remove foreign substances from the drive hole of the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

### 5.3 Sterilizing Rotor and Cap

Sterilize this rotor and the rotor cover in any following methods; autoclaving(121°C, 20minutes), gas sterilization method (ethylene oxide or formaldehyde), chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin), or ultraviolet rays (200-300nm) sterilization.

-  **WARNING:**
- Be sure to remove the rotor cover from the rotor before autoclaving, otherwise the rotor and the rotor cover may be deformed.
  - Never sterilize the rotor and the rotor cover by boiling, otherwise they will be brittle.

-  **CAUTION**
- After autoclaving, wait until the temperature in the autoclaving chamber is the same as the room temperature and then take out the rotor and the rotor cover.
  - Take out the rotor and the rotor cover carefully from the autoclaving chamber as they are wet and slippery.
  - Inspect the packing (for the cover) after sterilization. If you observe that the packing (for the cover) is damaged and is deteriorated, replace it with the new one.
  - Do not soak the rotor and the rotor cover in chemicals. Failure to do so can cause the rotor to be damaged due to the remaining liquid inside of the rotor.

### 5.4 Inspecting Rotor


Periodically check that the rotor is free from corrosion every 100 hours use. Especially be careful of the tube cavities of the rotor and the drive hole at the rotor bottom because the rotor can be brittle if these portions are corroded. If the rotor surface is discolored, dented or cracked, the rotor is corroded. Do not use such a corroded rotor and contact an authorized sales or service representative immediately for further inspection.

## 6. Packing List







Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	PART No. of Spare Parts (Qty)
①	Rotor	1	—
②	Packing (for the cover)	1*	S411818 (1)
③	1500PP Bottle (WM)	4	S310087A (4)
④	1500AL Cap	4	S310071B (2)
⑤	Cap Support	4	S412932A (4)
⑥	1000 Inner Cap (WM)2	4	S412933A (4)
⑦	O-Ring	4	S412935A (10)
⑧	Bottle Setter	1	S412942A (1)
⑨	Cap Opener	1	S204715 (1)
⑩	Bottle Vise	1	S309911A (1)
⑪	Rotor Set Handle	1	455650B (1)
⑫	Rotor Stand	1	S205816 (1) : For the rotor of which the red adapter is equipped at the bottom
⑬	Lubricant for Screw	1	84810601 (1)
⑭	Rotor Instruction Manual	1	S998279

※: the quantities of the parts that are mounted to the rotor

## 7. Spare Parts of the Rotor

Figure	Part No.	Name	Qty
	S411818	Packing (for the cover)	1
	S310087A	1500PP Bottle (WM)	4
	S310071B	1500AL Cap (The 1500AL Cap does not include the Cap Support, 1000 Inner Cap (WM)2 , and O-ring.)	2
	S412932A	Cap Support	4
	S412933A	1000 Inner Cap (WM)2	4
	S412935A	O-Ring	10
	S412942A	Bottle Setter	1
	S204715	Cap Opener	1
	S309911A	Bottle Vise	1
	455650B	Rotor Set Handle	1
	S205816 (*1)	Rotor Stand	1

(\*1): For the rotor of which the red adapter is equipped at the bottom

---

Figure	Part No.	Name	Qty
	84810601	Lubricant for Screw	1
	483719	Silicone Grease (Vacuum Grease)	1
	S998279	Rotor Instruction Manual	1

## 8. Decontamination

---

**⚠ WARNING:**

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

## 9. Rotor Retirement

---

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

## (Appendix) K Factors

---

- Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit (S).

Namely, the time taken for particles of 1S to settle from Rmin. to Rmax. is K factor.

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R \text{ max. }) - \ln(R \text{ min. })}{\omega^2} \cdot \frac{10^{13}}{3600}$$

$\omega$ : Angular velocity  $\omega = 2\pi/60 \times N$

N: Speed (rpm)

Rmin. : Distance from axis of rotation to top of solution (cm)

Rmax. : Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t (in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 8,500rpm of Model R9A2 rotor is 4313. The time taken for particles of 1000 S to be settled from Rmin. to Rmax. using this rotor is as follows:

$$t = \frac{4313}{1000} = 4.3 \text{ hours}$$

The time taken for settling of particles of 1000 S by the same rotor at a speed of 6,000 rpm is:

$$K = \left( \frac{8,500}{6,000} \right)^2 \times 4313 = 8656 \quad t = \frac{8656}{1000} = 8.7 \text{ hours}$$

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Copyright © 2021 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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