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(URL <https://www.himac-science.jp>)

Keep this rotor instruction manual and the centrifuge manual in the file.

Angle Rotors for Large Capacity Refrigerated Centrifuge

R7A

(Serial No. 4001 or later)

Before using the rotor, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Bottles	_____	3
3. How to use	_____	4
4. How to use the bottle vise	_____	7
5. Maintenance	_____	9
6. Packing list	_____	11
7. Spare parts	_____	12
8. Warranty of the rotor	_____	13
9. Decontamination	_____	14
10. Rotor retirement	_____	14
(Appendix) K Factors	_____	15

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.

These notes are defined as follows.

 **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.

 **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the surface of the cover-locking knob).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density or kinds of bottles.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that toxic or radioactive, or blood samples that pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Do not run the rotor without the rotor cover in position.
Tighten the cover-locking knob securely.
- Balance the tubes/bottles and the samples including the caps or adapters within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use tubes/bottles within their actual capacities.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge. The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature. Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/bottles if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.
Contact an authorized sales or service representative.

1. Specifications

1.1 Construction


Fig. 1-1 Construction of R7A rotor

1.2 Specifications

Applicable centrifuge:	CR7N/CR7and CR30NX*
Max. speed:	7,000 rpm
Max. RCF:	11,100 × g
Rotor No.:	57
Nominal rotor capacity:	1000mL × 6 bottles = 6000mL
Rotor Size:	φ 422 × 207.5 mm
Weight:	19.0 kg
Material:	Aluminum alloy

*: When you use this rotor in CR30NX centrifuges, refer to the "Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))".

Note: The centrifuge which you have been using may require replacement of the ROM (extra cost required) or modification of the centrifuge (extra cost required) to use this rotor. Contact an authorized sales or service representative for further information.

1.3 Cross-sectional view of rotor


Fig. 1-1 Cross-sectional view of R7A rotor

(*1) This bolt is not included in the rotor used for CR30NX, because it is not necessary to fix the rotor to the drive shaft in CR30NX.

1.4 Characteristic of rotor

Table 1-1 Characteristics

Rotor speed (rpm)	RCF (× g)		K factor*
	R = Rmin. (10.33cm)	R = Rmax. (20.34cm)	
1,000	115	227	171,622
2,000	462	910	42,906
3,000	1,040	2,050	19,069
4,000	1,850	3,640	10,726
5,000	2,890	5,690	6,865
6,000	4,160	8,190	4,767
6,500	4,880	9,610	4,062
7,000	5,660	11,100	3,502

*For details about K factor, refer to "(Appendix) K Factors"

Calculation formula

$$RCF = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)
N: speed(rpm)

1.5 Allowable speed of the rotor

 **WARNING :**
 Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic and kinds of tubes/bottles.
 Do not exceed the allowable speed of the rotor.

The maximum speed marked on the surface of the cover-locking knob is permitted under the conditions as follows.

-Average sample density : less than 1.2g/mL

Depending on sample density, kinds of bottles and centrifuge model, the rotor speed may be limited as follows. Do not exceed the allowable speed of the rotor.

(1) Allowable speed for sample density

 **WARNING :**
 When sample density exceed 1.2g/mL, calculate the allowable speed according to the following equation.

$$\text{Allowable speed (rpm)} = 7,000(\text{rpm}) \times \sqrt{\frac{1.2(\text{g/mL})}{\text{Average density of sample}(\text{g/mL})}}$$

(2) Temperature that is capable of cooling the sample depending on the rotor speed

Temperature that is capable of cooling the sample depends on the rotor speed (see table 1-2).

Table 1-2 Allowable rotor speeds and relative centrifugal force (RCF) for Temperature that is capable of cooling the sample

Item	Allowable speed (rpm)				Max.RCF (xg)			
	CR30NX	CR7N	CR7		CR30NX	CR7N	CR7	
Centrifuge								
Temperature	4°C	4°C	10°C	4°C	4°C	4°C	10°C	4°C
Rotor								
R7A	7,000	7,000	7,000	6,500	11,100	11,100	11,100	9,610

Note: The above numerical values are applicable when the room temperature is 25°C and the less.

2. Bottles

2.1 Applicable bottles

Table 2-1 Applicable bottles

Part No.	Name	Size (mm)	Actual capacity	Max. speed	Max. RCF
S307700A	1000 PP bottle (WM)	φ98 × 170	900mL	7,000rpm	11,100 × g
S308357A	1000 PC bottle (WM)	φ98 × 170	900mL	7,000rpm	11,100 × g

Refer to the attached Chemical resistance chart (Part No. S999313) for chemical resistance of each material.

2.2 Cleaning and sterilization of bottles and caps of bottles

Clean and sterilize bottles and caps of bottles according to Table 2-2.

● Autoclaving PP and PC bottles

Wash the bottles well. Put them on a heat-resistant bottle stand facing their mouths upward. Remove the cap of the bottle to prevent deformation or breakage during autoclaving. After autoclaving, wait until the temperature in the autoclaving chamber reduces to the room temperature, then take out the bottles.

	CAUTION: Do not place the cap on the mouth of the bottle. Otherwise the bottle may be deformed or damaged during autoclaving.
---	---

Table 2-2 Washing and sterilization of bottles

○: Applicable x: Inapplicable

Conditions		Material	PC	PP	
Washing	Washing solution	pH 5 or lower	x	x	
		pH 5 to 9	○	○	
		pH 9 or higher	x	○	
		Warm water (50°C or lower)	○	○	
	Ultrasonic washing	Neutral detergent (pH 7)	○	○	
Sterilization	High-pressure steam sterilization(Autoclaving)	115°C(0.7kg/cm ²), 30 minutes	○	○	
		121°C (1.0kg/cm ²), 20 minutes	○	○	
		126°C (1.4kg/cm ²), 15 minutes	x	x	
		Boiling	For 15 to 30 minutes	○	○
		Ultraviolet ray	200-300nm	x	x
		Gas	Ethylene oxide	x	○
	Formaldehyde		○	○	
		Chemical	Ethanol (70%)	x	○
	Hydrogen peroxide (3%)		○	○	
	Formalin (3%)		○	○	

2.3 Life expectancy of plastic bottle and cap of bottle

The life of plastic bottle and cap of bottle depends upon the sample, speed of rotor, temperature, etc. The standard life of plastic bottle when used with aqueous samples (of pH5 to pH9) is specified as follows. (in case they are used repeatedly for 1 hour at the maximum speed of rotor):

- Bottle and cap of bottle(PP,PC) 50 times

However, the standard life of the PC bottle is specified as follows depending upon the pre-treating conditions such as washing, sterilization, etc.:

Table 2-3 Standard Life of tube and the PC bottle

Sample	Washing and sterilization	Sterilization in gas or washing in hot water	Autoclaving at 121°C x 20min
Neutral(pH7)		50 times	10 times
Weak alkali(pH7 to pH9)		30 times	5 times

The standard life of plastic bottles as the above are an approximate guide when they are used repeatedly for 1 hour at the maximum speed of rotor. We do not warrant the life expectancies of bottles.

● Washing :

Wash the bottles with tap water or a dilute solution of neutral detergent. If the bottles stain heavily, soak them in a dilute solution of neutral detergent or wash them with a soft brush being careful not to scratch them. Then rinse the tubes/bottles in tap water and finally rinse them with distilled water. Dry the bottles in the air. Use a detergent of pH 5 to 9 because the PC bottles are susceptible to alkaline solutions. Avoid soaking the bottles in a dilute solution of neutral detergent for a long time. Otherwise, the bottles may be degraded in strength. For the caps of the bottles, remove them from the bottles. Then disassemble and wash them. And dry and store the caps in the same manner as the bottles.

3. How to use

3.1 Preparation of bottles

Check that the bottles are free from crack or deformation. Do not use remarkably deformed or cracked bottles.

⚠ WARNING :

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

⚠ CAUTION:

- Do not use any sample inapplicable to the bottles and the caps in use referring to the attached Chemical resistance chart. Otherwise, the bottles and the caps may be degraded.
- Chemical resistance and the strength of bottles varies with speed, temperature, and so on. Before using sample, fill bottles with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the specified sample capacity of the bottles, otherwise sample leakage or deformation of bottles may occur.
- Balance the sample within the allowable imbalance including the weight of the bottles and caps. Excessive imbalance can cause damage to the centrifuge and the rotor.
- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise the precipitation levels may be different by centrifugation and such operation may increase the level of imbalance. And it may result in damage to the centrifuge and the rotor.
- Run the rotor at the proper temperature (from 4 °C to 25°C). Otherwise the bottles may be damaged or deformed during operation. When running the rotor beyond this limit (under 4°C or above 25°C), test the bottles to ensure that there is no abnormality before actual operation.
- Tighten the cap of the bottle securely, otherwise the sample leakage during operation can cause the excessive imbalance, and it may result in damage to the centrifuge and the rotor.

(1) Inject the sample into the bottles and tighten the cap of the bottle securely. Do not exceed the specified sample capacity of the bottles, otherwise sample leakage may occur during operation.

(2) Balance the bottles that are symmetrically arranged. The allowable imbalance of this rotor is within 6g/bottles.

3.2 Preparation of the rotor

(1) Inspection of the rotor and the tube holder

- Check that the bottom of the tube holder and the drive hole of the rotor are not cracked and are not stained. And check that there is no foreign matter in the bottom of the tube holder and the drive hole of the rotor.
- Check that the drive pins of the drive hole are not deformed.
- Check that the screw portion of the rotor, the cover-locking knob, and the caps are not worn.

(2) Apply the lubricant for screw (Part No.84810601) to the screw portion of the rotor, the cover-locking knob, and the caps slightly.


3.2 Setting of bottles into rotor

⚠ CAUTION:


- Do not load only the one bottle to the rotor or do not load the bottles asymmetrically with respect to the drive shaft. Otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- Before using the tube holder, be sure to wipe off the water and the sample which remained into the tube holder. Otherwise it may be hard to remove the bottle from the tube holder.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, bottles, and sample are warmed or cooled to the desired temperature before operation.

Check that the tube cavities of the tube holder are free from foreign substances. Set the balanced bottles symmetrically with respect to the drive shaft as shown in the figure below.


Center of the drive shaft

Setting of two bottles

Setting of three bottles

Setting of four bottles


Fig. 3-1 Setting of bottles

3.4 Setting the cap to the tube holder

 **CAUTION:**

- Use the caps, whose color is silver, and the rotor, whose serial number is 4001 or later (serial number of the rotor is mentioned on the upper surface of the rotor), otherwise the caps or packing may come off during operation and it may cause the centrifuge and the rotor to be damaged.
- Do not set only the one cap to the tube holder or do not set the caps asymmetrically with respect to the drive shaft. Otherwise the centrifuge or the rotor may be damaged due to excessive imbalance operation.
- Be careful not to get your fingers caught between the caps when setting the cap.
- Be sure to set all caps and tighten them securely, otherwise samples may evaporate during operation.

- (1) Apply the lubricant for screw (Part No.84810601) to the screw portion of the caps.
- (2) Apply the silicone grease (vacuum grease) (the standard accessory of the centrifuge) to the packing and then fit it into the groove of the cap. Replace the packing with the new one if damaged or deteriorated.
- (3) Tighten the cap to the tube holder securely by the cap opener (Part No.S204715)(the standard accessory of this rotor). Tighten the cap to the tube holder until the caps touches the tube holder.

Fit the packing into the groove of the cap so that the protrusion of the packing faces downward

Fig. 3-2 Packing and cap

3.5 Mounting the rotor to the drive shaft

 **CAUTION:**

- Mount the rotor to the drive shaft not to overlap the drive pin of the rotor and the crown pin of the drive shaft. And then tighten the rotor lock bolt to fix the rotor to the drive shaft securely. If the rotor is not mounted properly, the rotor lock bolt can not be tightened sufficiently and the rotor may abnormally vibrate or come off and it may result in damage to the centrifuge and the rotor.

- (1) Mount the rotor without the rotor cover to the drive shaft. Then insert the rotor lock bolt into the center hole of the rotor and tighten the rotor lock bolt securely by the hexagonal wrench (the hexagonal wrench is the accessory of the rotor).

Above work is not necessary when you use this rotor in CR30NX.

3.6 Mounting the rotor cover

CAUTION:

- Be sure to mount the rotor cover to the rotor body. If the rotor without the rotor cover is rotated, it may result in the following inconveniences.
 - The loss of sample may be caused by the rotor temperature going up due to the friction between the rotor and the air.
 - The rotor speed may not reach the set speed. It results in the run stoppage of the centrifuge due to displaying the alarm code.
- Tighten the cover- locking knob securely before operation, otherwise the rotor cover may come off during operation and it may result in damage to the centrifuge and the rotor.

(1) Apply the lubricant for screw (Part No.84810601) to the screw portion of the cover-locking knob.

(2) Mount the rotor cover to the rotor body and tighten the cover-locking knob securely.

3.7 Operation

(1) Perform operation according to the instructions of the centrifuge instruction manual.

3.8 Taking out samples

CAUTION:

- Be careful not to get your fingers caught between the caps when removing the cap.
- Remove the bottle holding the rotor. If there is water or the sample in the tube holder, it may be hard to remove the bottle from the tube holder.

(1) After the operation, gently take out the rotor from the centrifuge and remove the rotor cover by turning the cover-locking knob.

(2) Remove the cap from the tube holder and then remove the bottles. Take out the sample in a proper manner. If it is hard to remove the cap or bottle from the tube holder, use the attached cap opener (Part No. S204715) (See Fig. 3-3). Mount small pin of the cap opener in the notch of the cap and then turn the cap opener counterclockwise. (When opening the cap of the bottle, use the large pins of the cap opener.)



Fig. 3-3 How to use the cap opener

4. How to use the bottle vise

- The tube vise is not included in this rotor. If necessary, please order the tube vise (Part No. S205037A).

4.1 Mounting the bottle vise

- (1) Check the thickness of the desk. If the distance between the bottom of the bottle vise and the handle (B) is shorter than the thickness of the desk, turn the handle (B) counterclockwise to make the distance longer between the bottom of the bottle vise and the handle (B).
- (2) Position the bottle vise on the desk and turn the handle (B) clockwise to mount the bottle vise on the desk.




4.2 Tightening the bottle

CAUTION:

- Do not tighten the handle (A) too tightly. Failure to do so might cause the bottle to be deformed.


- (1) Position the bottle on the center of the bottle vise.
- (2) If the bottle touches the fixing plate, turn the handle (A) counterclockwise to move the fixing plate.
- (3) Push the bottle against the fixing plate of the bottle vise and turn the handle (A) counterclockwise to move the fixing plate until the fixing plate touches the bottle lightly.
- (4) Tighten the handle (A) to the 90° to fix the bottle on the bottle vise.

4.3 Tightening the cap

Tighten the cap with the provided cap opener.

Mount the two pins of the cap opener in the holes of the cap and then turn the cap opener clockwise.

5. Maintenance

5.1 Corrosion resistance of rotor and cap

 **WARNING:**
Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor and the caps.

This rotor body and caps is made of aluminum alloy . The rotors made of aluminum alloy have high corrosion resistance and they are covered with an anodic oxidation coating. However, use of inapplicable chemicals can corrode the rotor body and caps and decrease the strength of them. Use samples that will not affect the rotor and caps referring to the attached chemical resistance chart.

5.2 Maintenance of rotor and cap

After the operation, properly maintain the rotor and cap to prevent corrosion that can cause the rotor and caps breakage.

●Normal maintenance

 **CAUTION:**
Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor and caps with tap water or a dilute solution of neutral detergent and rinse them out with distilled water. Wipe the rotor and caps with a soft cloth. Turn the rotor upside down with the rotor cover and caps removed to dry it well. Check that the rotor and caps is completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor and caps. Store the rotor and caps in a dry place.

●Maintenance after use of a corrosive sample

Wash the rotor and caps with tap water immediately after the operation. Then perform the normal maintenance.

●Maintenance when foreign substances are adhered to the rotor and caps

 **CAUTION:**
Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and caps can be discolored or corroded.

Soak the rotor and caps in warm water for one or two hours and wash the inside of the tube cavities and the drive hole at the rotor bottom with a soft brush to remove foreign substances. Do not forget to remove foreign substances adhered in the drive hole at the rotor bottom, otherwise the rotor and the drive shaft may be damaged.

5.3 Sterilizing rotor and cap

Sterilize this rotor and caps in any following methods; autoclaving(120 °C , 20minutes), gas sterilization method (ethylene oxide or formaldehyde) or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).

 **WARNING:** Be sure to remove the rotor cover and the caps from the rotor before autoclaving, otherwise the rotor, the rotor cover or the caps may be deformed.

 **WARNING:** Never sterilize the rotor, the rotor cover and the cap by boiling, otherwise they will be brittle.

 **CAUTION :** After autoclaving, wait until the temperature in the autoclaving chamber reduces to the room temperature, then take out the rotor, the rotor cover, and the caps.

 **CAUTION :** Take out the rotor , the rotor cover and the caps carefully from the autoclaving chamber as they are wet and slippery.

 **CAUTION :** Inspect the packing of the cap. If the packing is deteriorated, it may cause the damaged packing and rotor seal can not work properly when the sample leaks from the bottle.

 **CAUTION :** Do not dip the rotor in the formalin (3%) solution more than 2 hours.

5.4 Inspecting rotor

Periodically check that the rotor is free from corrosion every 100 hours use. Especially be careful of the tube cavities of the rotor and the drive hole at the rotor bottom because the rotor can be brittle if these portions are corroded. If the rotor surface is discolored, dented or cracked, the rotor is corroded. Do not use such a corroded rotor and contact an authorized sales or service representative immediately for further inspection.


6. Packing list

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.



No.	Name	Qty	Part No. of spare parts (Qty)
①	Rotor	1	—
②	Cap	6* ¹	S414217A (2)
③	Packing	6* ¹	S409741A (2)
④	Rotor lock bolt* ²	1* ¹	84551409 (1)
⑤	Cap opener	1	S204715 (1)
⑥	1000 cap (WM) ass'y	6	S410870A (2)
⑦	1000 bottle (WM) ass'y	6	S307700A (6)
⑧	Bottle setter	1	84550301 (1)
⑨	Hexagonal wrench* ²	1	84850702 (1)
⑩	Lubricant for screw	1	84810601 (1)
⑪	Rotor stand	1	S205816 (1) : For the rotor of which the red adapter is equipped at the bottom
⑫	Rotor instruction manual	1	S999892

*1: the quantities of the parts that are mounted to the rotor.

*2: Rotor for CR30NX does not include them.

7. Spare parts

	Part No.	Name	Qty
	S414217A	1000AL cap (WM) ass'y2	2
	S409741A	Packing ass'y	2
	84810601	Lubricant for screw	1
	483719	Silicone grease (vacuum grease)	1
	84551409	Rotor lock bolt (*2)	1
	S204715	Cap opener	1
	84550301	Bottle setter	1
	84850702	Hexagonal wrench (*2)	1
	S205816 (*1)	Rotor stand	1
	S999892	Rotor Instruction manual	1

(*1): For the rotor of which the red adapter is equipped at the bottom

(*2): They are not required for the rotor which is used in CR30NX.

8. Warranty of the rotor

● Warranty for Refrigerated Centrifuge Rotors

Eppendorf Himaс Technologies warrants its refrigerated centrifuge rotors under the terms and conditions set forth below.

1. Scope of Application

This warranty applies to:

- (1) failure or damage to the rotor due to rotor defects, in materials or workmanship, that are attributable to Eppendorf Himaс Technologies, and
- (2) failure or damage to an Eppendorf Himaс Technologies refrigerated centrifuge due to defects in a rotor covered by the warranty.

Note: This warranty covers the rotor or refrigerated centrifuge only, and Eppendorf Himaс Technologies shall not be liable for damage to samples or other items that may arise from failure of or damage to the rotor or refrigerated centrifuge.

2. Warranty Period

This warranty is valid for seven years from the date of delivery of the rotor by Eppendorf Himaс Technologies or an Eppendorf Himaс Technologies dealer.

3. Notes

- (1) This warranty extends only to the original buyer from Eppendorf Himaс Technologies or the original buyer from an Eppendorf Himaс Technologies dealer. Rotors whose ownership passes from the original buyer to a third person due to resale, transfer or the like are not covered by this warranty. Neither does this warranty cover cases where the ownership of a rotor is reacquired subsequently by the original buyer.
- (2) The warranty for the rotor proper (warranty period and terms) shall continue to apply in cases where rotor parts and the like are additionally purchased.
In such cases, the warranty period shall continue to be measured from the date of delivery of the rotor. If such parts are subject to a restricted warranty period, this warranty shall only cover failure or damage involving such parts that occurs both within the rotor warranty period and within the warranty period for the parts concerned.
- (3) This warranty does not cover any of the following, even during the warranty period:
 - ① Failure or damage due to abuse or misuse.
 - ② Failure or damage due to operation or maintenance in a manner contrary to the instructions in the operator's manuals for the rotor and refrigerated centrifuge.
 - ③ Failure or damage due to disassembly or modification without the permission of Eppendorf Himaс Technologies.
 - ④ Failure or damage due to use of a bucket, adapter, tube, bottle, or the like, other than those approved by Eppendorf Himaс Technologies for use with the product.
 - ⑤ Failure or damage due to disaster such as fire or earthquake.
 - ⑥ Failure or damage due to use with a refrigerated centrifuge not of Eppendorf Himaс Technologies manufacture without the permission of Eppendorf Himaс Technologies.
 - ⑦ Failure or damage due to use with a refrigerated centrifuge that has been modified without the permission of Eppendorf Himaс Technologies.
 - ⑧ Failure or damage due to a failed or damaged part (including a tube, bottle, adapter, etc) not submitted to Eppendorf Himaс Technologies.
 - ⑨ Replacement of consumable items.
 - ⑩ Parts subject to a restricted warranty period that has expired.

● Other Matters

Rotors that are delivered to Eppendorf Himaс Technologies for inspection or repair must be sterilized and cleaned beforehand.

Be aware that Eppendorf Himaс Technologies will sometimes decline to carry out inspection or repair if a rotor or other related equipment is or could be contaminated as a result of use without protection or under inadequate safety management in a radioisotope facility or a laboratory of the level P2 or higher.

9. Decontamination

 **WARNING:**

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

10. Rotor Retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure. Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

(Appendix) K Factors

- Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

Namely, the time taken for particles of 1S to settle from Rmin. to Rmax. is K factor.

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max.}) - \ln(R_{\min.})}{\omega^2} \cdot \frac{10^{13}}{3600}$$

ω : Angular velocity $\omega = 2\pi/60 \times N$

N : Speed (rpm)

Rmin. : Distance from axis of rotation to top of solution (cm)

Rmax. : Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t(in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, the K factor at the maximum speed of 7,000rpm of Model R7A rotor is 3,502. The time taken for particles of 3000 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{3,502}{3,000} = 1.2 \text{ hours}$$

The time taken for settling of particles of 3000 S by the same rotor at a speed of 5,000 rpm is:

$$K = \left[\frac{7,000}{5,000} \right]^2 \times 3,502 = 6,864 \quad t = \frac{6,864}{3,000} = 2.3 \text{ hours}$$

Eppendorf Himac Technologies Co., Ltd.

1060, Takeda, Hitachinaka City
Ibaraki Pref., 312-8502 Japan

URL: <https://www.himac-science.com>

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Copyright© 2020 Eppendorf Himaac Technologies Co., Ltd. All rights reserved.