eppendorf



Peak of Perfection

Made for a faster and easier daily routine: Eppendorf Liquid Handling Consumables

Our Benchmark: Your Requirements

How liquid handling consumables can influence your assay results

The reliability and consistency of your research results are our priority when developing consumables.

Chemicals like slip agents, plasticizers or biocides, used as manufacturing additives, can leach out of the plastic into the sample and substantially inhibit enzymatic assays and binding studies. As described in recent publications these chemicals, such as slip agents (e.g. oleamide), can slow down evaporation, skew absorbance readings and lead to erroneous DNA quantification. Some of these slip agents have also been shown to negatively affect the outcome of biological tests like enzyme activity or receptor-binding assays. Providing excellence for scientific experts ensures unaffected and reproducible results.

To achieve highest reliability and consistency for your experiments, Eppendorf has optimized materials and processes to minimize the risk of interference.

- > Eppendorf liquid handling consumables are made of highest quality, virgin polypropylene (PP) free of biocides, plasticizers and latex
- > Optimized, highly polished molds produced without the use of slip agents like oleamide, erucamide, stearamide
- > Used dyes are free of organic additives and heavy metals



»We have seen substantial inhibition of our enzyme assays by chemicals leaching from disposable plastic consumables. To obtain the best possible reproducibility we use consumables from manufacturers that can confirm the absence of critical manufacturing additives.«

Dr. Andrew Holt

Department of Pharmacology, University of Alberta, Canada



»We need to avoid that contaminants from the plastic material enter the sample and inhibit bacterial growth. The consumables that we use to analyse water samples should be of the highest purity to obtain reliable results.«

Karen Thomsen

Mikrobiologie-Zentrallabor, Hamburg Wasser GmbH, Germany



»Our DNA isolation protocols from both animal and plant material require grinding of tissue prior to and during the cell lysis process. As we work with high numbers of samples, breaking of consumables and subsequent sample loss can be critical. The excellent quality and stability of the 1 mL Deepwell Plates from Eppendorf convinced us as it improved the reliability of our process significantly.«

Dr. Paul Gooding

Plant Genomics Centre, Australian Genome Research Facility

Premium Quality Is Our Standard

Certified quality and purity

Continuous quality assurance throughout the entire production process – from the initial material to the finished product.

- > Manufactured from carefully selected, purest raw materials which comply with international purity criteria (FDA 21 CFR§177.1520 »Olefin Polymers«, FDA 21 CFR§178.2010 »Antioxidants and Stabilizers for Polymers«)
- > Fully automatic production under clean room conditions according to VDI 2083 class 6 and to U.S. Fed.Stand. 209D class 100.000: continuous bioburden and particle monitoring of the plant for impeccable production environment
- > SOP present for storage and every production step elimination of human error
- > Frequent production tool checks plus regular in-process quality and functional checks for consistent quality guarantee and absolute reproducibility
- > Full traceability for each product full control, each product to be traced back to material lot

- > Proof of compliance with requirements of standards, guidelines and regulations that apply to biological and industrial laboratories:
 - > Lot-specific purity certificates document testing by an independent lab for products with the purity grades
 - > Eppendorf Biopur®
 - > PCR clean
 - > Sterile
 - > PCR clean and Sterile
 - available at www.eppendorf.com/certificates
 - > General quality certificates as e.g.
 - > Certificates of Purity for PCR clean and Eppendorf Biopur products
 - > Certificates of Quality e.g. Statement on BSE/TSE
 - > Product specific certificates
 - > For trace metal
 - > For absence of surface active additives
 - > With special focus, e.g. filter efficiency for ep Dualfilter T.I.P.S.[®]
 - > Certificates verifying Quality Management System/ compliance with standards

Eppendorf Certificate eppendorf	Certificate of Purity	Eppendorf Certificate
Certificate of Purity - BIOPUR®	– Eppendorf PCR Clean & Sterile –	Eppendorr certificate
This package contains a high-quality conveniable manufactured under the Biopur Epismolof Pavity Standard, The Episendorf Biopur consumables are produced in a class 6	LADR GeldH Kecksnische Versogragszerem Brenen Besch Lossenzhler, Futernitze & Limavitanip(tk	Certificate of Quality
taccording to VDI 2053) and a class 100,000 faccording to U.S. Fed. Stand. 209 DI clean toom environment.	Lating Opt-Bot Thomas Veriget	Combitips advanced* Typical values for trace metal(s) The values in the table indicate typical values of trace metal concentrations which are obtained after involvating Combitips advanced with core, nitric acid for 10
For this product Eppendarf certifies the following: > Sterilie		minutes (see: Materials and Methods).
> store > human DNA-free > factorial DNA-free biopur	ep Dualifiter T.I.P.S.® 0.1 - 5 mL	As the indicated values were determined in a one-time measurement, they cannot be puramited for every lot of Combition advanced. Rather, they give an idea to what extent trace elements can be eluted from Combition advanced.
> Riverse free	■ 2023-06-28 STEMLE N	Texts Initial Material (1994)
> PCR inhibiter free > ATP-free	C Chater no.: 0000 277.565 Cdt. Anno. 002047047	11 Fair 1
These parameters are continuously monitored by an independent certified laboratory. Eggendert guarantees the conformity within the fullowing limits:	V C Berner Berner Berner ton Bernert ton	0.2 mi 2.00739 +0.00002 +0.00002 +0.0010 +0.0000 +0.0000 +0.00000 +0.00000 +0.00000 +0.00000 +0.00000 +0.00000
Sterlity in accontance with USP, Ph. Eur. 2.6.12		0.5 mi 0.5020 -0.0000 -0.0000 -0.0000 -0.0000 -0.0000 -0.0000 -0.0000 -0.0000 -0.0000
Peropersi (Endotoxini) < 0.001 EU/nit tlested according to Pic Eur. 2.6.14 (LAL text) Human DNA <2 pg		1.0 ee 0.001700 e0.00000 e0.00000 e0.00010 e0.0001 e0.0000 e0.00000 e0.00000 e0.00000 e0.00000 e0.00000
Bacterial DNA < 50 fg DNase not detectable (LOD: 1.0 x 10* Kunitz units)	* Nutrian EMA / Nutriana EMA * Genetiky / Santiky / San	2 5 mi 2 00107 to 00000 +0 00000 +0 00000 +0 00000 +0 00000 +0 00000 +0 00000 +0 00000
RNase ovit detectable (LOD: 1.0 x 10° Kurvitz units) PCR inhibition less than 10 targets amplifiable	Instructure fragmeneran Danas Rusadom Marter Calvebream Elevation accosting (2014) 2013 2019 Larar / Calvester - A for larar Rusa real forman all / Larar / Calvester - A forman all / * Desa	1 0 0 0 0 0 0000 0 0 00000 0 0 0000 0 0 0
ATP < 5.5 x 10 ⁻⁹ mg	Katrice / Institutes DNA Capacitor / POA Vortes Lind / Granzowsh not detected / robot randowshar sold COL Physics of Galaxies for errors ROD / Nyohwangsold tab you be	12 H4 E 0000 (4 0000) +0.0000 +0.0000 +0.0000 +0.0000 +0.0000 +0.0000 +0.0000 +0.0000 +0.0000
Lot specific certificates can be downloaded from the loternet at	* Status Biolizationerstatistic Control / Poli Vertes Noticel / Mathematication / Poli Vertes	28
www.appendorf.com/certificates.	Law / Georgeweit on Joshn recharacióne UDP, Ph. Eur. 2 6:13 0.007 / Hardingsment 1.2 x 10 ⁴ functionalita · · · · · · · · · · · · · · · · · · ·	58 est 1 00024* =0.0002 +0.0002 +0.0001 +0.0007 +0.0005 =0.00005 +0.00000 +0.00000
March, 2016 March 2016	* PCR Industries PCR Initializes Initial / Polyandrates: Initial	A CA C Co Ha Ma Ha Pa Sa
1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Level / Oversevent Level than 10 largets amplitude / (Ph. Go. 2.6.14) weaking at als 10 amplitude are used of the 2.6.140 Targets Level / Genevater 4.4.004 EE (BJL 1.E.1.3U.) INL	100 2 00000 200000 1000 1000 0000 0000
paid (by as high set	Yanging of the above described bit should conformity within the limits of detectors. The foll is retensed for vee.	
Barret Vors Dr. B. Schweiter	Die Printung des ober ausgewassenen Litte ergeb eine Einhaltung der Greinweite. Das Lati ist hierritt heipegeben.	
Policet Line Manager Vice President Quility Management B - Commissions Regulatory Affairs	15. St. 18 GPD ABOS JSBC The Willing Weight	Page 1 of 2
Yeor local dambets: www.spendarf.com/contact Sponorth 2 22311 Interburg Generativ generatingegeotaticsm.	The last reach well-well-well-show to be here beind. The certificate must not be opportunity without the toported of LAUK Gradies. To the statement on a manufaction and de Protoportunite. One Zertificat del descapamente richt ohne	
Egyanoloff, an Egyanolof lago of Egyanf ao regalement subsystem of Egyanolof Ad, Sarnane CS, George Fanna ao subsiste en enteringementen antipe Al como ference due systemis en accente European Unite the Proceeding Adu.	edottische Genetringung der LADR Getet vervalltigt werten	www.sppendorf.com
www.eppendorf.com	LACH Gedel BVZ Bremen, Beruch Laberanitel: & Unweitensive (Freekstr-Kanille: 22 (D-2005 Bremen) Gelmany Proze: +61-421 4201 500 (Faz: +69 421 4307:199 (mobile/activ/sign/@aboreantrum-bremen.de (move laboreantrum-bremen.de)	Contract Advances of the contract of the

> All certificates may be downloaded on www.eppendorf.com/certificates

A System You Can Rely On

Tailored to your application needs

Eppendorf has set industry standards in consumable purity levels. Building on the famous Eppendorf Quality, four additional purity grades tailor Eppendorf Consumables to various applicational needs – Sterile, PCR clean, Forensic DNA Grade and Biopur.

You have highest expectations on consistency and reproducibility? Make Eppendorf consumables your choice!

Paving the way for modern lab standards

All consumables are subject to internal process controls on function, tightness, precision, transparency, low wetting property and high chemical and thermal resistance. Additionally, all Sterile, PCR clean, Biopur, and Forensic DNA Grade products are tested by an independent, external laboratory for compliance. The lot-specific certificates issued by this laboratory may be downloaded here: www.eppendorf.com/certificates



- > DNase-free
- > RNase-free
- > PCR inhibitor-free
- > ATP-free

biopu

- > Pyrogen-free
- > Sterile

Example applications:

For highest purity demands in cell culture, nucleic acid analytics, or hygiene monitoring



	Sterile	PCR clean	PCR clean and Sterile	Forensic DNA Grade*	Biopur*
Lot testing (certified) for the following purity criteria					
Pyrogen-free (endotoxin-free)	•		•		•
Sterile (Ph. Eur./USP)	•		•		•
Human DNA-free	·	•	•	•	•
DNA-free (Human- and bacterial DNA)	·				•
DNase-free		•	•	•	•
RNase-free		•	•	•	•
PCR inhibitor-free		•	•	•	•
ATP-free					•
Methods (examples)	·				-
Bacteria and yeast culture			√		√
Cell and tissue culture	√		√		~~
Isolation and storage of DNA		 √√	√		√
Isolation and storage of RNA		√	√		11
DNA analysis (PCR, qPCR, restriction analysis,		<u></u>	√		√
hybridization, microarrays, sequencing)					
Mitochondrial DNA analysis				√√	√ √
Bacterial DNA analysis					√ √
RNA analysis		√	√	√	√ √
Application Areas (Examples)					
Molecular biology	√	 √√	√	 √√	√
Microbiology	√		√		√
Cell biology:					
> Stem cell research			\checkmark		$\checkmark\checkmark$
> Transgenic animals/plants					
Research:					
> Medical			\checkmark		$\checkmark\checkmark$
> Agriculture and aquaculture					
Quality control:					
> Food and beverage			√		~~
> Water supply					
> Environmental monitoring					
Forensic		√	√		√ √

✓ Recommended, ✓✓ Highly recommended

* Increased safety due to individually packaged/single-blistered products.

The Perfect Fit – epT.I.P.S.®

Each of your valuable samples deserve best treatment. See for yourself how Eppendorf pipette tips will save time and reduce costs.

With respect to material, fit, design and operating forces our pipette tips set new standards. The close environment of each sample should be adapted to its specific quality and purity needs. This can involve a specific purity level or the absence of certain substances, but also stability, reliability or geometry. The epT.I.P.S. pipette tips from Eppendorf are designed to cover all of the specific needs of your samples.

Our epT.I.P.S. – Eppendorf Totally Integrated Pipetting System – have been developed to work in perfect combination with our pipettes. This results in reduced tip attachment and ejection forces while maintaining a complete seal between tip and pipette. Additionally, the universal nose cone design allows the tips to be used with pipettes from other manufacturers. In the Eppendorf owned production facilities in the north of Germany, we can maintain the highest standards in the selection and processing of plastic materials. Ideal wetting properties, high transparency, and special certified purity levels are the visible expression of this production philosophy.

Every tip is manufactured to match specifically to its respective Eppendorf pipettes, thus ensuring the maximum precision and reliability you have come to expect and to rely on.

Each pipette tip together with the pipette forms a complete coordinated system. In our current Application Note 354 "The tip of the iceberg – How pipette tips influence results" we have compiled research results to provide you with comprehensive information about erroneous test results caused by using an uncoordinated pipetting system.

Learn more at: www.eppendorf.com/tip-quality

How can you recognize original Eppendorf pipette tips? The unmistakable feature of our tips is the "ring" of drop-shaped relief elements at the edge of the tip crown. Additionally, every tip has the "Eppendorf" name embossed on the edge.



epT.I.P.S.[®] Racks

- > Ideal for when absolute purity is needed available in Sterile/Pyrogen free and Eppendorf Biopur
- > Eppendorf Biopur[®] pipette tips provide maximum biological purity. Guaranteed PCR clean, sterile, pyrogen-free, ATP-free, bacterial DNA-free, they meet the most stringent requirements of the medical, pharmaceutical and foodstuffs industry as well as those of molecular biology and cell biology
- > Lot-specific certificates issued by an independent laboratory are available on the Internet at www.eppendorf.com/certificates
- > Packaged in racks of 96, 48 or 24 tips
- > Batch number and expiration date on each rack label



- epT.I.P.S.[®] Standard
- > Original, high-quality
 Eppendorf pipette tips
 packaged in resealable
 bags
- > Available in sizes from 10 µL to 10 mL
- > epT.I.P.S. Standard pipette tips can be autoclaved at 121 °C in 20 min.



epT.I.P.S.[®] Box/ epT.I.P.S.[®] Set > Contact-free transfer of

- trays to the working box > System optimized for
- use with multichannel pipettes
- > Color-coded trays for simple identification of tips and matching Eppendorf pipette
- > Reloads require less storage space than racks or boxes
- > Reload trays and epT.I.P.S. boxes are entirely autoclavable for future use at 121 °C



epT.I.P.S.® Reloads

- > Reduced waste as compared with disposable racks
- > The Reload system, depending on tip size, is
- packaged as either dualsided or in stack form > Reloads are available in
- two purity levels: Eppendorf Quality and PCR clean
- > Easy one handed reloading under clean bench



epT.I.P.S.[®] Singles

- > Individually wrapped pipette tips for the strictest lab environments when absolute purity must be guaranteed until each tip is ready for use in Eppendorf Biopur: guaranteed sterile, RNase-, DNA-, ATP- and pyrogen-free
- > Batch number and expiration date printed on each blisterpack
- > Continuous quality control of each batch by an independent laboratory – Batch-specific purity certificates available on www.eppendorf.com/ certificates

Long-Distance Tips

Extended length for reliable sampling

epT.I.P.S. L pipette tips enable you to reliably reach your sample while pipetting from and into conical tubes, high reagent bottles, narrow deep vessels, cell culture flasks or deepwell plates. The long and slim design of these elongated pipette tips gives you free access to your sample with reduced risk of touching the sides of tubes or wells. Cross-contamination during pipetting can be reduced to a minimum.

Product features and benefits

- > Safe sample access to Eppendorf 5 mL tubes, 15 mL conical tubes, test tubes, cell culture flasks, deepwell plates and other deep vessels
- > Highest protection for pipette and sample with ep Dualfilter T.I.P.S L
- > Available in the purity standards Eppendorf Quality, PCR clean, Eppendorf Biopur and as ep Dualfilter T.I.P.S. in PCR clean/Sterile
- $> 0.5 20 \mu$ L L, 46 mm for 0.2, 0.5, 1.5, and 2.0 mL tubes or plates
- > 50 1,250 μL L,103 mm for 1 mL deepwell plates and with multi-channel pipette
- > 0.2 5 mL L,175 mm for 15 mL, 50 mL conical tubes
- > 0.5 10 mL L, 243 mm for 75 cm² cell culture flasks,
 - 1 L reagent bottles, Erlenmeyer flasks, measuring cylinders



> epT.I.P.S. L are available in volumes of: 0.5–20 μL L 46 mm long, 50–1,250 μL L 103 mm long, 0.2–5 mL L 175 mm long and 0.5–10 mL L 243 mm long

ep Dualfilter T.I.P.S.®

Two filter layers are better than one

Eppendorf ep Dualfilter T.I.P.S. are the first filter tips with a two-phase filter for contamination protection. The unmistakable blue and white filter layers are made of flexible, hydrophobic material to fit perfectly in the tip cone and retain practically 100 % of all aerosols* and biomolecules. This unique filtering effect is achieved using various well-defined pore sizes in the two filter layers.

The white layer that faces the sample retains drops, splashes and aerosols. The blue layer forms a highly-effective double barrier that reliably binds biomolecules.

* An aerosol is a colloid of fine solid particles or liquid droplets, in air or another gas

The air flow rate through the filter is the same as with singlelayer filters – it simply makes your rapid pipetting tasks much safer. The ep Dualfilter T.I.P.S. hydrophobic filter surface also enables easier and complete recovery of samples.

Product features und benefits

- > Two filter layers provide dual protection
- > Dual protection for pipette and sample
- > Dual protection from aerosols and biomolecules
- > Free of PCR inhibitor additives
- > Eppendorf PCR clean, Sterile and pyrogen-free
- > Batch certified

Applications

- > DNA applications (e.g. PCR)
- > RNA applications (e.g. Gene expression analysis)
- > Protein applications (e.g. Antibody Research)
- > Cell Culture applications (e.g. Media)
- > Applications with radioactive substances

Dual protection against contamination right from the start.

ep Dualfilter T.I.P.S. are manufactured to the highest possible quality from pure, non-recycled materials under cleanroom conditions. They are characterized by defined flow dynamics, low wettability and high thermal stability.

ep Dualfilter T.I.P.S. are sterilized by electron beams and certified pyrogen-free and PCR clean (free from human DNA, DNase, RNase and PCR-inhibiting substances).



Aerosols are formed during the movement of liquids. Without a filter **1** the pipette is exposed to contamination by samples and aerosols. Conventional single-layer filters **2** do not fully block particles and molecules. Only ep Dualfilter T.I.P.S. **3** provide reliable protection even against the finest impurities.

Scientifically Proven

Excellent protection with ep Dualfilter T.I.P.S.®

For the introduction of the ep Dualfilter T.I.P.S.® filter tips in 2006, filter tips from various manufacturers were tested in a blind trial at the Fraunhofer Institute for Toxicology and Experimental Medicine (ITEM) in Hanover, Germany. Testing focused on the prevention of tip cone contamination by aerosols. Efficacy of the filters from various manufacturers to protect against salt aerosols and biomolecules (DNA fragments) was examined. Quantitative determination of DNA fragments was performed using real-time PCR.

The results indicate that ep Dualfilter T.I.P.S. filter tips are much better at protecting against aerosols than the other filter tips tested during the trial. With regard to particle permeability, ep Dualfilter T.I.P.S. are 55 to 677 times more effective. ep Dualfilter T.I.P.S. were 21 to 600 times more effective when protection against DNA fragments was examined for the first time. Both filters, ep Dualfilter T.I.P.S. and ep Dualfilter T.I.P.S. SealMax, are made of the same Polyethylene (PE) and have comparable pore size structure. Thus a further study with ep Dualfilter T.I.P.S. SealMax was conducted based on the standard EN 1822. This standard deals with filtration performance testing of filters as used, for example, for applications in clean room technology or pharmaceutical industry. The results of the filter efficiency testing done by a certified institute showed a minimum particle collection efficiency of the ep Dualfilter T.I.P.S. SealMax of 99.5 % with NaCl aerosol particle sizes of $0.05-0.5 \mu m$.



Performance of ep Dualfilter T.I.P.S. SealMax at MPPS (most penetrating particle size). Measurements performed according to EN 1822.

Filters retain different particle sizes with different performance. The MPPS displays the lowest performance. The labels on the right side of the graphic refer to different classes of MPPS within the standard EN 1822. All ep Dualfilter T.I.P.S. comply with class E12, whereas the most used ep Dualfilter T.I.P.S. variant, the 1,000 μ L tip, even complies with class H13. This means that it catches at least 99.98 % of the particles.

Don't Panic

Your pipette is safe with ep Dualfilter T.I.P.S. SealMax

Eppendorf's ep Dualfilter T.I.P.S. SealMax filter tips protect your pipette against both aerosols and liquids. When it comes to an accidental over-pipetting situation, the new violet layer of ep Dualfilter T.I.P.S. SealMax becomes a reliable barrier against sample liquid – no liquid will pass through the filter. Your pipette is safe at all times! Additionally, virtually 100 % of aerosols* and biomolecules are retained while the hydrophobic white layer protects against splashes and droplets. Further information on specific sample recovery features, PCR inhibition and aerosol protection are available on your local Eppendorf website.

* An aerosol is a colloid of fine solid particles or liquid droplets, in air or another gas.

Product features and benefits

- > Reliable blocking of liquid for comprehensive pipette protection
- > Maximum protection against aerosols and biomolecules for pipette and sample
- > No PCR-inhibition in case of sample contact
- > Available in PCR clean, Sterile (sterile, pyrogen-free)
- > Lot specific purity certificates

Extended applications

- > Applications with toxic or other dangerous sample material
- > Applications where limited sample amounts are available
- > Applications where sample preparation is especially time and cost consuming



Violet Layer

Forms a highly-effective barrier that blocks liquids and binds aerosols and biomolecules.

White Layer Designed to be hydrophobic to repel drops, splashes and to bind aerosols

epT.I.P.S.® LoRetention

Increased sensitivity with maximum reproducibility

Do you depend on your precious samples? Then don't leave any valuable residues in your tips when pipetting solutions that contain detergents! Based on a unique and innovative method of material modification, the new Eppendorf epT.I.P.S. LoRetention pipette tips allow you to transfer almost 100% of the liquid picked up (»pearl effect«) – for maximum recovery with unbeaten precision.

The pipetting of liquids containing detergents is ubiquitous in modern laboratory processes. The lower surface tension of detergents has a marked influence on the dosing properties of samples, thus making it noticeably more difficult to perform reproducible work with these liquids.

The »Pearl Effect«

The ultrahydrophobic extremely homogenous surface of the epT.I.P.S. LoRetention pipette tips is achieved through an innovative treatment at the molecular level – the Pearl Effect Technology.

The tips are uncoated, free from additives, and do not leach into the sample. Liquids roll off completely, so that only a tiny drop remains in the tip. In comparison, the adhesive layer of liquid in standard tips retains considerably more sample material.

For users of the epT.I.P.S. LoRetention pipette tips, this unique material property guarantees maximum sample recovery, improved reproducibility and higher sensitivity. Available in purity grade of PCR clean and Eppendorf Quality, also as PCR clean/Sterile ep Dualfilter T.I.P.S..



The »Pearl Effect« – Maximum recovery with epT.I.P.S. LoRetention Simple differentiation of racks and reloads through a new, clear lid with a white seal. The reusable seal supplied with the reloads can also be used to label existing epT.I.P.S. boxes. When open, the writing »LoRetention« can be seen from the inside.

The Facts Speak for Themselves

Maximum reproducibility in genomics epT.I.P.S.

LoRetention pipette tips are especially recommended for applications where highest precision of DNA/RNA analysis results are needed, such as for PCR and real-time PCR or NGS library preparation. For example expensive master mixes and enzyme solutions tend to adhere to the tips inner surface. Special treated low retention surfaces are created to repel detergent solutions to a maximum - for minimum loss of your valuable sample.

Examples of liquids with wetting effect:

- > Master mixes and NGS reagents
- > Enzyme solutions: restriction enzymes, ligation, DNase
- > DNA ladders for gel electrophoresis

epT.I.P.S. LoRetention perform considerably better in terms of precision and sample recovery than standard pipette tips as shown in fig. 1.



It is not only in molecular biology that high sensitivity detection methods require extreme reliability and reproducibility in pipetting. Also in protein analysis and purification the reagents and samples often contain detergents, like e.g. SDS-Page. By minimizing sample retention and improving reproducibility of pipetting, epT.I.P.S. LoRetention pipette tips and ep Dualfilter T.I.P.S. LoRetention filter tips are especially advantageous in proteomic applications.

Routine protein applications:

- > Isolation
- > Purification
- > Denaturation



Fig. 1: Comparison of residual liquid of epT.I.P.S. and epT.I.P.S. LoRetention with master mix for *real-time* PCR. Once the liquid had been dispensed, the residual liquid in the tips was determined. The measurements were repeated several times and the standard deviation was determined. epT.I.P.S. LoRetention resulted in the lowest residual liquid.

The facts speak for themselves – when compared to standard pipette tips, epT.I.P.S. LoRetention showed markedly better results in terms of precision and sample recovery, as displayed in fig. 1 and 2.



Fig. 2: Comparison of resistance of the low retention property when subjected to solvents typically used in proteomics. 200 μ L epT.I.P.S. and epT.I.P.S. LoRetention were treated with the solvents specified. An enzymatic buffer containing detergents was then pipetted and the residual moisture was determined. epT.I.P.S. LoRetention resulted in reproducible low levels of residual liquids.

Loading and Filling

Eppendorf GELoader®

Eppendorf's GELoader Tips have been designed to simplify the loading of samples onto polyacrylamide gels. These flexible, long and narrow tips prevent the gels from being damaged whilst allowing optimal handling of smallest volumes.

Product features and benefits

- > Special tip for gel electrophoresis
- > GELoader Tips and rack autoclavable (121°C, 20 min.)
- > Highest precision and accuracy when used with
- Eppendorf pipettes for 0.5 to 10 μL (gray control button)



Unmistakably Eppendorf – Thanks to the patented 3D design on the upper rim, Repeaters® are able to automatically and reliably detect the volume type of Combitips advanced®.

Eppendorf Microloader™

These unique tips are ideally suited to fill microcapillaries used for microinjection. The extremely long, fine and flexible tip provides also the ideal solution for all kind of applications in which additional immersion depth is needed while pipetting smallest volumes.

Product features and benefits

- > Extremely long, fine and flexible for filling of microcapillaries for microinjection
- > Rack package can be autoclaved by 121 °C
- > Ideal for recovering surplus solution from the capillary





Always ready to hand -

Combitips advanced are color coded for easy selection and secure connection – the Combitip rack also makes single-hand operation possible.



A perfect team – The Combitips advanced and ViscoTips[®] are optimized for all previous and new Repeater models – thus creating a perfect connection!

Time for a New Original



Eppendorf Combitips advanced®

The invention of the Eppendorf Repeater® manual dispenser and the Eppendorf Combitips® marks a milestone in the field of Liquid Handling! Executing long pipetting and dispensing series were thus made significantly simpler and faster. Thanks to its innovative sensor technology for automatic Combitip recognition, this dispensing system has become an indispensable tool for every laboratory. The increasingly high demands of modern laboratory work have also heightened the requirements for a highprecision dispensing tip. Our experts have thoroughly optimized the Combtips using sophisticated engineering.

The result: A revolutionary 360° evolution! Setting a completely new standard in dispensing systems.

The »Combitips« principle

- > Positive displacement principle (comparable to a syringe)
- > High-precision-dispensing regardless of the density and viscosity of the liquid
- > Prevents aerosol contamination with sealed piston for secure dispensing and provides protection from radioactive and toxic substances
- > Quick dispensing of long series with precise, repeated dispensing of identical volumes (in combination with the Repeater)
- > Individually color coded Quick identification of the desired Combitips facilitates the workflow
- > Variety and selection With 9 volume sizes (0.1 mL to 50 mL) and several purity levels you will always find the perfect Combitip for your application. The tips in Eppendorf Biopur, Sterile, PCR clean, and Forensic DNA Grade are individually blister-wrapped and feature an access tab which makes them easier to open, even with gloves

Technical Specifications Combitips advanced®

Combitips advanced	Min./max. volume	Increment/ step size	Max. dispensing volume	Test volume	Inaccuracy*1	Imprecision*
for Repeater [®] M4 and Cor						_
0.1 mL	, 	1μL	20 μL	2 μL	±1.6%	±3.0%
	20 μL	F.		20 μL	±1.0%	±2.0%
0.2 mL	2 μL	2 μL	40 μL	4 μL	±1.3%	±2.0%
	40 μL		·	40 μL	±0.8 %	±1.5%
0.5 mL	5 μL	5 μL	100 μL	10 μL	±0.9 %	±1.5%
	100 μL		·	100 μL	±0.8 %	±0.6 %
1 mL	10 μL	10 μL	200 μL	20 µL	±0.9 %	±0.9 %
	200 μL			_200 μL	±0.6 %	±0.4 %
2.5 mL	25 μL	25 μL	500 μL	50 μL	±0.8 %	±0.8%
	500 μL			500 μL	±0.5%	±0.3 %
5 mL	50 μL	50 μL	1,000 μL	100 μL	±0.6 %	±0.6 %
	1,000 μL			_1,000 μL	±0.5%	±0.25 %
10 mL	100 μL	100 μL	2,000 μL	200 μL	±0.5%	±0.6%
	2,000 μL			<u>2,000 μL</u>	<u>±0.5%</u>	<u>±0.25 %</u>
25 mL	250 μL	250 μL	5,000 μL	500 μL	±0.4 %	±0.6%
F.0. 1	<u>5,000 μL</u>			<u>5,000 μL</u>	$\pm 0.3\%$	<u>±0.25 %</u>
50 mL	500 μL	500 μL	10,000 μL	1,000 μL	±0.3%	±0.5%
for Repeater [®] stream/Xstr	<u>10,000 μL</u>			<u>10,000 μL</u>	±0.3 %	±0.3 %
				40		
0.1 mL	1 μL	0.1 μL	0.1 mL	10 μL	±1.6%	±2.5%
	100			50 μL	±1.0%	±1.5%
	<u>100 μL</u>	0.2 μL	0.2 mL	100 μL	$-\frac{\pm 1.0\%}{\pm 1.3\%}$	±0.5%
0.2 mL	2 μL	0.2 μL	0.2 ML	20 µL		±1.0 %
	200			100 μL	±1.0%	±1.0%
	200 μL			200 µL	<u>±1.0%</u>	±0.5 %
0.5 mL	5 μL	0.5 μL	0.5 mL	50 μL	±0.9%	±0.8%
	F00			250 μL	±0.9%	±0.5%
1 mL	<u>500 μL</u>			500 µL	$\pm 0.9\%$	<u>±0.3 %</u> ±0.55 %
ImL	10 μL	1 μL	1 mL	100 μL	±0.9%	
	1 1			500 μL	±0.6%	±0.3%
	<u>1 mL</u>	<u> </u>	2.5	1,000 μL	$\pm 0.6\%$	±0.2%
2.5 mL	25 μL	2.5 μL	2.5 mL	250 μL	±0.8%	±0.45 %
	2 E m.l			1,250 μL	±0.5%	±0.3%
5	2.5 mL			2,500 μL	$\pm 0.5\%$	$\pm 0.15\%$
5 mL	50 μL	5 μL	5 mL	500 μL	±0.8%	±0.35 % ±0.25 %
	E mal			2,500 μL	±0.5%	
10 mL	<u>5 mL</u>	10 μL	10 ml	5,000 μL	$\frac{\pm 0.5\%}{\pm 0.5\%}$	<u>+0.15 %</u> ±0.25 %
IU IIIL	100 μL	το με	10 mL	1,000 μL		
	10			5,000 μL	±0.4%	±0.25 %
	<u>10 mL</u>		25	10,000 μL	$\pm 0.4\%$	±0.15 %
25 mL	250 μL	25 μL	25 mL	2,500 μL	±0.3%	±0.35 %
	25			12,500 μL	±0.3 %	±0.25 %
	25 mL			25,000 μL	±0.3%	±0.15%
50 mL	500 μL	50 μL	50 mL	5,000 μL	±0.3%	±0.5%
	50 .			25,000 μL	±0.3%	±0.2%
	50 mL			50,000 μL	±0.3 %	±0.15 %

*1 The data for imprecision (random error) and inaccuracy (systematic error) according to EN ISO 8655 only apply when using original Eppendorf Combitips advanced.

Compatibility of Combitips advanced with standard laboratory tubes

Combitips advanced/ ViscoTip®	Eppendor	rf Safe-Lock	Tubes	Eppendorf Tubes [®] 5.0 mL	Conical t	tubes	Eppendorf	Deepwell Plate	es
	0.5 mL	1.5 mL	2.0 mL	5.0 mL	15 mL	50 mL	96/500 μL	96/1000 μL	96/2000 μL
0.2 mL	+	+	+	+	-	-	+	+	+
0.5 mL	+	+	+	+	-	-	+	+	+
1 mL	+	+	+	+	_	-	+	+	+
2.5 mL	+	+	+	+	+	+	+	+	+
5 mL	+	+	+	+	+	+	+	+	+
10 mL	+	+	+	_	_	+	+	+	+
25 mL	_	+	+	_	_	_	+	+	+
50 mL	-	+	+	-	-	-	+	+	+

ViscoTip®

Experience the new member of the Combitips advanced dispenser tip family. The ViscoTip is specifically designed and optimized for handling high viscosity liquids up to 14,000 mPa*s such as Glycerol 99.5%, Tween, oils, cremes, shampoos or honey. It significantly reduces operating forces while handling such liquids leading to enhanced ergonomics, increased working speed and longer charge life time of your Repeater battery.



Positive Displacement Tips



For Eppendorf Maxipettor®

The pipette tips for the Maxipettor are tailored to different vessels. For example, the Eppendorf Maxitips P are designed for aspirating 1 mL to 10 mL from beakers, and it pipettes according to the positive displacement principle. The Eppendorf Maxitip S for 2.5 mL to 10 mL forms a system with the Maxitip and a dispensing part. This system can be used for aspirating liquids from tall, narrow-neck vessels.

Product features and benefits

- > Maxitips P (fig. 1) for aspirating 1 mL to 10 mL from beakers using a positive displacement technique
- > Maxitips S (fig. 2): The dosing part, in combination with the Eppendorf Maxitip (fig. 3) forms a system for aspiration of between 2.5 mL and 10 mL liquid from high, narrowmouth vessels – air displacement principle
- > Valve for Maxitips ensures drip-free dispensing of liquids with a high vapour pressure
- > Maxitip is graduated for dispensing accurate volumes of nonaqueous liquids

Reliable Unattended Automation



epT.I.P.S.[®] Motion – Eppendorf Totally Integrated Pipetting System for Automation

For routine pipetting with the best precision in the industry it is important to have a perfect system of tips and automate. The epT.I.P.S. Motion have been developed to work in perfect combination with our ep*Motion*[®] platform.

The epT.I.P.S. Motion racks consist of two options for simple integration into the ep*Motion* liquid handling workstation: A single use box with color coded trays for easy volume identification and a TipHolder adapter to use the trays as a Reload-System. With color coded trays you can quickly and easily see the volume range of the tip and assure an even higher level of sample safety. The new Reload trays also offer a more environmentally friendly option and will be delivered in a PET blister packaging with a sealed cap assuring the same well established quality as all Eppendorf tips.

The optional TipHolder adapter (autoclavable) replaces the normal tip box component and reduces waste by up to 40 %. A perfect match for the most demanding user.

Product features and benefits

More safety

- > Each tip is checked for straightness before packaging
- > Color-coded trays for direct volume identification
- > Dispensing tool design assures an optimal fit each time
- > Optical sensor automatically identifies tip type Flexible handling
- > For volume ranges (0.2 µL to 1 mL) available with and w/o filter in multiple purity grades
- > Easily exchangeable, no additional labware files required
- > Available as SafeRacks for tip reuse with individual tip compartments
- Eco-friendly Reload System
- > Easy conversion with TipHolder adapter
- > TipHolder adapter is autoclavable
- > Reduces waste by up to 40 %





Volume range/epT.I.P.S epT.I.P.S.® LoRetentior		Standard	Reloads	Reloads	
(all tips shown are actu					
		Eppendorf Quality™	Eppendorf Quality™	PCR clean	
0.1–10 μL, 34 mm	dark gray	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
12		022492004	022491504 022493018 LoRetention	022491709022493010 LoRetention	
0.1–20 μL, 40 mm	medium gray	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
le		022492012	022491512	022491717	
0.5–20 μL L, 46 mm	light gray	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
9 <u>-</u>		022492021	022491521 022493020 LoRetention	022491725 022493012 LoRetention	
2–200 μL, 53 mm	yellow	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
		022492039	022491539 022493022 LoRetention	 022491733 022493014 LoRetention 	
20–300 μL, 55 mm	orange	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
		022492047	022491547	022491741	
50–1,000 μL, 71 mm	blue	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
		022492055	022491555 022493024 LoRetention	 022491750 022493016 LoRetention 	
50–1,250 μL, 76 mm	green	2 bags of 500 tips = 1,000 tips	10 trays of 96 tips = 960 tips	10 trays of 96 tips = 960 tips	
ATT:		022492063	022491563	022491768	

Box	Set	Singles	Racks
Eppendorf Quality™	Eppendorf Quality™	Eppendorf Biopur®	Eppendorf Sterile Eppendorf Biopur®
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips		
022491300	022491407 0030072251 LoRetention		
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips	100 tips, individually wrapped	Biopur 5 racks of 96 tips = 480 tips Sterile 10 racks of 96 tips = 960 tips
022491318	022491415	022491130	022491067 Biopur 022492250 Sterile
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips		
022491326	 022491423 0030072260 LoRetention 		
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips	100 tips, individually wrapped	Biopur 5 racks of 96 tips = 480 tips Sterile 10 racks of 96 tips = 960 tips
022491334	022491431 0030072278 LoRetention	<mark>-</mark> 022491148	022491083 Biopur 022492276 Sterile
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips		Biopur 5 racks of 96 tips = 480 tips Sterile 10 racks of 96 tips = 960 tips
022491342	0022491440		022491091 Biopur 022492284 Sterile
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips	100 tips, individually wrapped	Biopur 5 racks of 96 tips = 480 tips Sterile 10 racks of 96 tips = 960 tips
022491351	 022491458 0030072286 LoRetention 	022491156	 022491105 Biopur 022492292 Sterile
1 reusable box incl. 96 tips	1 reusable box incl. 5 trays of 96 tips		Biopur 5 racks of 96 tips = 480 tips Sterile 10 racks of 96 tips = 960 tips
022491369	022491466		022491113 Biopur
			ULL+72500 Stellie



Box	Set	Singles	Racks
Eppendorf Quality™	Eppendorf Quality™	Eppendorf Biopur [®] (sterile)	Eppendorf Biopur [®] (sterile)
1 reusable box incl. 96 tips ■ 022494016			5 racks of 96 tips = 480 tips ■ 022494014
1 reusable box incl. 48 tips	1 reusable box		5 racks of 48 tips = 240 tips
022491377	incl. 5 trays of 48 tips ■ 022491474		022491121
1 reusable box incl. 24 tips ■ 022491385			5 racks of 24 tips = 120 tips 00030075200 Biopur 022492314 Sterile
			5 racks of 24 tips = 120 tips ■ 0030075218
			5 racks of 24 tips = 120 tips ■ 022491164



Volume range/ep Dualfilter T.I.P.S.®, Racks ep Dualfilter T.I.P.S.[®] SealMax, ep Dualfilter T.I.P.S.® LoRetention (all tips shown are actual size) PCR clean, sterile incl. pyrogen-free **0.1–10 μL S,** 34 mm 10 racks of 96 tips = 960 tips dark gray 022491202 022493000 LoRetention 0030077806 SealMax **0.1–10 μL M,** 40 mm 10 racks of 96 tips = 960 tips medium gray 022491211 0030077768 Forensic DNA Grade **0.5–20 μL L,** 46 mm 10 racks of 96 tips = 960 tips light gray 022491229 022493002 LoRetention 0030077814 SealMax **2–20 μL,** 53 mm yellow 10 racks of 96 tips = 960 tips 022491270 0030077776 Forensic DNA Grade **2–100 μL,** 53 mm yellow 10 racks of 96 tips = 960 tips 022491237 022493006 LoRetention 1 0030077822 SealMax **2–200 μL,** 55 mm 10 racks of 96 tips = 960 tips yellow 022491296 0030077830 SealMax 0030077784 Forensic DNA Grade **20–300 μL,** 55 mm orange 10 racks of 96 tips = 960 tips 022491245 022493004 LoRetention 0030077849 SealMax **50–1,000 μL,** 76 mm blue 10 racks of 96 tips = 960 tips 022491253 022493008 LoRetention 0030077857 SealMax 0030077792 Forensic DNA Grade **50–1,250 μL L,** 103 mm **Δ** dark green 5 Racks of 96 Tips = 480 Tips 022494002



Special Tips

Volume range/pipette tip	GELoader®	uppendier
0.5–20 μL, 62 mm light gray	2 racks of 96 tips = 192 tips	
	022351656	
	Microloader	
0.5–20 μL, 100 mm ■ medium gray	2 racks of 96 tips = 192 tips	quality
	930001007	

Combitips advanced[®] and Accessories

Combitips advanced®	Color coding	Eppendorf Quality™ 100 pcs.	PCR clean 100 pcs. (4 reclos-	Forensic DNA Grade 100 pcs. (indivi-	Eppendorf Biopur [®] 100 pcs. (individually	Sterile 100 pcs (indivi-
		(4 bags × 25 pcs.)	able bags × 25 pcs.)	dually wrapped)	wrapped)	dually wrapped
0.1 mL	□ white	0030089405	0030089766		0030089618	0030089510
0.2 mL	light blue	0030089413	0030089774		0030089626	0030089529
0.5 mL	violet	0030089421	0030089782		0030089634	0030089537
1 mL	yellow	0030089430	0030089790	0030089855	0030089642	0030089545
2.5 mL	green	0030089448	0030089804	0030089863	0030089650	0030089553
5 mL	blue	0030089456	0030089812	0030089871	0030089669	0030089561
10 mL	erange	0030089464	0030089820		0030089677	0030089570
25 mL*	red	0030089472	0030089839		0030089685	0030089588
50 mL*	light gray	0030089480	0030089847		0030089693	0030089596
Eppendorf ViscoTips®						
10 mL	orange	0030089502				
Accessories						
Adapter advanced						
25 mL adapter (1 pc.)		r ed	0030089715			
50 mL adapter (1 pc.)		light gray	0030089723			
25 mL adapter (7 pcs.)		r ed				0030089731
50 mL adapter (7 pcs.)		light gray				0030089740
Combitip Rack (for			0030089758			
8 Combitips advanced [®] ,						
0.1 mL–10 mL)						
Combitips advanced [®]			0030089936			
Assortment pack (1 Co	•					
of each size, incl. adapt	ters)					
* 4 boxes of 25 pcs, each; each box c	ontaining one adapter.					

* 4 boxes of 25 pcs. each; each box containing one adapter.

Eppendorf Maxitip P and S for Maxipettor®

Description	Order no.
Maxipettor, single-channel, with continuous volume selection, 1 – 10 mL	022290002
Maxitip P, to remove liquid from smaller vessels 100 pieces	022291050
Maxitip S Starter kit to liquid from narrow-necked vessels and measuring flasks, consisting of 100 Maxitips, 10 dispensing parts, 10 valves	022291301
Cylinder with piston, 30 pieces	022291408
Maxitip S, graduated, 200 pieces	022291351
Maxitip S valve, 100 pieces	022291459

epMotion[®] Automated Pipetting System

Description

epT.I.P.S.® Motion pipette tips

Order no.

Automatic pipette tips in individual racks for use on the ep*Motion*[®]. The tip type and size is automatically recognized on the device. 96 epT.I.P.S.[®]/ rack, 10 racks per set. The refill racks from the reload products can be placed in an autoclavable TipHolder adapter. Three purity levels: Eppendorf Quality[™], PCR clean and sterile. Pipette tips available with or without a filter. Testing of production batches (certificates available).

Without filter

10 μ L, sterile, free of pyrogens, volume range 0,2-10 μ L, 10 x 96 tips in racks	0030015185
50 μL, sterile, free of pyrogens, volume range 1–50 μL, 10 × 96 tips in racks	0030015207
300 µL , sterile, free of pyrogens, volume range 20–300 µL, 10×96 tips in racks	0030015223
1,000 μL, sterile, free of pyrogens, volume range 40–1,000 μ L, 10 × 96 tips in racks	0030015240
50 μL, Eppendorf Quality [™] , volume range 1–50 μL, 10 × 96 tips in racks	0030014405
300 μL, Eppendorf Quality [™] , volume range 20–300 μL, 10 × 96 tips in racks	0030014448
1,000 μL, Eppendorf Quality [™] , volume range 40–1,000 μL, 10 × 96 tips in racks	0030014480
50 μL, Eppendorf Quality [™] , Reloads, volume range 1–50 μL, 24 × 96 tips	0030014421
300 μL, Eppendorf Quality™, Reloads, volume range 20–300 μL, 24 × 96 tips	0030014464
1,000 μL, Eppendorf Quality [™] , Reloads, volume range 40–1,000 μL, 24 × 96 tips	0030014502
SafeRacks, 50 μ L, Eppendorf Quality TM , volume range 1–50 μ L, 10 × 96 tips	0030014600
SafeRacks, 300 μ L, Eppendorf Quality TM , volume range 20–300 μ L, 10 × 96 tips	0030014626
SafeRacks, 1,000 μL, Eppendorf Quality™, volume range 40–1,000 μL, 10 × 96 tips	0030014642
With filter	
10 μ L, PCR clean, sterile, free of pyrogens, volume range 0,2-10 μ L, 10 x 96 tips in racks	0030015193
50 µL, PCR clean, sterile, free of pyrogens, volume range 1–50 µL, 10×96 tips in racks	0030015215
300 μL , PCR clean, sterile, free of pyrogens, volume range 20–300 μ L, 10 × 96 tips in racks	0030015231
1,000 µL, PCR clean, sterile, free of pyrogens, volume range 40–1,000 µL, 10×96 tips in racks	0030015258
50 μ L, PCR clean, volume range 1–50 μ L, 10 × 96 tips in racks	0030014413
300 μL , PCR clean, volume range 20–300 μ L, 10 × 96 tips in racks	0030014456
1,000 μL, PCR clean, volume range 40–1,000 μL, 10 \times 96 tips in r acks	0030014499
50 μL, PCR clean, Reloads, volume range 1–50 μL, 24 × 96 tips	0030014430
50 μ L, PCR clean, sterile, free of pyrogens, Reloads , volume range 1–50 μ L, 24 × 96 tips	0030014537
300 μL, PCR clean, Reloads, volume range 20–300 μL, 24 \times 96 tips	0030014472
300 μL , PCR clean, sterile, free of pyrogens, Reloads , volume range 20–300 μ L, 24 × 96 tips	0030014529
1,000 μL, PCR clean, Reloads, volume range 40–1,000 μ L, 24 × 96 tips	0030014510
SafeRacks, 50 μ L, PCR clean, volume range 1–50 μ L, 10 × 96 tips	0030014618
SafeRacks, 300 μL, PCR clean, volume range 20–300 μL, 10 × 96 tips	0030014634
SafeRacks, 1,000 μ L, PCR clean, volume range 40–1,000 μ L, 10 × 96 tips	0030014650

Reservoir rack modules

Are inserted in a Reservoir rack. They can be temperature controlled with a thermal module for heating and cooling.

PCR 0.2 mL, for 8 × 0.2 mL PCR tubes	960002601
Eppendorf Tubes [®] 5.0 mL, for 4 × 5 mL tubes	5075799340
Eppendorf Safe-Lock, for 4 × 0.5/1.5/2 mL tubes	960002620
Tubes Ø 12 mm, for 4 × Ø 12 mm tubes	960002630
Tubes Ø 16 mm, for 4 × Ø 16 mm tubes	960002640
15 mL conical tubes, for $4 \times \emptyset$ 17 mm tubes	960002650
50 mL conical tubes, for $2 \times \emptyset$ 29 mm tubes	960002660
10 mL reservoir, for use with reservoir rack, 5 × 10 large volume reservoirs, PCR clean	0030126521
30 mL reservoir, for use with reservoir rack, 5 × 10 large volume reservoirs, PCR clean	960051009
100 mL reservoir, for use with reservoir rack, 5 × 10 large volume reservoirs, PCR clean	960051017
195 mL low dead volume reservoir, single well with 96 cavities, PCR clean, 4 x 5 reservoirs	0030126556
400 mL reservoir, also for use in the epMotion® VAC, 10 pcs./set, made of PP	960002229

Eppendorf Pipette/epT.I.P.S.® Combinations

Eppendorf Research [®] plus	0.1 μL− 10 μL (S)	0.1 μL− 10 μL (M)	0.1 μL− 20 μL	0.5 μL– 20 μL L	2 μL– 20 μL	2 μL– 100 μL	2 μL– 200 μL	20 μL− 300 μL
	📕 dark gray	mediu	m gray	light gray	yellow	yellow	yellow	orange
ep T.I.P.S.®	√		√	√	√	√		√
ep T.I.P.S. [®] LoRetention	√			√			√	
ep Dualfilter T.I.P.S.®	√							
ep Dualfilter T.I.P.S. [®] LoRete								√
ep Dualfilter T.I.P.S. [®] SealMa	ax <u>√</u>		√					√
Fixed volume								
10 μL medium	gray 🗸	√	√	√				
20 μL 🗌 light gra	ау			√				
10 μLyellow					√	√	√	
20 μL yellow					√	√	√	√
25 μL, 50 μL, 📃 yellow 100 μL						√	~	√
200 μLyellow							√	√
200 μL, 250 μL, 📕 blue _500 μL, 1,000 μL								
Adjustable volume								
0.1 μL–2.5 μL 📕 dark gra	ıy √	√	√					
0.5 μL–10 μL 📃 medium	gray 🗸	√	√	√				
2 μL–20 μL 📃 light gra	ay 🗸	\checkmark	√	√				
2 μ–20 μL yellow					√	√	√	\checkmark
10 μL–100 μL 🧧 yellow					\checkmark	√	√	√
20 μL–200 μL yellow					\checkmark	\checkmark	√	√
30 μL–300 μL 🗧 orange					\checkmark	\checkmark	\checkmark	√
100 μL–1,000 μL 📃 blue								
0.5 mL–5 mL violet								
1 mL–10 mL turquois	ie							

Eppendorf Xplorer® Eppendorf Xplorer® plus	0.1 μL− 10 μL (S)	0.1 μL− 10 μL (M)	0.1 μL− 20 μL	0.5 μL– 20 μL L	2 μL– 20 μL	2 μL– 100 μL	2 μL– 200 μL	20 μL– 300 μL
	dark gray	📕 mediu	um gray	light gray	yellow	yellow	yellow	orange
ep T.I.P.S.®	√		√	√	√	√	√	√
ep T.I.P.S. [®] LoRetention	√						√	
ep Dualfilter T.I.P.S.®	√	~		√		√	√	√
ep Dualfilter T.I.P.S.® LoRetention	√					√		√
ep Dualfilter T.I.P.S. [®] SealMax	1		√			√	√	√
0.5 μL–10 μL 📃 medium gray	√	√	√	√				
5 μL–100 μL 📃 yellow					\checkmark	√	√	√
15 μL–300 μL – orange					~	\checkmark	\checkmark	√
50 μL–1000 μL 📃 blue								
50 μL–1200 μL 📕 green								
0.25 mL–5 mL violet								
0.5 mL–10 mL turquoise								
✓: Compatible, √: Limited volume								

✓ : Compatible, ✓ : Limited volume

50 μL- 1,000 μL √ √ √ √ √ √ √ √ √	50 μL− 1,250 μL green ✓	50 μL− 1,250 μL L dark green ✓	0.25 mL- 2.5 mL red √	0.1 mL− 5 mL ✓ ✓	0.2 mL- 5 mL L ✓	0.5 mL- 10 mL ■ turquoise ✓ ✓	0.5 mL- 10 mL L ■ turqouise ✓ ✓	GELoader Microloader light gray ✓
✓ 	✓ 	✓						
✓	√	✓		√	√			

50 μL− 1,000 μL ■ blue	50 μL− 1,250 μL green	50 μL– 1,250 μL L dark green	0.25 mL− 2.5 mL ■ red	0.1 mL− 5 mL ■ violet	0.2 mL− 5 mL L ■ violet	0.5 mL– 10 mL turquoise	0.5 mL– 10 mL L turqouise	GELoader Microloader light gray
√	√	√	√	√	√	√	√	√
✓								
✓		✓		✓	✓	✓	✓	
√								
√								
								√
√	√	✓						
	√	√						
				✓	√			
						✓	✓	

Eppendorf Pipette/epT.I.P.S.® Combinations

Eppendorf Reference	e [®] 2	0.1 μL– 10 μL (S)	0.1 μL– 10 μL (M)	0.1 μL– 20 μL	0.5 μL– 20 μL L	2 μL– 20 μL	2 μL– 100 μL	2 μL– 200 μL	20 μL– 300 μL
		dark gray	mediu		light gray	yellow	yellow	yellow	orange
ep T.I.P.S.®		√		√	√	1	 √		
ep T.I.P.S. [®] LoRetenti	on	√			√			√	
ep Dualfilter T.I.P.S.®		√	√		√	√	√	√	√
ep Dualfilter T.I.P.S.®	LoRetention	√			√		√		√
ep Dualfilter T.I.P.S.®	SealMax	√		√			√	√	√
Fixed volume									
1 μL, 2 μL 📃 da	ark gray	√	√	√					
5 μL, 10 μL 📃 m	nedium gray	√	√	√	√				
20 μL lig	ght gray			√	√				
10 μL 🧧 ye	ellow					√	√	√	
20 μL 🧧 ye	ellow							√	√
25 μL, 50 μL, <mark> </mark>	ellow						~	~	1
200 μL 🧧 ye	ellow							√	√
200 μL, 250 μL, bl 500 μL, 1,000 μL	lue								
2 mL, 2.5 mL 📕 re	ed								
Adjustable volume									
0.1 μL–2.5 μL 🔳 da	ark gray	√	√	√					
0.5 μL-10 μL 🛛 🕅 m	nedium gray	√	√	√	√				
2 μL–20 μL 📃 lig	ght gray	\checkmark	\checkmark	√	√				
2 μL-20 μL 🦳 ye	ellow					√	√	√	\checkmark
10 μL-100 μL 🧧 ye	ellow					~	√	√	\checkmark
20 μL-200 μL 🧧 ye	ellow					~	\checkmark		√
30 μL-300 μL 📃 οι	range					\checkmark	\checkmark	\checkmark	√
100 μL–1,000 μL 📃 bl	lue								
0.25 mL–2.5 mL 📕 re	ed								
	iolet								
1 mL–10 mL 📃 tu	urquoise								

✓ : Compatible, √ : Limited volume

50 μL- 1,000 μL blue √ √ √ √ √ √	50 μL− 1,250 μL green ✓	50 μL− 1,250 μL L dark green ✓	0.25 mL- 2.5 mL ■ red √	0.1 mL- 5 mL ✓ ✓ ✓	0.2 mL- 5 mL L ■ violet ✓ ✓	0.5 mL- 10 mL ■ turquoise √ 	0.5 mL− 10 mL L turqouise √ √	GELoader Microloader Iight gray ✓
								✓ ✓ ✓
		✓	✓					
		✓ 						

eppendorf



Make Your Lab a Better Place

The Eppendorf Liquid Handling Instruments Portfolio

As the inventor of the microliter system, we at Eppendorf have more than fifty years of experience in precise manual and automatic pipetting and transferring of the smallest quantities of liquids. The standards defined by Eppendorf quality assurance (for example, precision and accuracy guidelines) clearly surpass the standards required to obtain these results, and guarantee the reproducibility that our customers have come to depend on.

Our »PhysioCare Concept®« incorporates growing ergonomic demands into all that we do. The use of our liquid handling products has been proven to reduce physical strain at the workstation to a minimum. This has been perfected by the ergonomic design of our products and flawless fit of device to consumable (e.g., measurably reduced ejecting force when pipetting).

For more information visit www.eppendorf.com/pipetting





> Ergonomics in the lab: www.eppendorf.com/ergonomics

Your local distributor: www.eppendorf.com/contact Eppendorf North America, Inc. · Phone: 800-645-3050 info@eppendorf.com · www.eppendorfna.com Eppendorf Canada Ltd. · Phone: 800-263-8715 canada@eppendorf.com · www.eppendorf.ca

www.eppendorf.com

Eppendorf[®], the Eppendorf Brand Design, Biopur[®], ep Dualfilter T.I.P.S.[®], epT.I.P.S.[®], epMotion[®], Eppendorf Varitips[®], Eppendorf Reference[®], Eppendorf Research[®], Combitips advanced[®], GELoader[®], Combitips[®], Biomaster[®], Mastertip[®], PhysioCare Concept[®], Multipette[®], Varipette[®], ViscoTip[®], and Eppendorf Xplorer[®] are registered trademarks of Eppendorf AG, Germany. Eppendorf QualityTM, and Eppendorf MicroloaderTM are trademarks of Eppendorf AG, Germany. All rights reserved, including graphics and images. Copyright © 2018 by Eppendorf AG. Order No.: A03X X13 080/ENA1/6T/0818/MCP/K&R