New Brunswick™ TCA-3 Temperature Monitoring System

Manual de instrucciones

Copyright

Copyright © 2014 Eppendorf AG, Germany. No part of this publication may be reproduced without the prior permission of the copyright owner.

The company reserves the right to change information in this document without notice. Updates to information in this document reflect our commitment to continuing product development and improvement.

Trademarks

Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf AG, Germany.

New Brunswick™ is a trademark of Eppendorf AG, Germany.

VeriSign Secured® is a registered trademark of VeriSign, Inc., USA.

Microsoft® and Excel® are registered trademarks of Microsoft Corporation in the United States and other countries.

Trademarks are not marked in all cases with ™ or ® in this manual.

Eppendorf has attempted to identify the ownership of all trademarks from public records. Any omissions or errors are unintentional.

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The image shows a login page for New Brunswick, an Eppendorf company. The page is divided into two columns. The left column contains the login form with the following fields: 'Gateway Mac:', 'Login Name:', 'Password:', and a checkbox labeled 'Remember Gateway Mac'. Below these fields is a 'Login' button. At the bottom of the left column, there is a link: 'Click [HERE](#) to try the demo site.' The right column contains the slogan 'Monitoring made easy.' The entire login form area is circled in red.

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
The screenshot shows a web page with a blue header bar containing the text "Welcome to your On". Below the header, there is a section titled "Account Setup Required" in blue text. Underneath this title, the text reads: "To properly setup your account we need you to complete some basic". This is followed by the sentence: "In the following steps we will ask you to enter information for:". Below this, there is a numbered list with three items: "1. Primary Contact Info", "2. User Accounts Logins", and "3. The Default Alarm Contact(s)". At the bottom of this section, there is a bold note: "NOTE: This information is required so that we can contact you and securely configured.". At the bottom right of the page, there is a button with a blue border and a dotted background, containing the text "Next >>".

Welcome to your On

Account Setup Required

To properly setup your account we need you to complete some basic

In the following steps we will ask you to enter information for:

1. Primary Contact Info
2. User Accounts Logins
3. The Default Alarm Contact(s)

NOTE: This information is required so that we can contact you and securely configured.

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Account Se

Primary Contact

From time to time we may need to notify you of system updates, schedules, and other important information. Who is the primary contact for this monitoring system account?

Who is the primary contact for this monitoring system account?

Full Name*:	<input type="text"/>
Company*:	<input type="text"/>
Address Line1*:	<input type="text"/> <small>Street address, P.O. box, company name, c/o</small>
Address Line2:	<input type="text"/> <small>Apartment, suite, unit, building, floor, etc.</small>
City*:	<input type="text"/>
State/Province/Region*:	<input type="text"/>
ZIP/Postal Code*:	<input type="text"/>
Country*:	<input type="text" value="United States"/>
Timezone*:	<input type="text" value="(GMT-05:00) Eastern Standard Time (DST)"/>
Phone Number*:	<input type="text"/>
Email Address*:	<input type="text"/>
Repeat Email Address*:	<input type="text"/>

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Account Se

User Accounts

Your Admin and Reader accounts have been setup with default passwords, you must change these passwords, if you haven't already.

<p>Administrator</p> <p>Login Name*: <input type="text" value="Admin"/></p> <p>Password: <input type="text"/></p> <p>Repeat Password: <input type="text"/></p> <p>First Name: <input type="text"/></p> <p>Last Name: <input type="text"/></p>	<p>Read-Only Account</p> <p>Login Name*: <input type="text"/></p> <p>Password: <input type="text"/></p> <p>Repeat Password: <input type="text"/></p> <p>First Name: <input type="text"/></p> <p>Last Name: <input type="text"/></p>
--	--

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The screenshot shows the 'Account Set' interface. The main section is titled 'Default Alarm Contact'. It contains a message: 'To ensure system reliability, at least one of the following lists must'. Below this is an 'Edit Message' button. To the right, there is a 'Notify Type:' dropdown menu currently set to 'All at'. Below these are two sections: 'Default Email List' and 'Default Phon'. A modal dialog is open in the foreground, titled 'New Email:'. It has two input fields: 'New Email:' and 'Re Enter:'. At the bottom right of the dialog are 'cancel' and 'add' buttons. At the bottom of the main interface are '<< Prev' and 'Next >>' buttons.

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Account Set

Default Alarm Contact

To ensure system reliability, at least one of the following lists must

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
The screenshot shows the main dashboard of the New Brunswick software. It features a header with the company logo and navigation links. Below the header, there are sections for 'TRIGGERED ALARMS' and 'LATEST MEASUREMENTS'. The 'TRIGGERED ALARMS' section contains a table with columns for Alarm Name, POD, Measurement, Value, Alarm Time, and Chart. The 'LATEST MEASUREMENTS' section shows a table with columns for Measurement Name, Value, and status icons. Four red circles with numbers 1 through 4 highlight specific UI elements: 1. The top navigation menu (SETTINGS | SUPPORT | LOGOUT). 2. The 'SEND TO HISTORY' button. 3. The 'LATEST MEASUREMENTS' table. 4. The 'VIEW ALARM HISTORY' button.

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The screenshot shows the New Brunswick TCA-3 monitoring system interface. At the top, there is a logo for New Brunswick (an endpoint company) and navigation links for SETTINGS, SUPPORT, and LOGOUT. Below the logo, it says 'Welcome' and 'Time to Next Data Storage: 00:00:43'. The main section is titled 'TRIGGERED ALARMS' and contains a table with columns: ALARM NAME, POD, MEASUREMENT, VALUE, ALARM TIME, and CHART. Two alarms are listed: 'Low supply voltage' (GW-0014a002017d, Supply voltage (Battery), 2.28, 01:32:00 PM 01/25/10) and 'Power or Internet Outage' (N/A, N/A, 01:34:00 PM 01/25/10). Below the table are buttons for 'SEND TO HISTORY' and 'VIEW ALARM HISTORY'. A secondary section titled 'LATEST MEASUREMENTS' shows data for GW-0014a002017d as of 02/04/10 09:51:00 AM, with two RTD sensors: RTD #1 (6.58 C) and RTD #2 (-20.22 C). On the right side of the measurements table, there are three tabs: 'Gateway Data and Diagnostics', 'Alarms', and 'Preferences', with the 'Preferences' tab highlighted.

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
This screenshot shows the 'Gateway Preferences' dialog box for gateway GW-0014a002017d. The dialog has three tabs: 'General', 'Info', and 'View', with 'General' selected. It contains the following fields: 'Name' (GW-0014a002017d), 'Store data every:' (set to 1 Min), and 'Check for alarm every:' (set to 1 Min). There is a 'Notes:' text area below. At the bottom are 'OK', 'Cancel', and 'Apply' buttons. The dialog is overlaid on the main dashboard, which is partially visible in the background. A red circle highlights the 'Preferences' tab icon in the background dashboard.

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
The screenshot displays the New Brunswick TCA-3 web interface. At the top, there is a navigation bar with 'SETTINGS | SUPPORT | LOGOUT' and a 'Welcome' message. Below this, a table titled 'TRIGGERED ALARMS' lists two active alarms:

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
Low supply voltage	GW- 0014a002017d	Supply Voltage (Battery)	2.28	01:32:00 PM 01/25/10	
Power or Internet Outage	N/A	N/A		01:34:00 PM 01/25/10	

Below the table are controls for 'SEND TO HISTORY' and 'VIEW ALARM HISTORY'. A modal window titled 'GATEWAY COMMUNICATION ALARM' is open, showing configuration for 'GW- 0014a002017d'. The 'General' tab is active, showing:

- Alarm Enabled
- Alarm Description: Name: Power or Internet Outage
- Notify List: Default List (selected)

Buttons for 'OK', 'Cancel', and 'Apply' are at the bottom of the modal. In the background, a sidebar menu has 'Gateway Description and Configuration' highlighted with a red circle.

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


4.11 Alarmas de sensor

The screenshot shows the New Brunswick TCA-3 Temperature Monitoring System interface. At the top left is the New Brunswick logo. The top right has links for SETTINGS, SUPPORT, and LOGOUT. Below the logo, it says 'Welcome' and 'Time to Next Data Storage: 00:01:17'. The main section is titled 'TRIGGERED ALARMS' and contains a table with columns: ALARM NAME, POD, MEASUREMENT, VALUE, ALARM TIME, and CHART. Below the table, it says 'There are no triggered alarms'. There are also buttons for 'SEND TO HISTORY' and 'VIEW ALARM HISTORY'. Below this is the 'LATEST MEASUREMENTS' section, which shows a table of measurements for device 'GW- 0014a002016a' on '01/05/12 09:19:00 AM'. The table has columns for measurement type, value, and status icons. A red circle highlights the alarm icon for 'RTD #2'. To the right of the table are links for 'Generate Data Logs and Downloads', 'Alarms', and 'Preferences'.

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
There are no triggered alarms					

LATEST MEASUREMENTS					
GW- 0014a002016a				01/05/12 09:19:00 AM	
RTD #1		-245.76 C			
RTD #2		-237.61 C			
Supply Voltage (Battery)		4.05 Volts			
Supply Voltage (External)		4.73 Volts			
Thermocouple		1,177.51 C			

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[Add Notes](#)

Notify List: Default List

This will notify everyone on the "Default Notify List".
The list is located in the "Settings" menu on the main page.

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
The screenshot shows a software window titled "ALARM DEFINITION" for "GW- 0014a002017d - RTD #1". It has two tabs: "General" and "Trigger". The "Trigger" tab is selected. Inside the "Trigger" tab, there is a section for "Alarm Limits" with the following settings: "Sensor value*" set to "<" with a value of "20", "And/Or" set to "OR", and "Sensor value*" set to ">" with a value of "30". Below this is an "Advanced" section with "Trigger Mode" set to "Auto-Rearm" and a "Filter" section set to "3 occurrences in 5 consecutive samples". At the bottom of the dialog are "OK", "Cancel", and "Apply" buttons. A small note in the top right corner of the dialog says "* = required".

4.12 Gráfico de sensor

TRIGGERED ALARMS

ALARM NAME	POD	MEASUREMENT	VALUE	ALARM TIME	CHART
There are no triggered alarms					

Select: All, None, Inactive ■ Active Alarm ■ Inactive Alarm

LATEST MEASUREMENTS

Device	Time	RTD #1	RTD #2	Supply Voltage (Battery)	Supply Voltage (External)	Thermocouple
GW-0014a002016a	01/05/12 09:19:00 AM	-245.76 C	-237.61 C	4.05 Volts	4.73 Volts	1,177.51 C

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7 Certificados

Declaration of Conformity

The product named below fulfills the requirements of directives and standards listed. In the case of unauthorized modifications to the product or an unintended use this declaration becomes invalid.

Product name:

TCA-3
including accessories

Product type:

Remote logging temperature monitor

Relevant directives / standards:

2004/108/EC EN 61326-1
2006/95/EC EN 61010-1
2011/65/EU

Management Board

Portfolio Management

Date: January 06, 2014

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

MI288-9999-00

Eppendorf® and the Eppendorf Logo are registered trademarks of Eppendorf AG, Hamburg/Germany. All rights reserved incl. graphics and pictures. Copyright 2013 © by Eppendorf AG.

Certificados

TCA-3 Temperature Monitoring System
Español (ES)

Evaluate your manual

Give us your feedback.
www.eppendorf.com/manualfeedback