Register your instrument!
www.eppendorf.com/myeppendorf

MixMate®

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Copyright © 2022 Eppendorf SE, Germany. All rights reserved, including graphics and images. No part of this publication may be reproduced without the prior permission of the copyright owner.

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany.

MixMate® is a registered trademark of Eppendorf SE, Germany.

Registered trademarks and protected trademarks are not marked in all cases with ® or ™ in this manual.

U.S. Patents are listed on www.eppendorf.com/ip

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OSTRZEŻENIE! Electric shock due to damage to the device or mains/power cord.

- ▶ Only switch on the device if the device and the mains/power cord are undamaged.
- ▶ Only operate devices which have been installed or repaired properly.
- ▶ In case of danger, disconnect the device from the mains/power supply voltage. Disconnect the mains/power plug from the device or the earth/grounded socket. Use the disconnecting device intended for this purpose (e.g., the emergency switch in the laboratory).

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Touching parts under high voltage can cause an electric shock. Electric shocks cause injuries to the heart and respiratory paralysis.

- ▶ Ensure that the housing is closed and undamaged.
 - ▶ Do not remove the housing.
 - ▶ Make sure that no liquids can enter the device.
- Only authorized service staff may open the device.

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**OSTRZEŻENIE! Danger due to incorrect voltage supply.**

- ▶ Only connect the device to voltage sources which correspond with the electrical requirements on the name plate.
- ▶ Only use earth/grounded sockets with a protective earth (PE) conductor.
- ▶ Only use the mains/power cord supplied.

**UWAGA! Damage to the display due to mechanical pressure.**

- ▶ Do not apply mechanical pressure to the display.

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UWAGA! Damage to electronic components due to condensation.

Condensate may form in the device when it has been transported from a cool environment to a warmer environment.

- ▶ After installing the device, wait for at least 3 h. Only then connect the device to the mains/power line.

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
- ▶ Do not use any aggressive chemicals on the device or its accessories, such as strong and weak bases, strong acids, acetone, formaldehyde, halogenated hydrocarbons or phenol.
 - ▶ If the device has been contaminated by aggressive chemicals, clean it immediately using a mild cleaning agent.
-

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Select the location for the device according to the following criteria:

- Mains/power connection in accordance with the name plate
- Minimum distance to other devices and walls: 10 cm
- Do not place the device in a wet location.
- Resonance free table with horizontal even work surface
- The location must be well ventilated
- The location is protected against direct sunlight

The mains/power switch and the disconnecting device for the mains/power line must be easily accessible during operation (e.g., a residual current circuit breaker).

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OSTRZEŻENIE! Device fire hazard due to the ingress of liquid.

Liquid ingress can cause a fire due to a short circuit in the device.

- ▶ Do not allow any liquids to enter the inside of the housing.
- ▶ Only mix in sealed tubes and plates.
- ▶ If there was an ingress of liquid: switch off the device, pull the mains/power plug, and have the device cleaned by service technicians who have been authorized by Eppendorf.

UWAGA! Damage due to aggressive chemicals.

- ▶ Do not use any aggressive chemicals on the device or its accessories, such as strong and weak bases, strong acids, acetone, formaldehyde, halogenated hydrocarbons or phenol.
- ▶ If the device has been contaminated by aggressive chemicals, clean it immediately using a mild cleaning agent.

UWAGA! Corrosion due to aggressive cleaning agents and disinfectants.

- ▶ Do not use any corrosive cleaning agents, aggressive solvents or abrasive polishes.
- ▶ Do not incubate the accessories in aggressive cleaning agents or disinfectants for longer periods.

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ZAGROŻENIE! Electric shock due to the ingress of liquid.

- ▶ Switch off the device and disconnect it from the mains/power line before starting cleaning or disinfection.
- ▶ Do not allow any liquids to enter the inside of the housing.
- ▶ Do not spray clean or spray disinfect the housing.
- ▶ Do not reconnect the device to the mains/power line unless both the inside and outside of the device are completely dry.

**OSTRZEŻENIE! Device fire hazard due to the ingress of liquid.**

Liquid ingress can cause a fire due to a short circuit in the device.

- ▶ Do not allow any liquids to enter the inside of the housing.
- ▶ Only mix in sealed tubes and plates.
- ▶ If there was an ingress of liquid: switch off the device, pull the mains/power plug, and have the device cleaned by service technicians who have been authorized by Eppendorf.

**UWAGA! Damage from UV and other high-energy radiation.**

- ▶ Do not use UV, beta, gamma, or any other high-energy radiation for disinfection.
- ▶ Avoid storage in areas with high UV radiation.

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**OSTRZEŻENIE! Risk to health from contaminated device.**

1. Observe the information contained in the decontamination certificate. It is available as a PDF document on our webpage (www.eppendorf.com/decontamination).
 2. Decontaminate all parts to be shipped.
 3. Include the fully completed decontamination certificate in the shipment.
-

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Mains/power connection:	220 to 240 V $\pm 10\%$, 50 to 60 Hz 100 to 120 V $\pm 10\%$, 50 to 60 Hz
Power consumption:	40 W
Overvoltage category:	II

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Ambience	For indoor use only
Ambient temperature	2 °C – 40 °C
Relative humidity	10 % – 75 %
Atmospheric pressure	Use up to an altitude of 2000 m above MSL.
Pollution degree	2
Protection class	I

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Declaration of Conformity

The product named below fulfills the requirements of directives and standards listed. In the case of unauthorized modifications to the product or an unintended use this declaration becomes invalid. This declaration of conformity is issued under the sole responsibility of the manufacturer.

Product name:

MixMate®

Product type:

Mixer for micro test tubes and plates

Relevant directives / standards:

2014/35/EU: EN 61010-1, EN 61010-2-051

2014/30/EU: EN 61326-1

2011/65/EU: EN IEC 63000
(incl. (EU) 2015/863)

Further applied standards: IEC 61010-1, IEC 61010-2-051
UL 61010-1, UL 61010-2-051
CAN/CSA C22.2 No. 61010-1-12, CAN/CSA C22.2 No. 61010-2-051

Hamburg, April 04, 2022

Dr. Wilhelm Plüster
Management Board

Dr. Marlene Jentzsch
Senior Vice President
Business Division Separation & Instrumentation

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com

Eppendorf®, the Eppendorf Brand Design and MixMate® are registered trademarks of Eppendorf SE, Germany.
U.S. Design Patents are listed on www.eppendorf.com/ip.
All rights reserved, incl. graphics and images. Copyright ©2022 by Eppendorf SE.

www.eppendorf.com

ISO
9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf Certificate

Declaration of Conformity

Product: MixMate®
 Order No. International: 5353 000.014
 Order No. North America: 022674200

The device was tested in accordance with EN ISO 5349-1:2001 "Mechanical vibration – Measurement and assessment of human exposure to hand-transmitted vibration – Part 1"

Hand- arm vibration assessment of MixMate® in accordance with EN ISO 5349-1:2001

Tests were performed on the MixMate operated in the "touch vortex" operation mode (3,500 rpm) with 15 mL and 50 mL conical tubes. The total vibration value was measured, the daily vibration exposure and the exposure points calculated from it.

The exposure action value **EAV of 2.5 m/s² A(8)** is equivalent to 100 exposure points and is the level at and above which hand-arm-vibration management procedures must be initiated where employees are regularly exposed to this level of vibration exposure.

The exposure limit value **ELV of 5 m/s² A(8)** is equivalent to 400 exposure points and must not be exceeded on any work day. If it is, immediate measures must be taken to control vibration levels or reduce exposure times to limit daily vibration exposure to below the ELV.

Usage pattern: Each test consisted of a representative usage pattern equivalent to 5 successive touch vortex operations of 30 seconds with 5 second intervals in a total of 2 minutes 30 seconds "on" and 20 seconds "off". Tubes were filled with 2/3 of water.	Tube	Total vibration value	Daily exposure value	1 HOUR exposure points	Time - single tool use: For not exceeding	
		m/s ²	m/s ²		EAV	ELV
					2.5 m/s ² minutes	5 m/s ² minutes
	15 mL	5.2	1.8	54	111	444
	50 mL	9.4	3.3	177	34	136

Eppendorf Certificate

Exposure Action Value (EAV)

The table indicates that the operator using 15 mL conical tubes in touch vortex mode (3,500 rpm) with reported usage would reach the EAV in 111 minutes and using 50 mL conical tubes in 34 minutes. Based on a usage pattern of 30 seconds vortex/5 seconds gap for changing the tube, this means that approximately 190 × 15 mL or 58 × 50 mL conical tubes could produce exposures that reach the EAV. With a usage pattern of 15 seconds vortex/5 seconds gap, 333 × 15 mL tubes or 102 × 50 mL tubes may be vortexed without exceeding the EAV. With a usage pattern of 45 seconds vortex/5 seconds gap 133 × 15 mL tubes or 40 × 50 mL tubes may be vortexed.

Exposure Limit Value (ELV)

The table indicates that the operator using 15 mL conical tubes in touch vortex mode (3,500 rpm) with reported usage would reach the ELV in 444 minutes and using 50 mL conical tubes in 136 minutes. Based on a usage pattern of 30 seconds vortex/5 seconds gap, this means that approximately 761 × 15 mL conical tubes or 233 × 50 mL conical tubes could produce exposures that reach the ELV. Immediate action must thus be taken to limit exposures to below the ELV. With a usage pattern of 15 seconds vortex/5 seconds gap, 1332 × 15 mL tubes or 408 × 50 mL tubes may be vortexed until the ELV is reached. With a usage pattern of 45 seconds vortex/5 seconds gap 532 × 15 mL tubes or 163 × 50 mL tubes may be vortexed.

Vortex time (5 seconds gap)	Number of tubes to reach EAV/for not exceeding ELV			
	EAV		ELV	
	15 mL	50 mL	15 mL	50 mL
15 seconds	333	102	1332	408
30 seconds	190	58	761	233
45 seconds	133	40	532	163

Date: November 4, 2014

Page 2 of 2

0015 033.509-03

Dr. S. Scheeff
Global Product Manager
Sample Preparation Instruments

Dr. B. Schreiber
Vice President
Quality Management & Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO
13485
Certified

ISO
14001
Certified

Measurement of sound power level according to
DIN EN ISO 3744:1995-11 and DIN EN ISO 11204:1996-7

Allround mixer MixMate

Eppendorf AG, 22331 Hamburg, Germany

1. Machine:

- | | | |
|-----|-----------------------|------------------------|
| 1.1 | Type: | Allround mixer MixMate |
| 1.2 | Serial Number: | 5353 05003* |
| 1.3 | Year of construction: | 2008 |

2. Specifications:

See documents of the manufacturer

3. Measurement system:

- | | | |
|-----|-------------------------------|--------------|
| 3.1 | Calibrated sound level meter: | Norsonic 118 |
|-----|-------------------------------|--------------|

4. Measurement surface:

- | | | |
|-----|-----------------------|-------------------------------------|
| 4.1 | Measurement surface: | hemisphere with 10 measuring points |
| 4.2 | Measurement distance: | 1 m |

5. Measurement conditions:

- | | | |
|-----|-------------------------------------|------------------------------------|
| 5.1 | Environment: | Free field over a reflecting plane |
| 5.2 | Environmental correction K_{2A} : | 0.9 dB(A) |
| 5.3 | accuracy class | 2 |
| 5.4 | Operating conditions: | 1400 rpm |

- | | | |
|----|--|------------|
| 6. | <u>Measurement surface sound pressure level:</u> | 31.2 dB(A) |
| | <u>Sound power level:</u> | 39.0 dB(A) |

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Nr. 109SST127 Date of measurement: 25.04.2009

Dipl. - Ing. C. Michalke

*) Sample measurement for this device only.

Measurement of sound power level according to
DIN EN ISO 3744:1995-11 and DIN EN ISO 11204:1996-7

Allround mixer MixMate

Eppendorf AG, 22331 Hamburg, Germany

1. **Machine:**

- | | | |
|-----|-----------------------|------------------------|
| 1.1 | Type: | Allround mixer MixMate |
| 1.2 | Serial Number: | 5353 05003 * |
| 1.3 | Year of construction: | 2008 |

2. **Specifications:**

See documents of the manufacturer

3. **Measurement system:**

- | | | |
|-----|-------------------------------|--------------|
| 3.1 | Calibrated sound level meter: | Norsonic 118 |
|-----|-------------------------------|--------------|

4. **Measurement surface:**

- | | | |
|-----|-----------------------|-------------------------------------|
| 4.1 | Measurement surface: | hemisphere with 10 measuring points |
| 4.2 | Measurement distance: | 1 m |

5. **Measurement conditions:**

- | | | |
|-----|-------------------------------------|------------------------------------|
| 5.1 | Environment: | Free field over a reflecting plane |
| 5.2 | Environmental correction K_{2A} : | 0.9 dB(A) |
| 5.5 | accuracy class | 2 |
| 5.6 | Operating conditions: | 1650 rpm |

- | | | |
|----|---|------------|
| 6. | <u>Measurement surface sound pressure level:</u> | 33.1 dB(A) |
| | <u>Sound power level:</u> | 41.0 dB(A) |

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Nr. 109SST127 Date of measurement: 25.04.2009

Dipl. - Ing. C. Michalke

*) Sample measurement for this device only.

Measurement of sound power level according to
DIN EN ISO 3744:1995-11 and DIN EN ISO 11204:1996-7

Allround mixer MixMate

Eppendorf AG, 22331 Hamburg, Germany

1. Machine:

- | | | |
|-----|-----------------------|------------------------|
| 1.1 | Type: | Allround mixer MixMate |
| 1.2 | Serial Number: | 5353 05003* |
| 1.3 | Year of construction: | 2008 |

2. Specifications:

See documents of the manufacturer

3. Measurement system:

- | | | |
|-----|-------------------------------|--------------|
| 3.1 | Calibrated sound level meter: | Norsonic 118 |
|-----|-------------------------------|--------------|

4. Measurement surface:

- | | | |
|-----|-----------------------|-------------------------------------|
| 4.1 | Measurement surface: | hemisphere with 10 measuring points |
| 4.2 | Measurement distance: | 1 m |

5. Measurement conditions:

- | | | |
|-----|-------------------------------------|------------------------------------|
| 5.1 | Environment: | Free field over a reflecting plane |
| 5.2 | Environmental correction K_{2A} : | 0.9 dB(A) |
| 5.7 | accuracy class | 2 |
| 5.8 | Operating conditions: | 3000 rpm |

- | | | |
|----|--|------------|
| 6. | <u>Measurement surface sound pressure level:</u> | 46.1 dB(A) |
| | <u>Sound power level:</u> | 54.0 dB(A) |

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Nr. 109SST127 Date of measurement: 25.04.2009

Dipl. - Ing. C. Michalke

*) Sample measurement for this device only.

Evaluate Your Manual

Give us your feedback.
www.eppendorf.com/manualfeedback

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

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TÜV Nord Umweltschutz GmbH & Co. KG