Register your instrument!
www.eppendorf.com/myeppendorf

Top Buret M/H

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Copyright© 2022 Eppendorf SE, Germany. All rights reserved, including graphics and images. No part of this publication may be reproduced without the prior permission of the copyright owner.

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany.

Eppendorf trademarks and trademarks of third parties may appear in this manual. All trademarks are the property of their respective owners. The respective trademark name, representations and listed owners can be found here:
www.eppendorf.com/ip.

U.S. Patents and U.S. Design Patents are listed on www.eppendorf.com/ip.

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订购信息 (International)	描述
	Bottle thread adapter PP
4960 800.040	from GL 32 to GL 25,
4960 800.139	from GL 32 to GL 27
4960 800.058	from GL 32 to GL/S 28
4960 800.120	from GL 45 to GL 32
4960 800.155	from GL 45 to GL 38
4960 800.147	from GL 45 to S 40
4960 800.082	from GL 32 to NS 19/26
4960 800.090	from GL 32 to NS 24/29
4960 800.104	from GL 32 to NS 29/32

Eppendorf Certificate

Declaration of Conformity – China RoHS 2 for Eppendorf Top Buret® M/H

Eppendorf has made reasonable efforts to ensure that hazardous materials and substances may not be used in the Eppendorf Top Buret M/H.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB/T 26572 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

- > Lead (Pb): 0.1%
- > Mercury (Hg): 0.1%
- > Cadmium (Cd): 0.01%
- > Hexavalent Chromium (Cr(+VI)): 0.1%
- > Polybrominated Biphenyls (PBB): 0.1%
- > Polybrominated Diphenyl Ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not cause serious environmental pollution, serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Eppendorf Top Buret M/H is 40 years.

Products manufactured by Eppendorf may become components of other devices or can be used with other appliances. With these third-party products and devices in particular, please note the EFUP labeled on these products. Eppendorf will not take responsibility for the EFUP of those products and devices.

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This table is created according to SJ/T 11364-2014.

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Note: Batteries, glassware and accessories might not be content of the enclosed device and/or may have their own EFUP-marking and/or might be maintaining parts with changing EFUP-marking.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with Lead (Pb), Mercury (Hg), Cadmium (Cd), hexavalent Chromium (Cr(+VI)), polybrominated Biphenyls (PBB), and polybrominated Diphenyl Ethers (PBDE).

Date: Hamburg, July 21, 2022

Page 2 of 2

Monika Schneider
Vice President Global Quality Management
& Regulatory Affairs

Dr. Christian Eggert
Senior Vice President
Business Division Liquid Handling

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

Eppendorf, the Eppendorf Brand Design and Top Buret® are registered trademarks of Eppendorf SE, Germany. All rights reserved, including graphics and images. Copyright 2022 © by Eppendorf SE.

www.eppendorf.com

1812TopBuretCRoHS2-02

Declaration of Conformity

The product named below fulfills the requirements of directives and standards listed. In the case of unauthorized modifications to the product or an unintended use this declaration becomes invalid. This declaration of conformity is issued under the sole responsibility of the manufacturer.

Product name:

Eppendorf Top Buret M/H 4965

Product type:

Bottletop Buret

Relevant directives / standards:

2014/30/EU: EN 61326-1

2011/65/EU: EN IEC 63000
(incl. EU 2015/863)

Date: December 3rd, 2021

Dr. Wilhelm Plüster
Management Board

Dr. Christian Eggert
Senior Vice President
Business Division Liquid Handling

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com

Eppendorf® and the Eppendorf logo are registered trademarks of Eppendorf SE, Germany.
All rights reserved, incl. graphics and pictures. Copyright 2021 © by Eppendorf SE.

www.eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified

4965 900.999-06

Evaluate Your Manual

Give us your feedback.
www.eppendorf.com/manualfeedback

Your local distributor: www.eppendorf.com/contact
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com