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Literature

- [1] Fernández O, Izquierdo G, Fernández V, et al. Adipose-derived mesenchymal stem cells (AdMSC) for the treatment of secondary-progressive multiple sclerosis: A triple blinded, placebo controlled, randomized phase I/II safety and feasibility study. *PLoS One* 2018; 13(5):e0195891.
- [2] Baboolal TG, et al. Synovial fluid hyaluronan mediates MSC attachment to cartilage, a potential novel mechanism contributing to cartilage repair in osteoarthritis using knee joint distraction. *Ann Rheum Dis* 2016; 75:908–915.
- [3] Wassef MA, et al. Therapeutic efficacy of differentiated versus undifferentiated mesenchymal stem cells in experimental type I diabetes in rat. *Biochemistry and Biophysics Reports* 2016; 5:468–475.
- [4] Cheng L, Zhao W, Hill AF. Exosomes and their role in the intercellular trafficking of normal and disease associated prion proteins. *Mol. Aspects Med.* 2018; 60:62–68.
- [5] Joo HS, Suh JH, Lee HJ, Bang ES, Lee JM. Current Knowledge and Future Perspectives on Mesenchymal Stem Cell-Derived Exosomes as a New Therapeutic Agent. *Int J Mol Sci.* 2020; 21:727.
- [6] Chen CC, Liu L, Ma F, et al. Elucidation of exosome migration across the blood–brain barrier model in vitro. *Cell Mol. Bioeng.* 2016; 9(4):509–529.
- [7] Yin K, Wang S, Zhao RC. Exosomes from mesenchymal stem/stromal cells: a new therapeutic paradigm. *Biomark Res.* 2019; 7:8.
- [8] June SH, Jinkwan K. Mesenchymal Stem Cell-derived Exosomes: Applications in Cell-free Therapy. *Korean J Clin. Lab. Sci.* 2018; 50:391-398.
- [9] Parallel Fed-batch CHO Culture on SciVario® twin, the Flexible Controller for All Your Bioprocess Needs. *Application Note* 2020; App-432.
- [10] Siddiquee K, Sha M. Microcarriers-based Expansion of Adipose-Derived Mesenchymal Stem Cells in Shake Flasks. *Bioprocessing Journal* 2013; 12:11-17.
- [11] Teixeira FG, Panchalingam KM, Anjo SI, et al. Do hypoxia/normoxia culturing conditions change the neuroregulatory profile of Wharton Jelly mesenchymal stem cell secretome? *Stem Cell Res Ther.* 2015; 6:133.
- [12] Thippabhotla S, Zhong C, He M. 3D cell culture stimulates the secretion of in vivo like extracellular vesicles. *Sci Rep* 2019; 9:13012.

Ordering information

Description	Order no. international
SciVario® twin Fermenter/Bioreactor Control System, base unit, 100 – 240 V/50/60 Hz, for 2 vessels	7600100001
BioBLU® 1c Single-Use Vessel, cell culture, open pipe, 2 pitched-blade impellers, optical pH, sterile, 4 pieces	1386110500
Centrifuge MiniSpin® plus, non-refrigerated, with Rotor F-45-12-11, 230 V/50 – 60 Hz(CN)	5453000097
Centrifuge 5430 R, keypad, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 230 V/50 – 60 Hz (CN)	5428000090
Incubator CellXpert® C170i, inner door with 4 door segments, handle right side, 1 – 20 % oxygen control, 220 – 240 V/50 – 60 Hz (CN)	6731001044
Recirculation tube (from 3L decanter kit (M1287-119102))	M1287-908402
Adaptor, compression fitting, for Pg 13.5 port, (from 3L decanter kit (M1287-119102))	M1287-503001
Eppendorf Conical Tubes, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 50 mL	0030122178
Eppendorf Conical Tubes, sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 15 mL	0030122151

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · Barkhausenweg 1 · 22339 · Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com/bioprocess

www.eppendorf.com/SciVario