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1. Margiana R. et al., "Clinical application of mesenchymal stem cell in regenerative medicine: a narrative review" Stem Cell Research & Therapy, 2022, 13:366.
2. Abbaszadeh H. et al., "Regenerative potential of Warthon's jelly-derived mesenchymal stem cells: A new horizon of stem cell therapy" Journal of Cellular Physiology 2020, 235:12; 9230-9240.
3. Yuan X. et al., "Strategies for improving adipose-derived stem cells for tissue regeneration" Burns & Trauma 2022, 10, tkac028.
4. Rodriguez-Fuentes D.E. et al., "Mesenchymal stem cells current clinical applications: A systematic review" 2021, 52:93-101.
5. Olson TR. et al., "Peak MSC – are we there yet?" Front Med. 5:178. (5)
6. Martin C. et al., "Revisiting MSC expansion from critical quality attributes to critical culture process parameters" Process Biochem. 2017, 59, 231-243. (6)
7. Kusuma GD. et al., Effect of 2D and 3D culture microenvironments on mesenchymal stem cell-derived extracellular vesicles potencies. Front. Cell Dev. Biol. 10:819726. (6)
8. Egger D. et al., "Physiologic isolation and expansion of human mesenchymal stem/stromal cells for manufacturing of cell-based therapy products" Engineering in life Sciences 2021, 27;22(3-4):361-372. (8)
9. Lin Y.M. „Expansion in microcarriers-spinner cultures improves the chondrogenic potential of human early mesenchymal stromal cells" Cytotherapy 2026, 18;740-753. (7)
10. Salazar-Noratto G.E. et al., "Understanding and leveraging cell metabolism to enhance mesenchymal stem cell transplantation survival in tissue engineering and regenerative medicine applications. Stem Cells. 2020;38:22-33.
11. Folmes C.D.L. et al., "Metabolic plasticity in stem cell homeostasis and differentiation." Cell Stem Cell. 2012;11(5):596-606.
12. Zhu H. et al., "Inducible metabolic adaptation promotes mesenchymal stem cell therapy for ischemia." Arterioscler Thromb Vasc Biol. 2014; 34:870-876.
13. Yuan X. et al., "Metabolism in human mesenchymal stromal cells: a missing link between hMSC biomanufacturing and therapy." Front. Immunol. 2019;10:977.
14. Jossen V. et al., "Theoretical and practical issues that are relevant when scaling up hMSC microcarrier production processes" Stem Cell Intl. 2016:4760414
15. Tsai AC. et al., "Influence of microenvironment on mesenchymal stem cell therapeutic potency: from planar culture to microcarriers", Frontiers in Bioengineering and Biotechnology, 2020, 24;8;640.
16. Yin J. et al., "Manufacturing of primed mesenchymal stromal cells for therapy" Nat. Biomed.End. 2019 2:90-104.

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电子邮件: marketinfo@eppendorf.cn

www.eppendorf.com/bioprocess

Beckman Coulter® and Vi-CELL® are registered trademarks of Beckman Coulter, Inc., USA. Synthemax® and Corning® are registered trademarks of Corning Inc., USA. RoosterBio® is a registered trademark of RoosterBio Inc., USA. RoosterNourish™ and RoosterReplenish™ are trademarks of RoosterBio Inc., USA. CellTiter-Glo® and GloMax® are registered trademarks of Promega inc., USA. Glucose-Glo™, Lactate-Glo™, and Glutamine/Glutamate-Glo™ are trademarks of Promega inc., USA. Eppendorf®, the Eppendorf Brand Design, BioBLU®, CellXpert®, and ThermoMixer® are registered trademarks of Eppendorf SE, Germany. DASbox® and DASware® is a registered trademark of DASGIP Information and Process Technology GmbH, Germany. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE.

Eppendorf SE reserves the right to modify its products and services at any time. This application note is subject to change without notice. Although prepared to ensure accuracy, Eppendorf SE assumes no liability for errors, or for any damages resulting from the application or use of this information. Viewing the application note alone cannot as such provide for or replace reading and respecting the current version of the operating manual.