

# Borrelia burgdorferi

Multiporator/Eppendorf Eporator®

## Transformation Protocol

Protocol No. 4308 915.506 – 12/2001

<b>Microorganism</b>	<i>Borrelia burgdorferi</i> B31
<b>Cell type</b>	Bacteria, gram negative
<b>Molecules injected</b>	Linear DNA (in water)
<b>Growth medium</b>	Barbour-Stoenner-Kelly (BSK) II medium without gelatin
<b>Washing solution</b>	PBS, without divalent cations; EPS (272 mM sucrose, 15% glycerol)
<b>Electroporation solution</b>	EPS (272 mM sucrose, 15% glycerol)
<b>Outgrowth medium</b>	BSK II medium
<b>Cuvette</b>	2 mm gap width
<b>Reference</b>	Samuels D.S., et al • 1994 • Journal of Bacteriology 176 • 6045-6049

### Making electrocompetent cells:

1. Grow cells in 500 ml BSK II at 32 °C to a cell density of  $3-7 \times 10^7$  cells/ml.
2. Harvest by centrifugation (4,000 x g, 20 min., 4 °C). Wash twice in 60 ml cold PBS and pellet by centrifugation (3,000 x g, 10 min., 4 °C).
3. Wash three times in 20 ml cold EPS and pellet by centrifugation (2,000 x g, 10 min., 4 °C).
4. Resuspend the cells in 0.6 ml cold EPS.

### Electroporation of cells:

1. Add 0.3 to 1 µg DNA to 50 µl of electrocompetent cells on ice. Homogenize by gently mixing with pipette several times. Transfer mixture into a prechilled cuvette.
2. Wipe moisture from the cuvette and insert the cuvette into the device.
3. Electroporation:

<b>Mode</b>	Prokaryotes "O"
<b>Voltage (V)</b>	2,500 V
<b>Time constant (τ)</b>	5 ms

4. Immediately add 1 ml BSK II medium, transfer cells to a 15 ml tube and add another 9 ml of BSK II. Incubate for 24-48 hours at 32 °C with shaking.
5. Plate 0.1 ml of culture on solid BSK II medium. Pellet the remaining 9.9 ml (4,000 x g, 10 min.), resuspend in 1 ml of supernatant and plate on selective plates. Incubate for 14 days, at 32 °C in a humidified 5% CO<sub>2</sub> chamber.

### Expected results:

Transformation efficiency up to  $10^3$  transformants/µg of DNA.

**eppendorf**