

# Staphylococcus aureus

Multiporator/Eppendorf Eporator®

## Transformation Protocol

Protocol No. 4308 915.533 – 01/2002

<b>Microorganism</b>	<i>Staphylococcus aureus</i>
<b>Cell type</b>	Bacteria, gram positive
<b>Molecules injected</b>	Plasmid DNA (pC194)
<b>Growth medium</b>	Basic medium (BM: 1% peptone, 0.5% yeast extract, 0.1% glucose, 0.5% NaCl, 0.1% K <sub>2</sub> HPO <sub>4</sub> )
<b>Washing solution</b>	10% glycerol (in distilled water)
<b>Electroporation solution</b>	10% glycerol (in distilled water)
<b>Outgrowth medium</b>	SMMP50 medium (5.5 parts SMM buffer (1 M sucrose, 0.04 M maleic acid, 0.04 M MgCl <sub>2</sub> , pH 6.5), 4 parts 7% Panassy broth, 0.5 parts 10% BSA)
<b>Cuvette</b>	2 mm gap width

**Reference** Augustin J. and Götz F. • 1990 • FEMS Microbiology Letters 66 • 203-208

### Making electrocompetent cells:

1. Inoculate a flask of BM with a fresh overnight culture. Grow at 37 °C with shaking to an O.D.<sub>578</sub> of about 0.5 to 0.65.
2. Wash in one volume of 10% glycerol at room temperature, followed by washes with 1/2, 1/20 and 1/50 volumes.
3. Resuspend to final volume of  $1-5 \times 10^{10}$  cells/ml in 10% glycerol. Aliquote and freeze at -70 °C.

### Electroporation of cells:

1. Thaw cells for 5 min at room temperature. Add 1-2 µl plasmid DNA (100-500 ng) to 50 µl of electrocompetent cells. Homogenize by gently mixing with pipette several times. Incubate at room temperature for 30 minutes.
2. Transfer mixture into a cuvette and insert the cuvette into the device.
3. Electroporation:

<b>Mode</b>	Prokaryotes "O"
<b>Voltage (V)</b>	1,800 V
<b>Time constant (τ)</b>	5 ms

4. Add 1 ml SMMP50 medium to the cells and transfer to sterile tube. Incubate 90 minutes at 37 °C with moderate shaking.
5. Plate on selective BM plates. Incubate 20 hours at 37 °C.

### Expected results:

Transformation efficiency up to  $1.5 \times 10^5$  transformants/µg of DNA.

**eppendorf**