

Ustilago maydis

Multiporator/Eppendorf Eporator®

Transformation Protocol

Protocol No. 4308 915.536 – 01/2002

| | |
|---------------------------------|---|
| Microorganism | <i>Ustilago maydis</i> fbd11, protoplasts |
| Cell type | Phytopathogenic fungus, yeast-like growing sporidia (diploid) |
| Molecules injected | Plasmid DNA (pNEBUH, pSMUT) |
| Growth medium | YEPS (1% (w/v) yeast extract, 2% (w/v) bacto-tryptone, 2% (w/v) saccharose) |
| Washing solution | SCS (1 M sorbitol, 20 mM sodium citrate (pH 5.8)); 1 M sorbitol |
| Electroporation solution | 1 M sorbitol |
| Selection plates | 1.5% agar in YEPS with 1 M sorbitol in 94 mm diameter plastic petri dishes - 10 ml bottom layer containing 300 µl hygromycin B/ml agar - 10 ml top layer of agar (without hygromycin B) poured off 10 to 20 min. prior to use |
| Cuvette | 2 mm gap width |
| Reference | Robert Fischer • Department of Biology • Genetics • Philipps University of Marburg Karl-von-Frisch-Str. 1 • D-35032 Marburg Phone +49 6421 282 7080 • Fax +49 6421 282 8971 • e-mail: rfischer@mail.uni-marburg.de |

Preparation of protoplasts:

1. Inoculate 50 ml of YEPS with a primary culture of *Ustilago maydis* and grow cells overnight at 30 °C to a cell density of O.D.₆₀₀ of 0.5 to 0.7.
2. Harvest by centrifugation at 3,000 rpm for 5 min at room temperature.
3. Wash one time with 25 ml SCS.
4. Digest the cell walls by resuspending cells in 2 ml filter sterilized SCS containing 3 mg/ml Novozyme-234, incubate until 20 to 50% of protoplasts have been formed (observe by microscope), then add 25 ml 1 M ice-cold sorbitol.
5. Centrifuge for 7 min at 2,300 rpm at 4 °C.
6. Wash three times with 25 ml ice-cold 1 M sorbitol. Centrifuge for 7 min at 2,300 rpm at 4 °C.
7. Resuspend in 1 ml ice-cold 1 M sorbitol (density: 10⁷ to 10⁸ protoplasts/ml). Keep on ice.

Electroporation of protoplasts:

1. Add 1-2 µl pNEBUH DNA or 10-20 µl pSMUT DNA (1 µg/µl) to 400 µl of protoplasts. Homogenize by gently mixing with pipette several times. Transfer mixture into a prechilled cuvette.
2. Wipe moisture from the cuvette and insert the cuvette into the device.
3. Electroporation:

| | |
|--------------------------|-----------------|
| Mode | Prokaryotes "O" |
| Voltage (V) | 480 V |
| Time constant (τ) | 5 ms |

4. Incubation for 10 to 20 min. on ice is possible but not required.
5. Plate on selective plates (freshly poured top layer).

eppendorf

Ustilago maydis

Multiporator/Eppendorf Eporator®

Expected results:

The protoplast viability was not determined because of incomplete treatment with Novozyme-234: death of protoplasts may occur to certain percentage, but probably not of remaining cells with intact cell wall.

Transformation efficiency (transient) approx. 400 colonies/μg autonomously replicating vector pNEBUH.

Transformation efficiency (stable) approx. 1 colony/μg integration vector pSMUT.

Note:

The preparation of protoplasts is very critical to obtain a high transformation efficiency. Cell wall degradation by Novozyme-234 continues during washing steps in 1 M sorbitol causing burst of protoplasts.

eppendorf

Eppendorf AG · 22331 Hamburg · Germany · Phone +49 40-5 38 01-0 · Fax +49 40-5 38 01-556

e-mail: eppendorf@eppendorf.com · Internet: www.eppendorf.com