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## Literature

- [1] <https://cellculturedish.com/questions/two-questions-what-do-you-consider-typical-titer-for-aav-production-and-what-can-you-do-to-increase-it/>
- [2] Gebre MS, Brito LA, et al. Novel approaches for vaccine development. *Cell*. 2021; 184:1589-1603.
- [3] Pollard AJ, Bijker EM. A guide to vaccinology: from basic principles to new developments. *Nature Reviews Immunology*. 2021; 21:83-100.
- [4] Lal H, et al. Efficacy of an adjuvanted herpes zoster subunit vaccine in older adults. *N. Engl. J. Med.* 2015; 372:2087-2096.
- [5] Pardi N, Hogan MJ, Porter FW, Weissman D. mRNA vaccines — a new era in vaccinology. *Nature Reviews Drug Discovery*. 2018; 17:261-279.
- [6] Sasso E, Morena A, Zambrano N et al. New viral vectors for infectious diseases and cancer. *Seminars in Immunology*. 2020; 50:101430.
- [7] Mendonça SA, Lorincz R, Boucher P et al. Adenoviral vector vaccine platforms in the SARS-CoV-2 pandemic. *npj Vaccines*. 2021; 6:97.
- [8] Zabaleta N, Dai W, Bhatt U, et al. An AAV-based, room-temperature-stable, single-dose COVID-19 vaccine provides durable immunogenicity and protection in non-human primates. *Cell Host & Microbe*. 2021; 29:1437-1453.
- [9] Sha M. Vero Cell-based Vaccine Production: Cell lines, Media and Bioreactor Options. *Application Note*. 2021; App-23.
- [10] Xiaofeng (Kevin) Han, Sha M. Vero High-Density Vero Cell Perfusion Culture in BioBLU® 5p Single-Use Vessels. *Application Note*. 2017; App-359.
- [11] Xiaofeng (Kevin) Han, Sha M. Easy Perfusion for Anchorage-dependent Cell Culture using an Eppendorf Vessel equipped with Microcarrier Spin Filter. *Application Note*. 2020; App-414.
- [12] Rimmelzwaan GF, Baars M, Claas ECJ, Osterhaus ADME. Comparison of RNA hybridization, hemagglutination assay, titration of infectious virus and immunofluorescence as methods for monitoring influenza virus replication in vitro. *Journal of Virological Methods*. 1998; 74:57-66.
- [13] Voeten JTM, Brands R, Palache AM, et al. Characterization of high-growth reassortant influenza A viruses generated in MDCK cells cultured in serum-free medium. *Vaccine*. 1999; 17:1942-1950.
- [14] Pau MG, Ophorst C, Koldijk MH, Schouten G, Mehtali M, Uytdehaag F. The human cell line PER.C6 provides a new manufacturing system for the production of influenza vaccines. *Vaccine*. 2001; 19:2716-2721.
- [15] Kistner O, Barrett PN, Mundt W, Reiter M, Schober-Bendixen S, Dorner F. Development of a mammalian cell (Vero) derived candidate influenza virus vaccine. *Vaccine*. 1998; 16:960-968.
- [16] Youil R, Su Q, Toner TJ, et al. Comparative study of influenza virus replication in Vero and MDCK cell lines. *Journal of Virological Methods*. 2004; 120:23-31.
- [17] Vemula SV, Mittal SK. Production of adenovirus vectors and their use as a delivery system for influenza vaccines. *Expert Opinion on Biological Therapy*. 2010; 10:1469-1487.
- [18] Fuenmayor Garcés J, Gutiérrez Granados S, Gòdia Casablanco F, Cervera Gracia L. Cell line and plasmid influence on the efficiency of HEK293 cells transient transfection for the production of VLP-based vaccines. *New Biotechnology*. 2016; 33:S32-S33.
- [19] Escobar Ivirico JL, Sha M. HEK293 Suspension Cell Culture using Eppendorf BioFlo 320 with BioBLU 3c Single-Use Vessels. *Application Note*. 2021; App-447.
- [20] Fakhiri J, Landegger LD, Grimm D. Breaking the sound barrier: Towards next-generation AAV vectors for gene therapy of hearing disorders. *Hearing Research*. 2020; 15:10809.
- [21] Mendell JR, Al-Zaidy SA, Rodino-Klapac LR et al. Current Clinical Applications of In Vivo Gene Therapy with AAVs. *Molecular Therapy*. 2021; 29:464-488.
- [22] Xin KQ, Urabe M, Yang J, et al. A novel recombinant adeno-associated virus vaccine induces a long-term humoral immune response to human immunodeficiency virus. *Hum. Gene Ther.* 2001; 12:1047-1061.
- [23] Du L, He Y, Wang Y, et al. Recombinant adeno-associated virus expressing the receptor-binding domain of severe acute respiratory syndrome coronavirus S protein elicits neutralizing antibodies: implication for developing SARS vaccines. *Virology*. 2006; 353:6-16.
- [24] LeWitt P, Rezai AR, Leehey MA, et al. AAV2-GAD gene therapy for advanced Parkinson's disease: a double-blind, sham-surgery controlled, randomised trial. *Lancet Neurol*. 2011; 10:309-319.
- [25] Li C, Samulski RJ. Engineering adeno-associated virus vectors for gene therapy. *Nat Rev Genet*. 2020; 21:255-272.
- [26] Kuzmin DA, Shutova MV, Johnston NR, et al. The clinical landscape for AAV gene therapies. *Nature Reviews Drug Discovery*. 2021; 20:173-174.

- [27] Graham FL, Smiley J, Russell WC, Nairn R. Characteristics of a Human Cell Line Transformed by DNA from Human Adenovirus Type 5. J. gen. Virol. 1977; 36:59–7.
- [28] Louis N, Eveleigh C, Graham FL. Cloning and Sequencing of the Cellular– Viral Junctions from the Human Adenovirus Type 5 Transformed 293 Cell Line. Virology. 1997; 233:423–429.
- [29] Suttle A, Longworth J, Sha M. From Shaker to Pilot/Production Bioreactor: How Scale Up Assist Using the BioFlo® 720 Bioreactor Control System Can Help Your Antibody Production Workflow. Application Note. 2021; App-436.
- [30] Suttle A, Sha M. Parallel Fed-batch CHO Culture on SciVario® twin, the Flexible Controller for All Your Bioprocess Needs. Application Note. 2020; App-432.
- [31] Suttle A, Sha M. CHO Cell Culture in Eppendorf BioBLU® 10c Single-Use Vessels. Application Note. 2018; App-361.

#### Ordering information

| Description                                                                                                                                                                          | Order no.  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <b>BioBLU® 1c Single-Use Bioreactor</b> , cell culture, open pipe, 2 pitched-blade impellers, optical pH, sterile, 4 pieces                                                          | 1386110500 |
| <b>SciVario® twin Fermenter/Bioreactor Control System</b> , base unit, 100 – 240 V/50/60 Hz, for 2 vessels                                                                           | 7600100001 |
| <b>New Brunswick S41i</b> , 170 L, CO2 incubator shaker with inner shelf and touch screen control, stackable up to 2-fold, 1 (2. optional) shelves, 120 V/60 Hz, orbit 2.5 cm (1 in) | S41I120010 |
| <b>Centrifuge MiniSpin® plus</b> , non-refrigerated, with Rotor F-45-12-11, 120 V/50 – 60 Hz (US)                                                                                    | 022620207  |
| <b>Centrifuge 5430 R</b> , keypad, refrigerated, with Rotor FA-45-30-11 incl. rotor lid, 120 V/50 – 60 Hz (US)                                                                       | 022620601  |
| <b>Eppendorf Conical Tubes</b> , sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 50 mL                                                                   | 0030122178 |
| <b>Eppendorf Conical Tubes</b> , sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 15 mL                                                                   | 0030122151 |

**Your local distributor:** [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · Barkhausenweg 1 · 22339 · Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com/bioprocess](http://www.eppendorf.com/bioprocess)

AAVpro® is a registered trademark of Takara Bio Inc., Japan. FectoVIR® is a registered trademark of POLYPLUS TRANSFECTION Société Anonyme, France. ISM® and Mettler Toledo® are registered trademarks of Mettler-Toledo AG, Switzerland. Pluronic® and Life Technologies® are registered trademarks of Life Technologies Corp. USA. Cedex® is a registered trademark of Roche Diagnostics GmbH, Germany. BioBLU®, SciVario®, Eppendorf®, MiniSpin®, and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. Expi293F™, Orion Star™, and Thermo Scientific™ are trademarks of Thermo Fisher Scientific Inc., USA. All rights reserved, including graphics and images. Copyright © 2022 by Eppendorf SE.