## APPLICATION NOTE

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## Literature

- [1]MARDIROSIAN D, SHA M, RASHID K. Viral Vaccine Production [J].
- [2]BARRETT P N, MUNDT W, KISTNER O, et al. Vero cell platform in vaccine production: moving towards cell culture-based viral vaccines [J]. Expert review of vaccines, 2009, 8(5): 607-18.
- [3]KIESSLICH S, KAMEN A A. Vero cell upstream bioprocess development for the production of viral vectors and vaccines [J]. Biotechnology Advances, 2020, 44: 107608.
- [4]SHA M. Vero Cell-based Vaccine Production: Cell lines, Media and Bioreactor Options [J]. Application Note, 2021.
- [5]GE H. Microcarrier Cell Culture: Principles and Methods [J]. 1999.

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[www.eppendorf.com/bioprocess](http://www.eppendorf.com/bioprocess)

AAVpro® is a registered trademark of Takara Bio Inc., Japan. FectoVIR® is a registered trademark of POLYPLUS TRANSFECTION Société Anonyme, France. ISM® and Mettler Toledo® are registered trademarks of Mettler-Toledo AG, Switzerland. Pluronic® and Life Technologies® are registered trademarks of Life Technologies Corp. USA. Cedex® is a registered trademark of Roche Diagnostics GmbH, Germany. BioBLU®, SciVario®, Eppendorf®, MiniSpin®, and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany. Expi293F™, Orion Star™, and Thermo Scientific™ are trademarks of Thermo Fisher Scientific Inc., USA. All rights reserved, including graphics and images. Copyright © 2022 by Eppendorf SE.