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- [1] Smith GE, Summers MD, Fraser MJ. Production of human beta interferon in insect cells infected with a baculovirus expression vector. *Mol Cell Biol.* 1983 Dec;3(12):2156-65. doi: 10.1128/mcb.3.12.2156-2165.1983. PMID: 6318086; PMCID: PMC370086.
- [2] Pidre, M.L.; Arrías, P.N.; Amorós Morales, L.C.; Romanowski, V. The Magic Staff: A Comprehensive Overview of Baculovirus-Based Technologies Applied to Human and Animal Health. *Viruses* 2023, 15, 80. <https://doi.org/10.3390/v15010080>
- [3] Wickham, T.J., Davis, T., Granados, R.R. et al. Baculovirus defective interfering particles are responsible for variations in recombinant protein production as a function of multiplicity of infection. *Biotechnol Lett* 13, 483–488 (1991). <https://doi.org/10.1007/BF01049204>
- [4] Imasaki T, Wenzel S, Yamada K, Bryant ML, Takagi Y. Titer estimation for quality control (TEQC) method: A practical approach for optimal production of protein complexes using the baculovirus expression vector system. *PLoS One.* 2018 Apr 3;13(4):e0195356. doi: 10.1371/journal.pone.0195356. PMID: 29614134; PMCID: PMC5882171.
- [5] Fradin S. Development of a Scale-Down Model for rAAV Viral Vector Production Using a Sf9/BEV System. Eppendorf Application Note 303.

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电子邮件: marketinfo@eppendorf.cn

www.eppendorf.com

Abcam® is a registered trademark of Abcam PLC, United Kingdom. Mettler Toledo® is a registered trademark of Mettler-Toledo GmbH, Switzerland. ABM® is a registered trademark of Applied Biological Materials Inc., Canada. VI-CELL® is a registered trademark of Beckman Coulter, Inc., USA. Beckman Coulter® is a registered trademark of Beckman Coulter, Inc., USA. ThawSTAR® is a registered trademark of MedCision, LLC, USA. MedCision® is a registered trademark of Biocision, LLC, USA. Thermo Fisher Scientific® is a registered trademark of Thermo Fisher Scientific Inc., USA. Corning® is a registered trademark of Corning Inc., USA. Eppendorf®, the Eppendorf Brand Design, and BioBLU® are registered trademarks of Eppendorf SE, Germany. BioFlo® and Innova® are registered trademarks of Eppendorf, Inc., USA.

Eppendorf SE reserves the right to modify its products and services at any time. This application note is subject to change without notice. Although prepared to ensure accuracy, Eppendorf SE assumes no liability for errors, or for any damages resulting from the application or use of this information. Viewing the application note alone cannot as such provide for or replace reading and respecting the current version of the operating manual. All rights reserved, including graphics and images. Copyright © 2023 by Eppendorf SE.