

## BioFlo® 120 - Software Update

### Instructions:

Determine which software version is currently installed in your controller:

- > Turn on your control station
- > Select the Setup button (Unit View) from the bottom
- > Select the Software Tab
- > View the software version

If the installed software is lower than the latest version listed, perform the upgrade as listed below. You may want to print this page for easy access to the instructions while you work.

NOTE: Some recipes may not be compatible with the system after firmware upgrades. We recommend manually recording all recipe values prior to upgrading your firmware, as well as backing them up on an external USB memory stick before performing the upgrade.

1. Click the link from the Software downloads page for the corresponding control station
2. Choose a location on your hard-drive for the download (suggest "desktop")
3. After download is complete close the dialog box
4. Extract the downloaded .ZIP file to your hard drive (suggest "desktop")
5. Install a USB external memory device into the computer's USB drive connection
6. Go to the Extracted .EXE file with in extracted file on your hard drive and double click the file
7. Pop up menu comes up – Select Run
8. From the pop up menu - under "Unzip to folder" – click "Browse" and select USB drive (Removable Disk)
9. Select Unzip button and Okay
10. After download is complete select "Close" to close the dialog box
11. Remove the USB external memory device from the computer and insert it into your control station USB port
12. On your control station, select the 'Setup' tab
13. Select the 'Software' tab
14. Select the 'Refresh' button located in the 'Software' panel
15. In the 'Updates' list, select the file to be installed
16. Select the 'Update' button and the update routine will begin
17. After the update is complete – you will come to the main screen
18. Verify that the new version is running by selecting the 'Setup' button and then the 'System Settings' tab. Under "Software" confirm you have latest revision
19. The unit will now be running with the new version that was installed on it
20. The USB external memory device can be removed from the unit

Congratulations! You have installed the latest performance enhancements. Check the Eppendorf website regularly for new revisions.