

# BioFlo® 310 / 415 & CelliGen® 310 Initial Load Software

**NOTE: This update is only to be used when installing a NEW CPU board!**

Date/File/Version	Description
6-August-2013 <b>Update_M1297_Initial_L.01.exe</b> Version L.01	BioFlo 310 / 415 & CelliGen 310 software for initial installation.

## Instructions for loading software on a new computer board

On a new computer board it will be necessary to restore the operating system. This is done automatically; the only requirements will be to calibrate the touchscreen and then to update the software to latest version.

During this process the screen may be blank or may appear as if the system is not updating;  
**do not restart the unit! This process will take about 10 minutes to complete.**

1. Go to the Saved file on your hard drive and double click the file
2. Pop up menu comes up - Select Run
3. Select Unzip button and Okay
4. After download is complete select "Close" to close the dialog box
5. Go to folder 'local v234 d1'
6. Double click Setup.EXE
7. Click next on welcome Pop Up
8. Enter name and company
9. Select Directory
10. Select program folder
11. Click next to start copying files
12. When prompted to insert disk 2, Click Browse
13. Select 'BC LOCAL' folder from saved location
14. Select 'local v234 d2'
15. Click OK on next two screens
16. When prompted to insert disk 3, Click Browse
17. Select 'BC LOCAL' folder from saved location
18. Select 'local v234 d3'
19. Click OK on next two screens
20. Click Finish