Keep this rotor instruction manual and the centrifuge manual in the file.

Swing Rotors for Refrigerated Centrifuge

R3S

Before using these rotors, please carefully read this instruction manual and the centrifuge instruction manual for its efficient operation and for your safety.

Keep this instruction manual for your reference and refer to it as required.

Contents

1. Specifications	_____	1
2. Tubes/Bottles	_____	4
3. How to Use	_____	8
4. Maintenance After Operation	_____	11
5. Maintenance	_____	12
6. Packing List	_____	13
7. Decontamination	_____	14
8. Rotor Retirement	_____	14
(Appendix) K Factors	_____	15

Safety Reminder

- Centrifuge rotors rotating at high speed have considerable potential for damage to personal properties if used improperly.
For safe and proper use of this rotor, carefully read the centrifuge instruction manual and this rotor instruction manual before use and observe the instructions.

 **WARNING :** and  **CAUTION :** notes are used to call your attention in this manual to prevent personal injury or damage to the rotor and the centrifuge.
These notes are defined as follows.

-  **WARNING:** indicates a potentially hazardous situation which, if not avoided, could result in personal severe injury or possible death.
-  **CAUTION:** indicates a hazardous situation which, if not avoided, could result in personal injury or severe damage to the instrument.

WARNING

- Never use any material capable of producing flammable or explosive vapors.
- Your centrifuge and rotor are not designed to confine any sample particles dispersed due to leakage. Therefore, when using toxic or radioactive samples or pathogenic or infectious blood samples, make sure to prepare necessary safety measures at your own responsibility.
- Never exceed the maximum speed of the rotor (mentioned on the rotor body).
Always reduce rotor speed as instructed in this manual when rotor speed is limited due to sample density, kinds of microplates, etc.
- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor.
- If the centrifuge, rotor, or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the centrifuge, rotor, or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the centrifuge, rotor, or the accessory properly before requesting repairs from an authorized sales or service representative.
- It is your responsibility to sterilize and/or decontaminate the centrifuge, rotor, or parts properly before returning them to an authorized sales or service representative.

CAUTION

- Balance the tubes/bottles and the samples including the caps, adapters, or tube racks within the allowable imbalance of the rotor. Do not exceed the allowable imbalance.
- Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.
- Use the tubes/bottles within their actual capacities.
- Do not use tubes/bottles that have exceeded their life expectancy. Failure to do so could result in damage of tubes/bottles and the rotor and the centrifuge. The life expectancy of tubes/bottles depends on factors such as the characteristics of samples, speed of the rotor used, and temperature. Always check for deterioration and damage (cracks, deformation, and so on) on tubes/bottles before using them. Do not use the tubes/bottles if you find such a problem.
- Inspect and maintain the rotor after use. If abnormality is observed, do not use it.
Contact an authorized sales or service representative.

1. Specifications

1.1 Construction and Cross-sectional view of the rotor

The construction and cross-sectional view of the rotor are shown below.


Figure1-1 Construction of R3S rotor


Figure1-2 Cross-sectional view of R3S rotor
(Buckets [standard accessories])

1.2 Specifications

Items	Model
	R3S
Applicable centrifuge	For applicable centrifuges, refer to the “Applicable centrifuges (Rotors for high-speed refrigerated centrifuges (Part No. S998611))”. (Some centrifuges are on sale in only Japan. Visit our web site.)
Maximum speed	3,000 rpm*
Maximum RCF	1,940 xg
Nominal rotor capacity	2,000 mL (500 mL X 4)
Rotor No.	40
Number of buckets	4 pcs.
Weight	6.0 kg (including 500PP tubes) 6.9 kg (including tube racks)
Material of the rotor body	Aluminum alloy

*: Maximum speed depends on buckets (see “1.3 Rotor speed limitation depending on buckets”).

1.3 Rotor speed limitation depending on buckets

**WARNING :**

Do not exceed the maximum speed of the rotor. The rotor speed should be limited depending on sample characteristic and weight.

Do not exceed the allowable speed of the rotor.

Applicable buckets		Allowable speed Maximum RCF
Order no.	Name	
Standard accessories (4 pcs.)	Bucket (Material: Stainless steel)	3,000 rpm 1,940 xg
5721 244 317 (4 pcs.)	MP Bucket (Material: Aluminum alloy, Stainless steel, Chloroprene rubber)	3,000 rpm 1,670 xg
5721 244 316 (4 pcs.)	M Bucket (Material: Stainless steel)	2,000 rpm 700 xg

**CAUTION:**

Do not use tubes/bottles/adapters/tube racks and so on other than them specified in “2. Tubes/Bottles” of this manual. Failure to do so can cause them to become damaged during centrifugation.

1.4 Allowable speed of the rotor when using heavy samples

Limit the weight of contents (total weight of the tube rack, tubes, bottles, samples, etc.) in each bucket within the weight listed in Table 1-1 to operate this rotor at the maximum speed (see Section 1.3).

Table 1-1

Buckets	Allowable weight
Bucket (Standard accessories)	760 (g)
MP Bucket	500 (g)
M Bucket	300 (g)

CAUTION :

If the total weight of contents in the bucket exceeds the weight listed in Table 1-1, calculate the allowable speed according to the following equation and operate under the allowable speed. Operation over the allowable speed may damage the buckets and result in severe damage to the centrifuge.

- When using Bucket (standard accessories), see below.

$$\text{Allowable speed (rpm)} = 3,000 \text{ (rpm)} \times \sqrt{\frac{760 \text{ (g)}}{\text{Total weight of contents in the bucket (g)}}}$$

- When using MP Bucket, see below.

$$\text{Allowable speed (rpm)} = 3,000 \text{ (rpm)} \times \sqrt{\frac{500 \text{ (g)}}{\text{Total weight of contents in the bucket (g)}}}$$

- When using M Bucket, see below.

$$\text{Allowable speed (rpm)} = 2,000 \text{ (rpm)} \times \sqrt{\frac{300 \text{ (g)}}{\text{Total weight of contents in the bucket (g)}}}$$

1.5 Characteristic of rotor (Buckets [standard accessories])

Rotor speed (rpm)	Tube rack			500PP Tube		
	Maximum RCF (xg)		K factor ※	Maximum RCF (xg)		K factor ※
	Rmax. (18.6 cm)	Rmin. (9.7 cm)		Rmax. (19.3 cm)	Rmin. (9.7 cm)	
1,000	208	108	164,909	216	108	174,267
2,000	832	434	41,227	863	434	43,567
3,000	1,870	976	18,323	1,940	976	19,363

※For details about K factor, refer to “(Appendix) K Factors”

Calculation formula

$$\text{RCF} = 1.118 \times 10^{-5} \times R \times N^2 (\times g)$$

R: Rotating radius (cm)

N: speed (rpm)

2. Tubes/Bottles

2.1 Applicable tubes/bottles

 **CAUTION:**

- Do not use tubes/microplates which are not listed. If tubes/microplates which do not match for adapters or tube racks are used, tubes/microplates and adapters or tube racks may be damaged.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise, the sample leakage may occur during operation.
- Use the tubes/microplates that are same products. Otherwise, the tubes/microplates may be damaged due to their different form or the imbalance operation may occur. And it may result in the damage to the centrifuge and the rotor.
- Use the tubes/microplates at a temperature between 4°C and 25°C. Otherwise, tubes/microplates may be broken or deformed during operation.
- Before using tubes/microplates on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- Before operation, check that the tube/bottle does not touch the rotor body when the bucket is swung. If the tube/bottle touches the rotor body when the bucket is swung, the tube/ bottle and the rotor body might be damaged during operation.
- Use the 50mL self-standing culture tube only in the tube rack 29.5X5(B). The 50mL self-standing culture tube cannot be used in the tube rack/adaptor other than the tube rack 29.5X5(B).

- Microplate buckets, tube racks and tubes are not included in this rotor. Please order applicable ones selecting from the list on the following pages (designate the Order no.). Applicable tubes etc. are also listed on the following pages.
- When using tubes or microplates on the market, perform operation under the allowable speed or the allowable RCF specified by the manufacturer and do not use any tube or microplates whose allowable speed or RCF is not specified. Otherwise the tubes or microplates may be broken during operation. We will check the applicability of tubes or microplates upon request. Contact an authorized sales or service representative.
- Clean and sterilize tubes or microplates on the market according to the instructions of the manufacturer.
- Check the size of tubes before using the tubes (see the following list). When the tube length is longer than the tube length shown in the following list, you cannot use the tube in your adapters/tube racks.
- Maximum speeds listed are guidelines only. Because of variances in user methodologies, no guarantee of performance is expressed or implied.

Applicable tube racks and tubes for the buckets (standard accessories)

Tube rack			Tube				Max. speed (Max. RCF) (Max. radius)															
Order no.	Name	Nominal Capacity × Qty	Order no.	Name	Actual Capacity	Size (mm)																
_____	_____	_____	5721 411 064	500PP tube (4 pcs.)	500mL	(D) 103 × (L) 108	3,000rpm(1,940 × g) (Rmax.:19.3cm)															
5721 244 301	Tube rack BL (4 pcs.) (Material:*1)	10mLX120	5721 411 042	10 glass tube (50 pcs.)	8mL	(D)12 × (L) 105																
	(Marketed)*2		_____	*5	(D)12 × (L) 120																	
	<table><tr><td>d1</td><td>12</td></tr><tr><td>d2</td><td>14</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>55</td></tr></table>				d1	12	d2	14	L1	120	L2	55										
	d1		12																			
d2	14																					
L1	120																					
L2	55																					
5721 244 302	Tube rack GR (4 pcs.) (Material:*1)	7mLX96	(Marketed)*2	5-7 glass tube	*5	(D)14 × (L)120																
	<table><tr><td>d1</td><td>14</td></tr><tr><td>d2</td><td>16</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>54</td></tr></table>						d1	14	d2	16	L1	120	L2	54								
	d1						14															
	d2						16															
L1	120																					
L2	54																					
5721 244 303	Tube rack YE (4 pcs.) (Material:*1)	15mLX72	5721 411 050	15 glass tube (50 pcs.)	13mL	(D)16.5 × (L)105																
	(Marketed)*2						_____	*5	(D)16.5 × (L)120													
										<table><tr><td>d1</td><td>16.5</td></tr><tr><td>d2</td><td>19</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>54</td></tr></table>				d1	16.5	d2	19	L1	120	L2	54	
										d1	16.5											
d2	19																					
L1	120																					
L2	54																					
5721 244 304	Tube rack OR (4 pcs.) (Material:*1)	50mLX16	5721 411 059	50 glass tube (10 pcs.)	45mL	(D)35 × (L)100																
			(Marketed)*2	_____	*5	(D)36.5 × (L)120																
		10mLX16	5721 411 042	10 glass tube (50 pcs.)	8mL	(D)12 × (L) 105																
			(Marketed)*2	_____	*5	(D)12 × (L) 115																
			<table><tr><td></td><td>50mL</td><td>10mL</td></tr><tr><td>d1</td><td>36.5</td><td>12</td></tr><tr><td>d2</td><td>39</td><td>14</td></tr><tr><td>L1</td><td>120</td><td>120</td></tr><tr><td>L2</td><td>52</td><td>55</td></tr></table>					50mL	10mL	d1	36.5	12	d2	39	14	L1	120	120	L2	52	55	
				50mL	10mL																	
d1	36.5	12																				
d2	39	14																				
L1	120	120																				
L2	52	55																				
5721 244 305	Tube rack BL (B) (4 pcs.) (Material:*1)	10mLX120	5721 411 042	10 glass tube (50 pcs.)	8mL	(D)12 × (L) 105																
							(Marketed)*2	_____	*5	(D)12 × (L) 120												
											<table><tr><td>d1</td><td>12</td></tr><tr><td>d2</td><td>14</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>74</td></tr></table>				d1	12	d2	14	L1	120	L2	74
											d1	12										
d2	14																					
L1	120																					
L2	74																					
5721 244 306	Tube rack GR (B) (4 pcs.) (Material:*1)	7mLX96	(Marketed)*2	5-7 glass tube	*5	(D)14 × (L)120																
							<table><tr><td>d1</td><td>14</td></tr><tr><td>d2</td><td>16</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>75</td></tr></table>				d1	14	d2	16	L1	120	L2	75				
							d1	14														
							d2	16														
L1	120																					
L2	75																					
5721 244 307	Tube rack YE (B) (4 pcs.) (Material:*1)	15mLX72	5721 411 050	15 glass tube (50 pcs.)	13mL	(D)16.5 × (L)105																
							(Marketed)*2	_____	*5	(D)16.5 × (L)120												
											<table><tr><td>d1</td><td>16.5</td></tr><tr><td>d2</td><td>19</td></tr><tr><td>L1</td><td>120</td></tr><tr><td>L2</td><td>73</td></tr></table>				d1	16.5	d2	19	L1	120	L2	73
											d1	16.5										
d2	19																					
L1	120																					
L2	73																					
5721 244 308	Tube rack OR (B) (4 pcs.) (Material:*1)	50mLX16	5721 411 059	50 glass tube (10 pcs.)	45mL	(D)35 × (L)100																
			(Marketed)	_____	*5	(D)36.5 × (L)120																
		10mLX16	5721 411 042	10 glass tube (50 pcs.)	8mL	(D)12 × (L) 105																
			(Marketed)*2	_____	*5	(D)12 × (L) 115																
			<table><tr><td></td><td>50mL</td><td>10mL</td></tr><tr><td>d1</td><td>36.5</td><td>12</td></tr><tr><td>d2</td><td>39</td><td>14</td></tr><tr><td>L1</td><td>120</td><td>120</td></tr><tr><td>L2</td><td>71</td><td>71</td></tr></table>					50mL	10mL	d1	36.5	12	d2	39	14	L1	120	120	L2	71	71	
				50mL	10mL																	
d1	36.5	12																				
d2	39	14																				
L1	120	120																				
L2	71	71																				

<

*1: Polypropylene

*2: For glass tubes on the market, select proper glass tubes made of hard glass or PYREX® glass to prevent breakage of tubes due to centrifugal force. Before using glass tubes on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality. Remove the ground stopper from the glass tube before operation, otherwise the ground stopper and the glass tube might be damaged during operation.

*5: Actual capacities of tubes/bottles on the market may vary depending on the specifications. Refer to the manufacture's catalog.

Applicable tube racks and tubes for the buckets (standard accessories)

Tube rack			Tubes				Max. speed (Max. RCF) (Max. radius)
Order no.	Name	Nominal Capacity × Qty	Order no.	Name	Actual Capacity	Size (mm)	
5721 244 311	Tube rack 17X10(4 pcs.) (Material:*1)	15mLX40	(Marketed)	15 mL culture tube (conical bottom)	15mL	(D)17 × (L) 120	3,000rpm(1,870 × g) (Rmax.:18.6cm)
5721 244 314	50T7S Adapter (4 pcs.) (Material:*1)	50mLX28	(Marketed)	50 mL culture tube (conical bottom)	50mL	(D)30 × (L) 115	3,000rpm(1,870 × g) (Rmax.:18.6cm)
5721 244 315	Tube rack 29.5X5(B) (4 pcs.) (Material:*1)	50mLX 20	(Marketed)	50 mL self-standing culture tube 50 mL culture tube (conical bottom)	50mL	(D)29 × (L) 115	3,000rpm(1,830 × g) (Rmax.:18.2cm)
5721 244 312	Tube rack 13.2X18 (4 pcs.) (Material:*1)	4-7mLX72	(Marketed)	Blood collection tube	*6	(D)12.4-13.2 × (L)75-100	3,000rpm(1,810 × g) (Rmax.:1p8.0cm)
5721 244 313	Tube rack 15.5X18 (4 pcs.) (Material:*1)	6-10mLX72	(Marketed)	Blood collection tube	*6	(D)15.2-15.5 × (L)75-100	3,000rpm(1,810 × g) (Rmax.:18.0cm)
5722 112 015	Tube rack 45.5X2 (4 pcs.) (Material:*1)	100mLX8	(Marketed) *2	100mL centrifuge test tube*3	*5	(D)45 × (L)120	3,000rpm(1,810 × g) (Rmax.:18.0cm)
5721 244 309	Tube rack 61X1 (4 pcs.) (Material:*1)	225mLX4	(Marketed)	Falcon® 225mL conical tube*4	225mL	(D)61 × (L)137	3,000rpm(1,900 × g) (Rmax.:18.9cm)
		175mLX4		Falcon® 175mL conical tube*4	175mL	(D)61 × (L)118	
5721 244 310	Tube rack 62X1(4 pcs.) (Material:*1)	250mLX4	(Marketed)	250mL Nalgene™ bottle	*5	(D)61.8 × (L)127.7	3,000rpm(1,830 × g) (Rmax.:18.2cm)
Unit (mm)			Unit (mm)				

*1: Polypropylene

*2: For glass tubes on the market, select proper glass tubes made of hard glass or PYREX® glass to prevent breakage of tubes due to centrifugal force. Before using glass tubes on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality. Remove the ground stopper from the glass tube before operation, otherwise the ground stopper and the glass tube might be damaged during operation.

*3: Maker: NICHIDEN-RIKA GLASS CO., LTD.

*4: Be sure to use the cushion (Falcon® 352090) to prevent deformation.

*5: Actual capacities of tubes/bottles on the market may vary depending on the specifications. Refer to the manufacture's catalog.

*6: Actual capacities of blood collection tubes may vary depending on the applications. Refer to the manufacture's catalog.

Applicable microplates for the Microplate buckets

Bucket		Microplate				Max. speed (Max. RCF) (Max. radius)
Order no.	Name	Order no.	Name	Actual Capacity	Size (mm)	
5721 244 316	M bucket *7	(Marketed)*9	96-hole microplate	-	130 × 85 × 14	2,000rpm(700 × g) (Rmax.:15.7cm)
			120-hole microplate	-	130 × 92 × 12	
5721 244 317	MP bucket *8	(Marketed)*9	96-hole microplate	-	130 × 85 × 14	3,000rpm(1,670 × g) (Rmax.:16.6cm)
			96-hole deep well plate	-	130 × 85 × 55	

*7: The maximum height of microplates that can be held in a bucket is 32 mm.

*8: The maximum height of microplates that can be held in a bucket is 55 mm.

*9: Before using microplates on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.

3. How to Use

3.1 Preparation of the rotor

- (1) Mount the rotor into the drive shaft.

CAUTION:

Inspect the rotor before use. Refer to "5. Maintenance" for inspecting the rotor.

- (2) Check that there is no foreign matter in the buckets.

- (3) Engage the buckets to the rotor.

WARNING:

Be sure to set the four buckets of the same serial number. (The serial number is mentioned at the bottom of each bucket.)

Check that the buckets are securely engaged with the rotor pins.

CAUTION:

Before using the bucket, be sure to wipe off the water and the sample which remained into the bucket. Otherwise it may be hard to remove the tube racks/adapters from the bucket.

3.2 Preparation of the sample and the tubes/bottles

Inject the sample into the tubes/bottles.

WARNING:

- These rotors and the centrifuge are not explosion-proof. Never use explosive or flammable samples.
- For safety, there are limits on the use of bio-samples which require bio-isolation, such as pathogenic germs and DNA recombination, as well as RI substance in centrifuges. Perform strict safety controls when separating samples containing these substances.

CAUTION:

- Check the attached chemical resistance chart, and do not use any sample inapplicable to the rotor (including buckets), tube racks, and adapters. Using such a sample could corrode or deteriorate them.
- Chemical resistance and the strength of tubes/bottles varies with speed, temperature, and so on. Before using sample, fill tubes/bottles with water, buffer solution, etc. instead of sample and run them at the intended speed, temperature, and so on to ensure that there is no abnormality.
- Do not exceed the actual capacity specified in the manufacturer. Otherwise the sample leakage may occur during operation.
- Use the tubes/bottles that are same products. Otherwise, the tubes/bottles may be damaged due to their different form or the imbalance operation may occur, and it may result in the damage to the centrifuge and the rotor.
- Use the tubes/bottles at a temperature between 4°C and 25°C. Otherwise, tubes/bottles may be broken or deformed during operation.
- Before using tubes/bottles on the market, test them by filling with water instead of sample, and run them at the intended speed to ensure that there is no abnormality.
- It takes much time for this rotor to reach the desired temperature. We recommend that the rotor, buckets, adapters, tubes/bottles, and samples are warmed or cooled to the desired temperature before operation.

3.3 Balancing the microplates

The allowable imbalance of this rotor is 6g (Allowable imbalance: Difference in weight of contents (tube racks/adapters, tubes/bottles, and samples) in the symmetrically arranged buckets).


CAUTION:

- Use the samples that are same in compositions and same in density when placing them symmetrically. Otherwise, the precipitation levels may be different by centrifugation and there may be variations in position of center of gravity. And it may result in damage to the centrifuge and the rotor due to the imbalance operation.
- Balance between symmetrically arranged microplates filled with the same-density sample. Variations in sample density cause imbalance operation. And it may result in the damage to the centrifuge and the rotor.



Use of tubes with caps is recommended to prevent contamination of the sample with dust or water.

Check that there is no foreign matter in the tube racks/adapters then load the tubes. When loading tubes less than the number of tube cavities in the tube racks/adapters, load the tubes symmetrically with respect to the center of the tube racks/adapters (see below). Examples of correct setting and incorrect tube setting are shown below.

[Correct setting]

Arrange tubes symmetrically with regard to the handle.

[Incorrect setting]

- : Tubes are loaded.
○: Tubes are not loaded.

When it is impossible to load the tubes symmetrically for lack of a tube, prepare a dummy tube filled with the same-density solution and load the tubes symmetrically.


3.4 Mounting the rotor and the buckets

Mount the balanced buckets with respect to the center of the rotor. And Mount the rotor to the drive shaft.


WARNING:

- Be sure to set all the buckets, all the tube racks/adapters regardless of the presence or absence of sample. Operation with only two buckets is very dangerous because the rotor may be deformed and the buckets may come off. Never use any bucket or tube rack/adaptor that is not designated for this rotor by us.
- Be sure to set all the buckets of the same serial number (The serial number is shown on the bottom of the bucket). Otherwise, it may cause the imbalance operation and it may result in the damage to the centrifuge and the rotor.
- Securely hook the buckets on the rotor pins. Otherwise the buckets may come off while the rotor is running.

Correct setting

Incorrect setting

CAUTION:

- Do not load tube racks/adapters and tubes/bottles asymmetrically. Asymmetrically loading may cause imbalance operation and the centrifuge may be damaged.

Check that the buckets swing smoothly. If the swinging is not smooth, clean both the rotor pins and the pin holes of the buckets.

3.5 Operation

Run the rotor according to the procedures specified in the centrifuge instruction manual.

3.6 Taking out samples

- (1) After the operation, carefully remove the rotor from the centrifuge.
- (2) Remove the tube racks/adapters from the buckets. Then take out the separated sample from the tubes/bottles.

4. Maintenance After Operation

4.1 Corrosion resistance of the rotor

WARNING:

Check the attached chemical resistance chart (Part No. S999313), and do not use any sample inapplicable to the rotor, buckets, and tube racks/adapters.

See "1. Specifications" of this manual about materials of each part.

4.2 Maintenance of rotor

After the operation, properly maintain the rotor and buckets to prevent corrosion that can cause the rotor and buckets breakage.

● Normal maintenance

CAUTION:

Clean the inside of the drive hole (crown hole) of the rotor and the surface of the drive shaft (crown) of the centrifuge once a month. If the drive hole or the drive shaft is stained or any foreign matter is adhered, the rotor may be improperly installed and come off during operation.

Wash the rotor and buckets with tap water or a dilute solution of neutral detergent and rinse them out with distilled water. Wipe them with a soft cloth. Dry the rotor which the buckets and tube racks/adapters are removed from. Check that the rotor is completely dried and then put a light coat of silicone grease (vacuum grease) (standard accessory of the centrifuge) on the rotor. Store them in a dry place.

● Maintenance after use of a corrosive sample

Wash the rotor and buckets with tap water immediately after the operation. Then perform the normal maintenance.

WARNING:

Wash the rotor and buckets with tap water immediately after the operation. Then perform the normal maintenance.

● Maintenance when foreign substances are adhered to the rotor and buckets

CAUTION:

Use a neutral detergent having a pH between 5 and 9, otherwise the rotor and buckets can be discolored or corroded.

After operation, properly maintain the rotor to prevent corrosion. Adhesion of salt etc. may corrode the rotor body. Clean the rotor body periodically (once a month). If salt or any other foreign matter adheres to the rotor body, wipe it off the rotor body.

4.3 Maintenance of tube racks/adapters

Wash the tube racks/adapters with tap water or a dilute solution of neutral detergent having a pH between 5 and 9.

4.4 Sterilizing the rotor, buckets, and tube racks/adapters

- Sterilize this rotor and buckets in any following methods; ultraviolet rays sterilization, gas sterilization method (ethylene oxide or formaldehyde), or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).

WARNING:

- Never autoclave the rotor and buckets to avoid decrease in strength.
- Never sterilize rotor and buckets by boiling, otherwise they will be brittle.

- Sterilize tube racks/adapters in any following methods; gas sterilization method (ethylene oxide or formaldehyde), or chemical sterilization method (70% ethanol, 3% hydrogen peroxide, 3% formalin).


CAUTION:

- Never autoclave the tube racks/adapters to avoid deformation.
- Never sterilize the tube racks/adapters by boiling, otherwise they will be deformed.
- Do not dip the tube racks/adapters in the formalin (3%) solution more than 2 hours.


5. Maintenance

Periodically check every 100 hours or once a month whether the rotor and buckets are not corroded. Especially be careful of the inside of the bucket, the inside bottom of the bucket (see below), and the crown hole of the rotor bottom, because the rotor and the buckets can be brittle if these portions are corroded. If their surfaces are discolored, dented or cracked, they are corroded. Do not use such a corroded rotor or a bucket. Contact an authorized sales or service representative immediately for further inspection.

Check it.

Check it.

CAUTION:


If there is any abnormality such as corrosion or cracks, stop use of the rotor and the buckets and contact an authorized sales or service representative.

There is no need to apply silicone grease (vacuum grease) to the rotor pins because the rotor pins have a coating for smooth swinging of the buckets. Although this coating may be discolored or partially flake off as the rotor is repeatedly used, there is no problem in performance.

6. Packing List

Check the parts and accessories of the rotor when the rotor is delivered. Contact an authorized sales or service representative if there are a missing item and inferior goods.

No.	Name	Qty	Order no. of spare parts (Qty)
①	Rotor	1	—
②	Bucket	4	—
③	Rotor instruction manual	1	—

7. Decontamination

WARNING:

- If the centrifuge, rotor or an accessory is contaminated by samples that are toxic or radioactive, or blood samples that are pathogenic or infectious, be sure to decontaminate the item according to good laboratory procedures and methods.
- If there is a possibility that the rotor or an accessory is contaminated by samples that might impair human health (for example, samples that are toxic or radioactive, or blood samples that are pathogenic or infectious), it is your responsibility to sterilize or decontaminate the rotor or the accessory properly before requesting repairs from an authorized sales or service representative. Note that we cannot repair the centrifuge, rotor or the accessory unless sterilization or decontamination is completed.
- It is your responsibility to sterilize and/or decontaminate the rotor or parts properly before returning them to an authorized sales or service representative. In such cases, copy the attached decontamination sheet and fill out the copied sheet, then attach it to the item to be returned. We may ask you about the treatment for the rotor or parts if the decontamination is checked and judged as insufficient by us. It is your responsibility to bear the cost of sterilization or decontamination. Note that we cannot repair or inspect the rotor or the accessory unless sterilization or decontamination is completed.

8. Rotor Retirement

After many years of use, there will be inevitably some corrosion or stress corrosion. At some points, the combination of such damage and metal fatigue could make the rotor vulnerable to a failure.

Although a rotor may appear to be in a good condition, you should follow the rotor retirement recommendation shown below.

Rotor	Material	Retire After Years
Refrigerated centrifuge rotors	Aluminum alloy	15
	Stainless steel	

(Appendix) K Factors

• Calculation of separation time

The K factor can conveniently be used which makes it possible to estimate the settling time directly from Svedberg unit(S).

The K factor can be obtained from the formula show below.

$$K = \frac{\ln(R_{\max.}) - \ln(R_{\min.})}{\omega^2} \cdot \frac{10^{13}}{3,600}$$

ω : Angular velocity $\omega = 2\pi/60 \times N$

N: Speed (rpm)

Rmin.: Distance from axis of rotation to top of solution (cm)

Rmax.: Distance from axis of rotation to bottom of tube (cm)

Using the K factor, the settling time t (in hours) is expressed as follows:

$$t = \frac{K}{S}$$

For example, when using the 500PA tube, the K factor at the maximum speed of 3,000 rpm of Model R3S rotor is 19,363. The time taken for particles of 1,000 S to be settled from R min. to R max. using this rotor is as follows:

$$t = \frac{19,363}{1,000} = 19.4 \text{ hours}$$

The time taken for settling of particles of 1,000 S by the same rotor at a speed of 2,000 rpm is:

$$K = \left[\frac{3,000}{2,000} \right]^2 \times 19,363 = 43,567 \quad t = \frac{43,567}{1,000} = 43.6 \text{ hours}$$

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Copyright © 2023 Eppendorf Himac Technologies Co., Ltd. All rights reserved.

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