


VisionNize[®] Lab Suite

Overview of compatible Eppendorf devices

How Devices Can Be Connected



- 1 VisioNize box (hardware hub)
- 2 VisioNize sense gateway router
- 3 VisioNize sense environmental monitoring sensors
- 4 VisioNize sense external temperature sensor



VisioNize

Device Overview

Overview of compatible Eppendorf devices for VisioNize Lab Suite

This list contains Eppendorf devices which can be connected to the platform. Visit www.eppendorf.com/visionize for the most current compatibility list.

VisioNize touch enabled devices

> Eppendorf devices equipped with the VisioNize touch interface can be connected via ethernet cable:

Eppendorf ULT Freezers	Eppendorf CO₂ Incubators	Eppendorf Cyclers
CryoCube® F740i	CellXpert® C170i	Mastercycler® X40
CryoCube® F740hi		Mastercycler® X50a
CryoCube® F740iw	Eppendorf Shakers	Mastercycler® X50s
CryoCube® F740hiw	Innova® S44i*	Mastercycler® X50p
	<small>* only compatible after serial number S44IJH302546</small>	Mastercycler® X50h
	Eppendorf Centrifuge	Mastercycler® X50i
	Centrifuge 5910 Ri	Mastercycler® X50l
		Mastercycler® X50r
		Mastercycler® X50t

VisioNize box enabled

> Eppendorf devices which are connectable via VisioNize box to the VisioNize Lab Suite. The VisioNize box acts as a hardware hub and works as a translation system enabling data collection of the device data in VisioNize Lab Suite.

Eppendorf Centrifuges	Eppendorf Shakers	Eppendorf ULT Freezers
Centrifuge 5910 R	Innova® 40	CryoCube® F440n
Centrifuge 5424 R	Innova® 42	CryoCube® F440h
Centrifuge 5425	Innova® 43	CryoCube® F570
Centrifuge 5425 R	Innova® 44	CryoCube® F570h
Centrifuge 5427 R	New Brunswick™ S41i	CryoCube® F570hw
		CryoCube® FC660
Eppendorf CO₂ Incubators	Eppendorf Thermomixer	CryoCube® FC660h
Galaxy® 48R	Eppendorf ThermoMixer® C, F0.5, F1.5, F2.0, FP	CryoCube® F740
CellXpert® C170	Eppendorf ThermoStat™ C	CryoCube® F740h
	Eppendorf Cyclers	CryoCube® F101h
	Mastercycler® Nexus	Innova® U101
	Mastercycler® Nexus gradient	Innova® U535
	Mastercycler® Nexus flat	Innova® U570
	Mastercycler® Nexus X2	Innova® U725
	Mastercycler® Nexus GX2	Innova® U725-G
	Mastercycler® Nexus SX1	Innova® C760
	Mastercycler® Nexus GSX1	Innova® C585
		Premium U410
		Premium U570
		HEF U410

Your local distributor: www.eppendorf.com/contact

Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com