

# Decontamination Certificate

The fully completed decontamination declaration is a prerequisite for the acceptance and further processing of your service. If your shipment is not accompanied with this declaration, processing will be delayed until receipt of your decontamination certificate.

**Sender**

|                  |                          |         |
|------------------|--------------------------|---------|
| Company name:    |                          |         |
| Department name: |                          |         |
| Address          |                          |         |
| Zip code:        | City:                    | Country |
| Name:            | Phone:                   | Email:  |
| Client number:   | Contract / Quote number: |         |

Privacy information according to Articles 13 and 14 GDPR: <http://www.eppendorf.com/gdpr>

**Sent goods**

| Item | Instrument / Assembly / Product | Serial number / Batch number | Comments |
|------|---------------------------------|------------------------------|----------|
| 1    |                                 |                              |          |
| 2    |                                 |                              |          |
| 3    |                                 |                              |          |
| 4    |                                 |                              |          |
| 5    |                                 |                              |          |
| 6    |                                 |                              |          |
| 7    |                                 |                              |          |
| 8    |                                 |                              |          |

**Laboratory safety level**

Please specify the safety level of the laboratory where the device was used:

S1  S2  S3  S4

**I hereby confirm that:**

- The device(s) or product(s) mentioned above has/have not been in contact with hazardous substances.
- The device(s) or product(s) mentioned above has/have been in contact with hazardous substances.
- The above mentioned device(s) have not been in contact with the SARS-CoV-2 corona virus
- The above mentioned device(s) has/have been in contact with the SARS-CoV-2 coronavirus\*

\*Extract from "Recommendations of the RKI on hygiene measures in the context of the treatment and care of patients with an infection by SARS-CoV-2 (as of 1.4.2020)": Thermal disinfection methods should be used preferentially whenever possible. If this is not possible, disinfectants with proven, at least limited virucidal effectiveness should be used.

Description of the measures for the decontamination of the listed objects:

- We hereby declare that the objects to be serviced have been carefully cleaned and decontaminated prior to shipment. All objects (instruments and all other objects sent to service) are free of harmful chemicals, infectious or radioactive substances, as well as health-threatening proteins, DNA or bacteriological contamination. All listed devices are safe to handle and ship.

Signature of authorized person

Company / Department

City, Date

---

Important information: If dosing units (or parts there of) are autoclaved, they must be thoroughly dried before shipment. Internally moist dosing units cannot be calibrated accurately, as the residual moisture affects the volume. Please refer to your operating instructions for further information.