

Now available in
Chennai & Bengaluru



Pipette Calibration on Call

Mobile pipette calibration facility at your doorstep!



Actual calibration lab in the vehicle

Pipette Calibration is a very important part of any laboratory's routine. Calibration ensures that pipettes are always working according to specifications and keeping the pipette error-free while doing daily dispensing. To facilitate easier calibration and maintenance of your pipettes, Eppendorf is pleased to announce the availability of **"Mobile Pipette Calibration Facility"** in **Chennai & Bengaluru**. A fully equipped vehicle for this purpose is now available. Contact us to get the vehicle to your facility. The calibration will be conducted in the parking lot of your facility! What's more, the calibration certificate will also be issued instantly. No hassles to send the pipettes and wait for days to get the report!

Scope of pipette calibration

- > Initial inspection to determine pipette condition
- > Cleaning of pipette surface
- > Cleaning and lubrication of seals, O-rings and piston
- > Replacement of worn or defective parts, when necessary
- > 10 measurements per test volume (min, mid, max)
- > Adjustment according to ISO 8655
- > Calibration report provided
- > Calibration label pasted on pipettes

Ordering information

Description	Price (₹)
Fixed-volume pipette	600
Adjustable-volume pipette	750
8-channel pipette	1350
12-channel pipette	1600

Eppendorf India Limited
 Doshi Towers · 4th Floor · 156 Poonamalle High Road
 Kilpauk · Chennai 600 010
 Tel +91 44 42 111 314 / 341
 info@eppendorf.co.in · www.eppendorf.co.in

Eppendorf India Limited
 No.32, 2nd Floor · 1st Main (Temple) Road
 Between 10th & 11th Cross · Malleswaram
 Bengaluru 560 003
 Tel +91 80 41 535 870 / 871