









		epT.I.P.S.®		
Vessel Format	Type of tips	Eppendorf Xplorer® plus, Eppendorf Research® plus		
	Pipette	Eppendorf Xplorer® plus, Eppendorf Research® plus		
	No. of channels	4	6	8
	Volume (µL)	300 / 1,200	300 / 1,200	300 / 1,200
	Tip distance (mm)	9–33	9–20	9–14
	96 Wells (Tip distance 9 mm)	■	■	■
	48 Wells (Tip distance 13 mm)	■	■	■
	24 Wells (Tip distance 19 mm)	■	■	–
	12 Wells (Tip distance 26 mm)	■	–	–
	1.5 mL / 2.0 mL Tube (Tip distance 14, 18, 20 mm)	■	■	■
	Agarose gel	14, 18, 20 mm ■	14, 18, 20 mm –	14 mm –