


Promotion Eppendorf  
Conical Tubes 25 mL



# Space Matters!

Gain up to 33% storage space with Eppendorf Conical Tubes 25 mL

How would you benefit from shorter conical tubes? Nominated in the Scientist Choice Awards as one of the best new lab products of the year, our 25 mL tubes offer sample integrity in a new format that will optimize your workflow.

- > High g-Safe® centrifugation stability up to 17,000 x g
- > Lid tightness from -86 °C up to 100 °C
- > Reduce the risk of contamination and improve your sample recovery with this lower-depth / wide-opening format
- > At the same diameter as standard 50 mL conical tubes, you can easily integrate the 25 mL design into your laboratory periphery

Discover more!  
[www.eppendorf.com/tiptalk](http://www.eppendorf.com/tiptalk)



Learn how we achieve optimal sample and pellet visibility  
[www.eppendorf.com/25ml](http://www.eppendorf.com/25ml)

## Buy 3 and get the 4<sup>th</sup> for free

Order No.	Eppendorf Conical Tubes 25 mL with SnapTec® cap	
0030 118.405	Eppendorf Quality™, colorless, 200 tubes	Promotion Code: MR2020JULY
0030 118.413	PCR clean, colorless, 200 tubes	Promotion Code: MR2020JULY
0030 118.421	Sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 150 tubes	Promotion Code: MR2020JULY
Order No.	Eppendorf Conical Tubes 25 mL with screw cap	
0030 122.410	Eppendorf Quality™, colorless, 200 tubes	Promotion Code: MR2020JULY
0030 122.429	PCR clean, colorless, 200 tubes	Promotion Code: MR2020JULY
0030 122.437	Sterile, pyrogen-, DNase-, RNase-, human and bacterial DNA-free, colorless, 200 tubes	Promotion Code: MR2020JULY

## Save 25%

Order No.	Starter Packs Eppendorf Conical Tubes 25 mL	
0030 118.502	<b>Starter Pack</b> Eppendorf Conical Tubes 25 mL with SnapTec® cap, 4 Single Tube Stands, 6 adapters for rotors with bore for 50 mL conical tubes, PCR clean, colorless, 200 tubes	Promotion Code: 2020JULY25
0030 122.500	<b>Starter Pack</b> Eppendorf Conical Tubes 25 mL with screw cap, 4 Single Tube Stands, 6 adapters for rotors with bore for 50 mL conical tubes, PCR clean, colorless, 200 tubes	Promotion Code: 2020JULY25
Order No.	Eppendorf Storage Box	
0030 140.729	<b>Storage Box 3 × 3</b> , freezer box, for 9 tubes 25 mL, 2 pcs., height 89 mm, 3.5 in, polypropylene, for freezing to -86 °C, autoclavable, with lid and alphanumeric code	Promotion Code: 2020JULY25

## A Complete System

**Centrifugation:** Centrifuges with rotors for 50 mL conical tubes, compatible with two adapter variants for 25 mL tubes SnapTec cap (*order no. 5820 733.004*) and screw cap (*order no. 5820 734.000*)

**Heating and Mixing:** Eppendorf ThermoMixer® C and Eppendorf ThermoStat™ C: Adapter for 25 mL tubes ensures compatibility with Eppendorf SmartBlock™ 50 mL (*order no. 5365 000.109*)

**MixMate®:** Tube Holder 25/50 mL, for 4 × 25/50 mL conical tubes (*order no. 5353 040.156*)

**Single Tube Stand:** Ideally suited for weighing samples on the balance and for handling individual samples (*order no.: 0030 119.860*)

**Storage box made from polypropylene:** Lower height allows space-saving storage of up to 9 × 25 mL tubes in 3.5-inch Eppendorf Storage Boxes (*order no.: 0030 140.729*)

**Tube Rack 25 mL:** Eppendorf Racks for 50 mL conical tubes can be used with the new 25 mL tubes (*order no.: 0030 119.835*)



Click and buy:  
[www.eppendorf.com/consumable-promo-2020](http://www.eppendorf.com/consumable-promo-2020)

Contact us today to request your sample!  
[www.eppendorf.com/contact](http://www.eppendorf.com/contact)

No cash redemption. Errors and mistakes reserved. Promotion is valid in the Eppendorf eShop. Until 31<sup>st</sup> August 2020, you can use this code as many times as you like. Cannot be combined with any other offer(s). Shipping costs may apply. Not valid for Eppendorf's distributor partners.

Eppendorf®, the Eppendorf Brand Design, SnapTec® cap, g-Safe®, ThermoMixer® and MixMate® are registered trademarks of Eppendorf AG, Hamburg, Germany. Eppendorf Quality™ and ThermoStat™ are trademarks of Eppendorf AG, Hamburg, Germany. Patents are listed on [www.eppendorf.com/ip](http://www.eppendorf.com/ip). Order no. CONPRoJULYEUR2020v1. All rights reserved, including graphics and images. Copyright © 2020 by Eppendorf AG.