

epServices
for premium performance



Qualification Services

Be prepared for your audits!



IQ/OQ Services – GxP compliant



Supporting you in Qualification matters

Installation Qualification and Operational Qualification services assure that instrumentation is installed and running at premium performance.

Test series for process critical measured values are performed to ensure manufacturer specifications. IQ and OQ GxP compliant documentation including the relevant data is supplied to help you prepare and be ready for your quality and regulatory audits.

www.eppendorf.com/epServices

Operational Qualification according to instrument type

Installation Qualification & Operational Qualification services	IQ/OQ GxP	OQ GxP
Check of delivery status, and installation environment	■	–
Installation incl. configuration and functional checks	■	–
Detailed IQ report	■	–
Test of alarm functions	■	■
Verification of instrument parameters (e.g. temperature, speed)	■	■
Detailed OQ report and handling of deviation list	■	■
Signed approval by customer and service technician	■	■

Reduce your administrative workload!

Contact your local Eppendorf Service Representative for qualification services supporting the complex GxP and regulations via phone, e-mail, or our website contact form.

Further information: **www.eppendorf.com/epServices**

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · 22331 Hamburg · Germany

E-mail: eppendorf@eppendorf.com

www.eppendorf.com