




VisionNize[®] Lab Suite

Overview of compatible Eppendorf devices

How Devices Can Be Connected



- 1 VisioNize box (hardware hub)
- 2 VisioNize sense gateway router
- 3 VisioNize sense environmental monitoring sensors
- 4 VisioNize sense external temperature sensor



VisioNize

Device Overview

Overview of compatible Eppendorf devices for VisioNize Lab Suite

This list contains Eppendorf devices which can be connected to the platform. Visit www.eppendorf.com/visionize for the most current compatibility list.

VisioNize touch enabled devices

> Eppendorf devices equipped with the VisioNize touch interface can be connected via ethernet cable:

Eppendorf ULT Freezers	Eppendorf CO₂ Incubators	Eppendorf Cyclers
CryoCube® F740i	CellXpert® C170i	Mastercycler® X50a
CryoCube® F740hi		Mastercycler® X50s
CryoCube® F740iw	Eppendorf Shakers	Mastercycler® X50p
CryoCube® F740hiw	Innova® S44i*	Mastercycler® X50h
	<small>* only compatible after serial number S44IJH302546</small>	Mastercycler® X50i
		Mastercycler® X50l
Bioprocess Systems	Eppendorf Centrifuge	Mastercycler® X50r
SciVario® twin	Centrifuge 5910 Ri	Mastercycler® X50t

VisioNize box enabled

> Eppendorf devices which are connectable via VisioNize box to the VisioNize Lab Suite. The VisioNize box acts as a hardware hub and works as a translation system enabling data collection of the device data in VisioNize Lab Suite.

Eppendorf ULT Freezers	Eppendorf Shakers	Eppendorf Centrifuges
CryoCube® F440n	Innova® 40	Centrifuge 5910 R
CryoCube® F440h	Innova® 42	Centrifuge 5424 R
CryoCube® F570	Innova® 43	Centrifuge 5425
CryoCube® F570h	Innova® 44	Centrifuge 5425 R
CryoCube® F570hw	New Brunswick™ S41i	Centrifuge 5427 R
CryoCube® FC660		
CryoCube® FC660h	Eppendorf Thermomixer	Eppendorf CO₂ Incubators
CryoCube® F101h	Eppendorf ThermoMixer® C, F0.5, F1.5, F2.0, FP	Galaxy® 48R
Innova® U101	Eppendorf ThermoStat™ C	CellXpert® C170
Innova® U535		
Innova® U570	Eppendorf Cyclers	
Innova® U725	Mastercycler® Nexus	
Innova® U725-G	Mastercycler® Nexus gradient	
Innova® C760	Mastercycler® Nexus flat	
Innova® C585	Mastercycler® Nexus X2	
Premium U410	Mastercycler® Nexus GX2	
Premium U570	Mastercycler® Nexus SX1	
HEF U410	Mastercycler® Nexus GSX1	

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com