

**Instructions for use**  
Copper Chamber  
English (EN)

This document is a supplement to the operating manual for CellXpert® C170i and not a substitute. Therefore, also read the operating manual. The current version can be found on the Internet at [www.eppendorf.com/manuals](http://www.eppendorf.com/manuals).

## 1 Copper

This chamber has a copper surface.

Changes in color do not reflect degradation. The tempered color is the result of heat during the chamber's high-temperature disinfection cycles.

Copper effectively minimizes contamination. Over time, the surface develops a patina due to oxidation of the copper. This process is normal and increases the decontamination effect.

Do not polish the copper surface.