



Cell type: MCF-7

Performance of Eppendorf Cell Culture Consumables (TC treated)

Procedure: MCF-7 cells were seeded in appropriate culture medium (RPMI 1640, 10 % FBS). Cell Culture vessels were incubated at 37 °C with a humidified 5 % CO₂ atmosphere. Cell growth and morphology was monitored on day 1, 2 and 3 post-seeding.



Bars represent average count of two wells of an Eppendorf 6-Well Cell Culture Plate



Representative light microscope images (100X magnification)

The tissue culture (TC) treated surface of Eppendorf Cell Culture Consumables allows optimal attachment and homogeneous growth of **MCF-7** cells.

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · 22331 Hamburg · Germany
 eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com