



# Accelerate

Speed up your work with the Easypet® 3 electronic pipette controller and Eppendorf Serological Pipets



## Experience a New Dimension of Control

### **Experience a new dimension of electronic pipetting with the Eppendorf Easypet® 3.**

The Easypet 3 was designed with the latest technological developments and now, pipetting speed is controlled simply with the tips of your fingers allowing you to experience more sensitive control of the meniscus. The powerful motor delivers efficient performance to speed up large volume pipetting while also providing the utmost control for smaller volumes. The Easypet 3 has an excellent speed adjustment that allows for an optimized gravity drain option that can be used for easy dispensing of liquid without splashing.

Weighing only 160 grams, the lightweight, well-balanced Easypet 3 has been designed considering all ergonomic findings to sit comfortably in your hand, matching its shape to allow for fatigue-free pipetting. Easypet 3 can be used with all types of serological and volumetric pipets from volumes of 0.1–100 mL.

### **Additional product features:**

- > Smooth setting of pump speed for greater meniscus control
- > Lithium-polymer battery offers long runtime without the need for recharging
- > Battery status display helps ensure the Easypet 3 will not lose power mid-run
- > Operation during recharging is possible
- > Autoclavable pipet adapter for sterile applications
- > Quick release aspirating cone for easy exchange of membrane filters
- > Pressure equalization system releases aspirated vapors before entering the Easypet 3 to help eliminate corrosion



#### Speed control and precision

Intuitive and convenient speed adjustment simply with the tips of your fingers. The new valve technology enables more sensitive control of the meniscus resulting in increased precision.



#### Ergonomics

Lightweight, well-balanced and ergonomically improved design allows for stress-free work. Following our PhysioCare Concept®, we strive to improve your time at work.



#### Battery status display

Vibrant backlit LEDs provide optical feedback of the remaining battery life.



## Eppendorf Serological Pipets

Experience the quality and convenience in liquid handling you have trusted in for decades:

### Easypet® 3 and Eppendorf Serological Pipets.

Serological pipets from Eppendorf are designed to work in perfect harmony with your Easypet 3. The pipets offer precise graduations overlaid on the ultraclear material to facilitate quick and easy reading of the volume. The pipets are individually wrapped and packaged in a robust dispenser box to ensure consistent premium quality.

- > Clear and precise graduations for easy volume determination
- > Color-coding for easy identification of pipet sizes
- > Comprehensive packaging to ensure convenient storage and reliable usage
- > Ultraclear virgin polystyrene meets the requirement of USP Class VI
- > Sterility assurance level of  $10^{-6}$
- > Certified absence of detectable pyrogens, DNA, RNase & DNase, non-cytotoxic





»Easypet® 3—it has never been easier to combine speed, safety, precision and comfort.«



»When developing the Easypet 3 we always had our users in mind. By combining the speed adjustment control with our new valve technology which enables more sensitive control of the meniscus, we want to improve their daily work.«

**Tobias David**, Project Manager Product Development  
Eppendorf Liquid Handling

## Technical Data

Rechargeable battery type	Lithium-polymer
Capacity	1,100 mAh/3.7 V
Charging time	~3 h
Number of dispensings	~2,000 (with a 25 mL pipet)
Weight without rechargeable battery and pipet	134 g (0.295 lb)
<b>Power supply</b>	
Input	100 V–240 V AC, ±10 %, 50 Hz–60 Hz, 0.5 A
Output	5 V DC, 1.0 A

## Ordering Information

Description	Order No.
Easypet® 3 incl. power supply and Lithium-polymer rechargeable battery, wall mount, shelf stand (not available in USA), 2 membrane filters (unsterile) 0.45 µm	4430 000 018
<b>Accessories</b>	
Membrane filter, sterile, 0.45 µm, pack of 5	4421 601 009
Membrane filter, sterile, 0.2 µm, pack of 5	4430 606 005
Lithium-polymer rechargeable battery for Eppendorf Easypet® 3	4430 605 009
Wall mount	4430 604 002
Pipet Helper® mechanical pipette aid	4423 000 010
<b>Eppendorf Serological Pipets</b>	
<b>sterile, free of detect. pyrogens, DNA, RNase &amp; DNase, non-cytotoxic</b>	
■ 1 mL, yellow, 800 pcs. (4 bags × 200 pcs.)	0030 127 692
■ 2 mL, green, 600 pcs. (4 bags × 150 pcs.)	0030 127 706
■ 5 mL, blue, 400 pcs. (4 bags × 100 pcs.)	0030 127 714
■ 10 mL, orange, 400 pcs. (4 bags × 100 pcs.)	0030 127 722
■ 25 mL, red, 200 pcs. (4 bags × 50 pcs.)	0030 127 730
■ 50 mL, violet, 160 pcs. (4 bags × 40 pcs.)	0030 127 749

### Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf SE • 22339 Hamburg • Germany  
 eppendorf@eppendorf.com • www.eppendorf.com  
 Eppendorf North America, Inc. • Phone: 800-645-3050  
 info@eppendorf.com • www.eppendorfna.com  
 Eppendorf Canada Ltd. • Phone: 800-263-8715  
 canada@eppendorf.com • www.eppendorf.ca

[www.eppendorf.com/easypet](http://www.eppendorf.com/easypet)