

High Spec Centrifuge Now With Promotional Price



CE*

* The Centrifuge 5427 R is designed for separating liquid substance mixtures with different densities, in particular, for processing and analyzing samples from the human body in in-vitro diagnostic applications to ensure that the in-vitro diagnostic device can be used according to its intended purpose. This centrifuge including components is an in-vitro diagnostic device according to Directive 98/79/EC of the European Parliament and the Council dated October 27, 1998.

Double Capacity and Higher Speed

Centrifuge 5427 R

- > Speed up to 25,000 × g / 16,220rpm
- > Max. capacity: 48 x 1.5/2.0 mL
- > Small footprint to save precious space on the lab bench
- > Standard aerosol-tight rotor certified by Public Health Agency, Porton Down, UK for maximum user and sample safety
- > Autoclavable rotors
- > Eppendorf QuickLock® rotor lids for fast and ergonomic lid locking
- > Temperature control range from -10 °C to 40 °C
- > FastTemp® function optimized for fast 11 min pre-cooling
- > ECO shut-off can be programmed to engage after 1, 2, 4, or 8 hours of non-use for energy savings up to 88% and extended compressor life
- > Choice of 9 rotors for all your molecular application needs
- > Made in Germany

Save more by choosing what you really need:



Standard-throughput 24 x 1.5/2 mL rotor comes with aerosol-tight QuickLock® rotor lid. ¼ turn to close or open rotor lid to reduce stress off your wrist.
> Max. speed: 25,000 × g (16,220 rpm)



Outer rim to protect tube caps of spin column from shearing off

High-throughput 48 x 1.5/2 mL rotor
> Max. speed: 18,210 × g (12,700 rpm)
Dual functions:
> Safe centrifugation of up to 24 spin columns or,
> High capacity of up to 48 x 1.5/2 mL tubes



Medium-throughput 30 x 1.5/2 mL tubes
> Max. speed: 20,800 × g (14,000 rpm)
> Aerosol-tight QuickLock® rotor lid

Ordering Information

Package	Description	Speed	Capacity
1	Centrifuge 5427 R (EU-IVD), refrigerated, with Rotor FA-45-48-11	18,210 × g	48 x 1.5/2.0 mL
2	Centrifuge 5427 R (EU-IVD), refrigerated, with Rotor FA-45-30-11	20,800 × g	30 x 1.5/2.0 mL
3	Centrifuge 5427 R (EU-IVD), refrigerated with Rotor FA-45-24-11, aerosol-tight Eppendorf QuickLock® lid	25,001 × g	24 x 1.5/2 mL

For ordering information, please check our local website or contact your local sales representative

www.eppendorf.com/centrifugation

Click or scan for more Eppendorf promotions:

