

# WANTED

DEAD OR ELIMINATED

100  $\mu\text{m}$

75  $\mu\text{m}$

50  $\mu\text{m}$

25  $\mu\text{m}$



For the crimes of false and misleading results

## MYCOPLASMA

Description – About 0.1–0.3  $\mu\text{m}$  in diameter,  
invisible via brightfield microscopy.

Sometimes found attached  
to the cell membrane of eukaryotic cells

Testing by DNA staining, PCR, or biochemical  
detection necessary to catch this criminal