



# A New Spin on Clinical Trials

Your partnership for global support, sourcing and service - beyond centrifugation





# Your Strong Partner Beyond Centrifugation

Are you an Ancillary Supply Manager, Clinical Research Manager or Clinical Study Manager looking for a competent partner who will organize everything concerning the centrifuge for your clinical trials? A central point of contact, for the entire organization around the centrifuges such as global sourcing and global service?

The Eppendorf solution for your clinical trial are the highly versatile Centrifuges 5702 (air-cooled) and 5702 R (refrigerated) for the first spin of blood and urine samples on site. They accommodate all common blood, pediatric, urine and cell preparation tubes e.g. BD CPT 8ml Vacutainer (16x125 mm).

Choose Eppendorf as your partner for centrifuges in clinical trials. We are a global company with experience in equipping clinical trials sites with the right centrifuges for sample processing.





# Central Contact



## Simply contact your expert

A central point of contact for the worldwide coordination of your centrifuge solutions for your clinical trials. With excellent sales and service support for all price negotiations, global quotes including transport and services, e.g., repairs and IQ/OQ GxP - regardless of where study sites are located.

We assist you with your centrifuge product selection to assure your compliance with sponsors’ requirements. An expert is available to provide advice and assistance at every step of the process.

Contact us: [clinicaltrials@eppendorf.com](mailto:clinicaltrials@eppendorf.com)


# International Sourcing

## We simplify procurement and logistics

As a company that is active globally, Eppendorf has direct sales and service organizations in many countries and warehouses at strategically important locations around the globe. Benefit from this in your centrifuge sourcing process and save valuable time by simplifying logistics and customs clearance through regional procurement and drop-shipments from Eppendorf

warehouses to study sites. Thanks to our stock levels and our experienced fulfillment team, we are able to deliver the Centrifuge 5702 family quickly. As we know that there is often a need for action at short notice. Eppendorf has a global presence through its strong network of vendors, with thousands of sales partners. Contact your expert.

Service  
300 highly experienced  
service technicians



Global product availability  
•22 warehouses and  
distribution centers  
•29 subsidiaries  
•32 additional sales offices  
•Strong dealer network of a  
thousand dealers worldwide

# Global Service

## We reduce your administrative workload

The worldwide coordination of services to meet GxP compliance is time consuming. Are you concerned about downtimes while your clinical trial is ongoing? Your single point of contact takes over the organization and coordination of all service issues by managing your IQ/OQ, maintenance and warranty - independent of the location of the clinical trials site.

For customers working in regulated environments, we offer support in qualification matters to be prepared for audits. Through our tailored installation qualification and operational qualification services, as well as through

our GxP compliant documentation, we ensure the proper installation and optimum performance of your centrifuge. Test series of values critical to the process are performed to verify manufacturer specifications.

IQ/OQ GxP documentation will give you peace of mind due to reduce workload for your centrifuge qualification and in case of an audit your centrifuge documentation meets specifications.

We offer an extended warranty on parts, labor and travel expenses on repairs.

## IQ/OQ Services – GxP

Operational Qualification according to instrument type		
Installation Qualification & Operational Qualification services	IQ/OQ GxP	OQ GxP
Check of delivery status and installation environment	■	-
Installation incl. configuration and functional checks	■	-
Detailed IQ report	■	-
Test of alarm functions	■	■
Verification of instrument parameters (e.g. temperature, speed)	■	■
Detailed OQ report and handling of deviation list	■	■
Signed approval by customer and service technician	■	■



# Perfect Solution

> Low noise level  
5702 R (52 dB(A))

> New bucket for swing out rotor A-4-38 for extra-tall blood tubes (16x125 mm)

> Small footprint (1,264 cm²), 5702 R 38x58x26 cm (WxDxH)

> Various buckets for all common blood, urine, pediatric, and cell preparation tubes

> Automatic start after programmable timer has elapsed, e.g., for a standardized clotting time in serum isolation

## Versatile centrifuge for almost any clinical trial

Are you looking for a centrifuge which meets the requirements of sponsors for most applications? A centrifuge that is suitable for nearly all clinical trials, as well as for in-vitro diagnostics?

Our versatile, high quality, low downtime centrifuge meets most sponsor requirements and can therefore be used in almost any clinical trials - whether you offer your customer a purchase or rental model.

The 5702 family is the perfect solution for clinical trial sample collection sites for the initial spin of blood and urine samples. These versatile centrifuges accommodate all common blood, pediatric and urine tubes. The swing-out rotor offers significant advantages in terms of reduced sample contamination and higher


sample yields compared to a fixed-angle rotor. Importantly, this centrifuge accommodates the extra-tall (16x125 mm) cell preparation tubes e.g., the BD CPT 8 mL Vacutainer for PBMC isolation.

With a very compact footprint, this centrifuge fits on every lab bench, and its quiet operation benefits your work environment. User-friendly operation minimizes human error and allows for easy training of site staff, regardless of their level of experience with laboratory equipment.

Its capacity of up to 16 standard tubes (16x100 mm or 13x75 mm tubes) makes this centrifuge suitable for doctor's offices, but it is also ideal for processing samples from multiple patients in one run in larger sites such as hospitals, independent clinical sites or phase 1 units.

\*Intended use  
The Centrifuge 5702/5702 R/5702 RH is a non-automatic centrifuge for separating liquid substance mixtures from the human body and is specifically intended for use as an accessory with an in-vitro diagnostic device in order to facilitate the in-vitro diagnostic device to be used in accordance with its intended use. Eppendorf centrifuges are only intended for indoor use and for operation by trained and skilled personnel.

# Your Centrifugation Solution for Clinical Trials



Eppendorf offers you a reliable partnership to minimize your administrative effort around your centrifuge solutions for clinical studies. We build on a trusting

cooperation and offer you comprehensive support. As a global company, we can deliver on time to keep your clinical trials running.

Request your individual support:

> Contact us: [clinicaltrials@eppendorf.com](mailto:clinicaltrials@eppendorf.com)

**Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)**  
Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)

Eppendorf® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE. All rights reserved, including graphics and images.  
Order no. AA01 025 080/ENA/0,2T/0922/STEF · Carbon neutrally printed in Germany. Copyright © 2022 by Eppendorf SE.