

Centrifuge 5702 / 5702 R / 5702 RH

*CE



General Benefits

- > **Small footprint:** Medium capacity (40 x blood tubes) at a very compact footprint
- > **High versatility:** Spin a large variety of different blood and urine tubes including pediatric tubes (Microtainer® or MiniCollect®)
- > **Optimal conditions:** Optimum pellet formation due to swing-out rotor for contamination-free removal of the supernatant
- > **Strong sample safety:** Timer function and program buttons (Centrifuge 5702 R) for reproducible results

| Centrifuge 5702 | | Rotor A-4-38 | Round bucket | | | | | Rectangular bucket | | | 16 x 125 mm bucket | | |
|----------------------------------|-------------------------|-----------------------------|--------------|--------------|--------------|--------------|--------------|--------------------|--------------|--------------|--------------------|--------------|---|
| | | | | | | | | | | | | | |
| | | | 5702 737 003 | 5702 719 005 | 5702 735 000 | 5702 724 009 | 5702 723 002 | 5702 763 004 | 5702 710 008 | 5702 711 004 | 5702 712 000 | 5702 764 000 | |
| Blood sample | BD Vacutainer® | Microtainer® | 16* | 16* | - | - | - | 36* | 40* | - | - | - | |
| | | 1.8 – 4.8 mL (13 x 75 mm) | 16 | 16 | - | - | - | 36 | - | - | - | - | |
| | | 4.0 – 7.0 mL (13 x 100 mm) | 16 | 16 | - | - | - | - | 40 | 32 | - | - | |
| | | 8.0 – 11.0 mL (16 x 100 mm) | - | - | 16 | 16 | - | - | - | - | - | 24 | - |
| | | CPT™ 4.0 mL (13 x 100 mm) | 16 | 16 | - | - | - | - | 40 | - | - | - | - |
| | | CPT™ 8.0 mL (16 x 125 mm) | - | - | - | - | - | - | - | - | - | - | 4 |
| | Greiner® Vacuette | MiniCollect® | 16 | 16 | - | - | - | 36* | 40* | - | - | - | - |
| | | 0.9 – 4.0 mL (13 x 75 mm) | 16 | 16 | - | - | - | 36 | 40 | - | - | - | - |
| | | 3.5 – 6.0 mL (13 x 100 mm) | 16 | 16 | - | - | - | - | 40 | - | - | - | - |
| | | 6.0 – 10.0 mL (16 x 100 mm) | - | - | 16 | 16 | - | - | - | - | - | 24 | - |
| | Sarstedt Monovette® | 1.6 – 4.0 mL (ø 11/13 mm) | 16 | - | - | - | 8 | 36 | - | - | - | - | - |
| | | 4.0 – 5.0 mL (ø 11 mm) | 16 | - | - | - | - | 36 | - | - | - | - | - |
| | | 4.0 – 10.0 mL (ø 13 mm) | - | 16 | - | - | - | - | 20 | 32 | - | - | - |
| | | 4.0 – 10.0 mL (ø 15 mm) | - | - | 16 | 16 | - | - | - | - | - | 24 | - |
| Streck Blood Collection Tube | 5.0 mL (13 x 75 mm) | 16 | 16 | - | - | - | 36 | 40 | - | - | - | - | |
| | 10.0 mL (16 x 100 mm) | - | - | 16 | 16 | 8 | - | - | - | - | 24 | - | |

* Manufacturer adapter needed

Centrifuge 5702 / 5702 R / 5702 RH

| Centrifuge 5702 | | Rotor A-4-38 | Round bucket | | | | | Rectangular bucket | | | 16 x 125 mm bucket | |
|-----------------|-------------------------|---------------------------------------|--------------|--------------|--------------|--------------|--------------|--------------------|--------------|--------------|--------------------|--------------|
| | | | | | | | | | | | | |
| | | | 5702 737 003 | 5702 719 005 | 5702 735 000 | 5702 724 009 | 5702 723 002 | 5702 763 004 | 5702 710 008 | 5702 711 004 | 5702 712 000 | 5702 764 000 |
| Urine sample | BD Vacutainer® | 4.0 – 6.0 mL (ø 13 mm) | 20 | 16 | – | – | – | 36 | – | – | – | – |
| | | 8.0 – 11.0 mL (ø 16 mm) | – | – | 16 | 16 | – | – | – | – | 24 | – |
| | | 8.0 – 9.5 mL (conical, ø 16 mm) | – | – | – | – | 8 | – | – | – | 24 | – |
| | Greiner Vacuette® | 4.0 mL (ø 13 mm) | 16 | 16 | – | – | – | – | 40 | – | – | – |
| | | 6.0 – 10.0 mL (ø 16 mm) | – | – | 16 | 16 | – | – | – | – | 24 | – |
| | Sarstedt Monovette® | 3.2 – 6.0 mL (ø 13 mm) | 16 | 16 | – | – | – | – | 40 | – | – | – |
| | | 8.5 – 10.0 mL (ø 15 mm) | – | – | 16 | 16 | – | – | – | – | 32 | 24 |
| | Test Tube | 15.0 mL | – | – | – | – | 8 | – | – | – | 24 | – |

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com