

Enjoy Peace of Mind - Above and Beyond Sample Safety

Get connected to VisioNize® Lab Suite

Whether you are looking for a remote monitoring solution, a lab management system or ways to easily access device performance data, VisioNize Lab Suite has all the services to help you build a smart laboratory with enhanced sample safety, productivity and reliability.



For starters:

The Essential plan provides you with the basic functions to begin the digital management of your laboratory, starting with increasing your sample safety.

For more peace of mind:

Create device documentation at the press of a button whilst benefiting from enhanced notification options. Enjoy the enriched feature set in VisioNize Lab Suite for more convenience!

	For Free (3 connected devices included*)	For more peace of mind (Billed annually per connected device**)
	VisioNize Lab Suite Essentials	VisioNize Lab Sample Safety+
Real-Time Monitoring Monitor connected lab devices and track status, parameters, and performance history from anywhere at any time.	✓	✓
Instant Notifications Receive alarm, error or event notifications to address any issues quickly.	Email only	Email & SMS messaging
Non-Connectable Devices Manage an unlimited number of non-connectable devices by manual data entry, e.g. serial number and location.	✓	✓
Software Updates Receive remote device software updates for connected Eppendorf lab devices.	✓	✓
Task Management Plan and assign device-related tasks and get notified when due.	✓	✓
Document Management Quickly retrieve relevant device documentation, e.g. operating manuals and other SOPs.	✓	✓
Device Data Export Export device parameter, event data and any user interactions.	✓	✓
Lab Equipment Booking Book and plan the usage of your devices in one shared calendar to boost productivity and smooth execution of your experiments, e.g. PCR run.	✓	✓
Device Data Access Quickly retrieve all relevant device data and go back in time with the history view.	7 days	90 days
Audit Trail Track user interactions or alarm acknowledgements with time-stamped data in an audit trail.	—	✓
Response-Based Alarm Escalation Define escalation schemes to individual or multiple recipient groups based on presence in the laboratory.	—	✓

* Additional devices require device licenses which will be billed annually. Hardware and software upgrades are not included.

** All pricing information are available online at: www.eppendorf.com/visionize-subscription



Contact us for more information: www.eppendorf.com/visionize-subscription

